


Regardless of what
the day brings, our
new urostomy pouch
helps maximize
discretion and comfort.



Now Available in One Piece!
Premier
Urostomy Pouches

A Great Pouch Becomes Even Better

Introducing the improved Premier urostomy pouch. The updated pouch features several product innovations designed to deliver more comfort and greater control for your patients.

Enhanced Flow-Control Tap


The new tap is made with softer material and a streamlined shape. The goal of this design? To help increase comfort against the body, so your patients can do more, with less worry.

Barrier Options for Best Fit


The Premier urostomy pouch offers the options you need to promote optimal fit and skin health. Barriers feature the Flextend extended wear formulation, and you can choose flat or convex based on your patient's needs and body type.

New Multi-Chamber Design

The updated pouches now have multiple built-in chambers, which are designed to facilitate even distribution of urine. As urine flows into the pouch, it is dispersed among the chambers. This results in a more balanced collection of urine, and a lower pouch profile for your patients.



Multi Chamber,
Front View



Multi Chamber,
Top View, Sliced

Improved Adapter

The adapter's new side tab design helps facilitate ease of use when connecting any Premier one-piece urostomy pouch to a bedside drainage collection system.



New Enhancements plus the Features You Already Know and Trust

In addition to these pouch improvements, patients can continue to rely on the other benefits of Premier urostomy pouches including:

- Odour-barrier pouch film
- Transparent, ultra-clear, and beige options
- Soft, ComfortWear panels
- Belt tabs

The Goal is a Greater Sense of Control

No matter what the day presents, your patients deserve an uneventful experience with their pouch. These Hollister urostomy pouch enhancements are designed to deliver more opportunities to make that increased sense of control possible.



You asked, we responded!
Premier urostomy pouches with convex barrier are now offered in ultra-clear film to facilitate inspection of the stoma and pouch contents.



Choose your option.
Both flat and convex barrier options are available.

Also available with upgrades -
New Image Two Piece Urostomy Pouch

For more details, contact your Hollister representative.

Urostomy Pouches - Ordering Information

NEW! Premier One-Piece Urostomy Pouch

- Flexend barrier is durable and offers high resistance to erosion
- Tape border and belt tabs for additional security
- Soft, enhanced flow-control tap for better control and comfort
- Multi-chamber design for a slimmer pouch profile

Flat Barrier

Pouch Colour	Barrier Opening	Product No.
Ultra-Clear	Up to 2-1/2" (64 mm) ✂	84590
	Up to 2-1/2" (64 mm) ✂	84690
	3/4" (19 mm)	84697
Transparent	1" (25 mm)	84692
	1-1/4" (32 mm)	84698
	1-1/2" (38 mm)	84693
	1-3/4" (44 mm)	84699
	2" (51 mm)	84694
Beige	Up to 2-1/2" (64 mm) ✂	84490
	3/4" (19 mm)	84497
	1" (25 mm)	84492
	1-1/4" (32 mm)	84498
	1-1/2" (38 mm)	84493
	1-3/4" (44 mm)	84499
	2" (51 mm)	84494

Convex Barrier

Pouch Colour	Barrier Opening	Product No.
Ultra-Clear	Up to 1" (25 mm) ✂	84794
	Up to 1-1/2" (38 mm) ✂	84798
	Up to 2" (51 mm) ✂	8479111
	1/2" (13 mm)	84890
	5/8" (16 mm)	84891
	3/4" (19 mm)	84892
	7/8" (22 mm)	84893
	1" (25 mm)	84894
	1-1/8" (29 mm)	84895
	1-1/4" (32 mm)	84896
	1-3/8" (35 mm)	84897
	1-1/2" (38 mm)	84898
	1-3/4" (44 mm)	84899
Beige	Up to 1" (25 mm) ✂	84394
	Up to 1-1/2" (38 mm) ✂	84398
	Up to 2-1/8" (55 mm) ✂	8439111
	1/2" (13 mm)	84990
	5/8" (16 mm)	84991
	3/4" (19 mm)	84992
	7/8" (22 mm)	84993
	1" (25 mm)	84994
	1-1/8" (29 mm)	84995
	1-1/4" (32 mm)	84996
	1-3/8" (35 mm)	84997
	1-1/2" (38 mm)	84998
	1-3/4" (44 mm)	84999

New Image Two-Piece Urostomy Pouch

- Belt tabs for additional security
- Soft, enhanced flow-control tap for better control and comfort
- Multi-chamber design for a slimmer pouch profile

Pouch Colour	Flange Size	Pouch Length		Product No.
Ultra-Clear	1-3/4" (44 mm)	9" (23 cm)		18922
	2-1/4" (57 mm)	9" (23 cm)		18923
	2-3/4" (70 mm)	9" (23 cm)		18924
Transparent	1-3/4" (44 mm)	9" (23 cm)		18902
	2-1/4" (57 mm)	9" (23 cm)		18903
	2-3/4" (70 mm)	9" (23 cm)		18904
Beige	1-3/4" (44 mm)	9" (23 cm)		18912
	2-1/4" (57 mm)	9" (23 cm)		18913
	2-3/4" (70 mm)	9" (23 cm)		18914

✂ Indicates Cut-To-Fit product



Hollister Limited
95 Mary Street
Aurora, ON
L4G 1G3
1-800-263-7400

www.hollister.com

The Hollister logo, ComfortWear, Premier and New Image are trademarks of Hollister Incorporated.

© 2016 Hollister Limited.

912084-CA