

Making a Difference in the Journey of Life

2022 Product Catalogue



The fundamental values and principles
of our founder, John Dickinson Schneider,
form the foundation of our company today.

Our Mission

To make life more rewarding and dignified for the people who use our products and services.

Our Vision

To grow and prosper as an independent, employee-owned company and in the process become better human beings.

Our Heritage

The roots of Hollister Incorporated date back to 1921, when 23-year-old entrepreneur John Dickinson Schneider founded a small printing company in Chicago, which he named JDS Printer Craftsman. From the very beginning, John Schneider's business philosophy was "only first class is good enough," and the company became known for its quality products, quality service, and quality people.

Table of Contents

General Information	2	Accessories	
Ostomy Care		Adapt™ Barrier Rings	27
Choosing the Right Skin Barrier FIT	4	Adapt™	28
Choosing the Right Skin Barrier FORMULATION	5	Karaya	28
New Image™		m9 Odour Eliminators	29
Two-Piece Pouching System		Pouch Accessories	29
Flat Skin Barriers	8	Access Port	29
Soft Convex Skin Barriers	10	Wound Drainage Collectors	29
Convex Skin Barriers	10	With FlexWear Skin Barrier	29
Drainable Pouches	12	Skin Sheets	29
Closed Pouches	13	Irrigation Products	29
High Output Pouches	13	Stoma Caps	29
Urostomy Pouches	14	Continence Care	
Ostomy Kits	15	Infyna Chic™	32
Urostomy Kits	15	VaPro™ and VaPro™ F-Style	34
Premier™		Touch Free Hydrophilic Intermittent Catheters	
One-Piece Pouching System		VaPro™ Pocket and VaPro™ Plus Pocket	35
Flat Skin Barriers - Drainable Pouch	16	Touch Free Hydrophilic Intermittent Catheters	
High Output Pouches	17	Advance Plus™	
Flat Skin Barriers - Closed Pouch	17	Touch Free Intermittent Catheters	37
Soft Convex Skin Barriers - Drainable Pouch	18	Apogee™	
Convex Skin Barriers - Drainable Pouch	18	Intermittent Catheters	38
Flat Skin Barriers - Urostomy Pouch	21	Male External Catheters	
Soft Convex Skin Barriers - Urostomy Pouch	22	Hollister Extended and Everyday Wear	39
Convex Skin Barriers - Urostomy Pouch	22	InView Silicone	39
Karaya™		Urinary Collection Systems	
One-Piece Pouching System—		Leg Bags	40
Natural Skin Barriers		Bedside Drainage	40
Drainable Pouches	24	m9 Cleaner/Decrystallizer	40
Closed Pouches	24	External Collection Devices	
Pouchkins		Urinary Pouches	41
Pouchkins™		Fecal Collectors	41
Pediatric Pouching System		Critical Care	
One-Piece	26	Oral Endotracheal Tube Fasteners	44
Two-Piece	26	Tube Attachment Devices	44
		Numerical Index	45

General Information

Terms of Payment

Net 30 days from date of invoice. Electronic Fund Transfers (ETF) are possible. Please contact the customer order department.

Minimum Order

The minimum order requirement for a dealer is \$250.00. A charge of \$10.00 will be applied to any order below the minimum requirement.

Shipping

Prices for shipments to destinations in Canada are FOB buyer's address.

If the method and route of shipment are specified by the buyer, the buyer is required to pay the delivery charges, upon delivery of products to the carrier, all risk of loss and other incidents of ownership pass to the buyer. These terms may not be altered.

Pricing & Terms of Sale

GST/HST will be charged where applicable. Buyer agrees that they will not at anytime, directly or indirectly, sell or distribute, to customers which are located outside of Canada.

Return Goods Policy

All returns for credit require an authorization number from Hollister Limited. A Returns Goods Authorization (RGA) can be obtained by contacting a representative at Customer Service 1.800.263.7400 or by email canadacustomerservice@hollister.com. Written request for authorization is required.

Hollister will authorize returns under the following conditions:

1. A 15% restock fee will be applied to all returned goods over 30 days from the date of invoice.
2. Copy of invoice or purchase order number must be provided for the return request.
3. **Product returned by customers without an authorized RGA number will be disposed of by Hollister Limited and no credit will be issued.**

Hollister will NOT authorize the return of the following:

1. Goods over one (1) year from the date of invoice.
2. Partial boxes of merchandise.
3. Boxes in unsaleable condition (i.e markings, sticker for sale, sun fading, dented boxes due to improper handling/return packaging, sticker glue, etc.)

Damaged Product:

1. Customer must call Hollister Limited to receive an RGA number within 5 business days of receipt.
2. Receiving dock is responsible for documenting, on the bill of lading, when a shipment is received in damaged condition.
3. Customer must provide copy of invoice, invoice number or purchase order number for the return request.

4. Credit will **NOT** be issued if customer disposes of damaged/ defective product without approval by Hollister Limited.
5. If product cannot be returned due to extreme damage (leakage/broken) customer must provide pictures to Hollister Limited within 5 business days of receipt to gain approval to dispose and thereby obtain credit for the same.

Ostomy & Continence Retail Sample Policy

All sample requests at the Retail level are to be processed through Hollister Customer Service by calling 1.800.263.7400 or by email canadacustomerservice@hollister.com.

Retailer/client is to contact Customer Service and provide consumer information for sample processing. Samples will be shipped directly to consumers from Hollister Limited. Samples will **no longer** be available to ship directly to a retail setting.

Shortages

Customer must notify Hollister Limited within 5 business days of receipt on shortages. A shortage is a difference between what was physically received and the master pack list.

Customer Satisfaction Return Policy

All Customer Satisfaction Returns require authorization by your Hollister Territory Manager or a Customer Service Representative.

Unauthorized Customer Satisfaction Returns will not be processed.

Hollister will authorize Customer Satisfaction Returns under the following conditions:

1. The customer has been given a no-charge, full box replacement of Hollister product for their partial box return, and
2. A Customer Satisfaction Return form has been completed, and
3. The returned box has been given to, and the completed Customer Satisfaction Return form has been authorized by your Hollister Territory Manager or a Customer Service Representative.

Prior to Use

Prior to use, be sure to read the Instructions for Use for information regarding intended use, contraindications, warnings, precautions and instructions.

Latex Statement

The specifications for the Hollister Continence Care products listed below include natural rubber latex as a component; our suppliers add natural rubber latex as part of the production process. If a Hollister product(s) found in this catalog is not listed below, the specification for that product(s) and packaging does not include natural rubber latex as a component, and our component suppliers do not add natural rubber as part of their production process. Product # 9206, 9207, 9208, 9209, 9106, 9107, 9108, 9109

Ostomy Care

A woman with dark hair in a braid and glasses, wearing a white shirt with thin vertical stripes, is smiling and working with a piece of dough on a wooden cutting board. A laptop is open in front of her, and a glass jar with red beads is on the table. The background is a bright, modern kitchen with white shelves.

Healthy skin. Positive outcomes.

Healthy Skin. Positive Outcomes.

Choosing the Right Skin Barrier FIT

FIT is critical to get a good seal around the stoma to prevent leakage and promote comfort.

Barrier Examples



Flat Skin Barrier



Flat Skin Barrier
with Round or Oval
Convex Barrier Ring



Integrated Soft
Convex Barrier



Integrated Convex
Barrier



Integrated Convex
Barrier with Round
or Oval Convex
Barrier Ring

Possible Clinical Case Scenarios



Flat Skin Barrier

- Flat peristomal plane
- Adequate stoma protrusion
- Centered stoma opening



Flexible Convexity

- Stoma flush with skin
- Uneven peristomal plane
- Off-centered stoma opening



Soft Convexity

- Gentle pressure required to
peristomal plane
- Flush stoma
- Peristomal injury



Firm Convexity

- Creases or folds in peristomal plane
- Off-centered stoma opening
- Stoma retracted below
skin surface



Firm Convexity with
More Depth

- Recessed or flushed stoma
- Uneven abdominal contours
- Off-centered stoma opening

Prior to use, be sure to read the Instructions for Use for information regarding intended use, contraindications, warnings, precautions and instructions.

1. Colwell JC, Goldberg MT & Carmel JE (Eds). (2015). WOCN Society Core Curriculum Ostomy Management. Chapter 10. Philadelphia, PA: Wolters Kluwer
2. Hoeflok J et al. Use of Convexity in Ostomy Care: Results of an International Consensus Meeting. JWOCN, January/February 2017

Healthy Skin. Positive Outcomes.

Choosing the Right Skin Barrier FORMULATION

FORMULATION is the customized ingredients responsible for adhesion, erosion resistance and absorption that must be carefully formulated to create an environment where healthy skin can thrive.

Extended Wear Skin Barriers

CeraPlus

- Contains skin barrier characteristics (adhesion, absorption and erosion resistance) and is infused with ceramide—a naturally occurring lipid in the skin which protects against dryness
- Formulation designed to help protect the skin's natural moisture barrier and help maintain good peristomal skin health

Flextend

- Durable with resistance to erosion
- Swells when it comes in contact with liquid to help form a seal between the skin and the stoma
- Suitable for those who perspire heavily
- Strong adhesion and flexibility

FormaFlex

- Maintains its shape after being stretched or shaped
- Provides an easy-to-customize opening due to barrier material
- Easy to shape to individual stomas without the use of scissors

Standard Wear Skin Barriers

FlexWear

- Good balance of durability and gentleness providing moderate adhesion, and can handle frequent or infrequent pouch changes
- Resistant to fluid erosion
- Suitable for all skin types, depending on change frequency
- Designed to be worn more than one day

SoftFlex

- Appropriate choice for:
 - Those who have thin, fragile, or sensitive skin (such as those who have undergone chemo/radiation therapy)
 - People who desire frequent removal
- Designed to be gentle to the skin, which allows for frequent pouch changes

Karaya 5

- Well-suited for sensitive skin or in some cases where allergic reactions to synthetic barriers has occurred
- Mildly adhesive—needs a belt to secure pouching system

LEARN MORE



Follow the evidence to help ensure healthy skin

After surgery: Peristomal skin starts healthy



Top layers of the epidermis are removed with *every barrier change*.


It's unavoidable, and the damage isn't always visible.



Epidermal stripping is linked to increased levels
of transepidermal water loss (TEWL)^{1,2}

Healthy Skin:

Keeps moisture in and irritants out




Low TEWL
value = higher
skin integrity³

Ceramides are abundant in healthy skin.
They form a waterproof barrier to help keep
moisture in and bacteria out.

— = Ceramides

Damaged Skin:

Allows moisture out and irritants in



Higher TEWL
value = lower
skin integrity³

Damaged skin contains fewer ceramides.

💧 = Water droplet illustrates moisture loss from skin.

References:

1. McNichol L, Lund C, Rosen T, Gray M. Medical adhesives and patient safety: state of the science: consensus statements for the assessment, prevention, and treatment of adhesive-related skin injuries. J Wound Ostomy Continence Nurs. 2013;40(4):365-380.
2. Grove G, Houser T, Sibbald G, Salvadalena G. Measuring epidermal effects of ostomy skin barriers. Skin Res Technol. 2019;25:179-186. <https://doi.org/10.1111/srt.12630>
3. Nybaek H, Lophagen S, Karlsmark T, Bang Knudsen D, Jemec GBE. Stratum corneum integrity as a predictor for peristomal skin problems in ostomates. Br J Dermatol. 2010;162:357-361.

Healthy Skin. Positive Outcomes.

What Happens When Skin is Compromised

Some common signs of a Peristomal Skin Complication (PSC) include:
**ITCHING • REDNESS • PAIN • WEEPING SKIN
BLEEDING • BURNING**



PSCs are often classified as **Peristomal Medical Adhesive-Related Skin Injury (PMARSI)⁹** and **Peristomal Moisture-Associated Skin Damage (PMASD)¹**



PSCs associated with
**INCREASED COST
OF CARE⁶**



PSCs associated with
**LOWER QUALITY
OF LIFE⁷**

What Happens When You Choose CeraPlus Skin Barrier



RELATIVE REDUCTION
in the number of PSCs⁵



INCREASED SATISFACTION
with prevention of itching⁵



Significantly
**LOWER STOMA-RELATED
COST OF CARE⁵**



Significantly
**MORE LIKELY TO
RESOLVE PSCs⁵**



**TEWL IS
MINIMIZED^{2,8}**



Healthier peristomal skin =
HIGHER QUALITY OF LIFE⁷

You have a choice. Choose CeraPlus skin barrier and help ensure healthy peristomal skin for your patients from **DAY ONE.**

References:

- Colwell J, Pittman J, Raizman R, Salvadalena G. A Randomized Controlled Trial Determining Variances in Ostomy Skin Conditions (ADVOCATE). J Wound Ostomy Continence Nurs. 2018;45(1):37-42.
- Taneja, C., et al., Risk and Economic Burden of Peristomal Skin Complications Following Ostomy Surgery. J Wound Ostomy Continence Nurs. 2019. 46(2): 143-149
- Nichols T., Inglese G. The Burden of Peristomal Skin Complications on an Ostomy Population as Assessed by Health Utility and the Physical Component Summary of the SF-36v2 International Society for Pharmacoeconomics and Outcomes Research (ISPOR), Science Direct Value in Health; 2018; p. 92-93, Tables 3-4
- Nichols T, Houser T, Grove G. Comparing the skin stripping effects of three ostomy skin barriers infused with ceramide, honey or aloe; Journal of Stomal Therapy Australia 2019; 39(2):14-18
- Gray M, Colwell J, Doughty D, Goldberg M, Hoek J, Manson A, McNichol L, Rao S. (2013) Peristomal Moisture-Associated Skin Damage in Adults with Fecal Ostomies: A Comprehensive Review and Consensus. J Wound Ostomy Continence Nurs: 38(5): 541-53



New Image

Two-Piece Pouching System

Flat Skin Barriers

CeraPlus with Remois Technology*



Cut-to-fit

Product #	Flange Size	For Stomas Up To	Box of	
11202	1-3/4" (44 mm)	1-1/4" (32 mm)	5	With Tape Border
11203	2-1/4" (57 mm)	1-3/4" (44 mm)		
11204	2-3/4" (70 mm)	2-1/4" (57 mm)		
15102	1-3/4" (44 mm)	1-1/4" (32 mm)	5	Without Tape Border
15103	2-1/4" (57 mm)	1-3/4" (44 mm)		
15104	2-3/4" (70 mm)	2-1/4" (57 mm)		

CeraPlus with Remois Technology*

New!

Pre-sized

Product #	Flange Size	Barrier Opening	Box of	
112013	1-3/4" (44 mm)	7/8" (22 mm)	5	With Tape Border
112014	1-3/4" (44 mm)	1" (25 mm)		
112015	1-3/4" (44 mm)	1-1/8" (29 mm)		
112016	1-3/4" (44 mm)	1-1/4" (32 mm)		
112017	2-1/4" (57 mm)	1-3/8" (35 mm)		
112018	2-1/4" (57 mm)	1-1/2" (38 mm)	5	Without Tape Border
151113	1-3/4" (44 mm)	7/8" (22 mm)		
151114	1-3/4" (44 mm)	1" (25 mm)		
151115	1-3/4" (44 mm)	1-1/8" (29 mm)		
151116	1-3/4" (44 mm)	1-1/4" (32 mm)		
151117	2-1/4" (57 mm)	1-3/8" (35 mm)		
151118	2-1/4" (57 mm)	1-1/2" (38 mm)		

*Remois is a technology of Alcare Co., Ltd.

Flextend



Cut-to-fit

Product #	Flange Size	For Stomas Up To	Box of	
14602	1-3/4" (44 mm)	1-1/4" (32 mm)	5	With Tape Border
14603	2-1/4" (57 mm)	1-3/4" (44 mm)		
14604	2-3/4" (70 mm)	2-1/4" (57 mm)		
14606	4" (102 mm)	3-1/2" (89 mm)		
15602	1-3/4" (44 mm)	1-1/4" (32 mm)	5	Without Tape Border
15603	2-1/4" (57 mm)	1-3/4" (44 mm)		
15604	2-3/4" (70 mm)	2-1/4" (57 mm)		

Pre-sized

Product #	Flange Size	Barrier Opening	Box of	
14701	1-3/4" (44 mm)	5/8" (16 mm)	5	With Tape Border
14702	1-3/4" (44 mm)	3/4" (19 mm)		
14703	1-3/4" (44 mm)	7/8" (22 mm)		
14704	1-3/4" (44 mm)	1" (25 mm)		
14705	1-3/4" (44 mm)	1-1/8" (29 mm)		
14706	1-3/4" (44 mm)	1-1/4" (32 mm)		
14707	2-1/4" (57 mm)	1-3/8" (35 mm)		
14708	2-1/4" (57 mm)	1-1/2" (38 mm)		
14709	2-1/4" (57 mm)	1-3/4" (44 mm)		
16101	1-3/4" (44 mm)	5/8" (16 mm)	5	Without Tape Border
16102	1-3/4" (44 mm)	3/4" (19 mm)		
16103	1-3/4" (44 mm)	7/8" (22 mm)		
16104	1-3/4" (44 mm)	1" (25 mm)		
16105	1-3/4" (44 mm)	1-1/8" (29 mm)		
16106	1-3/4" (44 mm)	1-1/4" (32 mm)		
16107	2-1/4" (57 mm)	1-3/8" (35 mm)		
16108	2-1/4" (57 mm)	1-1/2" (38 mm)		

The CeraPlus Skin Barrier
delivers both Fit and Formulation

Ceramide Infused Barrier

- **Protects** the skin's natural moisture barrier
- **Decreases** TEWL (transepidermal water loss) from damaged or eroded skin
- **Maintains** peristomal skin health



New Image

Two-Piece Pouching System

Flat Skin Barriers

FormaFlex



Shape-to-fit

Product #	Flange Size	For Stomas Up To	Box of	
14102	1-3/4" (44 mm)	1-1/4" (32 mm)	5	With Tape Border
14103	2-1/4" (57 mm)	1-11/16" (43 mm)		
14104	2-3/4" (70 mm)	2-1/4" (57 mm)		

Colour Matching System

Two-Piece Pouching Systems

To help you select compatible *New Image* pouching system components, every box label is colour-coded.

1. Select the skin barrier flange size with the appropriate barrier stoma opening and note its colour code.
2. Select a suitable pouch with the desired features in the matching colour.

	1-3/4" (44 mm)
	2-1/4" (57 mm)
	2-3/4" (70 mm)
	4" (102 mm)

FlexWear



Cut-to-fit

Product #	Flange Size	For Stomas Up To	Box of	
14202	1-3/4" (44 mm)	1-1/4" (32 mm)	5	With Tape Border
14203	2-1/4" (57 mm)	1-3/4" (44 mm)		
14204	2-3/4" (70 mm)	2-1/4" (57 mm)		
14206	4" (102 mm)	3-1/2" (89 mm)		
15202	1-3/4" (44 mm)	1-1/4" (32 mm)	5	Without Tape Border
15203	2-1/4" (57 mm)	1-3/4" (44 mm)		
15204	2-3/4" (70 mm)	2-1/4" (57 mm)		


Pre-sized

Product #	Flange Size	Barrier Opening	Box of	
14301	1-3/4" (44 mm)	5/8" (16 mm)	5	With Tape Border
14302	1-3/4" (44 mm)	3/4" (19 mm)		
14303	1-3/4" (44 mm)	7/8" (22 mm)		
14304	1-3/4" (44 mm)	1" (25 mm)		
14305	1-3/4" (44 mm)	1-1/8" (29 mm)		
14306	1-3/4" (44 mm)	1-1/4" (32 mm)		
14307	2-1/4" (57 mm)	1-3/8" (35 mm)		
14308	2-1/4" (57 mm)	1-1/2" (38 mm)		
14309	2-1/4" (57 mm)	1-3/4" (44 mm)		
16401	1-3/4" (44 mm)	5/8" (16 mm)	5	Without Tape Border
16402	1-3/4" (44 mm)	3/4" (19 mm)		
16403	1-3/4" (44 mm)	7/8" (22 mm)		
16404	1-3/4" (44 mm)	1" (25 mm)		
16405	1-3/4" (44 mm)	1-1/8" (29 mm)		
16406	1-3/4" (44 mm)	1-1/4" (32 mm)		
16407	2-1/4" (57 mm)	1-3/8" (35 mm)		
16408	2-1/4" (57 mm)	1-1/2" (38 mm)		

Designed for Diverse Bodies

Adhesive (Tape) Border Barriers are an appropriate option for most patients, especially those with:

- Skin irregularities such as hernias, wrinkles, or scars
- Anatomical challenges due to obesity or cachexia
- Added need for a sense of security, or extra assurances to accommodate an active lifestyle
- Increased desire for comfort and flexibility



New Image

Two-Piece Pouching System

Soft Convex Skin Barriers

CeraPlus with Remois Technology*



Cut-to-fit

Product #	Flange Size	For Stomas Up To	Box of
11702	1-3/4" (44 mm)	1" (25 mm)	5 With Tape Border
11703	2-1/4" (57 mm)	1-1/2" (38 mm)	
11704	2-3/4" (70 mm)	2" (51 mm)	

Pre-sized

Product #	Flange Size	Barrier Opening	Box of
11903	1-3/4" (44 mm)	7/8" (22 mm)	5 With Tape Border
11904	1-3/4" (44 mm)	1" (25 mm)	
11905	2-1/4" (57 mm)	1-1/8" (29 mm)	
11906	2-1/4" (57 mm)	1-1/4" (32 mm)	

SCAN ME



Convex Skin Barriers

CeraPlus with Remois Technology*



Cut-to-fit

Product #	Flange Size	For Stomas Up To	Box of
11402	1-3/4" (44 mm)	1" (25 mm)	5 With Tape Border
11403	2-1/4" (57 mm)	1-1/2" (38 mm)	
11404	2-3/4" (70 mm)	2" (51 mm)	
15302	1-3/4" (44 mm)	1" (25 mm)	5 Without Tape Border
15303	2-1/4" (57 mm)	1-1/2" (38 mm)	
15304	2-3/4" (70 mm)	2" (51 mm)	


Pre-sized

Product #	Flange Size	Barrier Opening	Box of
11503	1-3/4" (44 mm)	7/8" (22 mm)	5 With Tape Border
11504	1-3/4" (44 mm)	1" (25 mm)	
11505	2-1/4" (57 mm)	1-1/8" (29 mm)	
11506	2-1/4" (57 mm)	1-1/4" (32 mm)	
153103	1-3/4" (44 mm)	7/8" (22 mm)	5 Without Tape Border
153104	1-3/4" (44 mm)	1" (25 mm)	
153105	2-1/4" (57 mm)	1-1/8" (29 mm)	
153106	2-1/4" (57 mm)	1-1/4" (32 mm)	

New!
New!
New!
New!

*Remois is a technology of Alcare Co., Ltd.

Soft Convex CeraPlus skin barrier is designed to help prevent leakage



New Image

Two-Piece Pouching System

Convex Skin Barriers

Flextend



Cut-to-fit

Product #	Flange Size	For Stomas Up To	Box of	
14802	1-3/4" (44 mm)	1" (25 mm)	5	With Tape Border
14803	2-1/4" (57 mm)	1-1/2" (38 mm)		
14804	2-3/4" (70 mm)	2" (51 mm)		
15802	1-3/4" (44 mm)	1" (25 mm)	5	Without Tape Border
15803	2-1/4" (57 mm)	1-1/2" (38 mm)		
15804	2-3/4" (70 mm)	2" (51 mm)		

Pre-sized

Product #	Flange Size	Barrier Opening	Box of	
13901	1-3/4" (44 mm)	5/8" (16 mm)	5	With Tape Border (Tapered)
13902	1-3/4" (44 mm)	3/4" (19 mm)		
13903	1-3/4" (44 mm)	7/8" (22 mm)		
13904	1-3/4" (44 mm)	1" (25 mm)		
13905	2-1/4" (57 mm)	1-1/8" (29 mm)		
13906	2-1/4" (57 mm)	1-1/4" (32 mm)		
13907	2-1/4" (57 mm)	1-3/8" (35 mm)		
13908	2-1/4" (57 mm)	1-1/2" (38 mm)		
13909	2-3/4" (70 mm)	1-5/8" (41 mm)		
13910	2-3/4" (70 mm)	1-3/4" (44 mm)		
13911	2-3/4" (70 mm)	2" (51 mm)		
14901	1-3/4" (44 mm)	5/8" (16 mm)	5	With Tape Border
14902	1-3/4" (44 mm)	3/4" (19 mm)		
14903	1-3/4" (44 mm)	7/8" (22 mm)		
14904	1-3/4" (44 mm)	1" (25 mm)		
14905	2-1/4" (57 mm)	1-1/8" (29 mm)		
14906	2-1/4" (57 mm)	1-1/4" (32 mm)		
14907	2-1/4" (57 mm)	1-3/8" (35 mm)		
14908	2-1/4" (57 mm)	1-1/2" (38 mm)		
14909	2-3/4" (70 mm)	1-5/8" (41 mm)		
14910	2-3/4" (70 mm)	1-3/4" (44 mm)		
14911	2-3/4" (70 mm)	2" (51 mm)		
15903	1-3/4" (44 mm)	7/8" (22 mm)	5	Without Tape Border
15904	1-3/4" (44 mm)	1" (25 mm)		
15905	2-1/4" (57 mm)	1-1/8" (29 mm)		
15906	2-1/4" (57 mm)	1-1/4" (32 mm)		

SCAN ME



FlexWear



Cut-to-fit

Product #	Flange Size	For Stomas Up To	Box of	
13402	1-3/4" (44 mm)	1" (25 mm)	5	With Tape Border (Tapered)
13403	2-1/4" (57 mm)	1-1/2" (38 mm)		
13404	2-3/4" (70 mm)	2" (51 mm)		
14402	1-3/4" (44 mm)	1" (25 mm)	5	With Tape Border
14403	2-1/4" (57 mm)	1-1/2" (38 mm)		
14404	2-3/4" (70 mm)	2" (51 mm)		
15402	1-3/4" (44 mm)	1" (25 mm)	5	Without Tape Border
15403	2-1/4" (57 mm)	1-1/2" (38 mm)		
15404	2-3/4" (70 mm)	2" (51 mm)		

Pre-sized

Product #	Flange Size	Barrier Opening	Box of	
13503	1-3/4" (44 mm)	7/8" (22 mm)	5	With Tape Border (Tapered)
13504	1-3/4" (44 mm)	1" (25 mm)		
13505	2-1/4" (57 mm)	1-1/8" (29 mm)		
13506	2-1/4" (57 mm)	1-1/4" (32 mm)		
14503	1-3/4" (44 mm)	7/8" (22 mm)	5	With Tape Border
14504	1-3/4" (44 mm)	1" (25 mm)		
14505	2-1/4" (57 mm)	1-1/8" (29 mm)		
14506	2-1/4" (57 mm)	1-1/4" (32 mm)		
15503	1-3/4" (44 mm)	7/8" (22 mm)	5	Without Tape Border
15504	1-3/4" (44 mm)	1" (25 mm)		
15505	2-1/4" (57 mm)	1-1/8" (29 mm)		
15506	2-1/4" (57 mm)	1-1/4" (32 mm)		

Choosing the Right Skin Barrier FIT

- FIT is critical to get a good seal around the stoma to prevent leakage and promote comfort.
- Convex barrier products can be utilized to help provide a secure and predictable seal.
- Product selection is made by matching stomal protrusion and peristomal contours to the appropriate barriers.
- The outward curving of the convex skin barrier is designed to interface with the immediate peristomal skin.

Refer to page 4 for additional information on *The Continuum of Convex Barrier Options*.

New Image

Two-Piece Pouching System

Drainable Pouches

Lock 'n Roll Closure System



Colour: Ultra-clear

Product #	Flange Size	Pouch Length	Filter	Box of
18002	1-3/4" (44 mm)	12" (30 cm)	No	10
18003	2-1/4" (57 mm)	12" (30 cm)		
18004	2-3/4" (70 mm)	12" (30 cm)		
18006	4" (102 mm)	12" (30 cm)		

Colour: Transparent

Product #	Flange Size	Pouch Length	Filter	Box of
18132	1-3/4" (44 mm)	12" (30 cm)	No	10
18133	2-1/4" (57 mm)	12" (30 cm)		
18134	2-3/4" (70 mm)	12" (30 cm)		
18192	1-3/4" (44 mm)	12" (30 cm)	AF300	10
18193	2-1/4" (57 mm)	12" (30 cm)		
18194	2-3/4" (70 mm)	12" (30 cm)		
18292	1-3/4" (44 mm)	7" (18 cm)	AF300	20
18293	2-1/4" (57 mm)	7" (18 cm)		
18294	2-3/4" (70 mm)	7" (18 cm)		

Colour: Beige

Product #	Flange Size	Pouch Length	Filter	Box of
18112	1-3/4" (44 mm)	12" (30 cm)	No	10
18113	2-1/4" (57 mm)	12" (30 cm)		
18114	2-3/4" (70 mm)	12" (30 cm)		
18182	1-3/4" (44 mm)	12" (30 cm)	AF300	10
18183	2-1/4" (57 mm)	12" (30 cm)		
18184	2-3/4" (70 mm)	12" (30 cm)		
18282	1-3/4" (44 mm)	7" (18 cm)	AF300	20
18283	2-1/4" (57 mm)	7" (18 cm)		
18284	2-3/4" (70 mm)	7" (18 cm)		

Clamp Closure



Colour: Ultra-clear

Product #	Flange Size	Pouch Length	Filter	Box of
18172	1-3/4" (44 mm)	12" (30 cm)	No	10
18173	2-1/4" (57 mm)	12" (30 cm)		
18174	2-3/4" (70 mm)	12" (30 cm)		
18176	4" (102 mm)	12" (30 cm)		

Colour: Transparent

Product #	Flange Size	Pouch Length	Filter	Box of
18102	1-3/4" (44 mm)	12" (30 cm)	No	10
18103	2-1/4" (57 mm)	12" (30 cm)		
18104	2-3/4" (70 mm)	12" (30 cm)		

Colour: Beige

Product #	Flange Size	Pouch Length	Filter	Box of
18122	1-3/4" (44 mm)	12" (30 cm)	No	10
18123	2-1/4" (57 mm)	12" (30 cm)		
18124	2-3/4" (70 mm)	12" (30 cm)		
18142	1-3/4" (44 mm)	12" (30 cm)	AF300	10
18143	2-1/4" (57 mm)	12" (30 cm)		
18144	2-3/4" (70 mm)	12" (30 cm)		
18202	1-3/4" (44 mm)	9" (23 cm)	No	10
18203	2-1/4" (57 mm)	9" (23 cm)		
18204	2-3/4" (70 mm)	9" (23 cm)		

Lock 'n Roll Microseal Closure Features:

- A secure closure with tiny microseals
- An easy-to-open-and-close seal
- The material used on the closure is water-resistant and can be securely closed even when it is wet, but will work best when dry.

New Image

Two-Piece Pouching System

Closed Pouches



Colour: Transparent

Product #	Flange Size	Pouch Length	Filter	Box of
18362	1-3/4" (44 mm)	9" (23 cm)	AF300	60
18363	2-1/4" (57 mm)	9" (23 cm)		
18364	2-3/4" (70 mm)	9" (23 cm)		

Colour: Beige QuietWear Fabric

Product #	Flange Size	Pouch Length	Filter	Box of
18332	1-3/4" (44 mm)	9" (23 cm)	No	30
18333	2-1/4" (57 mm)	9" (23 cm)		
18334	2-3/4" (70 mm)	9" (23 cm)		
18322	1-3/4" (44 mm)	9" (23 cm)	AF300	30
18323	2-1/4" (57 mm)	9" (23 cm)		
18324	2-3/4" (70 mm)	9" (23 cm)		

Colour: Beige

Product #	Flange Size	Pouch Length	Filter	Box of
18352	1-3/4" (44 mm)	7" (18 cm)	No	30
18353	2-1/4" (57 mm)	7" (18 cm)		
18354	2-3/4" (70 mm)	7" (18 cm)		
18372-30	1-3/4" (44 mm)	9" (23 cm)	AF300	30
18373-30	2-1/4" (57 mm)	9" (23 cm)		
18374-30	2-3/4" (70 mm)	9" (23 cm)		
18392	1-3/4" (44 mm)	7" (18 cm)	AF300	60
18393	2-1/4" (57 mm)	7" (18 cm)		
18394	2-3/4" (70 mm)	7" (18 cm)		

High Output Pouches



Soft Tap Closure System

Colour: Ultra-Clear

Product #	Flange Size	Pouch Length	Box of
18012	1-3/4" (44 mm)	12" (30 cm)	10
18013	2-1/4" (57 mm)	12" (30 cm)	
18014	2-3/4" (70 mm)	12" (30 cm)	
18016	4" (102 mm)	12" (30 cm)	

New Image

Two-Piece Pouching System

Urostomy Pouches

Enhanced Design



Colour: Ultra-Clear

Product #	Flange Size	Pouch Length	Box of
18922	1-3/4" (44 mm)	9" (23 cm)	10
18923	2-1/4" (57 mm)	9" (23 cm)	
18924	2-3/4" (70 mm)	9" (23 cm)	

Colour: Transparent

Product #	Flange Size	Pouch Length	Box of
18902	1-3/4" (44 mm)	9" (23 cm)	10
18903	2-1/4" (57 mm)	9" (23 cm)	
18904	2-3/4" (70 mm)	9" (23 cm)	

Colour: Beige

Product #	Flange Size	Pouch Length	Box of
18912	1-3/4" (44 mm)	9" (23 cm)	10
18913	2-1/4" (57 mm)	9" (23 cm)	
18914	2-3/4" (70 mm)	9" (23 cm)	

Original Design



Colour: Transparent

Product #	Flange Size	Pouch Length	Box of
18402	1-3/4" (44 mm)	9" (23 cm)	10
18403	2-1/4" (57 mm)	9" (23 cm)	
18404	2-3/4" (70 mm)	9" (23 cm)	

Colour: Beige

Product #	Flange Size	Pouch Length	Box of
18412	1-3/4" (44 mm)	9" (23 cm)	10
18413	2-1/4" (57 mm)	9" (23 cm)	
18414	2-3/4" (70 mm)	9" (23 cm)	

A Great Pouch Becomes Even Better

Enhanced Flow-Control Valve

The new valve is made with softer material and a streamlined shape.

New Multi-Chamber Design

The updated pouches now have multiple built-in chambers, which are designed to facilitate even distribution of urine.



Multi Chamber
Front View



Multi Chamber,
Top View, Sliced

SCAN ME





New Image

Two-Piece Pouching System

Ostomy Kits









Drainable

CeraPlus with Remois Technology*
Lock 'n Roll Closure System









Non-sterile				
Product #	Flange Size	For Stomas Up To	Pouch Length	Box of
19603	 2-1/4" (57 mm)	1-3/4" (44 mm)	12" (30 mm)	5
19604	 2-3/4" (70 mm)	2-1/4" (57 mm)	12" (30 mm)	

*Remois is a technology of Alcare Co., Ltd.

FlexWear
Lock 'n Roll Closure System

Non-sterile					
Product #		Flange Size	For Stomas Up To	Pouch Length	Box of
19002		1-3/4" (44 mm)	1-1/4" (32 mm)	12" (30 mm)	5
19003		2-1/4" (57 mm)	1-3/4" (44 mm)	12" (30 mm)	
19004		2-3/4" (70 mm)	2-1/4" (57 mm)	12" (30 mm)	
19006		4" (102 mm)	3-1/2" (89 mm)	12" (30 mm)	
Sterile					
Product #		Flange Size	For Stomas Up To	Pouch Length	Box of
19052		1-3/4" (44 mm)	1-1/4" (32 mm)	12" (30 mm)	5
19053		2-1/4" (57 mm)	1-3/4" (44 mm)	12" (30 mm)	
19054		2-3/4" (70 mm)	2-1/4" (57 mm)	12" (30 mm)	
19056		4" (102 mm)	3-1/2" (89 mm)	12" (30 mm)	

FlexWear
Clamp Closure System


Non-sterile					
Product #		Flange Size	For Stomas Up To	Pouch Length	Box of
19102		1-3/4" (44 mm)	1-1/4" (32 mm)	12" (30 mm)	5
19103		2-1/4" (57 mm)	1-3/4" (44 mm)	12" (30 mm)	
19104		2-3/4" (70 mm)	2-1/4" (57 mm)	12" (30 mm)	
19106		4" (102 mm)	3-1/2" (89 mm)	12" (30 mm)	
Sterile					
Product #		Flange Size	For Stomas Up To	Pouch Length	Box of
19152		1-3/4" (44 mm)	1-1/4" (32 mm)	12" (30 mm)	5
19153		2-1/4" (57 mm)	1-3/4" (44 mm)	12" (30 mm)	
19154		2-3/4" (70 mm)	2-1/4" (57 mm)	12" (30 mm)	
19156		4" (102 mm)	3-1/2" (89 mm)	12" (30 mm)	

Ostomy Kits Include:

- Skin barrier with tape border
- Ultra-clear drainable pouch
- Loop ostomy bridge with 4" flange







Urostomy Kits

CeraPlus with Remois Technology*

Non-sterile			
Product #	Flange Size	For Stomas Up To	Box of
19903	 2-1/4" (57 mm)	1-3/4" (44 mm)	5

*Remois is a technology of Alcare Co., Ltd.

Flextend

Non-sterile				
Product #		Flange Size	For Stomas Up To	Box of
19202		1-3/4" (44 mm)	1-1/4" (32 mm)	5
19203		2-1/4" (57 mm)	1-3/4" (44 mm)	
19204		2-3/4" (70 mm)	2-1/4" (57 mm)	
Sterile				
Product #		Flange Size	For Stomas Up To	Box of
19252		1-3/4" (44 mm)	1-1/4" (32 mm)	5
19253		2-1/4" (57 mm)	1-3/4" (44 mm)	
19254		2-3/4" (70 mm)	2-1/4" (57 mm)	

Urostomy Kits Include:

- Ultra-clear pouch
- Cut-to-fit Flextend Skin Barrier
- Drain tube adapter

Premier


One-Piece Pouching System


Flat Skin Barriers

Drainable Pouch

CeraPlus with Remois Technology*
Lock 'n Roll Closure System



Cut-to-fit With Tape Border 					
Product #	Colour	For Stomas Up To	Pouch Length	Filter	Box of
8931	Ultra-clear	2-1/2" (64 mm)	12" (30 cm)	No	10
8901	Beige with Viewing Option	2-1/2" (64 mm)	12" (30 cm)	AF300	


Cut-to-fit Coming Soon Without Tape Border 					
Product #	Colour	For Stomas Up To	Pouch Length	Filter	Box of
89011	Beige with Viewing Option	2-1/2" (64 mm)	12" (30 cm)	AF300	10

Pre-sized With Tape Border 					
Product #	Colour	Barrier Opening	Pouch Length	Filter	Box of
8925	Beige with Viewing Option	1" (25 mm)	12" (30 cm)	AF300	10
8930		1-3/16" (30 mm)	12" (30 cm)		
8935		1-3/8" (35 mm)	12" (30 cm)		
8940		1-9/16" (40 mm)	12" (30 cm)		

*Remois is a technology of Alcare Co., Ltd.




FlexWear
Lock 'n Roll Closure System

Cut-to-fit With Tape Border 					
Product #	Colour	For Stomas Up To	Pouch Length	Filter	Box of
8031	Ultra-clear	2-1/2" (64 mm)	12" (30 cm)	No	10
8081	Transparent	2-1/2" (64 mm)	12" (30 cm)	AF300	
8181	Beige	2-1/2" (64 mm)	12" (30 cm)	AF300	


FlexWear
Clamp Closure System




Cut-to-fit With Tape Border 					
Product #	Colour	For Stomas Up To	Pouch Length	Filter	Box of
8131	Ultra-clear	2-1/2" (64 mm)	12" (30 cm)	No	10

Flextend
Lock 'n Roll Closure System





Cut-to-fit With Tape Border 					
Product #	Colour	For Stomas Up To	Pouch Length	Filter	Box of
8531	Ultra-clear	2-1/2" (64 mm)	12" (30 cm)	No	10
8331	Transparent	2-1/2" (64 mm)	12" (30 cm)	AF300	
8551	Beige	2-1/2" (64 mm)	12" (30 cm)	No	

Pre-sized With Tape Border 				
Product #	Colour	Barrier Opening	Pouch Length	Box of
8538	Transparent	1-1/4" (32 mm)	12" (30 cm)	10
8552	Beige	1" (25 mm)	12" (30 cm)	10
8558		1-1/4" (32 mm)	12" (30 cm)	
8559		1-3/4" (44 mm)	12" (30 cm)	

Flextend
Clamp Closure System



Cut-to-fit With Tape Border 				
Product #	Colour	For Stomas Up To	Pouch Length	Box of
8651	Ultra-clear	2-1/2" (64 mm)	12" (30 cm)	10
8631	Transparent	2-1/2" (64 mm)	12" (30 cm)	10
8641		2-1/2" (64 mm)	9" (23 cm)	

Pre-sized With Tape Border 				
Product #	Colour	Barrier Opening	Pouch Length	Box of
8637	Transparent	3/4" (19 mm)	12" (30 cm)	10
8632		1" (25 mm)	12" (30 cm)	
8638		1-1/4" (32 mm)	12" (30 cm)	
8633		1-1/2" (38 mm)	12" (30 cm)	
8639		1-3/4" (44 mm)	12" (30 cm)	
8634		2" (51 mm)	12" (30 cm)	
8635		2-1/2" (64 mm)	12" (30 cm)	
8647	Beige	3/4" (19 mm)	9" (23 cm)	10
8642		1" (25 mm)	9" (23 cm)	
8648		1-1/4" (32 mm)	9" (23 cm)	
8643		1-1/2" (38 mm)	9" (23 cm)	
8649		1-3/4" (44 mm)	9" (23 cm)	
8644		2" (51 mm)	9" (23 cm)	

Premier


One-Piece Pouching System


Flat Skin Barriers


Drainable Pouch

SoftFlex
Lock 'n Roll Closure System



Cut-to-fit Without Tape Border 					
Product #	Colour	For Stomas Up To	Pouch Length	Filter	Box of
88300	Beige	2-1/8" (55 mm)	12" (30 cm)	AF300	10
88700		2-1/8" (55 mm)	7" (18 cm)		20


Cut-to-fit—Oval Skin Barrier Without Tape Border 					
Product #	Colour	For Stomas Up To	Pouch Length	Filter	Box of
88402	Transparent	2-1/2" x 3" (66 x 77 mm)	12" (30 cm)	AF300	10
88302	Beige	2-1/2" x 3" (66 x 77 mm)	12" (30 cm)		

Pre-sized Without Tape Border 					
Product #	Colour	Barrier Opening	Pouch Length	Filter	Box of
88325	Beige	1" (25 mm)	12" (30 cm)	AF300	10
88330		1-3/16" (30 mm)	12" (30 cm)		
88335		1-3/8" (35 mm)	12" (30 cm)		
88340		1-9/16" (40 mm)	12" (30 cm)		
88725		1" (25 mm)	7" (18 cm)	AF300	20
88730		1-3/16" (30 mm)	7" (18 cm)		
88735		1-3/8" (35 mm)	7" (18 cm)		
88740		1-9/16" (40 mm)	7" (18 cm)		

High Output Pouches

Flextend Design




Non-sterile Without Tape Border 			
Product #	Barrier Opening	Pouch Length	Box of
80070	2-3/4" (70 mm)	12" (30 cm)	10
80110	4-1/3" (110 mm)	12" (30 cm)	


Flat Skin Barriers


Closed Pouch

SoftFlex



Cut-to-fit Without Tape Border 					
Product #	Colour	For Stomas Up To	Pouch Length	Filter	Box of
82300	Beige	2-1/8" (55 mm)	9" (23 cm)	AF300	30
82100		2-1/8" (55 mm)	7" (18 cm)		

Cut-to-fit—Oval Skin Barrier Without Tape Border 					
Product #	Colour	For Stomas Up To	Pouch Length	Filter	Box of
82402	Transparent	2-1/2" x 3" (66 mm x 77 mm)	9" (23 cm)	AF300	30
82302	Beige	2-1/2" x 3" (66 mm x 77 mm)	9" (23 cm)		

Pre-sized Without Tape Border 					
Product #	Colour	Barrier Opening	Pouch Length	Filter	Box of
82325	Beige	1" (25 mm)	9" (23 cm)	AF300	30
82330		1-3/16" (30 mm)	9" (23 cm)		
82335		1-3/8" (35 mm)	9" (23 cm)		
82125		1" (25 mm)	7" (18 cm)		
82130		1-3/16" (30 mm)	7" (18 cm)		
82135		1-3/8" (35 mm)	7" (18 cm)		

Did you know?

The integrated AF300 Filter is waterproof, and designed to allow gas and air to escape, but keeps odour and output in the pouch. Liquid cannot penetrate the filter, therefore not requiring a protective sticker.

Premier


One-Piece Pouching System


Soft Convex Skin Barriers


Drainable Pouch

CeraPlus with Remois Technology*
Lock 'n Roll Closure System



Cut-to-fit					With Tape Border 
Product #		For Stomas Up To	Pouch Length	Filter	Box of
8954	Beige with Viewing Option	1" (25 mm)	12" (30 cm)	AF300	5
8958		1-1/2" (38 mm)	12" (30 cm)		
89511		2-1/8" (55 mm)	12" (30 cm)		

Cut-to-fit					Without Tape Border 
Product #		For Stomas Up To	Pouch Length	Filter	Box of
89518	Beige with Viewing Option	1-1/2" (38 mm)	12" (30 cm)	AF300	5
8195111		2-1/8" (55 mm)	12" (30 cm)		


Pre-sized					With Tape Border 
Product #		Barrier Opening	Pouch Length	Filter	Box of
8960	Beige with Viewing Option	3/4" (19 mm)	12" (30 cm)	AF300	5
8961		7/8" (22 mm)	12" (30 cm)		
8962		1" (25 mm)	12" (30 cm)		
8963		1-1/8" (29 mm)	12" (30 cm)		


*Remois is a technology of Alcare Co., Ltd.

Flextend

Lock 'n Roll Closure System



Cut-to-fit					With Tape Border 
Product #	Colour	Stoma Size	Pouch Length	Filter	Box of
8674	Beige with Viewing Option	5/8" - 1" (13 - 25 mm)	12" (30 cm)	AF300	5
8678		5/8" - 1-1/2" (13 - 38 mm)	12" (30 cm)		
86711		5/8" - 2-1/8" (13 - 55 mm)	12" (30 cm)		


Pre-sized					With Tape Border 
Product #	Colour	Barrier Opening	Pouch Length	Filter	Box of
8661	Beige with Viewing Option	7/8" (22 mm)	12" (30 cm)	AF300	5
8662		1" (25 mm)	12" (30 cm)		
8663		1-1/8" (29 mm)	12" (30 cm)		
8664		1-1/4" (32 mm)	12" (30 cm)		
8665		1-1/2" (38 mm)	12" (30 cm)		
8666		1-3/4" (44 mm)	12" (30 cm)		

Convex Skin Barriers

Drainable Pouch

CeraPlus with Remois Technology*
Lock 'n Roll Closure System



Cut-to-fit					With Tape Border 
Product #		For Stomas Up To	Pouch Length	Filter	Box of
8984	Ultra-clear	1" (25 mm)	12" (30 cm)	No	5
8988		1-1/2" (38 mm)	12" (30 cm)		
89811		2" (51 mm)	12" (30 cm)		
8914	Beige with Viewing Option	1" (25 mm)	12" (30 cm)	AF300	5
8918		1-1/2" (38 mm)	12" (30 cm)		
89111		2" (51 mm)	12" (30 cm)		

Pre-sized					With Tape Border 
Product #		Barrier Opening	Pouch Length	Filter	Box of
8990	Beige with Viewing Option	3/4" (19 mm)	12" (30 cm)	AF300	5
8991		7/8" (22 mm)	12" (30 cm)		
8992		1" (25 mm)	12" (30 cm)		
8993		1-1/8" (29 mm)	12" (30 cm)		

*Remois is a technology of Alcare Co., Ltd.

Soft Convex Skin Barriers

Soft convexity is designed to help provide gentle peristomal pressure for an individualized fit. It helps conform to accommodate uneven skin surfaces, and is a more flexible alternative to firm convexity.

- Lock 'n Roll closure
- Tape border
- Beige with viewing option
- AF300 filter
- Belt tabs



LEARN MORE



Premier

One-Piece Pouching System

Convex Skin Barriers

Drainable Pouch

Flexextend
Lock 'n Roll Closure System



Cut-to-fit With Tape Border

Colour: Transparent

Product #		For Stomas Up To	Pouch Length	Filter	Box of
Tapered Border	Non-Tapered Border				
8244	8524	1" (25 mm)	12" (30 cm)	No	5
8248	8528	1-1/2" (38 mm)	12" (30 cm)		
82411	85211	2" (51 mm)	12" (30 cm)		
8274	8584	1" (25 mm)	12" (30 cm)		
8278	8588	1-1/2" (38 mm)	12" (30 cm)	AF300	
82711	85811	2" (51 mm)	12" (30 cm)		

Colour: Beige

8214	8504	1" (25 mm)	12" (30 cm)	No	5
8218	8508	1-1/2" (38 mm)	12" (30 cm)		
82111	85011	2" (51 mm)	12" (30 cm)		

Pre-sized With Tape Border

Colour: Transparent

Product #		Barrier Opening	Pouch Length	Box of
Tapered Border	Non-Tapered Border			
8280	8510	3/4" (19 mm)	12" (30 cm)	5
8281	8511	7/8" (22 mm)	12" (30 cm)	
8282	8512	1" (25 mm)	12" (30 cm)	
8283	8513	1-1/8" (29 mm)	12" (30 cm)	
8284	8514	1-1/4" (32 mm)	12" (30 cm)	
8285	8515	1-3/8" (35 mm)	12" (30 cm)	
8286	85116	1-1/2" (38 mm)	12" (30 cm)	
8288	85118	1-3/4" (44 mm)	12" (30 cm)	
8289	85119	2" (51 mm)	12" (30 cm)	

Colour: Beige

8290	8590	3/4" (19 mm)	12" (30 cm)	5
8291	8591	7/8" (22 mm)	12" (30 cm)	
8292	8592	1" (25 mm)	12" (30 cm)	
8293	8593	1-1/8" (29 mm)	12" (30 cm)	
8294	8594	1-1/4" (32 mm)	12" (30 cm)	
8295	8595	1-3/8" (35 mm)	12" (30 cm)	
8296	8596	1-1/2" (38 mm)	12" (30 cm)	
8297	8597	1-5/8" (41 mm)	12" (30 cm)	

Flexextend
Clamp Closure System



Cut-to-fit With Tape Border

Colour: Transparent

Product #		For Stomas Up To	Pouch Length	Box of
Tapered Border	Non-Tapered Border			
82624	8624	1" (25 mm)	12" (30 cm)	5
82628	8628	1-1/2" (38 mm)	12" (30 cm)	
826211	86211	2" (51 mm)	12" (30 cm)	

Pre-sized With Tape Border

Colour: Transparent

Product #		Barrier Opening	Pouch Length	Box of
Tapered Border	Non-Tapered Border			
82610	8610	3/4" (19 mm)	12" (30 cm)	5
82611	8611	7/8" (22 mm)	12" (30 cm)	
82612	8612	1" (25 mm)	12" (30 cm)	
82613	8613	1-1/8" (29 mm)	12" (30 cm)	
82614	8614	1-1/4" (32 mm)	12" (30 cm)	
82615	8615	1-3/8" (35 mm)	12" (30 cm)	
82616	8616	1-1/2" (38 mm)	12" (30 cm)	
82617	8617	1-5/8" (41 mm)	12" (30 cm)	
82618	8618	1-3/4" (44 mm)	12" (30 cm)	
82619	8619	2" (51 mm)	12" (30 cm)	

Tapered Border is Designed for Enhanced Sense of Security

- Tapered barrier edges help prevent lift from moisture or motion, which can help achieve desired wear time
- Gentle on the skin while still giving the feeling of extra security
- Adhesive flexes and bends to move comfortably with the body
- A new flexible plastic release liner that covers both the skin barrier and the adhesive border to help improve ease-of-use



A great fit. A formulation designed to support healthy skin.



The CeraPlus Soft Convex One-Piece Urostomy Pouching System

Soft convex fit

Flexible design helps provide gentle pressure around the stoma for a better fit

Adhesive border barrier

Additional perimeter of thin adhesive provides extra peace of mind

Ceramide-infused barrier

CeraPlus skin barrier is infused with ceramide to help protect the skin's natural moisture barrier

Multi-chamber design

Multiple built-in chambers are designed to facilitate a balanced collection of urine, resulting in a discreet pouch profile

Soft, flow-control tap

Made with soft material and a streamlined shape, the tap helps provide comfort against the body and makes it easy to identify when valve is closed

Improved adapter

The new side-tab design helps promote ease of use when connecting to any Premier/Moderma Flex one-piece urostomy pouch to a bedside drainage collection system



Ostomy Care

Healthy skin. Positive outcomes.

Premier


One-Piece Pouching System

Flat Skin Barriers

Urostomy Pouch

CeraPlus with Remois Technology*
Enhanced Design



Cut-to-fit					With Tape Border 
Product #	Colour	For Stomas	Pouch Length	Filter	Box of
84200	Ultra-clear	2-1/2" (64 mm)	9" (23 cm)	No	10

*Remois is a technology of Alcare Co., Ltd.


The CeraPlus Skin Barrier
delivers both Fit and Formulation


Ceramide Infused Barrier

- **Protects** the skin's natural moisture barrier
- **Decreases** TEWL (transepidermal water loss) from damaged or eroded skin
- **Maintains** peristomal skin health

Flextend
Enhanced Design





Cut-to-fit					With Tape Border 
Product #	Colour	For Stomas Up To	Pouch Length	Box of	
84590	Ultra-clear	2-1/2" (64 mm)	9" (23 cm)	10	
84490	Beige	2-1/2" (64 mm)	9" (23 cm)		

Pre-sized					With Tape Border 
Product #	Colour	Barrier Opening	Pouch Length	Box of	
84697	Transparent	3/4" (19 mm)	9" (23 cm)	10	
84692		1" (25 mm)	9" (23 cm)		
84698		1-1/4" (32 mm)	9" (23 cm)		
84693		1-1/2" (38 mm)	9" (23 cm)		
84699		1-3/4" (44 mm)	9" (23 cm)		
84694		2" (51 mm)	9" (23 cm)		
84497	Beige	3/4" (19 mm)	9" (23 cm)	10	
84492		1" (25 mm)	9" (23 cm)		
84498		1-1/4" (32 mm)	9" (23 cm)		
84493		1-1/2" (38 mm)	9" (23 cm)		
84499		1-3/4" (44 mm)	9" (23 cm)		
84494		2" (51 mm)	9" (23 cm)		

Flextend
Original Design



Cut-to-fit					With Tape Border 
Product #	Colour	For Stomas Up To	Pouch Length	Box of	
8450	Ultra-clear	2-1/2" (64 mm)	9" (23 cm)	10	
8460	Transparent	2-1/2" (64 mm)	9" (23 cm)		
8440	Beige	2-1/2" (64 mm)	9" (23 cm)		

Pre-sized					With Tape Border 
Product #	Colour	Barrier Opening	Pouch Length	Box of	
8467	Transparent	3/4" (19 mm)	9" (23 cm)	10	
8462		1" (25 mm)	9" (23 cm)		
8468		1-1/4" (32 mm)	9" (23 cm)		
8463		1-1/2" (38 mm)	9" (23 cm)		
8469		1-3/4" (44 mm)	9" (23 cm)		
8464		2" (51 mm)	9" (23 cm)		
8447	Beige	3/4" (19 mm)	9" (23 cm)	10	
8442		1" (25 mm)	9" (23 cm)		
8448		1-1/4" (32 mm)	9" (23 cm)		
8443		1-1/2" (38 mm)	9" (23 cm)		
8449		1-3/4" (44 mm)	9" (23 cm)		


Premier


One-Piece Pouching System

Soft Convex Skin Barriers Urostomy Pouch

CeraPlus with Remois Technology*
Enhanced Design



Cut-to-fit				With Tape Border 
Product #	Colour	For Stomas Up To	Pouch Length	Box of
84134	Ultra-clear	1" (25 mm)	9" (23 cm)	5
84138		1-1/2" (38 mm)	9" (23 cm)	
841311		2-1/8" (55 mm)	9" (23 cm)	


Pre-sized				With Tape Border 
Product #	Colour	Barrier Opening	Pouch Length	Box of
8412	Beige with Viewing Option	3/4" (19 mm)	9" (23 cm)	5
8413		7/8" (22 mm)	9" (23 cm)	
8414		1" (25 mm)	9" (23 cm)	
8415		1-1/8" (29 mm)	9" (23 cm)	


*Remois is a technology of Alcare Co., Ltd.

Convex Skin Barriers Urostomy Pouch

CeraPlus with Remois Technology*
Enhanced Design



Cut-to-fit				With Tape Border 
Product #	Colour	For Stomas Up To	Pouch Length	Box of
84404	Ultra-clear	1" (25 mm)	9" (23 cm)	5
84408		1-1/2" (38 mm)	9" (23 cm)	
844011		2" (51 mm)	9" (23 cm)	

Pre-sized				With Tape Border 
Product #	Colour	Barrier Opening	Pouch Length	Box of
84292	Beige with Viewing Option	3/4" (19 mm)	9" (23 cm)	5
84293		7/8" (22 mm)	9" (23 cm)	
84294		1" (25 mm)	9" (23 cm)	
84295		1-1/8" (29 mm)	9" (23 cm)	

*Remois is a technology of Alcare Co., Ltd.

Premier Soft Convex CeraPlus Pouching Systems:

- Are designed to conform to the body without compromising fit
- Deliver both convex support and flexibility
- Protect the skin's natural moisture barrier and help maintain good peristomal skin health

LEARN MORE




A Great Pouch Becomes Even Better

Enhanced Flow-Control Valve

The new valve is made with softer material and a streamlined shape.

New Multi-Chamber Design

The updated pouches now have multiple built-in chambers, which are designed to facilitate even distribution of urine.



Multi Chamber
Front View



Multi Chamber,
Top View, Sliced

Premier


One-Piece Pouching System


Convex Skin Barriers

Urostomy Pouch

Flextend
Enhanced Design



Cut-to-fit			With Tape Border 	
Product #	Colour	For Stomas Up To	Pouch Length	Box of
84794	Ultra-clear	1" (25 mm)	9" (23 cm)	5
84798		1-1/2" (38 mm)	9" (23 cm)	
8479111		2" (51 mm)	9" (23 cm)	
84394	Beige	1" (25 mm)	9" (23 cm)	
84398		1-1/2" (38 mm)	9" (23 cm)	
8439111		2" (51 mm)	9" (23 cm)	



Pre-sized			With Tape Border 	
Product #	Colour	Barrier Opening	Pouch Length	Box of
84890	Ultra-clear	1/2" (13 mm)	9" (23 cm)	5
84891		5/8" (16 mm)	9" (23 cm)	
84892		3/4" (19 mm)	9" (23 cm)	
84893		7/8" (22 mm)	9" (23 cm)	
84894		1" (25 mm)	9" (23 cm)	
84895		1-1/8" (29 mm)	9" (23 cm)	
84896		1-1/4" (32 mm)	9" (23 cm)	
84897		1-3/8" (35 mm)	9" (23 cm)	
84898		1-1/2" (38 mm)	9" (23 cm)	
84899		1-3/4" (44 mm)	9" (23 cm)	
84990	Beige	1/2" (13 mm)	9" (23 cm)	5
84991		5/8" (16 mm)	9" (23 cm)	
84992		3/4" (19 mm)	9" (23 cm)	
84993		7/8" (22 mm)	9" (23 cm)	
84994		1" (25 mm)	9" (23 cm)	
84995		1-1/8" (29 mm)	9" (23 cm)	
84996		1-1/4" (32 mm)	9" (23 cm)	
84997		1-3/8" (35 mm)	9" (23 cm)	
84998		1-1/2" (38 mm)	9" (23 cm)	
84999		1-3/4" (44 mm)	9" (23 cm)	

LEARN MORE



Flextend
Original Design



Cut-to-fit			With Tape Border 	
Product #	Colour	For Stomas Up To	Pouch Length	Box of
8474	Transparent	1" (25 mm)	9" (23 cm)	5
8478		1-1/2" (38 mm)	9" (23 cm)	
84711		2" (51 mm)	9" (23 cm)	
Pre-sized			With Tape Border 	
Product #	Colour	Barrier Opening	Pouch Length	Box of
8480	Transparent	1/2" (13 mm)	9" (23 cm)	5
8481		5/8" (16 mm)	9" (23 cm)	
8482		3/4" (19 mm)	9" (23 cm)	
8483		7/8" (22 mm)	9" (23 cm)	
8484		1" (25 mm)	9" (23 cm)	
8485		1-1/8" (29 mm)	9" (23 cm)	
8486		1-1/4" (32 mm)	9" (23 cm)	
8487		1-3/8" (35 mm)	9" (23 cm)	
8488		1-1/2" (38 mm)	9" (23 cm)	
8489		1-3/4" (44 mm)	9" (23 cm)	
8490	Beige	1/2" (13 mm)	9" (23 cm)	5
8491		5/8" (16 mm)	9" (23 cm)	
8492		3/4" (19 mm)	9" (23 cm)	
8493		7/8" (22 mm)	9" (23 cm)	
8494		1" (25 mm)	9" (23 cm)	
8495		1-1/8" (29 mm)	9" (23 cm)	
8496		1-1/4" (32 mm)	9" (23 cm)	
8497		1-3/8" (35 mm)	9" (23 cm)	
8498		1-1/2" (38 mm)	9" (23 cm)	


Adhesive (Tape) Border Barriers

- Feature a thin layer of adhesive that surrounds the hydrocolloid barrier.
- Designed to be gentle on the skin, the adhesive moves with the body to easily conform to skin irregularities, creases and folds.
- Application is simple, with no additional accessories required to gain an added sense of protection. Go with the skin barrier that comes with built-in confidence.

Drainable Pouches


Clamp Closure System



Colour: Transparent				With Tape Border 
Product #	Stoma Size	Flange Size	Pouch Length	Box of
3232	5/8" (15 mm)	1" (25 mm)	9" (23 cm)	30
3238	7/8" (22 mm)	1-1/4" (32 mm)	9" (23 cm)	
3233	1-1/8" (29mm)	1-1/2" (38 mm)	9" (23 cm)	
3239	1-3/8" (35 mm)	1-3/4" (44 mm)	9" (23 cm)	
3234	1-1/2" (38 mm)	2" (51 mm)	9" (23 cm)	
3222	5/8" (15 mm)	1" (25 mm)	12" (30 cm)	30
3228	7/8" (22 mm)	1-1/4" (32 mm)	12" (30 cm)	
3223	1-1/8" (29mm)	1-1/2" (38 mm)	12" (30 cm)	
3229	1-3/8" (35 mm)	1-3/4" (44 mm)	12" (30 cm)	
3224	1-1/2" (38 mm)	2" (51 mm)	12" (30 cm)	
3225	2" (51 mm)	2-1/4" (64 mm)	12" (30 cm)	30
3226	2-1/2" (64 mm)	3" (76 mm)	12" (30 cm)	
3278	7/8" (22 mm)	1-1/4" (32 mm)	16" (40 cm)	
3273	1-1/8" (29mm)	1-1/2" (38 mm)	16" (40 cm)	
3279	1-3/8" (35 mm)	1-3/4" (44 mm)	16" (40 cm)	
3274	1-1/2" (38 mm)	2" (51 mm)	16" (40 cm)	30
3275	2" (51 mm)	2-1/4" (64 mm)	16" (40 cm)	
3276	2-1/2" (64 mm)	3" (76 mm)	16" (40 cm)	
3608	7/8" (22 mm)	1-1/4" (32 mm)	12" (30 cm)	10
3603	1-1/8" (29mm)	1-1/2" (38 mm)	12" (30 cm)	
3609	1-3/8" (35 mm)	1-3/4" (44 mm)	12" (30 cm)	
3604	1-1/2" (38 mm)	2" (51 mm)	12" (30 cm)	
3605	2" (51 mm)	2-1/4" (64 mm)	12" (30 cm)	

Closed Pouches



Colour: Transparent With Integrated Filter				With Tape Border 
Product #	Stoma Size	Flange Size	Pouch Length	Box of
3328	7/8" (22 mm)	1-1/4" (32 mm)	9" (23 cm)	30
3323	1-1/8" (29mm)	1-1/2" (38 mm)	9" (23 cm)	
3329	1-3/8" (35 mm)	1-3/4" (44 mm)	9" (23 cm)	
3324	1-1/2" (38 mm)	2" (51 mm)	9" (23 cm)	
3325	2" (51 mm)	2-1/4" (64 mm)	9" (23 cm)	
3326	2-1/2" (64 mm)	3" (76 mm)	9" (23 cm)	

Did you know?

- To select the appropriate *Karaya* pouch, determine the stoma size and select the matching Product Code
- The labeling on the *Karaya* box and pouch you have selected indicates the gasket size, NOT the stoma size
- The *Karaya* seal extends into the gasket opening, ensuring that the Product Code you have selected fits the stoma size you have measured

A photograph of a family of three outdoors. On the left, a man with short dark hair and a mustache, wearing a blue t-shirt, is looking towards the center. In the middle, a young child with curly brown hair, wearing a red polo shirt and a plaid jacket, is smiling broadly. On the right, a woman with long brown hair, wearing a white top, is looking at the child. The background is a bright, out-of-focus outdoor setting with greenery. The word "Pouchkins" is overlaid in red text in the upper left corner.

Pouchkins

With the Gentleness and Flexibility of SoftFlex


Pouchkins

Pediatric Pouching System

One-Piece SoftFlex

Premie/Newborn Pouches



Cut-to-fit				Without Tape Border	
Product #	Description	For Stomas Up To	Pouch Length	Box of	
Premie					
3777	Ultra-clear, no starter hole	5/8" (15 mm)	4-1/2" (11.5 cm)	15	
Newborn					
3778	Ultra-clear, no starter hole	1-3/8" (35 x 22 x 7/8" mm)	6-1/2" (16 cm)		

Drainable Pouches Lock 'n Roll Closure System



Cut-to-fit Without Tape Border				
Product #	Description	For Stomas Up To	Pouch Length	Box of
3795	Transparent, no starter hole	2" (51 mm)	7" (18 cm)	10
3796		1-1/2" (38 mm)	6-1/2" (16 cm)	

Urostomy Pouch



Cut-to-fit Without Tape Border				
Product #	Description	For Stomas Up To	Pouch Length	Box of
3797	Transparent, no starter hole	1-1/2" (38 mm)	6-1/2" (16 cm)	10

Pouchkins SoftFlex Skin Barrier:

- Designed to be gentle on sensitive infant skin
- No starter hole, which enables the opening to be placed where needed
- No tape by design
- Flexible and sized for babies and children
- None of the *Pouchkins* pediatric pouching systems are made with natural rubber latex
- Any 1-3/4" *New Image* pouch will fit this skin barrier

Two-Piece SoftFlex

Flat Skin Barrier



Cut-to-fit Without Tape Border				
Product #	Description	Flange Size	For Stomas Up To	Box of
3761	No starter hole	1-3/4" (44 mm)	1-1/4" (32 mm)	5

Drainable Pouches Lock 'n Roll Closure System



Colour: Transparent

Product #	Flange Size	Pouch Length	Filter	Box of
3799	1-3/4" (44 mm)	6-1/2" (16 cm)	No	10

Colour: Beige

18282	1-3/4" (44 mm)	7" (18 cm)	AF300	20
--------------	----------------	------------	-------	----

Urostomy Pouches



Enhanced Design

Product #	Colour	Flange Size	Pouch Length	Box of
18922	Ultra-clear	1-3/4" (44 mm)	9" (23 cm)	10
18912	Beige	1-3/4" (44 mm)	9" (23 cm)	10

Ostomy Belt



Product #	Length	Box of
3774	10 - 17" (25 - 43 cm)	1

Shadow Buddies dolls

- Exclusive to **Hollister**, *Shadow Buddies* dolls are special friends for children who have had or will have ostomy surgery
- Great teaching tool to help with education and care giving
- Huggable dolls with heart-shaped eyes and a stoma provide comfort and support



Shadow Buddies is a trademark of the Shadow Buddies Foundation


Accessories

Adapt CeraRing

Sealed with the confidence of ceramide

The security + confidence of Adapt meets the science of ceramide

- Available in Flat and Convex for each individual need
- Helps prevent leakage and protect skin—and may even extend barrier wear time
- Resists erosion—even when wet or exposed to discharge from the stoma
- Offers ultimate comfort and ease of use—alcohol-free rings are designed to be gentle for all skin types
- Extended wear
- Infused with ceramide — a natural component of the skin
- Help protect the skin's natural moisture barrier
- Help decrease transepidermal water loss (TEWL) from damaged or eroded skin



Adapt Barrier Rings

Flat CeraRing



Product #	Type	Outer Diameter	Thickness	Box of
8805	Standard	2" (48 mm)	4.5 mm	10
8815	Slim	2" (48 mm)	2.30 mm	

Convex CeraRing



Product #	Inner Diameter	Stretch Up To	Box of
89520	13/16" (20 mm)	1" (25 mm)	10
89530	1-3/16" (30 mm)	1-3/8" (35 mm)	
89540	1-9/16" (40mm)	1-3/4" (45 mm)	

Oval Convex CeraRing



Product #	Inner Diameter	Stretch Up To	Box of
89601	7/8" x 1-1/2" (22 x 38 mm)	1-1/8" x 1-3/4" (27 x 43 mm)	10
89602	1-3/16" x 1-7/8" (30 x 48 mm)	1-3/8" x 2-1/8" (35 x 53 mm)	
89603	1-1/2" x 2-3/16" (38 x 56 mm)	1-3/4" x 2-3/8" (43 x 61 mm)	

Flat Flextend



Product #	Type	Outer Diameter	Thickness	Box of
7805	Standard	2" (48 mm)	4.5 mm	10
7806	Standard	4" (102 mm)	2,30 mm	
7815	Slim	2" (48 mm)	2.30 mm	

Convex Flextend



Product #	Inner Diameter	Stretch Up To	Box of
79520	13/16" (20 mm)	1" (25 mm)	10
79530	1-3/16" (30 mm)	1-3/8" (35 mm)	
79540	1-9/16" (40 mm)	1-3/4" (45 mm)	









Oval Convex Flextend




Product #	Inner Diameter	Stretch Up To	Box of
79601	7/8" x 1-1/2" (22 x 38 mm)	1-1/8" x 1-3/4" (27 x 43 mm)	10
79602	1-3/16" x 1-7/8" (30 x 48 mm)	1-3/8" x 2-1/8" (35 x 53 mm)	
79603	1-1/2" x 2-3/16" (38 x 56 mm)	1-3/4" x 2-3/8" (43 x 61 mm)	

Accessories

Adapt


Product #		Description	Usage	Size		Box of
7760		Universal Remover Wipes <i>*featuring no-sting silicone-based, alcohol free formulation</i>	This skin friendly product helps remove adhesive residues.	50 Wipes per Carton		50
7737		Adapt Medical Adhesive Remover Spray <i>*featuring no-sting, silicone-based, alcohol free formulation</i>	360° spray technology to deliver controlled, quiet spraying from any angle, with minimal waste	Spray Can	1.7 fl oz (50 ml)	1
7906		Adapt Stoma Powder	Use to absorb moisture from broken skin around the stoma. Dust on. Brush off excess powder. Discontinue use once the condition is resolved.	Puff bottle	1 oz (28 g)	1
79402		Adapt Barrier Extenders	Helps provide extra security around the edges of skin barrier, reduces skin barrier edge lifting and may help increase skin barrier wear time. Provides extra adhesion to accommodate a peristomal hernia.	Each barrier extender covers a surface area of 48 cm ²		30
78500		Adapt Lubricating Deodorant	Neutralizes odours, eases emptying, and helps prevent pouch static and sticking. It is not a perfume or masking agent. Should be used in drainable and closed pouching systems only.	Bottle	8 oz (236 ml)	1
78501				Packets	.27 oz (8 ml)	50
79301		Adapt Paste	Used to fill or caulk uneven skin contours to create a flatter surface. It can help to prevent drainage from getting under the ostomy skin barrier.	Tube	.5 oz (14 g)	20
79300					2 oz (60 g)	1
79400		Adapt Barrier Strips with FlexTend Barrier (Alcohol Free)	Used to fill in or caulk uneven skin contours near the stoma to create a flatter surface and help prevent ostomy drainage from getting under the ostomy barrier.	Triangular Shape	5-1/8" (13 cm)	10
7917		Adapt Skin Protective Wipe <i>*featuring no-sting silicone-based, alcohol free formulation</i>	Protects skin during adhesive removal Not recommended for use with extended wear skin barriers	50 Wipes per Carton		50

Karaya




Product #		Description	Usage	Size		Box of
7905		Powder	Use to absorb moisture from broken skin around the stoma. Dust on. Brush off excess powder. Discontinue use once the condition is resolved.	Puff Bottle	2.5 oz (71 g)	12

Accessories


m9 Odour Eliminators

Product #	Description	Size	Box of
7715	 Drops	1 oz (30 ml)	12
7717		8 oz (240 ml)	6
7732	 Spray Unscented	2 oz (60 ml)	12
7733		8 oz (240 ml)	6
7735	Spray Apple Scent	8 oz (240 ml)	

Pouch Accessories

Product #	Description	Size	Box of
8770	 Drainable Pouch Clamp, Beige		20
7300	 Beige Colour Ostomy Belt	Medium 23" - 43" (58 - 109 cm)	10
7299		Large 34" - 65" (86 - 165 cm)	
7331	 Urostomy Drain Tube Adapter		10

Access Port


Product #	Description	Box of
9779	 Graduated nipple and punch allows introduction of catheter or drain through Wound Drainage Collector or pouch system	8

Wound Drainage Collectors


With FlexWear Skin Barrier

Product #	For Wounds Up To	Box of
9773	3" (7.6 cm)	10
9776	3.75" (9.5 cm)	10
9778	4" x 8" (10 cm x 20 cm)	3

Skin Sheets

Product #	Description	Size	Box of
8800	 Flexextend (Non-sterile)	4" x 4" (10 cm x 10 cm)	5
8801		8" x 8" (20 cm x 20 cm)	3
7800	Premium (Non-sterile)	4" x 4" (10 cm x 10 cm)	5
7801		8" x 8" (20 cm x 20 cm)	3
7700	HolliHesive (Non-sterile)	4" x 4" (10 cm x 10 cm)	5

Irrigation Products

Product #	Description	Information		Box of
7721	Cone Irrigator Kit	Includes: (1) Irrigator (1) Stoma Cone with Connector (1) Cleaning Brush		1
7723	Stoma Cone	With Connector Replacement Unit		10
7728	Irrigator Sleeve	With Attached Closure and Belt Tabs	2" (51 mm)	20
7724			3" (76 mm)	
18152	 New Image Two-Piece Irrigator Sleeve	Flange: 1-3/4" (44 mm)	Pouch: 30" (76 cm)	5
18153		Flange: 2-1/4" (57 mm)	Pouch: 30" (76 cm)	
18154		Flange: 2-3/4" (70 mm)	Pouch: 30" (76 cm)	
7740	Stoma Lubricant	Bottle	4 oz (118 ml)	10

Stoma Caps

with Microporous Tape

Pre-sized					With Tape Border
Product #	Colour	For Stomas Up To	Pouch Length	Filter	Box of
3184	White	2" (51 mm)	4-1/4" (11 cm)	Yes	30
3186		3" (76 mm)	4-1/4" (11 cm)		

SoftFlex Skin Barrier

Pre-sized					Without Tape Border
Product #	Colour	For Stomas Up To	Pouch Length	Filter	Box of
1796	Beige	1-5/16" (50 mm)	4" (10 cm)	Yes	30



Hollister Secure Start Services

Hollister Secure Start Services provides a lifetime of customized support for people in the ostomy and continence care communities.

We are proud to offer dedicated support for each and every patient along the continuum of care! Once enrolled, a Hollister Secure Start Services Coordinator will call to explain our services, answer any questions, and provide ongoing support.

A Hollister Secure Start Services Coordinator can assist with the following:

- Finding the right products
- Providing product information and condition-specific education
- Finding local resources like retailers and clinicians
- Product usage care tips
- Finding reimbursement for your catheter products



To learn more about Hollister Secure Start Services, for Ostomy, call us at **1.866.789.7574**, or email us at **securestartcanada@hollister.com**.

For Continence, call us at **1.800.263.7421**, or email us at **canadasscccontact@hollister.com**.

Or visit us at **www.hollister.com**.



Continence Care



People are at the center of everything we do.

Discretion by design

Introducing a discreet catheter with design that is truly inspired.

The Infyna Chic™ Intermittent Catheter was designed with input from clinicians to help provide a high level of discretion for women who use catheters. It is beautiful and easy to use. The Infyna Chic catheter can help a woman feel better about having to use one.

Infyna Chic catheter

Product #	Description	Size	Length	Box of
7008	Straight tip	8 Fr	5" (13 cm)	30
7010	Straight tip	10 Fr	5" (13 cm)	30
7012	Straight tip	12 Fr	5" (13 cm)	30
7014	Straight tip	14 Fr	5" (13 cm)	30



Features that help change her intermittent catheter experience

- Flip-top cap can be opened and closed with one hand
- Ready to use, hydrophilic catheter
- Just right stiffness helps support easy, touch-free insertion
- Catheter length helps give confidence that the bladder is fully drained
- Can attach a collection bag without removing the catheter from the case



Easy to open with just one hand



Easy to grip and remove from case



Does not leak when reclosed

LEARN MORE





Start Right
Here.

VaPro Pocket™ catheters offer the features that patients want.¹

Protection. Ease of Use. Discretion.

VaPro Pocket catheters provide 100% No Touch Protection with a protective tip and sleeve, and they are easy to use, with an improved protective tip that offers a more comfortable insertion. The pocket-size package design is more discreet—a quality users prefer.

Consider VaPro Pocket catheters right from the start.

Learn more by visiting [Hollister.com](https://www.hollister.com) or contact us at 1.800.263.7400.



Reference: 1. Hollister data on file.

Prior to use, be sure to read the Instructions for Use for information regarding Intended Use, Contraindications, Warnings, Precautions, and Instructions.

The Hollister logo, VaPro Plus Pocket, and VaPro Pocket are trademarks of Hollister Incorporated.
©2021 Hollister Incorporated. ALL-01469

CE
0050

VaPro Pocket
No Touch Intermittent Catheter

VaPro Plus Pocket
No Touch Intermittent Catheter



Continence Care



Touch Free Hydrophilic Intermittent Catheters

VaPro Catheters Provide 100% No Touch Protection with a Protective Tip and Sleeve:

VaPro catheters are hygienic by design to help protect patients from germs throughout the entire catheterization process.

Protective tip:

- helps shield the sterile catheter during insertion from bacteria located within the first 15 mm of the distal urethra
- helps reduce the risk of carrying bacteria into the urinary tract

Protective sleeve:

- allows for catheter to be gripped anywhere
- provides a barrier to the catheter to help keep germs away



VaPro™



Product #	Description	Size	Length	Box of
72124	Straight tip	12 Fr	16" (40 cm)	30
73124	Coudé tip		16" (40 cm)	
72142	Straight tip	14 Fr	8" (20 cm)	30
72144	Straight tip		16" (40 cm)	
73144	Coudé tip	16 Fr	16" (40 cm)	30
72164	Straight tip		16" (40 cm)	
73164	Coudé tip		16" (40 cm)	

VaPro Plus™



Product #	Description	Size	Length	Box of
74124	Straight tip	12 Fr	16" (40 cm)	30
74144	Straight tip	14 Fr	16" (40 cm)	30

VaPro Plus is a closed system with an integrated 1000 ml collection bag.

VaPro F-Style is a **FIRMER** Touch Free Hydrophilic Catheter.

VaPro™ F-Style



Product #	Description	Size	Length	Box of
7600124	Straight tip	12 Fr	16" (40 cm)	30
7600144	Straight tip	14 Fr	16" (40 cm)	30
7600164	Straight tip	16 Fr	16" (40 cm)	30

VaPro Plus™ F-Style



Product #	Description	Size	Length	Box of
7700124	Straight tip	12 Fr	16" (40 cm)	30
7700144	Straight tip	14 Fr	16" (40 cm)	30

VaPro Plus F-Style is a closed system with an integrated 1000 ml collection bag.

Improved! A 50%* smaller protective tip,
the same 100% No Touch Protection.



Original

New



Improved design features:

- **50%* smaller** and more ergonomically shaped protective tip, for a more comfortable insertion.
- **Transparent protective tip** to aid visualization
- **More discreet** packaging design
- **Less waste** as a result of new design

*As measured by volume of material used.



Easy-to open, easy-to-empty
collection bag for cleaner
handling

(applies to the VaPro Plus Pocket
catheters only)

VaPro Pocket




Product #	Description	Size	Length	Box of
70124	Straight tip	12 Fr	16" (40 cm)	30
70144	Straight tip	14 Fr	16" (40 cm)	30
70164	Straight tip	16 Fr	16" (40 cm)	30

VaPro Plus Pocket

Product #	Description	Size	Length	Box of
71084	Straight tip	8 Fr	16" (40 cm)	30
71104	Straight tip	10 Fr	16" (40 cm)	30
71122	Straight tip	12 Fr	8" (20 cm)	30
71124	Straight tip		16" (40 cm)	
71142	Straight tip	14 Fr	8" (20 cm)	30
71144	Straight tip		16" (40 cm)	
71164	Straight tip	16 Fr	16" (40 cm)	30

*VaPro Plus Pocket is a closed system with an integrated
1000 ml collection bag.*



New Study¹ Confirms Superior Protection*

Hollister asked an independent laboratory to conduct a study to see if the **superior protection** of VaPro™ No Touch Intermittent Catheters could be measured. The study compared VaPro catheters to three other brands.

*In comparison to other tested catheters.

Reference: 1. Hollister data on file, TR-00453, 2021.

Find Out What
Happened



Continence Care

Touch Free Intermittent Catheters

Closed System



Product #	Description	Size	Length	Box of
94064IC	Straight tip	6 Fr	16" (40 cm)	100
94084IC	Straight tip	8 Fr	16" (40 cm)	
94104IC	Straight tip	10 Fr	16" (40 cm)	
94124IC	Straight tip	12 Fr	16" (40 cm)	
95124IC	Coudé tip		16" (40 cm)	
94144IC	Straight tip	14 Fr	16" (40 cm)	
95144IC	Coudé tip		16" (40 cm)	
94164IC	Straight tip	16 Fr	16" (40 cm)	
95164IC	Coudé tip		16" (40 cm)	
94184IC	Straight tip	18 Fr	16" (40 cm)	

Closed System Kit



Product #	Description	Size	Length	Box of
96064IC	Straight tip	6 Fr	16" (40 cm)	100
96084IC	Straight tip	8 Fr	16" (40 cm)	
96104IC	Straight tip	10 Fr	16" (40 cm)	
96124IC	Straight tip	12 Fr	16" (40 cm)	
97124	Coudé tip		16" (40 cm)	
96144IC	Straight tip	14 Fr	16" (40 cm)	
97144	Coudé tip		16" (40 cm)	
96164IC	Straight tip	16 Fr	16" (40 cm)	
97164	Coudé tip		16" (40 cm)	
96184IC	Straight tip	18 Fr	16" (40 cm)	

Kit Includes:

- Advance Plus Closed System
- Gloves
- Benzalkonium Chloride antiseptic swab (BZK)
- Waterproof underpad

Advance Plus Touch Free Intermittent Catheters:

- Have a protective tip that helps shield the catheter from bacteria located in the first 15 mm of the distal urethra
- Have a user-regulated gel reservoir which allows for control of the catheter lubrication
- Smooth eyelets help provide trouble-free insertion and withdrawal

Intermittent Catheters



Product #	Former Product #	Description	Size	Length	Box of
10610	New	Straight tip	6 Fr	10" (25 cm)	30
10616	B1007	Straight tip	6 Fr	15" (38 cm)	
10806	B1016	Straight tip	8 Fr	6" (15 cm)	30
10810	New	Straight tip		10" (25 cm)	
10816	B1018	Straight tip		16" (40 cm)	
10826	New	Coudé	8 Fr	16" (40 cm)	
11006	B1033	Straight tip	10 Fr	6" (15 cm)	30
11010	New	Straight tip		10" (25 cm)	
11026	1035	Coudé, Olive tip		16" (40 cm)	
11016	B1034	Straight tip		16" (40 cm)	
11206	B1047	Straight tip	12 Fr	6" (15 cm)	30
11216	B1048	Straight tip		16" (40 cm)	
11226	B1050	Coudé tip		16" (40 cm)	
11406	B1062	Straight tip	14 Fr	6" (15 cm)	30
11446	B1063	Straight tip, No funnel		6" (15 cm)	
11416	B1064	Straight tip		16" (40 cm)	
B1061	B1061	Soft tip		16" (40 cm)	
11426	B1066	Coudé tip		16" (40 cm)	
B1065	B1065	Straight tip with curved packaging		16" (40 cm)	50
11616	B1082	Straight tip	16 Fr	16" (40 cm)	30
11626	B1086	Coudé tip		16" (40 cm)	
11816	B1096	Straight tip	18 Fr	16" (40 cm)	30
11826	1097	Coudé	18 Fr	16" (40 cm)	

Features of the Apogee Family of Intermittent Catheters:

- Two staggered smooth eyelets help provide trouble-free insertion
- Not made with natural rubber latex
- Designed for comfort
- Medical grade Polyvinyl Chloride (PVC)
- Non lubricated, clear catheter
- Straight and coudé tip options available

Apogee products were changed to a new number in 2021. We have kept the former product number for cross reference purposes only. This change does not affect the product, its performance or intended use in any way. If you have any questions regarding this change please contact your Hollister representative or Hollister's Customer Service Department at 1-800-263-7400.

Male External Catheters

Hollister

Latex



Extended Wear with Unique Inner Flap			Sheath length 2-3/4" (7 cm), adhesive length 1-3/5" (4 cm)
Product #	Description	Size	Box of
9206	Small/Pediatric	7/8" - 1" (22 - 25 mm)	30
9207	Medium	1" - 1-1/8" (26 - 30 mm)	
9208	Intermediate	1-1/8" - 1-1/4" (31 - 35 mm)	
9209	Large	1-3/8" - 1-1/2" (36 - 39 mm)	
Everyday Wear			Sheath length 2-3/4" (7 cm), adhesive length 1-3/5" (4 cm)
Product #	Description	Size	Box of
9106	Small/Pediatric	7/8" - 1" (22 - 25 mm)	30
9107	Medium	1" - 1-1/8" (26 - 30 mm)	
9108	Intermediate	1-1/8" - 1-1/4" (31 - 35 mm)	
9109	Large	1-3/8" - 1-1/2" (36 - 39 mm)	

Contains Natural Rubber Latex

Size Selection Guide

Use the Hollister Sizing Guide to determine the most appropriate size catheter

- Place sizing guide just behind the glans of the penis
- Select the size that most closely fits



InView

Silicone



Standard			Sheath length 2-3/4" (7 cm), adhesive length 1-3/5" (4 cm)
Product #	Description	Size	Box of
97525	Small	1" (25 mm)	30
97529	Medium	1-1/8" (29 mm)	
97532	Intermediate	1-1/4" (32 mm)	
97536	Large	1-1/2" (36 mm)	
97541	Extra Large	1-5/8" (41 mm)	100
97525-100	Small	1" (25 mm)	
97529-100	Medium	1-1/8" (29 mm)	
97532-100	Intermediate	1-1/4" (32 mm)	
97536-100	Large	1-1/2" (36 mm)	
97541-100	Extra Large	1-5/8" (41 mm)	

Special			Sheath length 1-3/5" (4 cm), adhesive length 1-3/5" (4 cm)
Product #	Description	Size	Box of
97425	Small	1" (25 mm)	30
97429	Medium	1-1/8" (29 mm)	
97432	Intermediate	1-1/4" (32 mm)	
97436	Large	1-1/2" (36 mm)	
97441	Extra Large	1-5/8" (41 mm)	100
97425-100	Small	1" (25 mm)	
97429-100	Medium	1-1/8" (29 mm)	
97432-100	Intermediate	1-1/4" (32 mm)	
97436-100	Large	1-1/2" (36 mm)	
97441-100	Extra Large	1-5/8" (41 mm)	

Extra			Sheath length 2-3/4" (7 cm), adhesive length 2-3/4" (7 cm)
Product #	Description	Size	Box of
97625	Small	1" (25 mm)	30
97629	Medium	1-1/8" (29 mm)	
97632	Intermediate	1-1/4" (32 mm)	
97636	Large	1-1/2" (36 mm)	
97641	Extra Large	1-5/8" (41 mm)	100
97625-100	Small	1" (25 mm)	
97629-100	Medium	1-1/8" (29 mm)	
97632-100	Intermediate	1-1/4" (32 mm)	
97636-100	Large	1-1/2" (36 mm)	
97641-100	Extra Large	1-5/8" (41 mm)	

Not Made with Natural Rubber Latex

Urinary Collection Systems

Leg Bags

Leg Bag



Product #	Description	Capacity	Box of
9814	Medium	18 oz (540 ml)	10
9805	Large	30 oz (900 ml)	

Combination Pack



Product #	Description	Capacity	Box of
9348	Medium	18 oz (540 ml)	12 packs
9349	Large	30 oz (900 ml)	

Extension Tubing with Connector



Product #	Description	Length	Box of
9345	Non-sterile	18" (45 cm)	10
9346	Sterile	18" (45 cm)	

Vented Combination Pack



Product #	Description	Capacity	Box of
9645	Medium	18 oz (540 ml)	10 packs
9655	Large	30 oz (900 ml)	

Fabric Strap



Product #	Description	Capacity	Box of
9342	Small	15" (38 cm)	10 pairs
9343	Medium	23" (58 cm)	

Bedside Drainage



Product #	Description	Capacity	Box of
9839	2000 ml Collection System	68 oz (2000 ml)	20

Hollister Leg Bags Are:

- Sterile
- Not Made with Natural Rubber Latex
- Designed with an anti-reflux valve

m9 Cleaner/Decrystallizer



Concentrate & Wash Bottle

Product #	Description	Size	Box of
7736	1 Bottle of Concentrate 1 Washing Bottle	16 oz (480 ml)	12

External Collection Devices

Urinary Pouches

Female

SoftFlex Skin Barrier



Product #	Description	Box of
9840	Length 7-1/2" (19 cm) Width 4-1/2" (11cm)	10

Male

Flextend Skin Barrier

Thin barrier for flexibility, no additional tape needed



Product #	Description	Box of
9873	Length 7-1/2" (19 cm) Width 4-1/2" (11cm)	10

SoftFlex Skin Barrier



Product #	Description	Box of
9811	Length 7-1/2" (19 cm) Width 4-1/2" (11cm)	10

Hollister's Urinary Pouches:

- Not Made with Natural Rubber Latex
- Have a cut-to-fit skin barrier
- Have a odour-barrier pouch film

Fecal Collectors

Pouches

Flextend Skin Barrier



Product #	Description	Pouch Length	Box of
9823	Medium, Ultra-clear	10" (25 cm)	10
9880	Large, Transparent	12" (30 cm)	

Transparent Pouches

SoftFlex Skin Barrier



Product #	Description	Pouch Length	Box of
9822	Medium	10" (25 cm)	10
9821	Large	12" (30 cm)	



Now
Online

Our NEW Clinical Education Platform is Now Open



Clinical Research and Support



Continued Learning



Patient Education



For Clinicians.
By Clinicians.



Hollister®

Continence Care

A man and a woman are seen from behind, sitting closely together on a dark wooden bench. The man is wearing a light-colored, textured sweater, and the woman has curly blonde hair and is wearing a grey top. They are looking out over a scenic view that includes a body of water and a lush, green forested hillside under a bright, slightly hazy sky. The overall mood is peaceful and contemplative.

Critical Care

Dedicated to excellence.

Oral Endotracheal Tube Fasteners

AnchorFast



Product #	Description	Box of
9799	<ul style="list-style-type: none"> ET tube wrap with one-click security clamp locks the ET tube in place Gliding tube shuttle allows the ET tube to be easily repositioned from side to side while the tube remains secure FlexTend Skin Barrier For tubes sized 5.0-10.0mm inner diameter 	12

AnchorFast SlimFit



Product #	Description	Box of
9787	<ul style="list-style-type: none"> ET tube wrap with one-click security Gliding tube shuttle allows for side-to-side repositioning of the ET tube Thumb pad with increased surface area for finger or thumb placement when opening Designed to accommodate smaller faces Secures ET tubes ranging in size 5.0-8.0 mm 	12

AnchorFast Guard Select



Product #	Description	Box of
9700	<ul style="list-style-type: none"> ET tube wrap with one-click security clamp locks the ET tube in place Gliding tube shuttle allows the ET tube to be easily repositioned from side to side while the tube remains secure Integrated tube protection sleeve designed to help prevent occlusion Secures oral ET tubes ranging in size from 5.0 to 9.0 mm inner diameter and ET tubes with integrated subglottic suction lumen sizes 6.0 to 8.0 mm inner diameter 	12

AnchorFast Replacement Neck Band Straps



Product #	Description	Box of
9788	Replacement Neck Band Straps	4

Tube Attachment Devices

Horizontal Tube (HTAD)



Product #	Description	Box of
9781	Single use, sterile	5

Vertical Tube (VTAD)



Product #	Description	Box of
9782	Single use, sterile	5

Drain Tube (DTAD)



Product #	Description	Box of
9780	Single use, sterile	5

Feeding Tube (FTAD)



Product #	Description	Box of
9786	<ul style="list-style-type: none"> FlexWear skin barrier conforms to patient's nose Single use, sterile 	20

Suction Tube (STAD)



Product #	Description	Box of
9785	<ul style="list-style-type: none"> FlexWear skin barrier conforms to patient's nose Single use, sterile 	20

Numerical Index

Product #	Page	Product #	Page	Product #	Page	Product #	Page	Product #	Page	Product #	Page
1035	38	7014	32	8290	19	8496	23	8635	16	9106	39
1097	38	7299	29	8291	19	8497	23	8637	16	9107	39
1796	29	7300	29	8292	19	8498	23	8638	16	9108	39
3184	29	7331	29	8293	19	8504	19	8639	16	9109	39
3186	29	7700	29	8294	19	8508	19	8641	16	9206	39
3222	24	7715	29	8295	19	8510	19	8642	16	9207	39
3223	24	7717	29	8296	19	8511	19	8643	16	9208	39
3224	24	7721	29	8297	19	8512	19	8644	16	9209	39
3225	24	7723	29	8331	16	8513	19	8647	16	9342	40
3226	24	7724	29	8412	22	8514	19	8648	16	9343	40
3228	24	7728	29	8413	22	8515	19	8649	16	9345	40
3229	24	7732	29	8414	22	8524	19	8651	16	9346	40
3232	24	7733	29	8415	22	8528	19	8661	18	9348	40
3233	24	7735	29	8440	21	8531	16	8662	18	9349	40
3234	24	7736	40	8442	21	8538	16	8663	18	9645	40
3238	24	7737	28	8443	21	8551	16	8664	18	9655	40
3239	24	7740	29	8447	21	8552	16	8665	18	9700	44
3273	24	7760	28	8448	21	8558	16	8666	18	9773	29
3274	24	7800	29	8449	21	8559	16	8674	18	9776	29
3275	24	7801	29	8450	21	8584	19	8678	18	9778	29
3276	24	7805	27	8460	21	8588	19	8770	29	9779	29
3278	24	7806	27	8462	21	8590	19	8800	29	9780	44
3279	24	7815	27	8463	21	8591	19	8801	29	9781	44
3323	24	7905	28	8464	21	8592	19	8805	27	9782	44
3324	24	7906	28	8467	21	8593	19	8815	27	9785	44
3325	24	7917	28	8468	21	8594	19	8901	16	9786	44
3326	24	8031	16	8469	21	8595	19	8914	18	9787	44
3328	24	8081	16	8474	23	8596	19	8918	18	9788	44
3329	24	8131	16	8478	23	8597	19	8925	16	9799	44
3603	24	8181	16	8480	23	8610	19	8930	16	9805	40
3604	24	8214	19	8481	23	8611	19	8931	16	9811	41
3605	24	8218	19	8482	23	8612	19	8935	16	9814	40
3608	24	8244	19	8483	23	8613	19	8940	16	9821	41
3609	24	8248	19	8484	23	8614	19	8954	18	9822	41
3761	26	8274	19	8485	23	8615	19	8958	18	9823	41
3774	26	8278	19	8486	23	8616	19	8960	18	9839	40
3777	26	8280	19	8487	23	8617	19	8961	18	9840	41
3778	26	8281	19	8488	23	8618	19	8962	18	9873	41
3795	26	8282	19	8489	23	8619	19	8963	18	9880	41
3796	26	8283	19	8490	23	8624	19	8984	18	10610	38
3797	26	8284	19	8491	23	8628	19	8988	18	10616	38
3799	26	8285	19	8492	23	8631	16	8990	18	10806	38
7008	32	8286	19	8493	23	8632	16	8991	18	10810	38
7010	32	8288	19	8494	23	8633	16	8992	18	10816	38
7012	32	8289	19	8495	23	8634	16	8993	18	10826	38

Product #	Page	Product #	Page	Product #	Page	Product #	Page	Product #	Page	Product #	Page
11006	38	13907	11	14902	11	16402	9	18204	12	19004	15
11010	38	13908	11	14903	11	16403	9	18282	12	19006	15
11016	38	13909	11	14904	11	16404	9	18282	26	19052	15
11026	38	13910	11	14905	11	16405	9	18283	12	19053	15
11202	8	13911	11	14906	11	16406	9	18284	12	19054	15
11203	8	14102	9	14907	11	16407	9	18292	12	19056	15
11204	8	14103	9	14908	11	16408	9	18293	12	19102	15
11206	38	14104	9	14909	11	18002	12	18294	12	19103	15
11216	38	14202	9	14910	11	18003	12	18322	13	19104	15
11226	38	14203	9	14911	11	18004	12	18323	13	19106	15
11402	10	14204	9	15102	8	18006	12	18324	13	19152	15
11403	10	14206	9	15103	8	18012	13	18332	13	19153	15
11404	10	14301	9	15104	8	18013	13	18333	13	19154	15
11406	38	14302	9	15202	9	18014	13	18334	13	19156	15
11416	38	14303	9	15203	9	18016	13	18352	13	19202	15
11426	38	14304	9	15204	9	18102	12	18353	13	19203	15
11446	38	14305	9	15302	10	18103	12	18354	13	19204	15
11503	10	14306	9	15303	10	18104	12	18362	13	19252	15
11504	10	14307	9	15304	10	18112	12	18363	13	19253	15
11505	10	14308	9	15402	11	18113	12	18364	13	19254	15
11506	10	14309	9	15403	11	18114	12	18372-30	13	19603	15
11616	38	14402	11	15404	11	18122	12	18373-30	13	19604	15
11626	38	14403	11	15503	11	18123	12	18374-30	13	19903	15
11702	10	14404	11	15504	11	18124	12	18392	13	70124	35
11703	10	14503	11	15505	11	18132	12	18393	13	70144	35
11704	10	14504	11	15506	11	18133	12	18394	13	70164	35
11816	38	14505	11	15602	8	18134	12	18402	14	71084	35
11826	38	14506	11	15603	8	18142	12	18403	14	71104	35
11903	10	14602	8	15604	8	18143	12	18404	14	71122	35
11904	10	14603	8	15802	11	18144	12	18412	14	71124	35
11905	10	14604	8	15803	11	18152	29	18413	14	71142	35
11906	10	14606	8	15804	11	18153	29	18414	14	71144	35
13402	11	14701	8	15903	11	18154	29	18902	14	71164	35
13403	11	14702	8	15904	11	18172	12	18903	14	72124	34
13404	11	14703	8	15905	11	18173	12	18904	14	72142	34
13503	11	14704	8	15906	11	18174	12	18912	14	72144	34
13504	11	14705	8	16101	8	18176	12	18912	26	72164	34
13505	11	14706	8	16102	8	18182	12	18913	14	73124	34
13506	11	14707	8	16103	8	18183	12	18914	14	73144	34
13901	11	14708	8	16104	8	18184	12	18922	14	73164	34
13902	11	14709	8	16105	8	18192	12	18922	26	74124	34
13903	11	14802	11	16106	8	18193	12	18923	14	74144	34
13904	11	14803	11	16107	8	18194	12	18924	14	78500	28
13905	11	14804	11	16108	8	18202	12	19002	15	78501	28
13906	11	14901	11	16401	9	18203	12	19003	15	79300	28

Numerical Index


Product #	Page	Product #	Page	Product #	Page	Product #	Page	Product #	Page
79301	28	84394	23	85118	19	96144IC	37	151114	8
79400	28	84398	23	85119	19	96164IC	37	151115	8
79402	28	84404	22	85211	19	96184IC	37	151116	8
79520	27	84408	22	85811	19	97124	37	151117	8
79530	27	84490	21	86211	19	97144	37	151118	8
79540	27	84492	21	86711	18	97164	37	153103	10
79601	27	84493	21	88300	17	97425	39	153104	10
79602	27	84494	21	88302	17	97425-100	39	153105	10
79603	27	84497	21	88325	17	97429	39	153106	10
80070	17	84498	21	88330	17	97429-100	39	826211	19
80110	17	84499	21	88335	17	97432	39	841311	22
82100	17	84590	21	88340	17	97432-100	39	844011	22
82111	19	84692	21	88402	17	97436	39	7600124	34
82125	17	84693	21	88700	17	97436-100	39	7600144	34
82130	17	84694	21	88725	17	97441	39	7600164	34
82135	17	84697	21	88730	17	97441-100	39	7700124	34
82300	17	84698	21	88735	17	97525	39	7700144	34
82302	17	84699	21	88740	17	97525-100	39	8195111	18
82325	17	84711	23	89011	16	97529	39	8439111	23
82330	17	84794	23	89111	18	97529-100	39	8479111	23
82335	17	84798	23	89511	18	97532	39	B1007	38
82402	17	84890	23	89518	18	97532-100	39	B1016	38
82411	19	84891	23	89520	27	97536	39	B1018	38
82610	19	84892	23	89530	27	97536-100	39	B1033	38
82611	19	84893	23	89540	27	97541	39	B1034	38
82612	19	84894	23	89601	27	97541-100	39	B1047	38
82613	19	84895	23	89602	27	97625	39	B1048	38
82614	19	84896	23	89603	27	97625-100	39	B1050	38
82615	19	84897	23	89811	18	97629	39	B1061	38
82616	19	84898	23	94064IC	37	97629-100	39	B1061	38
82617	19	84899	23	94084IC	37	97632	39	B1062	38
82618	19	84990	23	94104IC	37	97632-100	39	B1063	38
82619	19	84991	23	94124IC	37	97636	39	B1064	38
82624	19	84992	23	94144IC	37	97636-100	39	B1065	38
82628	19	84993	23	94164IC	37	97641	39	B1065	38
82711	19	84994	23	94184IC	37	97641-100	39	B1066	38
84134	22	84995	23	95124IC	37	112013	8	B1082	38
84138	22	84996	23	95144IC	37	112014	8	B1086	38
84200	21	84997	23	95164IC	37	112015	8	B1096	38
84292	22	84998	23	96064IC	37	112016	8		
84293	22	84999	23	96084IC	37	112017	8		
84294	22	85011	19	96104IC	37	112018	8		
84295	22	85116	19	96124IC	37	151113	8		

www.hollister.com

The Hollister logo, Adapt, Advance, Advance Plus, Anchorfast, Anchorfast Slim Fit, Anchorfast Guard Select, Apogee, CeraPlus, Extended Wear, Everyday, Flexextend, FlexWear, FormaFlex, Infyna Chic, InView, Karaya 5, Lock 'n Roll, m9, New Image, Pouchkins, Premier, QuietWear, Secure Start, SoftFlex, VaPro, VaPro F-Style, VaPro Plus, VaPro Pocket and VaPro Plus Pocket are trademarks of Hollister Incorporated.

All other trademarks and copyrights are the property of their respective owners.

©2021 Hollister Incorporated



Contact Information

Customer Service

Toll Free: 1.800.263.7400

Fax

Toll Free: 1.800.432.8846

Email

canadacustomerservice@hollister.com

Business Hours

Monday – Friday; 8:30 am - 5:00 pm (EST) Excluding Holidays

Hollister Limited

Aurora, ON

Hollister Secure Start Services

Hollister Secure Start Services provides a lifetime of customized support for people in the ostomy and continence care communities.

We are proud to offer dedicated support for each and every patient along the continuum of care. Once enrolled, a Hollister Secure Start Services Coordinator will call to explain our services, answer any questions, and provide ongoing support.

A Hollister Secure Start Services Coordinator can assist with the following:

- Finding the right products
- Providing product information and condition-specific education
- Finding local resources like retailers and clinicians
- Product usage care tips
- Find reimbursement for your catheter products

Ostomy Care:

1.866.789.7574, or email us at
securestartcanada@hollister.com

Continence Care:

1.800.263.7421, or email us at
canadassccccontact@hollister.com



ENROLL HERE

Our NEW Clinical Education Platform is Now Open

Visit our NEW Clinical
Education Platform now.

Scan the QR code or visit
www.hollister.com/ostomyclinicaleducation

 SCAN ME



©2021 Hollister Incorporated

Ostomy Care
Healthy skin. Positive outcomes.