

# NEW!



## VaPro Plus

Touch Free Hydrophilic Intermittent Catheter

### 100% No Touch Protection

#### Protective Sleeve

promotes easy handling while also helping to keep bacteria away

#### Protective Tip


helps bypass bacteria in the first 15 mm of the distal urethra

#### Integrated Collection Bag

no receptacle needed

#### Ring Cap

helps keep protective tip clean and protected after package is opened



#### Simple insertion and removal

- Ready to use with pure catheter hydration technology
- Designed to be evenly lubricated for easier insertion
- Helps provide ease of insertion and withdrawal with smooth eyelets

#### Touch Free

- The catheter can be gripped anywhere along the protective sleeve
- The protective tip helps shield the sterile catheter, during insertion and removal, from bacteria located within the first 15 mm of the distal urethra

Description	Length	Product Code	Quantity per box
12 Fr Straight	16" (40 cm)	74124	30
14 Fr Straight	16" (40 cm)	74144	30

**Continence Care**  
People First



**Hollister**

95 Mary Street  
Aurora, ON L4G 1G3  
1-800-263-7400  
[www.hollister.com](http://www.hollister.com)