

Maintaining the
Right Environment
Makes all the
Difference



The CeraPlus Skin Barrier
*with Remois Technology**

Now Available on Premier One-Piece Pouching Systems!



Peristomal Skin Deserves Better

The importance of peristomal skin health in ostomy care can't be overstated. The condition of this skin greatly affects how well the skin barrier can adhere, as well as the general comfort and well-being of a patient living with a stoma.

Unfortunately, peristomal skin problems are too common. About half of those living with a stoma will report peristomal skin issues at some point in their lives.

Many even think that skin problems are a normal part of living with a stoma.

We disagree. Patients do not need to suffer, and our goal is to help maintain healthy peristomal skin from the start.

The CeraPlus Skin Barrier

Now Available on Premier One-Piece Pouching Systems!

The CeraPlus skin barrier is infused with ceramide, a naturally occurring lipid in the skin, which protects against dryness, and features a formulation to help:

- Protect the skin's natural moisture barrier
- Maintain good peristomal skin health
- Decrease transepidermal water loss (TEWL) from damaged or eroded skin

And now the CeraPlus skin barrier is available on Premier one-piece pouching systems. Your patients will continue to enjoy exclusive features such as the AF300 filter, and Lock 'n Roll microseal closure.

The CeraPlus skin barrier continues to be available as a two-piece product with tape border, to be paired with the New Image pouching system.


To continue to build on our two-piece CeraPlus portfolio, we are also introducing new two-piece products without tape border!



The Ceramide Advantage

Ceramide is a natural component of human skin, consisting of waxy lipid molecules. Ceramide links cells of the outer skin (epidermis) together to form a waterproof, protective barrier.

Ceramides help prevent water loss that can lead to skin damage and dryness. Many skin moisturizing products contain ceramide as a component to help protect the skin against loss of its natural moisture.



Ceramide:

- Is a natural component of the human skin
- Links the cells of the stratum corneum together to form a waterproof, protective barrier
- Helps prevent water loss that can lead to skin damage and dryness

*84% of responding clinicians reported they were likely to continue to use and recommend the CeraPlus skin barrier for their patients.***

Our new CeraPlus skin barrier is designed to support healthy skin from the start—so your patients maintain good skin health long after they leave your care.

Experience the CeraPlus skin barrier difference.

Contact your Hollister representative for details and a trial sample.

** Results from a 2014 product evaluation (N=111) of the CeraPlus skin barrier in Canada and The Netherlands.

The CeraPlus Skin Barrier

with Remois Technology*

The CeraPlus skin barrier is now available on **Premier** one-piece pouching systems. **Premier** pouching systems integrate exclusive features such as the AF300 filter, and Lock 'n Roll Microseal closure.

The CeraPlus skin barrier continues to be available as a two-piece product with tape border, to be paired with the New Image pouching system. To continue to build on the two-piece CeraPlus portfolio, we are also introducing new two-piece products without tape border!

Now available on both
Premier One-Piece
pouching systems and
New Image Two-Piece
skin barriers.

NEW! Premier CeraPlus One-Piece Pouching Systems

- Drainable pouch with Lock 'n Roll Microseal closure
- Tape border

Flat Barrier

Pouch Colour	Barrier Opening		Product No.
Ultra-Clear	Up to 2 1/2" (64 mm)		8931
	Up to 2 1/2" (64 mm)		8901
Beige with Viewing Option and AF300 filter	1" (25 mm)		8925
	1 3/16" (30 mm)		8930
	1 3/8" (35 mm)		8935
	1 9/16" (40 mm)		8940

Convex Barrier

Pouch Colour	Barrier Opening		Product No.
Ultra-Clear	Up to 1" (25 mm)		8984
	Up to 1 1/2" (38 mm)		8988
	Up to 2" (51 mm)		89811
Beige with Viewing Option & AF300 filter	Up to 1" (25 mm)		8914
	Up to 1 1/2" (28 mm)		8918
	Up to 2" (51 mm)		89111
	3/4" (19 mm)		8990
	7/8" (22 mm)		8991
	1" (25 mm)		8992
1 1/8" (29 mm)		8993	

Soft Convex Barrier

Pouch Colour	Barrier Opening		Product No.
Beige with Viewing Option & AF300 filter	Up to 1" (25 mm)		8954
	Up to 1 1/2" (38 mm)		8958
	Up to 2 1/8" (55 mm)		89511
	3/4" (19 mm)		8960
	7/8" (22 mm)		8961
	1" (25 mm)		8962
1 1/8" (29 mm)		8963	

New Image CeraPlus Two-Piece Skin Barriers

- Tape Border

Flat Barrier

Flange Size	Barrier Opening		Product No.
1 3/4" (44 mm)	Up to 1 1/4" (32 mm)		11202
2 1/4" (57 mm)	Up to 1 3/4" (44 mm)		11203
2 3/4" (70 mm)	Up to 2 1/4" (57 mm)		11204

Convex Barrier

Flange Size	Barrier Opening		Product No.
1 3/4" (44 mm)	Up to 1" (25 mm)		11402
2 1/4" (57 mm)	Up to 1 1/2" (38 mm)		11403
2 3/4" (70 mm)	Up to 2" (51 mm)		11404
1 3/4" (44 mm)	7/8" (22 mm)		11503
1 3/4" (44 mm)	1" (25 mm)		11504
2 1/4" (57 mm)	1 1/8" (29 mm)		11505
2 1/4" (57 mm)	1 1/4" (32 mm)		11506

NEW! New Image CeraPlus Two-Piece Skin Barriers

- Without tape border

Flat Barrier

Flange Size	Barrier Opening		Product No.
1 3/4" (44 mm)	Up to 1 1/4" (32 mm)		15102
2 1/4" (57 mm)	Up to 1 3/4" (44 mm)		15103
2 3/4" (70 mm)	Up to 2 1/4" (57 mm)		15104

Convex Barrier

Flange Size	Barrier Opening		Product No.
1 3/4" (44 mm)	Up to 1" (25 mm)		15302
2 1/4" (57 mm)	Up to 1 1/2" (38 mm)		15303
2 3/4" (70 mm)	Up to 2" (51 mm)		15304

indicates Cut-To-Fit Barrier