A photograph of a young Black couple sitting on a light-colored sofa. The woman, on the left, has curly hair and is wearing a light grey t-shirt and patterned shorts. She is laughing and has her hands raised in the air. The man, on the right, has a beard and is wearing a white t-shirt and dark pants. He is also laughing, holding a tablet in his left hand, and has his right hand raised. The background is a bright, minimalist interior with white walls and framed artwork.

Making a Difference in the Journey of Life

2021 Product Catalogue

A 100-Year Legacy...

The fundamental values and principles of our founder, John Dickinson Schneider, form the foundation of our company today.

Our Mission

To make life more rewarding and dignified for the people who use our products and services.

Our Vision

To grow and prosper as an independent, employee-owned company and in the process become better human beings.

Our Heritage

The roots of Hollister Incorporated date back to 1921, when 23-year-old entrepreneur John Dickinson Schneider founded a small printing company in Chicago, which he named JDS Printer Craftsman. From the very beginning, John Schneider's business philosophy was "only first class is good enough," and the company became known for its quality products, quality service, and quality people.



Table of Contents

| | | | |
|--|----|---|----|
| General Information | 2 | Accessories | |
| Ostomy Care | | Adapt™ Barrier Rings | 26 |
| Choosing the Right Skin Barrier FIT | 4 | Adapt™ Accessories | 27 |
| Choosing the Right Skin Barrier FORMULATION | 5 | m9 Odour Eliminators | 28 |
| New Image™ | | Pouch Accessories | 28 |
| Two-Piece Pouching System | | Karaya Accessories | 28 |
| Flat Skin Barriers | 7 | Skin Sheets | 28 |
| Soft Convex Skin Barriers | 9 | Loop Ostomy Bridge | 28 |
| Convex Skin Barriers | 9 | Irrigation Products | 28 |
| Drainable Pouches | 11 | Wound Drainage Collectors | 29 |
| Closed Pouches | 12 | Access Port | 29 |
| High Output Pouches | 12 | Stoma Caps | 29 |
| Urostomy Pouches | 13 | with Microporous Tape | 29 |
| Ostomy Kits | 14 | Continence Care | |
| Urostomy Kits | 14 | Infyna Chic™ | 32 |
| Premier™ | | VaPro™ | |
| One-Piece Pouching System | | Touch Free Hydrophilic Intermittent Catheters | 34 |
| Flat Skin Barriers - Drainable Pouch | 17 | VaPro | |
| High Output Pouches | 18 | Improved VaPro Pocket and VaPro Plus Pocket | 35 |
| Flat Skin Barriers - Closed Pouch | 18 | Look for the new VaPro Pocket and VaPro Plus Pocket catheters | 35 |
| Soft Convex Skin Barriers - Drainable Pouch | 19 | Advance Plus™ | |
| Convex Skin Barriers - Drainable Pouch | 19 | Touch Free Intermittent Catheters | 36 |
| Flat Skin Barriers - Urostomy Pouch | 21 | Apogee™ | |
| Soft Convex Skin Barriers - Urostomy Pouch | 22 | Apogee Family of Intermittent Catheters | 38 |
| Convex Skin Barriers - Urostomy Pouch | 22 | Male External Catheters | |
| Karaya™ | | Hollister Extended and Everyday Wear | 39 |
| One-Piece Pouching System—Natural Skin Barriers | | InView Silicone | 39 |
| Drainable Pouches | 24 | Urinary Collection Systems | |
| Closed Pouches | 24 | Leg Bags | 40 |
| Urostomy Pouches | 24 | Bedside Drainage | 40 |
| Pouchkins™ | | m9 Cleaner/Decrystallizer | 40 |
| Pediatric Pouching System | | External Collection Devices | |
| One-Piece | 25 | Urinary Pouches | 41 |
| Two-Piece | 25 | Fecal Collectors | 41 |
| | | Critical Care | |
| | | Oral Endotracheal Tube Fasteners | 44 |
| | | Tube Attachment Devices | 44 |
| | | Numerical Index | 45 |

General Information

Terms of Payment

Net 30 days from date of invoice. Electronic Fund Transfers (ETF) are possible. Please contact the customer order department.

Minimum Order

The minimum order requirement for a dealer is \$250.00. A charge of \$10.00 will be applied to any order below the minimum requirement.

Shipping

Prices for shipments to destinations in Canada are FOB buyer's address.

If the method and route of shipment are specified by the buyer, the buyer is required to pay the delivery charges, upon delivery of products to the carrier, all risk of loss and other incidents of ownership pass to the buyer. These terms may not be altered.

Pricing & Terms of Sale

GST/HST will be charged where applicable. Buyer agrees that they will not at anytime, directly or indirectly, sell or distribute, to customers which are located outside of Canada.

Return Goods Policy

All returns for credit require an authorization number from Hollister Limited. A Returns Goods Authorization (RGA) can be obtained by contacting a representative at Customer Service 1.800.263.7400 or by email canadacustomerservice@hollister.com. Written request for authorization is required.

Hollister will authorize returns under the following conditions:

1. A 15% restock fee will be applied to all returned goods over 30 days from the date of invoice.
2. Copy of invoice or purchase order number must be provided for the return request.
3. **Product returned by customers without an authorized RGA number will be disposed of by Hollister Limited and no credit will be issued.**

Hollister will NOT authorize the return of the following:

1. Goods over one (1) year from the date of invoice.
2. Partial boxes of merchandise.
3. Boxes in unsaleable condition (i.e markings, sticker for sale, sun fading, dented boxes due to improper handling/return packaging, sticker glue, etc.)

Damaged Product:

1. Customer must call Hollister Limited to receive an RGA number within 5 business days of receipt.
2. Receiving dock is responsible for documenting, on the bill of lading, when a shipment is received in damaged condition.
3. Customer must provide copy of invoice, invoice number or purchase order number for the return request.

4. Credit will **NOT** be issued if customer disposes of damaged/ defective product without approval by Hollister Limited.
5. If product cannot be returned due to extreme damage (leakage/broken) customer must provide pictures to Hollister Limited within 5 business days of receipt to gain approval to dispose and thereby obtain credit for the same.

Ostomy & Continence Retail Sample Policy

All sample requests at the Retail level are to be processed through Hollister Customer Service by calling 1.800.263.7400 or by email canadacustomerservice@hollister.com.

Retailer/client is to contact Customer Service and provide consumer information for sample processing. Samples will be shipped directly to consumers from Hollister Limited. Samples will **no longer** be available to ship directly to a retail setting.

Shortages

Customer must notify Hollister Limited within 5 business days of receipt on shortages. A shortage is a difference between what was physically received and the master pack list.

Customer Satisfaction Return Policy

All Customer Satisfaction Returns require authorization by your Hollister Territory Manager or a Customer Service Representative.

Unauthorized Customer Satisfaction Returns will not be processed.

Hollister will authorize Customer Satisfaction Returns under the following conditions:

1. The customer has been given a no-charge, full box replacement of Hollister product for their partial box return, and
2. A Customer Satisfaction Return form has been completed, and
3. The returned box has been given to, and the completed Customer Satisfaction Return form has been authorized by your Hollister Territory Manager or a Customer Service Representative.

Prior to Use

Prior to use, be sure to read the Instructions for Use for information regarding intended use, contraindications, warnings, precautions and instructions.

Latex Statement

The specifications for the Hollister Continence Care products listed below include natural rubber latex as a component; our suppliers add natural rubber latex as part of the production process. If a Hollister product(s) found in this catalog is not listed below, the specification for that product(s) and packaging does not include natural rubber latex as a component, and our component suppliers do not add natural rubber as part of their production process. Product # 9206, 9207, 9208, 9209, 9106, 9107, 9108, 9109

Ostomy Care



Healthy skin. Positive outcomes.

Healthy Skin. Positive Outcomes.

Choosing the Right Skin Barrier FIT

FIT is critical to get a good seal around the stoma to prevent leakage and promote comfort.

Barrier Examples



Flat Skin Barrier



Flat Skin Barrier
with Round or Oval
Convex Barrier Ring



Integrated Soft
Convex Barrier



Integrated Convex
Barrier



Integrated Convex
Barrier with Round
or Oval Convex
Barrier Ring

Possible Clinical Case Scenarios



Flat Skin Barrier

- Flat peristomal plane
- Adequate stoma protrusion
- Centered stoma opening



Flexible Convexity

- Stoma flush with skin
- Uneven peristomal plane
- Off-centered stoma opening



Soft Convexity

- Gentle pressure required to
peristomal plane
- Flush stoma
- Peristomal injury



Firm Convexity

- Creases or folds in peristomal plane
- Off-centered stoma opening
- Stoma retracted below
skin surface



Firm Convexity with
More Depth

- Recessed or flushed stoma
- Uneven abdominal contours
- Off-centered stoma opening

Prior to use, be sure to read the Instructions for Use for information regarding intended use, contraindications, warnings, precautions and instructions.

1. Colwell JC, Goldberg MT & Carmel JE (Eds). (2015). WOCN Society Core Curriculum Ostomy Management. Chapter 10. Philadelphia, PA: Wolters Kluwer
2. Hoeflok J et al. Use of Convexity in Ostomy Care: Results of an International Consensus Meeting. JWOCN, January/February 2017

Healthy Skin. Positive Outcomes.

Choosing the Right Skin Barrier FORMULATION

FORMULATION is the customized ingredients responsible for adhesion, erosion resistance and absorption that must be carefully formulated to create an environment where healthy skin can thrive.

Extended Wear Skin Barriers

CeraPlus

- Contains skin barrier characteristics (adhesion, absorption and erosion resistance) and is infused with ceramide—a naturally occurring lipid in the skin which protects against dryness
- Formulation designed to help protect the skin's natural moisture barrier and help maintain good peristomal skin health

Flexextend

- Durable with resistance to erosion
- Swells when it comes in contact with liquid to help form a seal between the skin and the stoma
- Suitable for those who perspire heavily
- Strong adhesion and flexibility

FormaFlex

- Maintains its shape after being stretched or shaped
- Provides an easy-to-customize opening due to barrier material
- Easy to shape to individual stomas without the use of scissors

Standard Wear Skin Barriers

FlexWear

- Good balance of durability and gentleness providing moderate adhesion, and can handle frequent or infrequent pouch changes
- Resistant to fluid erosion
- Suitable for all skin types, depending on change frequency
- Designed to be worn more than one day

SoftFlex

- Appropriate choice for:
 - Those who have thin, fragile, or sensitive skin (such as those who have undergone chemo/radiation therapy)
 - People who desire frequent removal
- Designed to be gentle to the skin, which allows for frequent pouch changes

Karaya 5

- Well-suited for sensitive skin or in some cases where allergic reactions to synthetic barriers has occurred
- Mildly adhesive—needs a belt to secure pouching system

CeraPlus™ skin barrier
with Remois Technology*

vs. SenSura Mio

The Results are Clear: CeraPlus Skin Barrier Helps Ensure Healthy Skin

In the Measuring Epidermal Effects of Ostomy Skin Barriers study (Grove et al, 2019), the CeraPlus skin barrier was significantly more effective at maintaining skin health than the SenSura Mio barrier.



**CeraPlus
Skin Barrier**

**SenSura
Mio Barrier**

| | | |
|---|---|---|
| Infused with ceramide, the skin's natural protection against water loss | ✓ | ✗ |
| Observed to be better at minimizing skin irritation | ✓ | ✗ |
| Observed to be better at managing Transepidermal water loss (TEWL) | ✓ | ✗ |
| Observed to be better at protecting skin against damage from repeated barrier removal | ✓ | ✗ |
| Observed to be better at maintaining skin health | ✓ | ✗ |



**CeraPlus Skin Barrier Minimizes Potential Impact
of Peristomal Skin Complications Right from the Start!**

For more information, please contact your Hollister Representative

LEARN MORE



Grove G, Houser T, Sibbald G, Salvadana G. Measuring epidermal effects of ostomy skin barriers. *Skin Res Technol*. 2019;25:179-186. <https://doi.org/10.1111/srt.12630>

The Hollister logo, CeraPlus and "Healthy skin. Positive Outcomes." are trademarks of Hollister Incorporated.

All other trademarks and copyrights are the property of their respective owners.

©2019 Hollister Incorporated

CE Prior to use, be sure to read the Instructions for Use for information regarding Intended Use, Contraindications, Warnings, Precautions, and Instructions. Not all products are CE marked.



*Remois is a technology of Alcare Co., Ltd.

New Image

Two-Piece Pouching System

Flat Skin Barriers

CeraPlus with Remois Technology*



Cut-to-fit

| Product # | Flange Size | For Stomas Up To | Box of | |
|--------------|----------------|------------------|--------|---------------------|
| 11202 | 1-3/4" (44 mm) | 1-1/4" (32 mm) | 5 | With Tape Border |
| 11203 | 2-1/4" (57 mm) | 1-3/4" (44 mm) | | |
| 11204 | 2-3/4" (70 mm) | 2-1/4" (57 mm) | | |
| 15102 | 1-3/4" (44 mm) | 1-1/4" (32 mm) | 5 | Without Tape Border |
| 15103 | 2-1/4" (57 mm) | 1-3/4" (44 mm) | | |
| 15104 | 2-3/4" (70 mm) | 2-1/4" (57 mm) | | |

*Remois is a technology of Alcare Co., Ltd.

Coming Soon

New Flat Pre-Sized CeraPlus

with Remois Technology*

*Remois is a technology of Alcare Co., Ltd.

Flextend



Cut-to-fit

| Product # | Flange Size | For Stomas Up To | Box of | |
|--------------|----------------|------------------|--------|---------------------|
| 14602 | 1-3/4" (44 mm) | 1-1/4" (32 mm) | 5 | With Tape Border |
| 14603 | 2-1/4" (57 mm) | 1-3/4" (44 mm) | | |
| 14604 | 2-3/4" (70 mm) | 2-1/4" (57 mm) | | |
| 14606 | 4" (102 mm) | 3-1/2" (89 mm) | 5 | Without Tape Border |
| 15602 | 1-3/4" (44 mm) | 1-1/4" (32 mm) | | |
| 15603 | 2-1/4" (57 mm) | 1-3/4" (44 mm) | | |
| 15604 | 2-3/4" (70 mm) | 2-1/4" (57 mm) | | |


Pre-sized

| Product # | Flange Size | Barrier Opening | Box of | |
|--------------|----------------|-----------------|--------|---------------------|
| 14701 | 1-3/4" (44 mm) | 5/8" (16 mm) | 5 | With Tape Border |
| 14702 | 1-3/4" (44 mm) | 3/4" (19 mm) | | |
| 14703 | 1-3/4" (44 mm) | 7/8" (22 mm) | | |
| 14704 | 1-3/4" (44 mm) | 1" (25 mm) | | |
| 14705 | 1-3/4" (44 mm) | 1-1/8" (29 mm) | | |
| 14706 | 1-3/4" (44 mm) | 1-1/4" (32 mm) | | |
| 14707 | 2-1/4" (57 mm) | 1-3/8" (35 mm) | | |
| 14708 | 2-1/4" (57 mm) | 1-1/2" (38 mm) | | |
| 14709 | 2-1/4" (57 mm) | 1-3/4" (44 mm) | | |
| 16101 | 1-3/4" (44 mm) | 5/8" (16 mm) | 5 | Without Tape Border |
| 16102 | 1-3/4" (44 mm) | 3/4" (19 mm) | | |
| 16103 | 1-3/4" (44 mm) | 7/8" (22 mm) | | |
| 16104 | 1-3/4" (44 mm) | 1" (25 mm) | | |
| 16105 | 1-3/4" (44 mm) | 1-1/8" (29 mm) | | |
| 16106 | 1-3/4" (44 mm) | 1-1/4" (32 mm) | | |
| 16107 | 2-1/4" (57 mm) | 1-3/8" (35 mm) | | |
| 16108 | 2-1/4" (57 mm) | 1-1/2" (38 mm) | | |

Designed for Diverse Bodies

Adhesive (Tape) Border Barriers are an appropriate option for most patients, especially those with:

- Skin irregularities such as hernias, wrinkles, or scars
- Anatomical challenges due to obesity or cachexia
- Added need for a sense of security, or extra assurances to accommodate an active lifestyle
- Increased desire for comfort and flexibility



New Image

Two-Piece Pouching System

Flat Skin Barriers

FormaFlex



| Product # | Flange Size | For Stomas Up To | Box of | |
|--------------|----------------|------------------|--------|------------------|
| 14102 | 1-3/4" (44 mm) | 1-1/4" (32 mm) | 5 | With Tape Border |
| 14103 | 2-1/4" (57 mm) | 1-11/16" (43 mm) | | |
| 14104 | 2-3/4" (70 mm) | 2-1/4" (57 mm) | | |

Colour Matching System

Two-Piece Pouching Systems

To help you select compatible *New Image* pouching system components, every box label is colour-coded.

1. Select the skin barrier flange size with the appropriate barrier stoma opening and note its colour code.
2. Select a suitable pouch with the desired features in the matching colour.

| | |
|--|----------------|
| | 1-3/4" (44 mm) |
| | 2-1/4" (57 mm) |
| | 2-3/4" (70 mm) |
| | 4" (102 mm) |

FlexWear



Cut-to-fit

| Product # | Flange Size | For Stomas Up To | Box of | |
|--------------|----------------|------------------|--------|---------------------|
| 14202 | 1-3/4" (44 mm) | 1-1/4" (32 mm) | 5 | With Tape Border |
| 14203 | 2-1/4" (57 mm) | 1-3/4" (44 mm) | | |
| 14204 | 2-3/4" (70 mm) | 2-1/4" (57 mm) | | |
| 14206 | 4" (102 mm) | 3-1/2" (89 mm) | | |
| 15202 | 1-3/4" (44 mm) | 1-1/4" (32 mm) | 5 | Without Tape Border |
| 15203 | 2-1/4" (57 mm) | 1-3/4" (44 mm) | | |
| 15204 | 2-3/4" (70 mm) | 2-1/4" (57 mm) | | |

Pre-sized

| Product # | Flange Size | Barrier Opening | Box of | |
|--------------|----------------|-----------------|--------|---------------------|
| 14301 | 1-3/4" (44 mm) | 5/8" (16 mm) | 5 | With Tape Border |
| 14302 | 1-3/4" (44 mm) | 3/4" (19 mm) | | |
| 14303 | 1-3/4" (44 mm) | 7/8" (22 mm) | | |
| 14304 | 1-3/4" (44 mm) | 1" (25 mm) | | |
| 14305 | 1-3/4" (44 mm) | 1-1/8" (29 mm) | | |
| 14306 | 1-3/4" (44 mm) | 1-1/4" (32 mm) | | |
| 14307 | 2-1/4" (57 mm) | 1-3/8" (35 mm) | | |
| 14308 | 2-1/4" (57 mm) | 1-1/2" (38 mm) | | |
| 14309 | 2-1/4" (57 mm) | 1-3/4" (44 mm) | | |
| 16401 | 1-3/4" (44 mm) | 5/8" (16 mm) | 5 | Without Tape Border |
| 16402 | 1-3/4" (44 mm) | 3/4" (19 mm) | | |
| 16403 | 1-3/4" (44 mm) | 7/8" (22 mm) | | |
| 16404 | 1-3/4" (44 mm) | 1" (25 mm) | | |
| 16405 | 1-3/4" (44 mm) | 1-1/8" (29 mm) | | |
| 16406 | 1-3/4" (44 mm) | 1-1/4" (32 mm) | | |
| 16407 | 2-1/4" (57 mm) | 1-3/8" (35 mm) | | |
| 16408 | 2-1/4" (57 mm) | 1-1/2" (38 mm) | | |

New Image

Two-Piece Pouching System

Soft Convex Skin Barriers **New!**

CeraPlus with Remois Technology*



Cut-to-fit

| Product # | Flange Size | For Stomas Up To | Box of |
|--------------|----------------|------------------|-----------------------|
| 11702 | 1-3/4" (44 mm) | 1" (25 mm) | 5 With Tape Border |
| 11703 | 2-1/4" (57 mm) | 1-1/2" (38 mm) | |
| 11704 | 2-3/4" (70 mm) | 2" (51 mm) | |

Pre-sized

| Product # | Flange Size | Barrier Opening | Box of |
|--------------|----------------|-----------------|-----------------------|
| 11903 | 1-3/4" (44 mm) | 7/8" (22 mm) | 5 With Tape Border |
| 11904 | 1-3/4" (44 mm) | 1" (25 mm) | |
| 11905 | 2-1/4" (57 mm) | 1-1/8" (29 mm) | |
| 11906 | 2-1/4" (57 mm) | 1-1/4" (32 mm) | |

Soft Convexity



Flexible design may help provide gentle pressure around the stoma to help obtain the right fit.

Convex Skin Barriers

CeraPlus with Remois Technology*



Cut-to-fit

| Product # | Flange Size | For Stomas Up To | Box of |
|--------------|----------------|------------------|--------------------------|
| 11402 | 1-3/4" (44 mm) | 1" (25 mm) | 5 With Tape Border |
| 11403 | 2-1/4" (57 mm) | 1-1/2" (38 mm) | |
| 11404 | 2-3/4" (70 mm) | 2" (51 mm) | |
| 15302 | 1-3/4" (44 mm) | 1" (25 mm) | 5 Without Tape Border |
| 15303 | 2-1/4" (57 mm) | 1-1/2" (38 mm) | |
| 15304 | 2-3/4" (70 mm) | 2" (51 mm) | |

Pre-sized

| Product # | Flange Size | Barrier Opening | Box of |
|--------------|----------------|-----------------|-----------------------|
| 11503 | 1-3/4" (44 mm) | 7/8" (22 mm) | 5 With Tape Border |
| 11504 | 1-3/4" (44 mm) | 1" (25 mm) | |
| 11505 | 2-1/4" (57 mm) | 1-1/8" (29 mm) | |
| 11506 | 2-1/4" (57 mm) | 1-1/4" (32 mm) | |

*Remois is a technology of Alcare Co., Ltd.

The CeraPlus Skin Barrier delivers both Fit and Formulation

Ceramide Infused Barrier

- **Protects** the skin's natural moisture barrier
- **Decreases** TEWL (transepidermal water loss) from damaged or eroded skin
- **Maintains** peristomal skin health



New Image

Two-Piece Pouching System

Convex Skin Barriers

Flextend



Cut-to-fit

| Product # | Flange Size | For Stomas Up To | Box of | |
|--------------|----------------|------------------|--------|---------------------|
| 14802 | 1-3/4" (44 mm) | 1" (25 mm) | 5 | With Tape Border |
| 14803 | 2-1/4" (57 mm) | 1-1/2" (38 mm) | | |
| 14804 | 2-3/4" (70 mm) | 2" (51 mm) | | |
| 15802 | 1-3/4" (44 mm) | 1" (25 mm) | 5 | Without Tape Border |
| 15803 | 2-1/4" (57 mm) | 1-1/2" (38 mm) | | |
| 15804 | 2-3/4" (70 mm) | 2" (51 mm) | | |

Pre-sized

| Product # | Flange Size | Barrier Opening | Box of | |
|--------------|----------------|-----------------|--------|----------------------------|
| 13901 | 1-3/4" (44 mm) | 5/8" (16 mm) | 5 | With Tape Border (Tapered) |
| 13902 | 1-3/4" (44 mm) | 3/4" (19 mm) | | |
| 13903 | 1-3/4" (44 mm) | 7/8" (22 mm) | | |
| 13904 | 1-3/4" (44 mm) | 1" (25 mm) | | |
| 13905 | 2-1/4" (57 mm) | 1-1/8" (29 mm) | | |
| 13906 | 2-1/4" (57 mm) | 1-1/4" (32 mm) | | |
| 13907 | 2-1/4" (57 mm) | 1-3/8" (35 mm) | | |
| 13908 | 2-1/4" (57 mm) | 1-1/2" (38 mm) | | |
| 13909 | 2-3/4" (70 mm) | 1-5/8" (41 mm) | | |
| 13910 | 2-3/4" (70 mm) | 1-3/4" (44 mm) | | |
| 13911 | 2-3/4" (70 mm) | 2" (51 mm) | | |
| 14901 | 1-3/4" (44 mm) | 5/8" (16 mm) | 5 | With Tape Border |
| 14902 | 1-3/4" (44 mm) | 3/4" (19 mm) | | |
| 14903 | 1-3/4" (44 mm) | 7/8" (22 mm) | | |
| 14904 | 1-3/4" (44 mm) | 1" (25 mm) | | |
| 14905 | 2-1/4" (57 mm) | 1-1/8" (29 mm) | | |
| 14906 | 2-1/4" (57 mm) | 1-1/4" (32 mm) | | |
| 14907 | 2-1/4" (57 mm) | 1-3/8" (35 mm) | | |
| 14908 | 2-1/4" (57 mm) | 1-1/2" (38 mm) | | |
| 14909 | 2-3/4" (70 mm) | 1-5/8" (41 mm) | | |
| 14910 | 2-3/4" (70 mm) | 1-3/4" (44 mm) | | |
| 14911 | 2-3/4" (70 mm) | 2" (51 mm) | | |
| 15903 | 1-3/4" (44 mm) | 7/8" (22 mm) | 5 | Without Tape Border |
| 15904 | 1-3/4" (44 mm) | 1" (25 mm) | | |
| 15905 | 2-1/4" (57 mm) | 1-1/8" (29 mm) | | |
| 15906 | 2-1/4" (57 mm) | 1-1/4" (32 mm) | | |

SCAN ME



FlexWear



Cut-to-fit

| Product # | Flange Size | For Stomas Up To | Box of | |
|--------------|----------------|------------------|--------|----------------------------|
| 13402 | 1-3/4" (44 mm) | 1" (25 mm) | 5 | With Tape Border (Tapered) |
| 13403 | 2-1/4" (57 mm) | 1-1/2" (38 mm) | | |
| 13404 | 2-3/4" (70 mm) | 2" (51 mm) | | |
| 14402 | 1-3/4" (44 mm) | 1" (25 mm) | 5 | With Tape Border |
| 14403 | 2-1/4" (57 mm) | 1-1/2" (38 mm) | | |
| 14404 | 2-3/4" (70 mm) | 2" (51 mm) | | |
| 15402 | 1-3/4" (44 mm) | 1" (25 mm) | 5 | Without Tape Border |
| 15403 | 2-1/4" (57 mm) | 1-1/2" (38 mm) | | |
| 15404 | 2-3/4" (70 mm) | 2" (51 mm) | | |

Pre-sized

| Product # | Flange Size | Barrier Opening | Box of | |
|--------------|----------------|-----------------|--------|----------------------------|
| 13503 | 1-3/4" (44 mm) | 7/8" (22 mm) | 5 | With Tape Border (Tapered) |
| 13504 | 1-3/4" (44 mm) | 1" (25 mm) | | |
| 13505 | 2-1/4" (57 mm) | 1-1/8" (29 mm) | | |
| 13506 | 2-1/4" (57 mm) | 1-1/4" (32 mm) | | |
| 14503 | 1-3/4" (44 mm) | 7/8" (22 mm) | 5 | With Tape Border |
| 14504 | 1-3/4" (44 mm) | 1" (25 mm) | | |
| 14505 | 2-1/4" (57 mm) | 1-1/8" (29 mm) | | |
| 14506 | 2-1/4" (57 mm) | 1-1/4" (32 mm) | | |
| 15503 | 1-3/4" (44 mm) | 7/8" (22 mm) | 5 | Without Tape Border |
| 15504 | 1-3/4" (44 mm) | 1" (25 mm) | | |
| 15505 | 2-1/4" (57 mm) | 1-1/8" (29 mm) | | |
| 15506 | 2-1/4" (57 mm) | 1-1/4" (32 mm) | | |

Choosing the Right Skin Barrier FIT

- FIT is critical to get a good seal around the stoma to prevent leakage and promote comfort.
- Convex barrier products can be utilized to help provide a secure and predictable seal.
- Product selection is made by matching stomal protrusion and peristomal contours to the appropriate barriers.
- The outward curving of the convex skin barrier is designed to interface with the immediate peristomal skin.

Refer to page 4 for additional information on The Continuum of Convex Barrier Options.

New Image

Two-Piece Pouching System

Drainable Pouches

Lock 'n Roll Closure System



Colour: Ultra-clear

| Product # | Flange Size | Pouch Length | Filter | Box of |
|--------------|----------------|--------------|--------|--------|
| 18002 | 1-3/4" (44 mm) | 12" (30 cm) | No | 10 |
| 18003 | 2-1/4" (57 mm) | 12" (30 cm) | | |
| 18004 | 2-3/4" (70 mm) | 12" (30 cm) | | |
| 18006 | 4" (102 mm) | 12" (30 cm) | | |

Colour: Transparent

| Product # | Flange Size | Pouch Length | Filter | Box of |
|--------------|----------------|--------------|--------|--------|
| 18132 | 1-3/4" (44 mm) | 12" (30 cm) | No | 10 |
| 18133 | 2-1/4" (57 mm) | 12" (30 cm) | | |
| 18134 | 2-3/4" (70 mm) | 12" (30 cm) | | |
| 18192 | 1-3/4" (44 mm) | 12" (30 cm) | AF300 | 10 |
| 18193 | 2-1/4" (57 mm) | 12" (30 cm) | | |
| 18194 | 2-3/4" (70 mm) | 12" (30 cm) | | |
| 18292 | 1-3/4" (44 mm) | 7" (18 cm) | AF300 | 20 |
| 18293 | 2-1/4" (57 mm) | 7" (18 cm) | | |
| 18294 | 2-3/4" (70 mm) | 7" (18 cm) | | |

Colour: Beige

| Product # | Flange Size | Pouch Length | Filter | Box of |
|--------------|----------------|--------------|--------|--------|
| 18112 | 1-3/4" (44 mm) | 12" (30 cm) | No | 10 |
| 18113 | 2-1/4" (57 mm) | 12" (30 cm) | | |
| 18114 | 2-3/4" (70 mm) | 12" (30 cm) | | |
| 18182 | 1-3/4" (44 mm) | 12" (30 cm) | AF300 | 10 |
| 18183 | 2-1/4" (57 mm) | 12" (30 cm) | | |
| 18184 | 2-3/4" (70 mm) | 12" (30 cm) | | |
| 18282 | 1-3/4" (44 mm) | 7" (18 cm) | AF300 | 20 |
| 18283 | 2-1/4" (57 mm) | 7" (18 cm) | | |
| 18284 | 2-3/4" (70 mm) | 7" (18 cm) | | |

Lock 'n Roll Microseal Closure Features:

- A secure closure with tiny microseals
- An easy-to-open-and-close seal
- The material used on the closure is water-resistant and can be securely closed even when it is wet, but will work best when dry.

Clamp Closure



Colour: Ultra-clear

| Product # | Flange Size | Pouch Length | Filter | Box of |
|--------------|----------------|--------------|--------|--------|
| 18172 | 1-3/4" (44 mm) | 12" (30 cm) | No | 10 |
| 18173 | 2-1/4" (57 mm) | 12" (30 cm) | | |
| 18174 | 2-3/4" (70 mm) | 12" (30 cm) | | |
| 18176 | 4" (102 mm) | 12" (30 cm) | | |

Colour: Transparent

| Product # | Flange Size | Pouch Length | Filter | Box of |
|--------------|----------------|--------------|--------|--------|
| 18102 | 1-3/4" (44 mm) | 12" (30 cm) | No | 10 |
| 18103 | 2-1/4" (57 mm) | 12" (30 cm) | | |
| 18104 | 2-3/4" (70 mm) | 12" (30 cm) | | |

Colour: Beige

| Product # | Flange Size | Pouch Length | Filter | Box of |
|--------------|----------------|--------------|--------|--------|
| 18122 | 1-3/4" (44 mm) | 12" (30 cm) | No | 10 |
| 18123 | 2-1/4" (57 mm) | 12" (30 cm) | | |
| 18124 | 2-3/4" (70 mm) | 12" (30 cm) | | |
| 18142 | 1-3/4" (44 mm) | 12" (30 cm) | AF300 | 10 |
| 18143 | 2-1/4" (57 mm) | 12" (30 cm) | | |
| 18144 | 2-3/4" (70 mm) | 12" (30 cm) | | |
| 18202 | 1-3/4" (44 mm) | 9" (23 cm) | No | 10 |
| 18203 | 2-1/4" (57 mm) | 9" (23 cm) | | |
| 18204 | 2-3/4" (70 mm) | 9" (23 cm) | | |

New Image

Two-Piece Pouching System

Closed Pouches



Colour: Transparent

| Product # | Flange Size | Pouch Length | Filter | Box of |
|--------------|----------------|--------------|--------|--------|
| 18362 | 1-3/4" (44 mm) | 9" (23 cm) | AF300 | 60 |
| 18363 | 2-1/4" (57 mm) | 9" (23 cm) | | |
| 18364 | 2-3/4" (70 mm) | 9" (23 cm) | | |

Colour: Beige QuietWear Fabric

| Product # | Flange Size | Pouch Length | Filter | Box of |
|--------------|----------------|--------------|--------|--------|
| 18332 | 1-3/4" (44 mm) | 9" (23 cm) | No | 30 |
| 18333 | 2-1/4" (57 mm) | 9" (23 cm) | | |
| 18334 | 2-3/4" (70 mm) | 9" (23 cm) | | |
| 18322 | 1-3/4" (44 mm) | 9" (23 cm) | AF300 | 30 |
| 18323 | 2-1/4" (57 mm) | 9" (23 cm) | | |
| 18324 | 2-3/4" (70 mm) | 9" (23 cm) | | |

Colour: Beige

| Product # | Flange Size | Pouch Length | Filter | Box of |
|-----------------|----------------|--------------|--------|--------|
| 18352 | 1-3/4" (44 mm) | 7" (18 cm) | No | 30 |
| 18353 | 2-1/4" (57 mm) | 7" (18 cm) | | |
| 18354 | 2-3/4" (70 mm) | 7" (18 cm) | | |
| 18372-30 | 1-3/4" (44 mm) | 9" (23 cm) | AF300 | 30 |
| 18373-30 | 2-1/4" (57 mm) | 9" (23 cm) | | |
| 18374-30 | 2-3/4" (70 mm) | 9" (23 cm) | | |
| 18392 | 1-3/4" (44 mm) | 7" (18 cm) | AF300 | 60 |
| 18393 | 2-1/4" (57 mm) | 7" (18 cm) | | |
| 18394 | 2-3/4" (70 mm) | 7" (18 cm) | | |

High Output Pouches



Soft Tap Closure System

Colour: Ultra-Clear

| Product # | Flange Size | Pouch Length | Box of |
|--------------|----------------|--------------|--------|
| 18012 | 1-3/4" (44 mm) | 12" (30 cm) | 10 |
| 18013 | 2-1/4" (57 mm) | 12" (30 cm) | |
| 18014 | 2-3/4" (70 mm) | 12" (30 cm) | |
| 18016 | 4" (102 mm) | 12" (30 cm) | |

Coming Soon

New Enhanced Design High Output Pouches

- ✓ Clear and Beige
- ✓ Three Flange Sizes

New Image

Two-Piece Pouching System

Urostomy Pouches

Enhanced Design



Colour: Ultra-Clear

| Product # | Flange Size | Pouch Length | Box of |
|--------------|----------------|--------------|--------|
| 18922 | 1-3/4" (44 mm) | 9" (23 cm) | 10 |
| 18923 | 2-1/4" (57 mm) | 9" (23 cm) | |
| 18924 | 2-3/4" (70 mm) | 9" (23 cm) | |

Colour: Transparent

| Product # | Flange Size | Pouch Length | Box of |
|--------------|----------------|--------------|--------|
| 18902 | 1-3/4" (44 mm) | 9" (23 cm) | 10 |
| 18903 | 2-1/4" (57 mm) | 9" (23 cm) | |
| 18904 | 2-3/4" (70 mm) | 9" (23 cm) | |

Colour: Beige

| Product # | Flange Size | Pouch Length | Box of |
|--------------|----------------|--------------|--------|
| 18912 | 1-3/4" (44 mm) | 9" (23 cm) | 10 |
| 18913 | 2-1/4" (57 mm) | 9" (23 cm) | |
| 18914 | 2-3/4" (70 mm) | 9" (23 cm) | |

Original Design



Colour: Ultra-Clear

| Product # | Flange Size | Pouch Length | Box of |
|--------------|----------------|--------------|--------|
| 18422 | 1-3/4" (44 mm) | 9" (23 cm) | 10 |
| 18423 | 2-1/4" (57 mm) | 9" (23 cm) | |
| 18424 | 2-3/4" (70 mm) | 9" (23 cm) | |

Colour: Transparent

| Product # | Flange Size | Pouch Length | Box of |
|--------------|----------------|--------------|--------|
| 18402 | 1-3/4" (44 mm) | 9" (23 cm) | 10 |
| 18403 | 2-1/4" (57 mm) | 9" (23 cm) | |
| 18404 | 2-3/4" (70 mm) | 9" (23 cm) | |

Colour: Beige

| Product # | Flange Size | Pouch Length | Box of |
|--------------|----------------|--------------|--------|
| 18412 | 1-3/4" (44 mm) | 9" (23 cm) | 10 |
| 18413 | 2-1/4" (57 mm) | 9" (23 cm) | |
| 18414 | 2-3/4" (70 mm) | 9" (23 cm) | |

A Great Pouch Becomes Even Better

Enhanced Flow-Control Valve

The new valve is made with softer material and a streamlined shape.

New Multi-Chamber Design

The updated pouches now have multiple built-in chambers, which are designed to facilitate even distribution of urine.



Multi Chamber
Front View



Multi Chamber,
Top View, Sliced

SCAN ME



New Image

Two-Piece Pouching System

Ostomy Kits

Drainable









CeraPlus Barrier

Lock 'n Roll Closure System

| Non-sterile | | | | |
|--------------|----------------|------------------|--------------|--------|
| Product # | Flange Size | For Stomas Up To | Pouch Length | Box of |
| 19602 | 1-3/4" (44 mm) | 1-1/4" (32 mm) | 12" (30 mm) | 5 |
| 19603 | 2-1/4" (57 mm) | 1-3/4" (44 mm) | 12" (30 mm) | |
| 19604 | 2-3/4" (70 mm) | 2-1/4" (57 mm) | 12" (30 mm) | |

FlexWear Barrier

Lock 'n Roll Closure System

| Non-sterile | | | | | |
|-------------|---|----------------|------------------|--------------|--------|
| Product # | | Flange Size | For Stomas Up To | Pouch Length | Box of |
| 19002 |  | 1-3/4" (44 mm) | 1-1/4" (32 mm) | 12" (30 mm) | 5 |
| 19003 |  | 2-1/4" (57 mm) | 1-3/4" (44 mm) | 12" (30 mm) | |
| 19004 |  | 2-3/4" (70 mm) | 2-1/4" (57 mm) | 12" (30 mm) | |
| 19006 |  | 4" (102 mm) | 3-1/2" (89 mm) | 12" (30 mm) | |
| Sterile | | | | | |
| Product # | | Flange Size | For Stomas Up To | Pouch Length | Box of |
| 19052 |  | 1-3/4" (44 mm) | 1-1/4" (32 mm) | 12" (30 mm) | 5 |
| 19053 |  | 2-1/4" (57 mm) | 1-3/4" (44 mm) | 12" (30 mm) | |
| 19054 |  | 2-3/4" (70 mm) | 2-1/4" (57 mm) | 12" (30 mm) | |
| 19056 |  | 4" (102 mm) | 3-1/2" (89 mm) | 12" (30 mm) | |

Urostomy Kits

CeraPlus Barrier

| Non-sterile | | | |
|--------------|----------------|------------------|--------|
| Product # | Flange Size | For Stomas Up To | Box of |
| 19902 | 1-3/4" (44 mm) | 1-1/4" (32 mm) | 5 |
| 19903 | 2-1/4" (57 mm) | 1-3/4" (44 mm) | |
| 19904 | 2-3/4" (70 mm) | 2-1/4" (57 mm) | |

Urostomy Kits Include:

- Ultra-clear pouch
- Cut-to-fit Flexextend Skin Barrier
- Drain tube adapter









CeraPlus Barrier

Clamp Closure System

| Non-sterile | | | | |
|--------------|----------------|------------------|--------------|--------|
| Product # | Flange Size | For Stomas Up To | Pouch Length | Box of |
| 19802 | 1-3/4" (44 mm) | 1-1/4" (32 mm) | 12" (30 mm) | 5 |
| 19803 | 2-1/4" (57 mm) | 1-3/4" (44 mm) | 12" (30 mm) | |
| 19804 | 2-3/4" (70 mm) | 2-1/4" (57 mm) | 12" (30 mm) | |

FlexWear Barrier







Clamp Closure System

| Non-sterile | | | | | |
|-------------|---|----------------|------------------|--------------|--------|
| Product # | | Flange Size | For Stomas Up To | Pouch Length | Box of |
| 19102 |  | 1-3/4" (44 mm) | 1-1/4" (32 mm) | 12" (30 mm) | 5 |
| 19103 |  | 2-1/4" (57 mm) | 1-3/4" (44 mm) | 12" (30 mm) | |
| 19104 |  | 2-3/4" (70 mm) | 2-1/4" (57 mm) | 12" (30 mm) | |
| 19106 |  | 4" (102 mm) | 3-1/2" (89 mm) | 12" (30 mm) | |
| Sterile | | | | | |
| Product # | | Flange Size | For Stomas Up To | Pouch Length | Box of |
| 19152 |  | 1-3/4" (44 mm) | 1-1/4" (32 mm) | 12" (30 mm) | 5 |
| 19153 |  | 2-1/4" (57 mm) | 1-3/4" (44 mm) | 12" (30 mm) | |
| 19154 |  | 2-3/4" (70 mm) | 2-1/4" (57 mm) | 12" (30 mm) | |
| 19156 |  | 4" (102 mm) | 3-1/2" (89 mm) | 12" (30 mm) | |

Ostomy Kits Include:

- Skin barrier with tape border
- Ultra-clear drainable pouch
- Loop ostomy bridge with 4" flange

Flexextend Barrier

| Non-sterile | | | | |
|-------------|---|----------------|------------------|--------|
| Product # | | Flange Size | For Stomas Up To | Box of |
| 19202 |  | 1-3/4" (44 mm) | 1-1/4" (32 mm) | 5 |
| 19203 |  | 2-1/4" (57 mm) | 1-3/4" (44 mm) | |
| 19204 |  | 2-3/4" (70 mm) | 2-1/4" (57 mm) | |
| Sterile | | | | |
| Product # | | Flange Size | For Stomas Up To | Box of |
| 19252 |  | 1-3/4" (44 mm) | 1-1/4" (32 mm) | 5 |
| 19253 |  | 2-1/4" (57 mm) | 1-3/4" (44 mm) | |
| 19254 |  | 2-3/4" (70 mm) | 2-1/4" (57 mm) | |

LEARN MORE



Follow the evidence to help ensure healthy skin

After surgery: Peristomal skin starts healthy



Top layers of the epidermis are removed with *every barrier change*.


It's unavoidable, and the damage isn't always visible.



Epidermal stripping is linked to increased levels
of transepidermal water loss (TEWL)^{1,2}

Healthy Skin:

Keeps moisture in and irritants out




Ceramides are abundant in healthy skin.
They form a waterproof barrier to help keep
moisture in and bacteria out.

— = Ceramides

Damaged Skin:

Allows moisture out and irritants in



Damaged skin contains fewer ceramides.

— = Water droplet illustrates moisture loss from skin.

References:

1. McNichol L, Lund C, Rosen T, Gray M. Medical adhesives and patient safety: state of the science: consensus statements for the assessment, prevention, and treatment of adhesive-related skin injuries. J Wound Ostomy Continence Nurs. 2013;40(4):365-380.
2. Grove G, Houser T, Sibbald G, Salvadalena G. Measuring epidermal effects of ostomy skin barriers. Skin Res Technol. 2019;25:179-186. <https://doi.org/10.1111/srt.12630>
3. Nybaek H, Lophagen S, Karlsmark T, Bang Knudsen D, Jemec GBE. Stratum corneum integrity as a predictor for peristomal skin problems in ostomates. Br J Dermatol. 2010;162:357-361.

Healthy Skin. Positive Outcomes.

What Happens When Skin is Compromised

Some common signs of a Peristomal Skin Complication (PSC) include:
**ITCHING • REDNESS • PAIN • WEEPING SKIN
BLEEDING • BURNING**



PSCs are often classified as **Peristomal Medical Adhesive-Related Skin Injury (PMARSI)⁹** and **Peristomal Moisture-Associated Skin Damage (PMASD)¹**



PSCs associated with
**INCREASED COST
OF CARE⁶**



PSCs associated with
**LOWER QUALITY
OF LIFE⁷**

What Happens When You Choose CeraPlus Skin Barrier



RELATIVE REDUCTION
in the number of PSCs⁵



INCREASED SATISFACTION
with prevention of itching⁵



Significantly
**LOWER STOMA-RELATED
COST OF CARE⁵**



Significantly
**MORE LIKELY TO
RESOLVE PSCs⁵**



**TEWL IS
MINIMIZED^{2,8}**



Healthier peristomal skin =
HIGHER QUALITY OF LIFE⁷

You have a choice. Choose CeraPlus skin barrier and help ensure healthy peristomal skin for your patients from **DAY ONE.**

References:

- Colwell J, Pittman J, Raizman R, Salvadalena G. A Randomized Controlled Trial Determining Variances in Ostomy Skin Conditions (ADVOCATE). J Wound Ostomy Continence Nurs. 2018;45(1):37-42.
- Taneja, C., et al., Risk and Economic Burden of Peristomal Skin Complications Following Ostomy Surgery. J Wound Ostomy Continence Nurs. 2019. 46(2): 143-149
- Nichols T, Inglese G. The Burden of Peristomal Skin Complications on an Ostomy Population as Assessed by Health Utility and the Physical Component Summary of the SF-36v2 International Society for Pharmacoeconomics and Outcomes Research (ISPOR), Science Direct Value in Health; 2018; p. 92-93, Tables 3-4
- Nichols T, Houser T, Grove G. Comparing the skin stripping effects of three ostomy skin barriers infused with ceramide, honey or aloe; Journal of Stomal Therapy Australia 2019; 39(2):14-18
- Gray M, Colwell J, Doughty D, Goldberg M, Hoek J, Manson A, McNichol L, Rao S. (2013) Peristomal Moisture-Associated Skin Damage in Adults with Fecal Ostomies: A Comprehensive Review and Consensus. J Wound Ostomy Continence Nurs: 38(5): 541-53



Premier


One-Piece Pouching System


Flat Skin Barriers

Drainable Pouch

CeraPlus with Remois Technology*
Lock 'n Roll Closure System




| Cut-to-fit With Tape Border  | | | | | |
|---|---------------------------|------------------|--------------|--------|--------|
| Product # | Colour | For Stomas Up To | Pouch Length | Filter | Box of |
| 8931 | Ultra-clear | 2-1/2" (64 mm) | 12" (30 cm) | No | 10 |
| 8901 | Beige with Viewing Option | 2-1/2" (64 mm) | 12" (30 cm) | AF300 | |


| Pre-sized With Tape Border  | | | | | |
|--|---------------------------|-----------------|--------------|--------|--------|
| Product # | Colour | Barrier Opening | Pouch Length | Filter | Box of |
| 8925 | Beige with Viewing Option | 1" (25 mm) | 12" (30 cm) | AF300 | 10 |
| 8930 | | 1-3/16" (30 mm) | 12" (30 cm) | | |
| 8935 | | 1-3/8" (35 mm) | 12" (30 cm) | | |
| 8940 | | 1-9/16" (40 mm) | 12" (30 cm) | | |

*Remois is a technology of Alcare Co., Ltd.

Flextend
Lock 'n Roll Closure System




| Cut-to-fit With Tape Border  | | | | | |
|---|-------------|------------------|--------------|--------|--------|
| Product # | Colour | For Stomas Up To | Pouch Length | Filter | Box of |
| 8531 | Ultra-clear | 2-1/2" (64 mm) | 12" (30 cm) | No | 10 |
| 8331 | Transparent | 2-1/2" (64 mm) | 12" (30 cm) | AF300 | |
| 8551 | Beige | 2-1/2" (64 mm) | 12" (30 cm) | No | |

| Pre-sized With Tape Border  | | | | |
|--|-------------|-----------------|--------------|--------|
| Product # | Colour | Barrier Opening | Pouch Length | Box of |
| 8538 | Transparent | 1-1/4" (32 mm) | 12" (30 cm) | 10 |
| 8552 | Beige | 1" (25 mm) | 12" (30 cm) | 10 |
| 8558 | | 1-1/4" (32 mm) | 12" (30 cm) | |
| 8559 | | 1-3/4" (44 mm) | 12" (30 cm) | |


FlexWear
Lock 'n Roll Closure System



| Cut-to-fit With Tape Border  | | | | | |
|---|-------------|------------------|--------------|--------|--------|
| Product # | Colour | For Stomas Up To | Pouch Length | Filter | Box of |
| 8031 | Ultra-clear | 2-1/2" (64 mm) | 12" (30 cm) | No | 10 |
| 8081 | Transparent | 2-1/2" (64 mm) | 12" (30 cm) | AF300 | |
| 8181 | Beige | 2-1/2" (64 mm) | 12" (30 cm) | AF300 | |


FlexWear
Clamp Closure System




| Cut-to-fit With Tape Border  | | | | | |
|---|-------------|------------------|--------------|--------|--------|
| Product # | Colour | For Stomas Up To | Pouch Length | Filter | Box of |
| 8131 | Ultra-clear | 2-1/2" (64 mm) | 12" (30 cm) | No | 10 |

Flextend
Clamp Closure System



| Cut-to-fit With Tape Border  | | | | |
|---|-------------|------------------|--------------|--------|
| Product # | Colour | For Stomas Up To | Pouch Length | Box of |
| 8651 | Ultra-clear | 2-1/2" (64 mm) | 12" (30 cm) | 10 |
| 8631 | Transparent | 2-1/2" (64 mm) | 12" (30 cm) | 10 |
| 8641 | | 2-1/2" (64 mm) | 9" (23 cm) | |

| Pre-sized With Tape Border  | | | | |
|--|-------------|-----------------|--------------|--------|
| Product # | Colour | Barrier Opening | Pouch Length | Box of |
| 8637 | Transparent | 3/4" (19 mm) | 12" (30 cm) | 10 |
| 8632 | | 1" (25 mm) | 12" (30 cm) | |
| 8638 | | 1-1/4" (32 mm) | 12" (30 cm) | |
| 8633 | | 1-1/2" (38 mm) | 12" (30 cm) | |
| 8639 | | 1-3/4" (44 mm) | 12" (30 cm) | |
| 8634 | | 2" (51 mm) | 12" (30 cm) | |
| 8635 | | 2-1/2" (64 mm) | 12" (30 cm) | |
| 8647 | Beige | 3/4" (19 mm) | 9" (23 cm) | 10 |
| 8642 | | 1" (25 mm) | 9" (23 cm) | |
| 8648 | | 1-1/4" (32 mm) | 9" (23 cm) | |
| 8643 | | 1-1/2" (38 mm) | 9" (23 cm) | |
| 8649 | | 1-3/4" (44 mm) | 9" (23 cm) | |
| 8644 | | 2" (51 mm) | 9" (23 cm) | |

Premier

One-Piece Pouching System

Flat Skin Barriers

Drainable Pouch

SoftFlex
Lock 'n Roll Closure System



| Cut-to-fit Without Tape Border | | | | | |
|--------------------------------|--------|------------------|--------------|--------|--------|
| Product # | Colour | For Stomas Up To | Pouch Length | Filter | Box of |
| 88300 | Beige | 2-1/8" (55 mm) | 12" (30 cm) | AF300 | 10 |
| 88700 | | 2-1/8" (55 mm) | 7" (18 cm) | | 20 |

| Cut-to-fit—Oval Skin Barrier Without Tape Border | | | | | |
|--|-------------|--------------------------|--------------|--------|--------|
| Product # | Colour | For Stomas Up To | Pouch Length | Filter | Box of |
| 88402 | Transparent | 2-1/2" x 3" (66 x 77 mm) | 12" (30 cm) | AF300 | 10 |
| 88302 | Beige | 2-1/2" x 3" (66 x 77 mm) | 12" (30 cm) | | |

| Pre-sized Without Tape Border | | | | | |
|-------------------------------|--------|-----------------|--------------|--------|--------|
| Product # | Colour | Barrier Opening | Pouch Length | Filter | Box of |
| 88325 | Beige | 1" (25 mm) | 12" (30 cm) | AF300 | 10 |
| 88330 | | 1-3/16" (30 mm) | 12" (30 cm) | | |
| 88335 | | 1-3/8" (35 mm) | 12" (30 cm) | | |
| 88340 | | 1-9/16" (40 mm) | 12" (30 cm) | | |
| 88725 | | 1" (25 mm) | 7" (18 cm) | AF300 | 20 |
| 88730 | | 1-3/16" (30 mm) | 7" (18 cm) | | |
| 88735 | | 1-3/8" (35 mm) | 7" (18 cm) | | |
| 88740 | | 1-9/16" (40 mm) | 7" (18 cm) | | |

Flat Skin Barriers

Closed Pouch

SoftFlex



| Cut-to-fit Without Tape Border | | | | | |
|--------------------------------|--------|------------------|--------------|--------|--------|
| Product # | Colour | For Stomas Up To | Pouch Length | Filter | Box of |
| 82300 | Beige | 2-1/8" (55 mm) | 9" (23 cm) | AF300 | 30 |
| 82100 | | 2-1/8" (55 mm) | 7" (18 cm) | | |

| Cut-to-fit—Oval Skin Barrier Without Tape Border | | | | | |
|--|-------------|-----------------------------|--------------|--------|--------|
| Product # | Colour | For Stomas Up To | Pouch Length | Filter | Box of |
| 82402 | Transparent | 2-1/2" x 3" (66 mm x 77 mm) | 9" (23 cm) | AF300 | 30 |
| 82302 | Beige | 2-1/2" x 3" (66 mm x 77 mm) | 9" (23 cm) | | |

| Pre-sized Without Tape Border | | | | | |
|-------------------------------|--------|-----------------|--------------|--------|--------|
| Product # | Colour | Barrier Opening | Pouch Length | Filter | Box of |
| 82325 | Beige | 1" (25 mm) | 9" (23 cm) | AF300 | 30 |
| 82330 | | 1-3/16" (30 mm) | 9" (23 cm) | | |
| 82335 | | 1-3/8" (35 mm) | 9" (23 cm) | | |
| 82125 | | 1" (25 mm) | 7" (18 cm) | | |
| 82130 | | 1-3/16" (30 mm) | 7" (18 cm) | | |
| 82135 | | 1-3/8" (35 mm) | 7" (18 cm) | | |

High Output Pouches

Flextend



| Non-sterile Without Tape Border | | | |
|---------------------------------|-----------------|--------------|--------|
| Product # | Barrier Opening | Pouch Length | Box of |
| 80070 | 2-3/4" (70 mm) | 12" (30 cm) | 10 |
| 80110 | 4-1/3" (110 mm) | 12" (30 cm) | |

Coming Soon

New Enhanced Design High Output Pouches

- ✓ Flextend Skin Barrier
- ✓ Flat and Soft Convex
- ✓ Cut-to-fit and Pre-sized options

Premier

One-Piece Pouching System

Soft Convex Skin Barriers Drainable Pouch

CeraPlus with Remois Technology*
Lock 'n Roll Closure System



| Cut-to-fit | | | With Tape Border | | |
|------------|---------------------------|------------------|------------------|--------|--------|
| Product # | | For Stomas Up To | Pouch Length | Filter | Box of |
| 8954 | Beige with Viewing Option | 1" (25 mm) | 12" (30 cm) | AF300 | 5 |
| 8958 | | 1-1/2" (38 mm) | 12" (30 cm) | | |
| 89511 | | 2-1/8" (55 mm) | 12" (30 cm) | | |

| Pre-sized | | | With Tape Border | | |
|-----------|---------------------------|-----------------|------------------|--------|--------|
| Product # | | Barrier Opening | Pouch Length | Filter | Box of |
| 8960 | Beige with Viewing Option | 3/4" (19 mm) | 12" (30 cm) | AF300 | 5 |
| 8961 | | 7/8" (22 mm) | 12" (30 cm) | | |
| 8962 | | 1" (25 mm) | 12" (30 cm) | | |
| 8963 | | 1-1/8" (29 mm) | 12" (30 cm) | | |

*Remois is a technology of Alcare Co., Ltd.

Flextend Lock 'n Roll Closure System



| Cut-to-fit | | | | With Tape Border | |
|------------|---------------------------|----------------------------|--------------|------------------|--------|
| Product # | Colour | Stoma Size | Pouch Length | Filter | Box of |
| 8674 | Beige with Viewing Option | 5/8" - 1" (13 - 25 mm) | 12" (30 cm) | AF300 | 5 |
| 8678 | | 5/8" - 1-1/2" (13 - 38 mm) | 12" (30 cm) | | |
| 86711 | | 5/8" - 2-1/8" (13 - 55 mm) | 12" (30 cm) | | |

| Pre-sized | | | | With Tape Border | |
|-----------|---------------------------|-----------------|--------------|------------------|--------|
| Product # | Colour | Barrier Opening | Pouch Length | Filter | Box of |
| 8661 | Beige with Viewing Option | 7/8" (22 mm) | 12" (30 cm) | AF300 | 5 |
| 8662 | | 1" (25 mm) | 12" (30 cm) | | |
| 8663 | | 1-1/8" (29 mm) | 12" (30 cm) | | |
| 8664 | | 1-1/4" (32 mm) | 12" (30 cm) | | |
| 8665 | | 1-1/2" (38 mm) | 12" (30 cm) | | |
| 8666 | | 1-3/4" (44 mm) | 12" (30 cm) | | |

Convex Skin Barriers Drainable Pouch

CeraPlus with Remois Technology*
Lock 'n Roll Closure System



| Cut-to-fit | | | | | With Tape Border |
|------------|---------------------------|------------------|--------------|--------|------------------|
| Product # | | For Stomas Up To | Pouch Length | Filter | Box of |
| 8984 | Ultra-clear | 1" (25 mm) | 12" (30 cm) | No | 5 |
| 8988 | | 1-1/2" (38 mm) | 12" (30 cm) | | |
| 89811 | | 2" (51 mm) | 12" (30 cm) | | |
| 8914 | Beige with Viewing Option | 1" (25 mm) | 12" (30 cm) | AF300 | 5 |
| 8918 | | 1-1/2" (38 mm) | 12" (30 cm) | | |
| 89111 | | 2" (51 mm) | 12" (30 cm) | | |

| Pre-sized | | | | | With Tape Border |
|-----------|---------------------------|-----------------|--------------|--------|------------------|
| Product # | | Barrier Opening | Pouch Length | Filter | Box of |
| 8990 | Beige with Viewing Option | 3/4" (19 mm) | 12" (30 cm) | AF300 | 5 |
| 8991 | | 7/8" (22 mm) | 12" (30 cm) | | |
| 8992 | | 1" (25 mm) | 12" (30 cm) | | |
| 8993 | | 1-1/8" (29 mm) | 12" (30 cm) | | |

*Remois is a technology of Alcare Co., Ltd.

Soft Convex Skin Barriers

Soft convexity is designed to help provide gentle peristomal pressure for an individualized fit. It helps conform to accommodate uneven skin surfaces, and is a more flexible alternative to firm convexity.

- Lock 'n Roll closure
- Tape border
- Beige with viewing option
- AF300 filter
- Belt tabs



LEARN MORE



Premier

One-Piece Pouching System

Convex Skin Barriers

Drainable Pouch

Flextend
Lock 'n Roll Closure System



Cut-to-fit With Tape Border 

Colour: Transparent

| Product # | | For Stomas Up To | Pouch Length | Filter | Box of |
|----------------|--------------------|------------------|--------------|--------|--------|
| Tapered Border | Non-Tapered Border | | | | |
| 8244 | 8524 | 1" (25 mm) | 12" (30 cm) | No | 5 |
| 8248 | 8528 | 1-1/2" (38 mm) | 12" (30 cm) | | |
| 82411 | 85211 | 2" (51 mm) | 12" (30 cm) | | |
| 8274 | 8584 | 1" (25 mm) | 12" (30 cm) | | |
| 8278 | 8588 | 1-1/2" (38 mm) | 12" (30 cm) | AF300 | |
| 82711 | 85811 | 2" (51 mm) | 12" (30 cm) | | |

Colour: Beige

| | | | | | |
|--------------|--------------|----------------|-------------|----|---|
| 8214 | 8504 | 1" (25 mm) | 12" (30 cm) | No | 5 |
| 8218 | 8508 | 1-1/2" (38 mm) | 12" (30 cm) | | |
| 82111 | 85011 | 2" (51 mm) | 12" (30 cm) | | |

Pre-sized With Tape Border 

Colour: Transparent

| Product # | | Barrier Opening | Pouch Length | Box of |
|----------------|--------------------|-----------------|--------------|--------|
| Tapered Border | Non-Tapered Border | | | |
| 8280 | 8510 | 3/4" (19 mm) | 12" (30 cm) | 5 |
| 8281 | 8511 | 7/8" (22 mm) | 12" (30 cm) | |
| 8282 | 8512 | 1" (25 mm) | 12" (30 cm) | |
| 8283 | 8513 | 1-1/8" (29 mm) | 12" (30 cm) | |
| 8284 | 8514 | 1-1/4" (32 mm) | 12" (30 cm) | |
| 8285 | 8515 | 1-3/8" (35 mm) | 12" (30 cm) | |
| 8286 | 85116 | 1-1/2" (38 mm) | 12" (30 cm) | |
| 8288 | 85118 | 1-3/4" (44 mm) | 12" (30 cm) | |
| 8289 | 85119 | 2" (51 mm) | 12" (30 cm) | |

Colour: Beige

| | | | | |
|-------------|-------------|----------------|-------------|---|
| 8290 | 8590 | 3/4" (19 mm) | 12" (30 cm) | 5 |
| 8291 | 8591 | 7/8" (22 mm) | 12" (30 cm) | |
| 8292 | 8592 | 1" (25 mm) | 12" (30 cm) | |
| 8293 | 8593 | 1-1/8" (29 mm) | 12" (30 cm) | |
| 8294 | 8594 | 1-1/4" (32 mm) | 12" (30 cm) | |
| 8295 | 8595 | 1-3/8" (35 mm) | 12" (30 cm) | |
| 8296 | 8596 | 1-1/2" (38 mm) | 12" (30 cm) | |
| 8297 | 8597 | 1-5/8" (41 mm) | 12" (30 cm) | |


Flextend
Clamp Closure System



Cut-to-fit With Tape Border 

Colour: Transparent

| Product # | | For Stomas Up To | Pouch Length | Box of |
|----------------|--------------------|------------------|--------------|--------|
| Tapered Border | Non-Tapered Border | | | |
| 82624 | 8624 | 1" (25 mm) | 12" (30 cm) | 5 |
| 82628 | 8628 | 1-1/2" (38 mm) | 12" (30 cm) | |
| 826211 | 86211 | 2" (51 mm) | 12" (30 cm) | |

Pre-sized With Tape Border 

Colour: Transparent

| Product # | | Barrier Opening | Pouch Length | Box of |
|----------------|--------------------|-----------------|--------------|--------|
| Tapered Border | Non-Tapered Border | | | |
| 82610 | 8610 | 3/4" (19 mm) | 12" (30 cm) | 5 |
| 82611 | 8611 | 7/8" (22 mm) | 12" (30 cm) | |
| 82612 | 8612 | 1" (25 mm) | 12" (30 cm) | |
| 82613 | 8613 | 1-1/8" (29 mm) | 12" (30 cm) | |
| 82614 | 8614 | 1-1/4" (32 mm) | 12" (30 cm) | |
| 82615 | 8615 | 1-3/8" (35 mm) | 12" (30 cm) | |
| 82616 | 8616 | 1-1/2" (38 mm) | 12" (30 cm) | |
| 82617 | 8617 | 1-5/8" (41 mm) | 12" (30 cm) | |
| 82618 | 8618 | 1-3/4" (44 mm) | 12" (30 cm) | |
| 82619 | 8619 | 2" (51 mm) | 12" (30 cm) | |

Tapered Border is Designed for Enhanced Sense of Security

- Tapered barrier edges help prevent lift from moisture or motion, which can help achieve desired wear time
- Gentle on the skin while still giving the feeling of extra security
- Adhesive flexes and bends to move comfortably with the body
- A new flexible plastic release liner that covers both the skin barrier and the adhesive border to help improve ease-of-use



Premier


One-Piece Pouching System

Flat Skin Barriers

Urostomy Pouch

CeraPlus with Remois Technology*
Enhanced Design



| Cut-to-fit | | | | | With Tape Border  |
|--------------|-------------|----------------|--------------|--------|--|
| Product # | Colour | For Stomas | Pouch Length | Filter | Box of |
| 84200 | Ultra-clear | 2-1/2" (64 mm) | 9" (23 cm) | No | 10 |

*Remois is a technology of Alcare Co., Ltd.


The CeraPlus Skin Barrier
delivers both Fit and Formulation


Ceramide Infused Barrier

- **Protects** the skin's natural moisture barrier
- **Decreases** TEWL (transepidermal water loss) from damaged or eroded skin
- **Maintains** peristomal skin health

Flextend
Enhanced Design





| Cut-to-fit | | | | | With Tape Border  |
|--------------|-------------|------------------|--------------|--------|--|
| Product # | Colour | For Stomas Up To | Pouch Length | Box of | |
| 84590 | Ultra-clear | 2-1/2" (64 mm) | 9" (23 cm) | 10 | |
| 84490 | Beige | 2-1/2" (64 mm) | 9" (23 cm) | | |

| Pre-sized | | | | | With Tape Border  |
|--------------|-------------|-----------------|--------------|--------|--|
| Product # | Colour | Barrier Opening | Pouch Length | Box of | |
| 84697 | Transparent | 3/4" (19 mm) | 9" (23 cm) | 10 | |
| 84692 | | 1" (25 mm) | 9" (23 cm) | | |
| 84698 | | 1-1/4" (32 mm) | 9" (23 cm) | | |
| 84693 | | 1-1/2" (38 mm) | 9" (23 cm) | | |
| 84699 | | 1-3/4" (44 mm) | 9" (23 cm) | | |
| 84694 | | 2" (51 mm) | 9" (23 cm) | | |
| 84497 | Beige | 3/4" (19 mm) | 9" (23 cm) | 10 | |
| 84492 | | 1" (25 mm) | 9" (23 cm) | | |
| 84498 | | 1-1/4" (32 mm) | 9" (23 cm) | | |
| 84493 | | 1-1/2" (38 mm) | 9" (23 cm) | | |
| 84499 | | 1-3/4" (44 mm) | 9" (23 cm) | | |
| 84494 | | 2" (51 mm) | 9" (23 cm) | | |

Flextend
Original Design



| Cut-to-fit | | | | | With Tape Border  |
|-------------|-------------|------------------|--------------|--------|--|
| Product # | Colour | For Stomas Up To | Pouch Length | Box of | |
| 8450 | Ultra-clear | 2-1/2" (64 mm) | 9" (23 cm) | 10 | |
| 8460 | Transparent | 2-1/2" (64 mm) | 9" (23 cm) | | |
| 8440 | Beige | 2-1/2" (64 mm) | 9" (23 cm) | | |

| Pre-sized | | | | | With Tape Border  |
|-------------|-------------|-----------------|--------------|--------|--|
| Product # | Colour | Barrier Opening | Pouch Length | Box of | |
| 8467 | Transparent | 3/4" (19 mm) | 9" (23 cm) | 10 | |
| 8462 | | 1" (25 mm) | 9" (23 cm) | | |
| 8468 | | 1-1/4" (32 mm) | 9" (23 cm) | | |
| 8463 | | 1-1/2" (38 mm) | 9" (23 cm) | | |
| 8469 | | 1-3/4" (44 mm) | 9" (23 cm) | | |
| 8464 | | 2" (51 mm) | 9" (23 cm) | | |
| 8447 | Beige | 3/4" (19 mm) | 9" (23 cm) | 10 | |
| 8442 | | 1" (25 mm) | 9" (23 cm) | | |
| 8448 | | 1-1/4" (32 mm) | 9" (23 cm) | | |
| 8443 | | 1-1/2" (38 mm) | 9" (23 cm) | | |
| 8449 | | 1-3/4" (44 mm) | 9" (23 cm) | | |


Premier


One-Piece Pouching System

Soft Convex Skin Barriers Urostomy Pouch

CeraPlus with Remois Technology*
Enhanced Design



| Cut-to-fit | | | | With Tape Border  |
|---------------|-------------|------------------|--------------|--|
| Product # | Colour | For Stomas Up To | Pouch Length | Box of |
| 84134 | Ultra-clear | 1" (25 mm) | 9" (23 cm) | 5 |
| 84138 | | 1-1/2" (38 mm) | 9" (23 cm) | |
| 841311 | | 2-1/8" (55 mm) | 9" (23 cm) | |


| Pre-sized | | | | With Tape Border  |
|-------------|---------------------------|-----------------|--------------|--|
| Product # | Colour | Barrier Opening | Pouch Length | Box of |
| 8412 | Beige with Viewing Option | 3/4" (19 mm) | 9" (23 cm) | 5 |
| 8413 | | 7/8" (22 mm) | 9" (23 cm) | |
| 8414 | | 1" (25 mm) | 9" (23 cm) | |
| 8415 | | 1-1/8" (29 mm) | 9" (23 cm) | |

*Remois is a technology of Alcare Co., Ltd.

Convex Skin Barriers Urostomy Pouch

CeraPlus with Remois Technology*
Enhanced Design



| Cut-to-fit | | | | With Tape Border  |
|---------------|-------------|------------------|--------------|--|
| Product # | Colour | For Stomas Up To | Pouch Length | Box of |
| 84404 | Ultra-clear | 1" (25 mm) | 9" (23 cm) | 5 |
| 84408 | | 1-1/2" (38 mm) | 9" (23 cm) | |
| 844011 | | 2" (51 mm) | 9" (23 cm) | |

| Pre-sized | | | | |
|--------------|---------------------------|-----------------|--------------|--------|
| Product # | Colour | Barrier Opening | Pouch Length | Box of |
| 84292 | Beige with Viewing Option | 3/4" (19 mm) | 9" (23 cm) | 5 |
| 84293 | | 7/8" (22 mm) | 9" (23 cm) | |
| 84294 | | 1" (25 mm) | 9" (23 cm) | |
| 84295 | | 1-1/8" (29 mm) | 9" (23 cm) | |

*Remois is a technology of Alcare Co., Ltd.

Premier Soft Convex CeraPlus Pouching Systems:

- Are designed to conform to the body without compromising fit
- Deliver both convex support and flexibility
- Protect the skin's natural moisture barrier and help maintain good peristomal skin health



LEARN MORE




A Great Pouch Becomes Even Better

Enhanced Flow-Control Valve

The new valve is made with softer material and a streamlined shape.

New Multi-Chamber Design

The updated pouches now have multiple built-in chambers, which are designed to facilitate even distribution of urine.



Multi Chamber
Front View



Multi Chamber,
Top View, Sliced

Premier


One-Piece Pouching System


Convex Skin Barriers

Urostomy Pouch

Flextend
Enhanced Design



| Cut-to-fit | | | With Tape Border  | |
|----------------|-------------|------------------|--|--------|
| Product # | Colour | For Stomas Up To | Pouch Length | Box of |
| 84794 | Ultra-clear | 1" (25 mm) | 9" (23 cm) | 5 |
| 84798 | | 1-1/2" (38 mm) | 9" (23 cm) | |
| 8479111 | | 2" (51 mm) | 9" (23 cm) | |
| 84394 | Beige | 1" (25 mm) | 9" (23 cm) | |
| 84398 | | 1-1/2" (38 mm) | 9" (23 cm) | |
| 8439111 | | 2" (51 mm) | 9" (23 cm) | |


| Pre-sized | | | With Tape Border  | |
|--------------|-------------|-----------------|--|--------|
| Product # | Colour | Barrier Opening | Pouch Length | Box of |
| 84890 | Ultra-clear | 1/2" (13 mm) | 9" (23 cm) | 5 |
| 84891 | | 5/8" (16 mm) | 9" (23 cm) | |
| 84892 | | 3/4" (19 mm) | 9" (23 cm) | |
| 84893 | | 7/8" (22 mm) | 9" (23 cm) | |
| 84894 | | 1" (25 mm) | 9" (23 cm) | |
| 84895 | | 1-1/8" (29 mm) | 9" (23 cm) | |
| 84896 | | 1-1/4" (32 mm) | 9" (23 cm) | |
| 84897 | | 1-3/8" (35 mm) | 9" (23 cm) | |
| 84898 | | 1-1/2" (38 mm) | 9" (23 cm) | |
| 84899 | | 1-3/4" (44 mm) | 9" (23 cm) | |
| 84990 | Beige | 1/2" (13 mm) | 9" (23 cm) | 5 |
| 84991 | | 5/8" (16 mm) | 9" (23 cm) | |
| 84992 | | 3/4" (19 mm) | 9" (23 cm) | |
| 84993 | | 7/8" (22 mm) | 9" (23 cm) | |
| 84994 | | 1" (25 mm) | 9" (23 cm) | |
| 84995 | | 1-1/8" (29 mm) | 9" (23 cm) | |
| 84996 | | 1-1/4" (32 mm) | 9" (23 cm) | |
| 84997 | | 1-3/8" (35 mm) | 9" (23 cm) | |
| 84998 | | 1-1/2" (38 mm) | 9" (23 cm) | |
| 84999 | | 1-3/4" (44 mm) | 9" (23 cm) | |


LEARN MORE



Flextend
Original Design



| Cut-to-fit | | | With Tape Border  | |
|--------------|-------------|------------------|--|--------|
| Product # | Colour | For Stomas Up To | Pouch Length | Box of |
| 8474 | Transparent | 1" (25 mm) | 9" (23 cm) | 5 |
| 8478 | | 1-1/2" (38 mm) | 9" (23 cm) | |
| 84711 | | 2" (51 mm) | 9" (23 cm) | |

| Pre-sized | | | With Tape Border  | |
|-------------|-------------|-----------------|--|--------|
| Product # | Colour | Barrier Opening | Pouch Length | Box of |
| 8480 | Transparent | 1/2" (13 mm) | 9" (23 cm) | 5 |
| 8481 | | 5/8" (16 mm) | 9" (23 cm) | |
| 8482 | | 3/4" (19 mm) | 9" (23 cm) | |
| 8483 | | 7/8" (22 mm) | 9" (23 cm) | |
| 8484 | | 1" (25 mm) | 9" (23 cm) | |
| 8485 | | 1-1/8" (29 mm) | 9" (23 cm) | |
| 8486 | | 1-1/4" (32 mm) | 9" (23 cm) | |
| 8487 | | 1-3/8" (35 mm) | 9" (23 cm) | |
| 8488 | | 1-1/2" (38 mm) | 9" (23 cm) | |
| 8489 | | 1-3/4" (44 mm) | 9" (23 cm) | |
| 8490 | Beige | 1/2" (13 mm) | 9" (23 cm) | 5 |
| 8491 | | 5/8" (16 mm) | 9" (23 cm) | |
| 8492 | | 3/4" (19 mm) | 9" (23 cm) | |
| 8493 | | 7/8" (22 mm) | 9" (23 cm) | |
| 8494 | | 1" (25 mm) | 9" (23 cm) | |
| 8495 | | 1-1/8" (29 mm) | 9" (23 cm) | |
| 8496 | | 1-1/4" (32 mm) | 9" (23 cm) | |
| 8497 | | 1-3/8" (35 mm) | 9" (23 cm) | |
| 8498 | | 1-1/2" (38 mm) | 9" (23 cm) | |

Adhesive (Tape) Border Barriers

- Feature a thin layer of adhesive that surrounds the hydrocolloid barrier.
- Designed to be gentle on the skin, the adhesive moves with the body to easily conform to skin irregularities, creases and folds.
- Application is simple, with no additional accessories required to gain an added sense of protection. Go with the skin barrier that comes with built-in confidence.

Drainable Pouches

Clamp Closure System



Colour: Transparent With Tape Border 

| Product # | Stoma Size | Pouch Length | Box of |
|-------------|----------------|--------------|--------|
| 3232 | 5/8" (15 mm) | 9" (23 cm) | 30 |
| 3238 | 7/8" (22 mm) | 9" (23 cm) | |
| 3233 | 1-1/8" (29mm) | 9" (23 cm) | |
| 3239 | 1-3/8" (35 mm) | 9" (23 cm) | |
| 3234 | 1-1/2" (38 mm) | 9" (23 cm) | |
| 3222 | 5/8" (15 mm) | 12" (30 cm) | 30 |
| 3228 | 7/8" (22 mm) | 12" (30 cm) | |
| 3223 | 1-1/8" (29mm) | 12" (30 cm) | |
| 3229 | 1-3/8" (35 mm) | 12" (30 cm) | |
| 3224 | 1-1/2" (38 mm) | 12" (30 cm) | |
| 3225 | 2" (51 mm) | 12" (30 cm) | 30 |
| 3226 | 2-1/2" (64 mm) | 12" (30 cm) | |
| 3278 | 7/8" (22 mm) | 16" (40 cm) | |
| 3273 | 1-1/8" (29mm) | 16" (40 cm) | |
| 3279 | 1-3/8" (35 mm) | 16" (40 cm) | |
| 3274 | 1-1/2" (38 mm) | 16" (40 cm) | 10 |
| 3275 | 2" (51 mm) | 16" (40 cm) | |
| 3276 | 2-1/2" (64 mm) | 16" (40 cm) | |
| 3608 | 7/8" (22 mm) | 12" (30 cm) | |
| 3603 | 1-1/8" (29mm) | 12" (30 cm) | |
| 3609 | 1-3/8" (35 mm) | 12" (30 cm) | |
| 3604 | 1-1/2" (38 mm) | 12" (30 cm) | |
| 3605 | 2" (51 mm) | 12" (30 cm) | |

Colour: Beige With Tape Border 

| Product # | Stoma Size | Pouch Length | Box of |
|-------------|----------------|--------------|--------|
| 3112 | 5/8" (15 mm) | 12" (30 cm) | 30 |
| 3119 | 1-3/8" (35 mm) | 12" (30 cm) | |

Colour: Transparent With Replacement Filter With Tape Border 

| Product # | Stoma Size | Pouch Length | Box of |
|-------------|--------------|--------------|--------|
| 3668 | 7/8" (22 mm) | 12" (30 cm) | 10 |

Closed Pouches



Colour: Transparent Without Tape Border 

| Product # | Stoma Size | Pouch Length | Box of |
|-------------|----------------|--------------|--------|
| 7164 | 1-1/2" (38 mm) | 9" (23 cm) | 30 |

Colour: Transparent With Integrated Filter With Tape Border 

| Product # | Stoma Size | Pouch Length | Box of |
|-------------|----------------|--------------|--------|
| 3328 | 7/8" (22 mm) | 9" (23 cm) | 30 |
| 3323 | 1-1/8" (29mm) | 9" (23 cm) | |
| 3329 | 1-3/8" (35 mm) | 9" (23 cm) | |
| 3324 | 1-1/2" (38 mm) | 9" (23 cm) | |
| 3325 | 2" (51 mm) | 9" (23 cm) | |
| 3326 | 2-1/2" (64 mm) | 9" (23 cm) | |

Urostomy Pouches



Colour: Transparent With Tape Border 

| Product # | Stoma Size | Pouch Length | Box of |
|-------------|----------------|--------------|--------|
| 1432 | 5/8" (15 mm) | 9" (23 cm) | 10 |
| 1438 | 7/8" (22 mm) | 9" (23 cm) | |
| 1433 | 1-1/8" (29mm) | 9" (23 cm) | |
| 1439 | 1-3/8" (35 mm) | 9" (23 cm) | |
| 1434 | 1-1/2" (38 mm) | 9" (23 cm) | |

- To select the appropriate *Karaya* pouch, determine the stoma size and select the matching Product Code
- The labeling on the *Karaya* box and pouch you have selected indicates the gasket size, NOT the stoma size
- The *Karaya* seal extends into the gasket opening, ensuring that the Product Code you have selected fits the stoma size you have measured


Pouchkins

Pediatric Pouching System

One-Piece SoftFlex

Premie/Newborn Pouches



| Cut-to-fit | | | | Without Tape Border |  |
|------------|---------------------------------|-------------------------------|------------------|---------------------|---|
| Product # | Description | For Stomas Up To | Pouch Length | Box of | |
| Premie | | | | | |
| 3777 | Ultra-clear, no starter hole | 5/8" (15 mm) | 4-1/2" (11.5 cm) | 15 | |
| Newborn | | | | | |
| 3778 | Ultra-clear, no starter hole | 1-3/8" (35 x 22 x 7/8" mm) | 6-1/2" (16 cm) | | |

Drainable Pouches Lock 'n Roll Closure System



| Cut-to-fit Without Tape Border | | | | |
|--------------------------------|------------------------------|------------------|----------------|--------|
| Product # | Description | For Stomas Up To | Pouch Length | Box of |
| 3795 | Transparent, no starter hole | 2" (51 mm) | 7" (18 cm) | 10 |
| 3796 | | 1-1/2" (38 mm) | 6-1/2" (16 cm) | |

Urostomy Pouch



| Cut-to-fit Without Tape Border | | | | |
|--------------------------------|------------------------------|------------------|----------------|--------|
| Product # | Description | For Stomas Up To | Pouch Length | Box of |
| 3797 | Transparent, no starter hole | 1-1/2" (38 mm) | 6-1/2" (16 cm) | 10 |

Pouchkins SoftFlex Skin Barrier:

- Designed to be gentle on sensitive infant skin
- No starter hole, which enables the opening to be placed where needed
- No tape by design
- Flexible and sized for babies and children
- None of the *Pouchkins* pediatric pouching systems are made with natural rubber latex
- Any 1-3/4" *New Image* pouch will fit this skin barrier

Two-Piece SoftFlex

Flat Skin Barrier



| Cut-to-fit Without Tape Border | | | | |
|--------------------------------|-----------------|----------------|------------------|--------|
| Product # | Description | Flange Size | For Stomas Up To | Box of |
| 3761 | No starter hole | 1-3/4" (44 mm) | 1-1/4" (32 mm) | 5 |

Drainable Pouches Lock 'n Roll Closure System



Colour: Transparent

| Product # | Flange Size | Pouch Length | Filter | Box of |
|--------------|----------------|----------------|--------|--------|
| 3799 | 1-3/4" (44 mm) | 6-1/2" (16 cm) | No | 10 |
| 18292 | 1-3/4" (44 mm) | 7" (18 cm) | AF300 | 20 |

Urostomy Pouches



| Product # | Flange Size | Pouch Length | Box of |
|--------------|----------------|--------------|--------|
| 18402 | 1-3/4" (44 mm) | 9" (23 cm) | 10 |

Enhanced Design

| Product # | Colour | Flange Size | Pouch Length | Box of |
|--------------|-------------|----------------|--------------|--------|
| 18922 | Ultra-clear | 1-3/4" (44 mm) | 9" (23 cm) | 10 |
| 18912 | Beige | 1-3/4" (44 mm) | 9" (23 cm) | 10 |

Ostomy Belt



| Product # | Length | Box of |
|-------------|-----------------------|--------|
| 3774 | 10 - 17" (25 - 43 cm) | 1 |

Adapt Barrier Rings



Flat CeraRing



| Product # | Type | Outer Diameter | Thickness | Box of |
|-------------|----------|----------------|-----------|--------|
| 8805 | Standard | 2" (48 mm) | 4.5 mm | 10 |
| 8815 | Slim | 2" (48 mm) | 2.30 mm | |

Convex CeraRing



| Product # | Inner Diameter | Stretch Up To | Box of |
|--------------|-----------------|----------------|--------|
| 89520 | 13/16" (20 mm) | 1" (25 mm) | 10 |
| 89530 | 1-3/16" (30 mm) | 1-3/8" (35 mm) | |
| 89540 | 1-9/16" (40mm) | 1-3/4" (45 mm) | |

Oval Convex CeraRing



| Product # | Inner Diameter | Stretch Up To | Box of |
|--------------|----------------------------------|---------------------------------|--------|
| 89601 | 7/8" x 1-1/2" (22 x 38 mm) | 1-1/8" x 1-3/4" (27 x 43 mm) | 10 |
| 89602 | 1-3/16" x 1-7/8" (30 x 48 mm) | 1-3/8" x 2-1/8" (35 x 53 mm) | |
| 89603 | 1-1/2" x 2-3/16" (38 x 56 mm) | 1-3/4" x 2-3/8" (43 x 61 mm) | |

Flat Flextend



| Product # | Type | Outer Diameter | Thickness | Box of |
|-------------|----------|----------------|-----------|--------|
| 7805 | Standard | 2" (48 mm) | 4.5 mm | 10 |
| 7806 | Standard | 4" (102 mm) | 4.5 mm | |
| 7815 | Slim | 2" (48 mm) | 2.30 mm | |

Convex Flextend



| Product # | Inner Diameter | Stretch Up To | Box of |
|--------------|-----------------|----------------|--------|
| 79520 | 13/16" (20 mm) | 1" (25 mm) | 10 |
| 79530 | 1-3/16" (30 mm) | 1-3/8" (35 mm) | |
| 79540 | 1-9/16" (40 mm) | 1-3/4" (45 mm) | |

Oval Convex Flextend




| Product # | Inner Diameter | Stretch Up To | Box of |
|--------------|----------------------------------|---------------------------------|--------|
| 79601 | 7/8" x 1-1/2" (22 x 38 mm) | 1-1/8" x 1-3/4" (27 x 43 mm) | 10 |
| 79602 | 1-3/16" x 1-7/8" (30 x 48 mm) | 1-3/8" x 2-1/8" (35 x 53 mm) | |
| 79603 | 1-1/2" x 2-3/16" (38 x 56 mm) | 1-3/4" x 2-3/8" (43 x 61 mm) | |

Adapt CeraRing

Sealed with the confidence of ceramide










The security + confidence of Adapt meets the science of ceramide

- Available in Flat and Convex for each individual need
- Helps prevent leakage and protect skin—and may even extend barrier wear time
- Resists erosion—even when wet or exposed to discharge from the stoma
- Offers ultimate comfort and ease of use—alcohol-free rings are designed to be gentle for all skin types
- Extended wear
- Infused with ceramide — a natural component of the skin
- Help protect the skin's natural moisture barrier
- Help decrease transepidermal water loss (TEWL) from damaged or eroded skin




Accessories

Adapt Accessories





| Product # | | Description | Usage | Size | | Box of |
|-----------|---|--|--|---|-------------------|--------|
| 7760 |  | Universal Remover Wipes <i>*featuring no-sting silicone-based, alcohol free formulation</i> | Silicone-based and alcohol-free, this skin friendly product helps remove adhesive residues. | 50 Wipes per Carton | | 50 |
| 7737 |  | Adapt Medical Adhesive Remover Spray <i>*featuring no-sting, silicone-based, alcohol free formulation</i> | 360° spray technology to deliver controlled, quiet spraying from any angle, with minimal waste | Spray Can | 1.7 fl oz (50 ml) | 1 |
| 7906 |  | Adapt Stoma Powder | Use to absorb moisture from broken skin around the stoma. Dust on. Brush off excess powder. Stop using when skin heals. | Puff bottle | 1 oz (28 g) | 1 |
| 79402 |  | Adapt Barrier Extenders | Helps provide extra security around the edges of skin barrier, reduces skin barrier edge lifting and may help increase skin barrier wear time. Provides extra adhesion to accommodate a peristomal hernia. | Each barrier extender covers a surface area of 48 cm ² | | 30 |
| 78500 |  | Adapt Lubricating Deodorant | Neutralizes odours, eases emptying, and helps prevent pouch static and sticking. It is not a perfume or masking agent. Should be used in drainable and closed pouching systems only. | Bottle | 8 oz (236 ml) | 1 |
| 78501 | | | | Packets | .27 oz (8 ml) | 50 |
| 79301 |  | Adapt Paste | Used to fill or caulk uneven skin contours to create a flatter surface. It can help to prevent drainage from getting under the ostomy skin barrier. | Tube | .5 oz (14 g) | 20 |
| 79300 | | | | | 2 oz (60 g) | 1 |
| 79400 |  | Adapt Barrier Strips with FlexTend Barrier (Alcohol Free) | Used to fill in or caulk uneven skin contours near the stoma to create a flatter surface and help prevent ostomy drainage from getting under the ostomy barrier. | Triangular Shape | 5-1/8" (13 cm) | 10 |
| 7730 |  | Medical Adhesive Spray | Improves adhesive contact between skin and barrier Not recommended for use with extended wear barriers | Spray Can <i>Contains no CFCs</i> | 3.8 oz (112 ml) | 1 |
| 7917 |  | Adapt Skin Protective Wipe <i>*featuring no-sting silicone-based, alcohol free formulation</i> | Protects skin during adhesive removal Not recommended for use with extended wear skin barriers | 50 Wipes per Carton | | 50 |

Accessories

m9 Odour Eliminators

| Product # | | Description | Size | Box of |
|-------------|---|-------------------|---------------|--------|
| 7715 |  | Drops | 1 oz (30 ml) | 12 |
| 7717 | | | 8 oz (240 ml) | 6 |
| 7732 |  | Spray Unscented | 2 oz (60 ml) | 12 |
| 7733 | | | 8 oz (240 ml) | 6 |
| 7735 | | Spray Apple Scent | 8 oz (240 ml) | |

Pouch Accessories

| Product # | | Description | Size | Box of |
|-------------|---|---|--|--------|
| 8770 |  | Drainable Pouch Clamp, Beige | | 20 |
| 7300 |  | Beige Colour Ostomy Belt | Medium 23"- 43" (58 - 109 cm) | 10 |
| 7299 | | | Large 34"- 65" (86 - 165 cm) | |
| 7766 |  | Replacement Filter | 20 filters and 2 plugs/bag (product series 366x, 383x, and 389x) | 5 bags |
| 7331 |  | Enhanced Design Urostomy Drain Tube Adapter | | 10 |

Karaya Accessories

| Product # | | Description | Size | Box of |
|-------------|---|-------------|---------------------------|--------|
| 7910 |  | Paste | Tube 4.5 oz (128 g) | 12 |
| 7905 | | | Puff Bottle 2.5 oz (71 g) | |

Skin Sheets



| Product # | Description | Size | Box of |
|-------------|---------------------------|-------------------------|--------|
| 8800 | Flexextend (Non-sterile) | 4" x 4" (10 cm x 10 cm) | 5 |
| 8801 | | 8" x 8" (20 cm x 20 cm) | 3 |
| 7800 | Premium (Non-sterile) | 4" x 4" (10 cm x 10 cm) | 5 |
| 7801 | | 8" x 8" (20 cm x 20 cm) | 3 |
| 7700 | HolliHesive (Non-sterile) | 4" x 4" (10 cm x 10 cm) | 5 |

Skin barrier sheets are designed to be the most resistant to discharge and is the barrier of choice for many urostomies and ileostomies. It may be cut for customized use in stoma care. Available in extended wear (Flexextend) and standard wear (Premium & HolliHesive).

Loop Ostomy Bridge



| Product # | Description | Length | Box of |
|-------------|-------------|----------------|--------|
| 7767 | Sterile | 3-1/2" (89 mm) | 20 |

Irrigation Products



| Product # | Description | Information | | Box of |
|-----------|--------------------------------------|---|--------------------|--------|
| 7721 | Cone Irrigator Kit | Includes: (1) Irrigator (1) Stoma Cone with Connector (1) Cleaning Brush | | 1 |
| 7723 | Stoma Cone | With Connector Replacement Unit | | 10 |
| 7728 | Irrigator Sleeve | With Attached Closure and Belt Tabs | 2" (51 mm) | 20 |
| 7724 | | | 3" (76 mm) | |
| 18152 | New Image Two-Piece Irrigator Sleeve | Flange: 1-3/4" (44 mm) | Pouch: 30" (76 cm) | 5 |
| 18153 | | Flange: 2-1/4" (57 mm) | Pouch: 30" (76 cm) | |
| 18154 | | Flange: 2-3/4" (70 mm) | Pouch: 30" (76 cm) | |
| 7740 | Stoma Lubricant | Bottle | 4 oz (118 ml) | 10 |

Accessories

Wound Drainage Collectors

With FlexWear Skin Barrier



| Product # | For Wounds Up To | | Box of |
|-------------|-------------------------|-------------|--------|
| 9701 | 3" (7.6 cm) | Sterile | 5 |
| 9702 | 3.75" (9.5 cm) | | 5 |
| 9703 | 4" x 8" (10 cm x 20 cm) | | 3 |
| 9773 | 3" (7.6 cm) | Non-sterile | 10 |
| 9776 | 3.75" (9.5 cm) | | 10 |
| 9778 | 4" x 8" (10 cm x 20 cm) | | 3 |

With Microporous Tape

| Product # | For Wounds Up To | | Box of |
|-------------|-------------------------|---------|--------|
| 9772 | 3" (7.6 cm) | Sterile | 10 |
| 9775 | 3.75" (9.5 cm) | | 10 |
| 9777 | 4" x 8" (10 cm x 20 cm) | | 3 |

Wound Drainage Collectors

For management of draining wounds, fistulas, or draining tube sites

- Available with and without FlexWear skin barriers
- Sterile and non-sterile drainable pouches
- Odour-barrier film
- Resealable drainage cap and large access window
- Without starter hole

Access Port




| Product # | Description | | Box of |
|-------------|--|---------|--------|
| 9779 | Graduated nipple and punch allows introduction of catheter or drain through Wound Drainage Collector or pouch system | Sterile | 8 |

Stoma Caps


with Microporous Tape



| Pre-sized | | | | | With Tape Border  |
|-------------|--------|------------------|----------------|--------|--|
| Product # | Colour | For Stomas Up To | Pouch Length | Filter | Box of |
| 3184 | White | 2" (51 mm) | 4-1/4" (11 cm) | Yes | 30 |
| 3186 | | 3" (76 mm) | 4-1/4" (11 cm) | | |

SoftFlex Skin Barrier



| Pre-sized | | | | | Without Tape Border  |
|-------------|--------|------------------|--------------|--------|---|
| Product # | Colour | For Stomas Up To | Pouch Length | Filter | Box of |
| 1796 | Beige | 1-5/16" (50 mm) | 4" (10 cm) | Yes | 30 |



Hollister Secure Start Services

Hollister Secure Start Services provides a lifetime of customized support for people in the ostomy and continence care communities.

We are proud to offer dedicated support for each and every patient along the continuum of care! Once enrolled, a Hollister Secure Start Services Coordinator will call to explain our services, answer any questions, and provide ongoing support.

A Hollister Secure Start Services Coordinator can assist with the following:

- Finding the right products
- Providing product information and condition-specific education
- Finding local resources like retailers and clinicians
- Product usage care tips



To learn more about Hollister Secure Start Services, for Ostomy, call us at **1.866.789.7574**, or email us at **securestartcanada@hollister.com**.

For Continence, call us at **1.800.263.7421**, or email us at **canadasscccontact@hollister.com**.

Or visit us at **www.hollister.com**.



Continence Care

A woman with dark hair tied back is seated in a blue wheelchair, positioned at a white desk. She is holding a paintbrush and working on a canvas mounted on a wooden easel. The desk is cluttered with various art supplies, including paint tubes and brushes. A large, adjustable desk lamp with a silver shade is positioned over the desk. In the background, a large window offers a view of a city skyline. To the left of the desk, there is a dark-colored sofa and a small potted plant on a windowsill. The floor is covered with a red, shaggy rug. The overall scene is brightly lit by natural light from the window.

People are at the center of everything we do.

Discretion by design

Introducing a discreet catheter with design that is truly inspired.

The Infyna Chic™ Intermittent Catheter was designed with input from clinicians to help provide a high level of discretion for women who use catheters. It is beautiful and easy to use. The Infyna Chic catheter can help a woman feel better about having to use one.

Infyna Chic catheter

| Product # | Description | Size | Length | Box of |
|-------------|--------------|-------|------------|--------|
| 7008 | Straight tip | 8 Fr | 5" (13 cm) | 30 |
| 7010 | Straight tip | 10 Fr | 5" (13 cm) | 30 |
| 7012 | Straight tip | 12 Fr | 5" (13 cm) | 30 |
| 7014 | Straight tip | 14 Fr | 5" (13 cm) | 30 |



Features that help change her intermittent catheter experience

- Flip-top cap can be opened and closed with one hand
- Ready to use, hydrophilic catheter
- Just right stiffness helps support easy, touch-free insertion
- Catheter length helps give confidence that the bladder is fully drained
- Can attach a collection bag without removing the catheter from the case



Easy to open with just one hand



Easy to grip and remove from case



Does not leak when reclosed

LEARN MORE





Another reason to choose VaPro™ No Touch F-Style Intermittent Catheter

The VaPro F-Style catheter, a firmer-style catheter, offers patients the freedom to choose more firmness and the extra protection that the portfolio of VaPro intermittent catheters delivers.

VaPro™
No Touch Intermittent Catheter



Continence Care



Touch Free Hydrophilic Intermittent Catheters

VaPro Catheters Provide 100% No Touch Protection with a Protective Tip and Sleeve:

VaPro catheters are hygienic by design to help protect patients from germs throughout the entire catheterization process.

Protective tip:

- helps shield the sterile catheter during insertion from bacteria located within the first 15 mm of the distal urethra
- helps reduce the risk of carrying bacteria into the urinary tract

Protective sleeve:

- allows for catheter to be gripped anywhere
- provides a barrier to the catheter to help keep germs away



VaPro™



| Product # | Description | Size | Length | Box of |
|--------------|--------------|-------|-------------|--------|
| 72124 | Straight tip | 12 Fr | 16" (40 cm) | 30 |
| 73124 | Coudé tip | | 16" (40 cm) | |
| 72142 | Straight tip | 14 Fr | 8" (20 cm) | 30 |
| 72144 | Straight tip | | 16" (40 cm) | |
| 73144 | Coudé tip | 16 Fr | 16" (40 cm) | 30 |
| 72164 | Straight tip | | 16" (40 cm) | |
| 73164 | Coudé tip | | 16" (40 cm) | |

VaPro Plus™



| Product # | Description | Size | Length | Box of |
|--------------|--------------|-------|-------------|--------|
| 74124 | Straight tip | 12 Fr | 16" (40 cm) | 30 |
| 74144 | Straight tip | 14 Fr | 16" (40 cm) | 30 |

VaPro Plus is a closed system with an integrated 1000 ml collection bag.

VaPro F-Style is a **FIRMER** Touch Free Hydrophilic Catheter.

VaPro™ F-Style



| Product # | Description | Size | Length | Box of |
|----------------|--------------|-------|-------------|--------|
| 7600124 | Straight tip | 12 Fr | 16" (40 cm) | 30 |
| 7600144 | Straight tip | 14 Fr | 16" (40 cm) | 30 |
| 7600164 | Straight tip | 16 Fr | 16" (40 cm) | 30 |

VaPro Plus™ F-Style



| Product # | Description | Size | Length | Box of |
|----------------|--------------|-------|-------------|--------|
| 7700124 | Straight tip | 12 Fr | 16" (40 cm) | 30 |
| 7700144 | Straight tip | 14 Fr | 16" (40 cm) | 30 |

VaPro Plus F-Style is a closed system with an integrated 1000 ml collection bag.

Improved! A 50%* smaller protective tip,
the same 100% No Touch Protection.



Original

New



New Artwork!



New design features:

- **50%* smaller** and more ergonomically shaped protective tip, for a more comfortable insertion.
- **Transparent protective tip** to aid visualization
- **More discreet** packaging design
- **Less waste** as a result of new design

*As measured by volume of material used.



Easy-to open, easy-to-empty
collection bag for cleaner
handling

(applies to the VaPro Plus Pocket
catheters only)

VaPro Pocket



| Product # | Description | Size | Length | Box of |
|--------------|--------------|-------|-------------|--------|
| 70124 | Straight tip | 12 Fr | 16" (40 cm) | 30 |
| 70144 | Straight tip | 14 Fr | 16" (40 cm) | 30 |
| 70164 | Straight tip | 16 Fr | 16" (40 cm) | 30 |

VaPro Plus Pocket

| Product # | Description | Size | Length | Box of |
|--------------|--------------|-------|-------------|--------|
| 71084 | Straight tip | 8 Fr | 16" (40 cm) | 30 |
| 71104 | Straight tip | 10 Fr | 16" (40 cm) | 30 |
| 71122 | Straight tip | 12 Fr | 8" (20 cm) | 30 |
| 71124 | Straight tip | | 16" (40 cm) | |
| 71142 | Straight tip | 14 Fr | 8" (20 cm) | 30 |
| 71144 | Straight tip | | 16" (40 cm) | |
| 71164 | Straight tip | 16 Fr | 16" (40 cm) | 30 |

*VaPro Plus Pocket is a closed system with an integrated
1000 ml collection bag.*

Touch Free Intermittent Catheters

Closed System



| Product # | Description | Size | Length | Box of |
|----------------|--------------|-------|-------------|--------|
| 94064IC | Straight tip | 6 Fr | 16" (40 cm) | 100 |
| 94084IC | Straight tip | 8 Fr | 16" (40 cm) | |
| 94104IC | Straight tip | 10 Fr | 16" (40 cm) | |
| 94124IC | Straight tip | 12 Fr | 16" (40 cm) | |
| 95124IC | Coudé tip | | 16" (40 cm) | |
| 94144IC | Straight tip | 14 Fr | 16" (40 cm) | |
| 95144IC | Coudé tip | | 16" (40 cm) | |
| 94164IC | Straight tip | 16 Fr | 16" (40 cm) | |
| 95164IC | Coudé tip | | 16" (40 cm) | |
| 94184IC | Straight tip | 18 Fr | 16" (40 cm) | |

Closed System Kit




| Product # | Description | Size | Length | Box of |
|----------------|--------------|-------|-------------|--------|
| 96064IC | Straight tip | 6 Fr | 16" (40 cm) | 100 |
| 96084IC | Straight tip | 8 Fr | 16" (40 cm) | |
| 96104IC | Straight tip | 10 Fr | 16" (40 cm) | |
| 96124IC | Straight tip | 12 Fr | 16" (40 cm) | |
| 97124 | Coudé tip | | 16" (40 cm) | |
| 96144IC | Straight tip | 14 Fr | 16" (40 cm) | |
| 97144 | Coudé tip | | 16" (40 cm) | |
| 96164IC | Straight tip | 16 Fr | 16" (40 cm) | |
| 97164 | Coudé tip | | 16" (40 cm) | |
| 96184IC | Straight tip | 18 Fr | 16" (40 cm) | |

Kit Includes:

- Advance Plus Closed System
- Gloves
- Benzalkonium Chloride antiseptic swab (BZK)
- Waterproof underpad

Advance Plus Touch Free Intermittent Catheters:

- Have a protective tip that helps shield the catheter from bacteria located in the first 15 mm of the distal urethra
- Have a user-regulated gel reservoir which allows for control of the catheter lubrication
- Smooth eyelets help provide trouble-free insertion and withdrawal

A photograph of a young man with a beard, wearing a blue t-shirt, blue jeans, and a blue baseball cap, sitting in a wheelchair. He is looking off to the side with a slight smile. A black and white dog is sitting next to him, looking up at him. The background is a blurred outdoor setting with trees and grass.

*"It's not possible for anyone to live life
and never have any problems. But for every
problem, there is a solution. A person just
has to go out and find it."*

Trusted Quality. Reliable Performance. Cost Effective.

Apogee IC intermittent catheters feature smooth eyelets that are designed to aid comfortable insertion and removal. Non lubricated and not made with natural rubber latex.



Smooth eyelets aid comfortable insertion and removal

Apogee IC
Intermittent Catheter

Apogee Family of Intermittent Catheters



| New Product # | Existing Product # | Description | Size | Length | Box of |
|------------------|--------------------|------------------------------------|-------|-------------|--------|
| 10616 | B1007 | Straight tip | 6 Fr | 15" (38 cm) | 30 |
| 10806 | B1016 | Straight tip | 8 Fr | 6" (15 cm) | 30 |
| 10816 | B1018 | Straight tip | | 16" (40 cm) | |
| 10826 | New | Coudé | 8 Fr | 16" (40 cm) | |
| 11006 | B1033 | Straight tip | 10 Fr | 6" (15 cm) | 30 |
| 11026 | 1035 | Coudé, Olive tip | | 16" (40 cm) | |
| 11016 | B1034 | Straight tip | | 16" (40 cm) | |
| 11206 | B1047 | Straight tip | 12 Fr | 6" (15 cm) | 30 |
| 11216 | B1048 | Straight tip | | 16" (40 cm) | |
| 11226 | B1050 | Coudé tip | 14 Fr | 16" (40 cm) | 30 |
| 11406 | B1062 | Straight tip | | 6" (15 cm) | |
| 11446 | B1063 | Straight tip, No funnel | | 6" (15 cm) | |
| 11416 | B1064 | Straight tip | | 16" (40 cm) | |
| no change | B1061 | Soft tip | | 16" (40 cm) | |
| 11426 | B1066 | Coudé tip | 16 Fr | 16" (40 cm) | 50 |
| no change | B1065 | Straight tip with curved packaging | | 16" (40 cm) | |
| 11616 | B1082 | Straight tip | 16 Fr | 16" (40 cm) | 30 |
| 11626 | B1086 | Coudé tip | | 16" (40 cm) | |
| 11816 | B1096 | Straight tip | 18 Fr | 16" (40 cm) | 30 |
| 11826 | New | Coudé | 18 Fr | 16" (40 cm) | |

New!

10" (25 cm) Apogee Intermittent Catheters

| New Product # | Description | Size | Length | Box of |
|---------------|--------------|-------|-------------|--------|
| 10610 | Straight tip | 6 Fr | 10" (25 cm) | 30 |
| 10810 | Straight tip | 8 Fr | 10" (25 cm) | |
| 11010 | Straight tip | 10 Fr | 10" (25 cm) | |

*Apogee products will be changing to a new number throughout 2021. This change does not affect the product, its performance or intended use in any way. If you have any questions regarding this change please contact your Hollister representative or Hollister's Customer Service Department at **1-800-263-7400**.

Features of the Apogee Family of Intermittent Catheters:

- Two staggered smooth eyelets help provide trouble-free insertion
- Not made with natural rubber latex
- Designed for comfort
- Medical grade Polyvinyl Chloride (PVC)
- Non lubricated, clear catheter
- Straight and coudé tip options available

Male External Catheters

Hollister

Latex



| Extended Wear with Unique Inner Flap | | | Sheath length 2-3/4" (7 cm), adhesive length 1-3/5" (4 cm) |
|---|-----------------|------------------------------|---|
| Product # | Description | Size | Box of |
| 9206 | Small/Pediatric | 7/8" - 1" (22 - 25 mm) | 30 |
| 9207 | Medium | 1" - 1-1/8" (26 - 30 mm) | |
| 9208 | Intermediate | 1-1/8" - 1-1/4" (31 - 35 mm) | |
| 9209 | Large | 1-3/8" - 1-1/2" (36 - 39 mm) | |
| Everyday Wear | | | Sheath length 2-3/4" (7 cm), adhesive length 1-3/5" (4 cm) |
| Product # | Description | Size | Box of |
| 9106 | Small/Pediatric | 7/8" - 1" (22 - 25 mm) | 30 |
| 9107 | Medium | 1" - 1-1/8" (26 - 30 mm) | |
| 9108 | Intermediate | 1-1/8" - 1-1/4" (31 - 35 mm) | |
| 9109 | Large | 1-3/8" - 1-1/2" (36 - 39 mm) | |

Contains Natural Rubber Latex

Size Selection Guide

Use the Hollister Sizing Guide to determine the most appropriate size catheter

- Place sizing guide just behind the glans of the penis
- Select the size that most closely fits



InView

Silicone



| Standard | | | Sheath length 2-3/4" (7 cm), adhesive length 1-3/5" (4 cm) |
|------------------|--------------|----------------|--|
| Product # | Description | Size | Box of |
| 97525 | Small | 1" (25 mm) | 30 |
| 97529 | Medium | 1-1/8" (29 mm) | |
| 97532 | Intermediate | 1-1/4" (32 mm) | |
| 97536 | Large | 1-1/2" (36 mm) | |
| 97541 | Extra Large | 1-5/8" (41 mm) | 100 |
| 97525-100 | Small | 1" (25 mm) | |
| 97529-100 | Medium | 1-1/8" (29 mm) | |
| 97532-100 | Intermediate | 1-1/4" (32 mm) | |
| 97536-100 | Large | 1-1/2" (36 mm) | |
| 97541-100 | Extra Large | 1-5/8" (41 mm) | |

| Special | | | Sheath length 1-3/5" (4 cm), adhesive length 1-3/5" (4 cm) |
|------------------|--------------|----------------|--|
| Product # | Description | Size | Box of |
| 97425 | Small | 1" (25 mm) | 30 |
| 97429 | Medium | 1-1/8" (29 mm) | |
| 97432 | Intermediate | 1-1/4" (32 mm) | |
| 97436 | Large | 1-1/2" (36 mm) | |
| 97441 | Extra Large | 1-5/8" (41 mm) | 100 |
| 97425-100 | Small | 1" (25 mm) | |
| 97429-100 | Medium | 1-1/8" (29 mm) | |
| 97432-100 | Intermediate | 1-1/4" (32 mm) | |
| 97436-100 | Large | 1-1/2" (36 mm) | |
| 97441-100 | Extra Large | 1-5/8" (41 mm) | |

| Extra | | | Sheath length 2-3/4" (7 cm), adhesive length 2-3/4" (7 cm) |
|------------------|--------------|----------------|--|
| Product # | Description | Size | Box of |
| 97625 | Small | 1" (25 mm) | 30 |
| 97629 | Medium | 1-1/8" (29 mm) | |
| 97632 | Intermediate | 1-1/4" (32 mm) | |
| 97636 | Large | 1-1/2" (36 mm) | |
| 97641 | Extra Large | 1-5/8" (41 mm) | 100 |
| 97625-100 | Small | 1" (25 mm) | |
| 97629-100 | Medium | 1-1/8" (29 mm) | |
| 97632-100 | Intermediate | 1-1/4" (32 mm) | |
| 97636-100 | Large | 1-1/2" (36 mm) | |
| 97641-100 | Extra Large | 1-5/8" (41 mm) | |

Not Made with Natural Rubber Latex

Urinary Collection Systems

Leg Bags

Leg Bag



| Product # | Description | Capacity | Box of |
|-------------|-------------|----------------|--------|
| 9814 | Medium | 18 oz (540 ml) | 10 |
| 9805 | Large | 30 oz (900 ml) | |

Combination Pack



| Product # | Description | Capacity | Box of |
|-------------|-------------|----------------|----------|
| 9348 | Medium | 18 oz (540 ml) | 12 packs |
| 9349 | Large | 30 oz (900 ml) | |

Extension Tubing with Connector




| Product # | Description | Length | Box of |
|-------------|-------------|-------------|--------|
| 9345 | Non-sterile | 18" (45 cm) | 10 |
| 9346 | Sterile | 18" (45 cm) | |

Vented Combination Pack



| Product # | Description | Capacity | Box of |
|-------------|-------------|----------------|----------|
| 9645 | Medium | 18 oz (540 ml) | 10 packs |
| 9655 | Large | 30 oz (900 ml) | |

Fabric Strap



| Product # | Description | Capacity | Box of |
|-------------|-------------|-------------|----------|
| 9342 | Small | 15" (38 cm) | 10 pairs |
| 9343 | Medium | 23" (58 cm) | |

Bedside Drainage



| Product # | Description | Capacity | Box of |
|-------------|---------------------------|-----------------|--------|
| 9839 | 2000 ml Collection System | 68 oz (2000 ml) | 20 |

Hollister Leg Bags Are:

- Sterile
- Not Made with Natural Rubber Latex
- Designed with an anti-reflux valve

m9 Cleaner/Decrystallizer



Concentrate & Wash Bottle

| Product # | Description | Size | Box of |
|-------------|---|----------------|--------|
| 7736 | 1 Bottle of Concentrate 1 Washing Bottle | 16 oz (480 ml) | 12 |

External Collection Devices

Urinary Pouches

Female

SoftFlex Skin Barrier

Adapt Paste, tape strips and Skin Gel Wipe included



| Product # | Description | Box of |
|-------------|---|--------|
| 9840 | Length 7-1/2" (19 cm) Width 4-1/2" (11cm) | 10 |

Male

Flextend Skin Barrier

Thin barrier for flexibility, no additional tape needed



| Product # | Description | Box of |
|-------------|---|--------|
| 9873 | Length 7-1/2" (19 cm) Width 4-1/2" (11cm) | 10 |

SoftFlex Skin Barrier

Tape strips and Skin Gel Wipe included



| Product # | Description | Box of |
|-------------|---|--------|
| 9811 | Length 7-1/2" (19 cm) Width 4-1/2" (11cm) | 10 |

Hollister's Urinary Pouches:

- Not Made with Natural Rubber Latex
- Have a cut-to-fit skin barrier
- Have a odour-barrier pouch film

Fecal Collectors

Pouches

Flextend Skin Barrier



| Product # | Description | Pouch Length | Box of |
|-------------|---------------------|--------------|--------|
| 9823 | Medium, Ultra-clear | 10" (25 cm) | 10 |
| 9880 | Large, Transparent | 12" (30 cm) | |

Transparent Pouches

SoftFlex Skin Barrier



| Product # | Description | Pouch Length | Box of |
|-------------|-------------|--------------|--------|
| 9822 | Medium | 10" (25 cm) | 10 |
| 9821 | Large | 12" (30 cm) | |



With 100% No Touch Protection,
whatever you touch
won't touch VaPro catheters.



VaPro

Touch Free Hydrophilic
Intermittent Catheter



Continence Care



Critical Care

Dedicated to excellence.

Oral Endotracheal Tube Fasteners

AnchorFast



| Product # | Description | Box of |
|-------------|--|--------|
| 9799 | <ul style="list-style-type: none"> ET tube wrap with one-click security clamp locks the ET tube in place Gliding tube shuttle allows the ET tube to be easily repositioned from side to side while the tube remains secure FlexTend Skin Barrier For tubes sized 5.0-10.0mm inner diameter | 12 |

AnchorFast SlimFit



| Product # | Description | Box of |
|-------------|---|--------|
| 9787 | <ul style="list-style-type: none"> ET tube wrap with one-click security Gliding tube shuttle allows for side-to-side repositioning of the ET tube Thumb pad with increased surface area for finger or thumb placement when opening Designed to accommodate smaller faces Secures ET tubes ranging in size 5.0-8.0 mm | 12 |

AnchorFast Guard Select



| Product # | Description | Box of |
|-------------|--|--------|
| 9700 | <ul style="list-style-type: none"> ET tube wrap with one-click security clamp locks the ET tube in place Gliding tube shuttle allows the ET tube to be easily repositioned from side to side while the tube remains secure Integrated tube protection sleeve designed to help prevent occlusion Secures oral ET tubes ranging in size from 5.0 to 9.0 mm inner diameter and ET tubes with integrated subglottic suction lumen sizes 6.0 to 8.0 mm inner diameter | 12 |

AnchorFast Replacement Neck Band Straps



| Product # | Description | Box of |
|-------------|------------------------------|--------|
| 9788 | Replacement Neck Band Straps | 4 |

Tube Attachment Devices

Horizontal Tube (HTAD)



| Product # | Description | Box of |
|-------------|---------------------|--------|
| 9781 | Single use, sterile | 5 |

Vertical Tube (VTAD)



| Product # | Description | Box of |
|-------------|---------------------|--------|
| 9782 | Single use, sterile | 5 |

Drain Tube (DTAD)



| Product # | Description | Box of |
|-------------|---------------------|--------|
| 9780 | Single use, sterile | 5 |

Feeding Tube (FTAD)



| Product # | Description | Box of |
|-------------|---|--------|
| 9786 | <ul style="list-style-type: none"> FlexWear skin barrier conforms to patient's nose Single use, sterile | 20 |

Suction Tube (STAD)



| Product # | Description | Box of |
|-------------|---|--------|
| 9785 | <ul style="list-style-type: none"> FlexWear skin barrier conforms to patient's nose Single use, sterile | 20 |

Numerical Index

| Product # | Page | Product # | Page | Product # | Page | Product # | Page | Product # | Page | Product # | Page |
|-----------|------|-----------|------|-----------|------|-----------|------|-----------|------|-----------|------|
| 1035 | 38 | 3795 | 25 | 8248 | 20 | 8484 | 23 | 8614 | 20 | 8954 | 19 |
| 1432 | 24 | 3796 | 25 | 8274 | 20 | 8485 | 23 | 8615 | 20 | 8958 | 19 |
| 1433 | 24 | 3797 | 25 | 8278 | 20 | 8486 | 23 | 8616 | 20 | 8960 | 19 |
| 1434 | 24 | 3799 | 25 | 8280 | 20 | 8487 | 23 | 8617 | 20 | 8961 | 19 |
| 1438 | 24 | 7008 | 32 | 8281 | 20 | 8488 | 23 | 8618 | 20 | 8962 | 19 |
| 1439 | 24 | 7010 | 32 | 8282 | 20 | 8489 | 23 | 8619 | 20 | 8963 | 19 |
| 1796 | 29 | 7012 | 32 | 8283 | 20 | 8490 | 23 | 8624 | 20 | 8984 | 19 |
| 3112 | 24 | 7014 | 32 | 8284 | 20 | 8491 | 23 | 8628 | 20 | 8988 | 19 |
| 3119 | 24 | 7164 | 24 | 8285 | 20 | 8492 | 23 | 8631 | 17 | 8990 | 19 |
| 3184 | 29 | 7299 | 28 | 8286 | 20 | 8493 | 23 | 8632 | 17 | 8991 | 19 |
| 3186 | 29 | 7300 | 28 | 8288 | 20 | 8494 | 23 | 8633 | 17 | 8992 | 19 |
| 3222 | 24 | 7331 | 28 | 8289 | 20 | 8495 | 23 | 8634 | 17 | 8993 | 19 |
| 3223 | 24 | 7700 | 28 | 8290 | 20 | 8496 | 23 | 8635 | 17 | 9106 | 39 |
| 3224 | 24 | 7715 | 28 | 8291 | 20 | 8497 | 23 | 8637 | 17 | 9107 | 39 |
| 3225 | 24 | 7717 | 28 | 8292 | 20 | 8498 | 23 | 8638 | 17 | 9108 | 39 |
| 3226 | 24 | 7721 | 28 | 8293 | 20 | 8504 | 20 | 8639 | 17 | 9109 | 39 |
| 3228 | 24 | 7723 | 28 | 8294 | 20 | 8508 | 20 | 8641 | 17 | 9206 | 39 |
| 3229 | 24 | 7724 | 28 | 8295 | 20 | 8510 | 20 | 8642 | 17 | 9207 | 39 |
| 3232 | 24 | 7728 | 28 | 8296 | 20 | 8511 | 20 | 8643 | 17 | 9208 | 39 |
| 3233 | 24 | 7730 | 27 | 8297 | 20 | 8512 | 20 | 8644 | 17 | 9209 | 39 |
| 3234 | 24 | 7732 | 28 | 8331 | 17 | 8513 | 20 | 8647 | 17 | 9342 | 40 |
| 3238 | 24 | 7733 | 28 | 8412 | 22 | 8514 | 20 | 8648 | 17 | 9343 | 40 |
| 3239 | 24 | 7735 | 28 | 8413 | 22 | 8515 | 20 | 8649 | 17 | 9345 | 40 |
| 3273 | 24 | 7736 | 40 | 8414 | 22 | 8524 | 20 | 8651 | 17 | 9346 | 40 |
| 3274 | 24 | 7737 | 27 | 8415 | 22 | 8528 | 20 | 8661 | 19 | 9348 | 40 |
| 3275 | 24 | 7740 | 28 | 8440 | 21 | 8531 | 17 | 8662 | 19 | 9349 | 40 |
| 3276 | 24 | 7760 | 27 | 8442 | 21 | 8538 | 17 | 8663 | 19 | 9645 | 40 |
| 3278 | 24 | 7766 | 28 | 8443 | 21 | 8551 | 17 | 8664 | 19 | 9655 | 40 |
| 3279 | 24 | 7767 | 28 | 8447 | 21 | 8552 | 17 | 8665 | 19 | 9700 | 44 |
| 3323 | 24 | 7800 | 28 | 8448 | 21 | 8558 | 17 | 8666 | 19 | 9701 | 29 |
| 3324 | 24 | 7801 | 28 | 8449 | 21 | 8559 | 17 | 8674 | 19 | 9702 | 29 |
| 3325 | 24 | 7805 | 26 | 8450 | 21 | 8584 | 20 | 8678 | 19 | 9703 | 29 |
| 3326 | 24 | 7806 | 26 | 8460 | 21 | 8588 | 20 | 8770 | 28 | 9772 | 29 |
| 3328 | 24 | 7815 | 26 | 8462 | 21 | 8590 | 20 | 8800 | 28 | 9773 | 29 |
| 3329 | 24 | 7905 | 28 | 8463 | 21 | 8591 | 20 | 8801 | 28 | 9775 | 29 |
| 3603 | 24 | 7906 | 27 | 8464 | 21 | 8592 | 20 | 8805 | 26 | 9776 | 29 |
| 3604 | 24 | 7910 | 28 | 8467 | 21 | 8593 | 20 | 8815 | 26 | 9777 | 29 |
| 3605 | 24 | 7917 | 27 | 8468 | 21 | 8594 | 20 | 8901 | 17 | 9778 | 29 |
| 3608 | 24 | 8031 | 17 | 8469 | 21 | 8595 | 20 | 8914 | 19 | 9779 | 29 |
| 3609 | 24 | 8081 | 17 | 8474 | 23 | 8596 | 20 | 8918 | 19 | 9780 | 44 |
| 3668 | 24 | 8131 | 17 | 8478 | 23 | 8597 | 20 | 8925 | 17 | 9781 | 44 |
| 3761 | 25 | 8181 | 17 | 8480 | 23 | 8610 | 20 | 8930 | 17 | 9782 | 44 |
| 3774 | 25 | 8214 | 20 | 8481 | 23 | 8611 | 20 | 8931 | 17 | 9785 | 44 |
| 3777 | 25 | 8218 | 20 | 8482 | 23 | 8612 | 20 | 8935 | 17 | 9786 | 44 |
| 3778 | 25 | 8244 | 20 | 8483 | 23 | 8613 | 20 | 8940 | 17 | 9787 | 44 |

| Product # | Page | Product # | Page | Product # | Page | Product # | Page | Product # | Page | Product # | Page |
|-----------|------|-----------|------|-----------|------|-----------|------|-----------|------|-----------|------|
| 9788 | 44 | 11826 | 38 | 14506 | 10 | 15603 | 7 | 18142 | 11 | 18402 | 25 |
| 9799 | 44 | 11903 | 9 | 14602 | 7 | 15604 | 7 | 18143 | 11 | 18403 | 13 |
| 9805 | 40 | 11904 | 9 | 14603 | 7 | 15802 | 10 | 18144 | 11 | 18404 | 13 |
| 9811 | 41 | 11905 | 9 | 14604 | 7 | 15803 | 10 | 18152 | 28 | 18412 | 13 |
| 9814 | 40 | 11906 | 9 | 14606 | 7 | 15804 | 10 | 18153 | 28 | 18413 | 13 |
| 9821 | 41 | 13402 | 10 | 14701 | 7 | 15903 | 10 | 18154 | 28 | 18414 | 13 |
| 9822 | 41 | 13403 | 10 | 14702 | 7 | 15904 | 10 | 18172 | 11 | 18422 | 13 |
| 9823 | 41 | 13404 | 10 | 14703 | 7 | 15905 | 10 | 18173 | 11 | 18423 | 13 |
| 9839 | 40 | 13503 | 10 | 14704 | 7 | 15906 | 10 | 18174 | 11 | 18424 | 13 |
| 9840 | 41 | 13504 | 10 | 14705 | 7 | 16101 | 7 | 18176 | 11 | 18902 | 13 |
| 9873 | 41 | 13505 | 10 | 14706 | 7 | 16102 | 7 | 18182 | 11 | 18903 | 13 |
| 9880 | 41 | 13506 | 10 | 14707 | 7 | 16103 | 7 | 18183 | 11 | 18904 | 13 |
| 10610 | 38 | 13901 | 10 | 14708 | 7 | 16104 | 7 | 18184 | 11 | 18912 | 13 |
| 10616 | 38 | 13902 | 10 | 14709 | 7 | 16105 | 7 | 18192 | 11 | 18912 | 25 |
| 10806 | 38 | 13903 | 10 | 14802 | 10 | 16106 | 7 | 18193 | 11 | 18913 | 13 |
| 10810 | 38 | 13904 | 10 | 14803 | 10 | 16107 | 7 | 18194 | 11 | 18914 | 13 |
| 10816 | 38 | 13905 | 10 | 14804 | 10 | 16108 | 7 | 18202 | 11 | 18922 | 13 |
| 10826 | 38 | 13906 | 10 | 14901 | 10 | 16401 | 8 | 18203 | 11 | 18922 | 25 |
| 11006 | 38 | 13907 | 10 | 14902 | 10 | 16402 | 8 | 18204 | 11 | 18923 | 13 |
| 11010 | 38 | 13908 | 10 | 14903 | 10 | 16403 | 8 | 18282 | 11 | 18924 | 13 |
| 11016 | 38 | 13909 | 10 | 14904 | 10 | 16404 | 8 | 18283 | 11 | 19002 | 14 |
| 11026 | 38 | 13910 | 10 | 14905 | 10 | 16405 | 8 | 18284 | 11 | 19003 | 14 |
| 11202 | 7 | 13911 | 10 | 14906 | 10 | 16406 | 8 | 18292 | 11 | 19004 | 14 |
| 11203 | 7 | 14102 | 8 | 14907 | 10 | 16407 | 8 | 18292 | 25 | 19006 | 14 |
| 11204 | 7 | 14103 | 8 | 14908 | 10 | 16408 | 8 | 18293 | 11 | 19052 | 14 |
| 11206 | 38 | 14104 | 8 | 14909 | 10 | 18002 | 11 | 18294 | 11 | 19053 | 14 |
| 11216 | 38 | 14202 | 8 | 14910 | 10 | 18003 | 11 | 18322 | 12 | 19054 | 14 |
| 11226 | 38 | 14203 | 8 | 14911 | 10 | 18004 | 11 | 18323 | 12 | 19056 | 14 |
| 11402 | 9 | 14204 | 8 | 15102 | 7 | 18006 | 11 | 18324 | 12 | 19102 | 14 |
| 11403 | 9 | 14206 | 8 | 15103 | 7 | 18012 | 12 | 18332 | 12 | 19103 | 14 |
| 11404 | 9 | 14301 | 8 | 15104 | 7 | 18013 | 12 | 18333 | 12 | 19104 | 14 |
| 11406 | 38 | 14302 | 8 | 15202 | 8 | 18014 | 12 | 18334 | 12 | 19106 | 14 |
| 11416 | 38 | 14303 | 8 | 15203 | 8 | 18016 | 12 | 18352 | 12 | 19152 | 14 |
| 11426 | 38 | 14304 | 8 | 15204 | 8 | 18102 | 11 | 18353 | 12 | 19153 | 14 |
| 11446 | 38 | 14305 | 8 | 15302 | 9 | 18103 | 11 | 18354 | 12 | 19154 | 14 |
| 11503 | 9 | 14306 | 8 | 15303 | 9 | 18104 | 11 | 18362 | 12 | 19156 | 14 |
| 11504 | 9 | 14307 | 8 | 15304 | 9 | 18112 | 11 | 18363 | 12 | 19202 | 14 |
| 11505 | 9 | 14308 | 8 | 15402 | 10 | 18113 | 11 | 18364 | 12 | 19203 | 14 |
| 11506 | 9 | 14309 | 8 | 15403 | 10 | 18114 | 11 | 18372-30 | 12 | 19204 | 14 |
| 11616 | 38 | 14402 | 10 | 15404 | 10 | 18122 | 11 | 18373-30 | 12 | 19252 | 14 |
| 11626 | 38 | 14403 | 10 | 15503 | 10 | 18123 | 11 | 18374-30 | 12 | 19253 | 14 |
| 11702 | 9 | 14404 | 10 | 15504 | 10 | 18124 | 11 | 18392 | 12 | 19254 | 14 |
| 11703 | 9 | 14503 | 10 | 15505 | 10 | 18132 | 11 | 18393 | 12 | 19602 | 14 |
| 11704 | 9 | 14504 | 10 | 15506 | 10 | 18133 | 11 | 18394 | 12 | 19603 | 14 |
| 11816 | 38 | 14505 | 10 | 15602 | 7 | 18134 | 11 | 18402 | 13 | 19604 | 14 |

Numerical Index

| Product # | Page | Product # | Page | Product # | Page | Product # | Page | Product # | Page |
|-----------|------|-----------|------|-----------|------|-----------|------|-----------|------|
| 19802 | 14 | 82300 | 18 | 84699 | 21 | 89111 | 19 | 97532-100 | 39 |
| 19803 | 14 | 82302 | 18 | 84711 | 23 | 89511 | 19 | 97536 | 39 |
| 19804 | 14 | 82325 | 18 | 84794 | 23 | 89520 | 26 | 97536-100 | 39 |
| 19902 | 14 | 82330 | 18 | 84798 | 23 | 89530 | 26 | 97541 | 39 |
| 19903 | 14 | 82335 | 18 | 84890 | 23 | 89540 | 26 | 97541-100 | 39 |
| 19904 | 14 | 82402 | 18 | 84891 | 23 | 89601 | 26 | 97625 | 39 |
| 70124 | 35 | 82411 | 20 | 84892 | 23 | 89602 | 26 | 97625-100 | 39 |
| 70144 | 35 | 82610 | 20 | 84893 | 23 | 89603 | 26 | 97629 | 39 |
| 70164 | 35 | 82611 | 20 | 84894 | 23 | 89811 | 19 | 97629-100 | 39 |
| 71084 | 35 | 82612 | 20 | 84895 | 23 | 94064IC | 36 | 97632 | 39 |
| 71104 | 35 | 82613 | 20 | 84896 | 23 | 94084IC | 36 | 97632-100 | 39 |
| 71122 | 35 | 82614 | 20 | 84897 | 23 | 94104IC | 36 | 97636 | 39 |
| 71124 | 35 | 82615 | 20 | 84898 | 23 | 94124IC | 36 | 97636-100 | 39 |
| 71142 | 35 | 82616 | 20 | 84899 | 23 | 94144IC | 36 | 97641 | 39 |
| 71144 | 35 | 82617 | 20 | 84990 | 23 | 94164IC | 36 | 97641-100 | 39 |
| 71164 | 35 | 82618 | 20 | 84991 | 23 | 94184IC | 36 | 826211 | 20 |
| 72124 | 34 | 82619 | 20 | 84992 | 23 | 95124IC | 36 | 841311 | 22 |
| 72142 | 34 | 82624 | 20 | 84993 | 23 | 95144IC | 36 | 844011 | 22 |
| 72144 | 34 | 82628 | 20 | 84994 | 23 | 95164IC | 36 | 7600124 | 34 |
| 72164 | 34 | 82711 | 20 | 84995 | 23 | 96064IC | 36 | 7600144 | 34 |
| 73124 | 34 | 84134 | 22 | 84996 | 23 | 96084IC | 36 | 7600164 | 34 |
| 73144 | 34 | 84138 | 22 | 84997 | 23 | 96104IC | 36 | 7700124 | 34 |
| 73164 | 34 | 84200 | 21 | 84998 | 23 | 96124IC | 36 | 7700144 | 34 |
| 74124 | 34 | 84292 | 22 | 84999 | 23 | 96144IC | 36 | 8439111 | 23 |
| 74144 | 34 | 84293 | 22 | 85011 | 20 | 96164IC | 36 | 8479111 | 23 |
| 78500 | 27 | 84294 | 22 | 85116 | 20 | 96184IC | 36 | B1007 | 38 |
| 78501 | 27 | 84295 | 22 | 85118 | 20 | 97124 | 36 | B1016 | 38 |
| 79300 | 27 | 84394 | 23 | 85119 | 20 | 97144 | 36 | B1018 | 38 |
| 79301 | 27 | 84398 | 23 | 85211 | 20 | 97164 | 36 | B1033 | 38 |
| 79400 | 27 | 84404 | 22 | 85811 | 20 | 97425 | 39 | B1034 | 38 |
| 79402 | 27 | 84408 | 22 | 86211 | 20 | 97425-100 | 39 | B1047 | 38 |
| 79520 | 26 | 84490 | 21 | 86711 | 19 | 97429 | 39 | B1048 | 38 |
| 79530 | 26 | 84492 | 21 | 88300 | 18 | 97429-100 | 39 | B1050 | 38 |
| 79540 | 26 | 84493 | 21 | 88302 | 18 | 97432 | 39 | B1061 | 38 |
| 79601 | 26 | 84494 | 21 | 88325 | 18 | 97432-100 | 39 | B1062 | 38 |
| 79602 | 26 | 84497 | 21 | 88330 | 18 | 97436 | 39 | B1063 | 38 |
| 79603 | 26 | 84498 | 21 | 88335 | 18 | 97436-100 | 39 | B1064 | 38 |
| 80070 | 18 | 84499 | 21 | 88340 | 18 | 97441 | 39 | B1065 | 38 |
| 80110 | 18 | 84590 | 21 | 88402 | 18 | 97441-100 | 39 | B1066 | 38 |
| 82100 | 18 | 84692 | 21 | 88700 | 18 | 97525 | 39 | B1082 | 38 |
| 82111 | 20 | 84693 | 21 | 88725 | 18 | 97525-100 | 39 | B1086 | 38 |
| 82125 | 18 | 84694 | 21 | 88730 | 18 | 97529 | 39 | B1096 | 38 |
| 82130 | 18 | 84697 | 21 | 88735 | 18 | 97529-100 | 39 | | |
| 82135 | 18 | 84698 | 21 | 88740 | 18 | 97532 | 39 | | |

www.hollister.com

The Hollister logo, Adapt, Advance, Advance Plus, Anchorfast, Anchorfast Slim Fit, Anchorfast Guard Select, Apogee, CeraPlus, Extended Wear, Everyday, Flexextend, FlexWear, FormaFlex, Infyna Chic, InView, Karaya 5, Lock 'n Roll, m9, New Image, Pouchkins, Premier, QuietWear, Secure Start, SoftFlex, VaPro, VaPro F-Style, VaPro Plus, VaPro Pocket and VaPro Plus Pocket are trademarks of Hollister Incorporated.

All other trademarks and copyrights are the property of their respective owners.

©2020 Hollister Incorporated



Contact Information

Customer Service

Toll Free: 1.800.263.7400

Local: 905.727.4344

Fax

Toll Free: 1.800.432.8846

Local: 905.727.1614

Email

canadacustomerservice@hollister.com

Business Hours

Monday – Friday; 8:30 am - 5:00 pm (EST) Excluding Holidays

Hollister Limited

Aurora, ON

Hollister Secure Start Services

Hollister Secure Start Services provides a lifetime of customized support for people in the ostomy and continence care communities.

We are proud to offer dedicated support for each and every patient along the continuum of care. Once enrolled, a Hollister Secure Start Services Coordinator will call to explain our services, answer any questions, and provide ongoing support.

A Hollister Secure Start Services Coordinator can assist with the following:

- Finding the right products
- Providing product information and condition-specific education
- Finding local resources like retailers and clinicians
- Product usage care tips
- Find reimbursement for your catheter products

Ostomy Care:

1.866.789.7574, or email us at
securestartcanada@hollister.com


Continence Care:

1.800.263.7421, or email us at
canadassccccontact@hollister.com



ENROLL HERE

We're here for you, no matter where 'here' is.



We know that the way clinicians connect with patients has seen much change in recent times. Through it all, Hollister remains committed to our mission of making life more rewarding and dignified for the people who use our products and services. We've kept our customer support programs accessible and approachable— even if remotely—and we're committed to continuing our focus on maintaining the supply of essential Ostomy and Continence Care products our customers need. You can rely on us, so they can rely on you.

SCAN ME



©2020 Hollister Incorporated