

Making a Difference in the Journey of Life



Sara Levitt
Hollister End-User and Ambassador,
Montreal, QC

2026 Product Catalogue



The fundamental values and principles
of our founder, John Dickinson Schneider,
form the foundation of our company today.

Our Mission

To make life more rewarding and
dignified for the people who use
our products and services.

Our Vision

To grow and prosper as an
independent, employee-owned
company and in the process
become better human beings.

Our Heritage

The roots of Hollister Incorporated date back to 1921, when 23-year-old entrepreneur John Dickinson Schneider founded a small printing company in Chicago, which he named JDS Printer Craftsman. From the very beginning, John Schneider's business philosophy was "only first class is good enough," and the company became known for its quality products, quality service, and quality people.

Table of Contents

Ostomy Care	5	Pouchkins™	37
CeraPlus™ Ostomy Products.....	9	Pediatric Pouching System.....	38
Two-Piece Barriers	10	Continence Care	39
Two-Piece Barriers + Kits	11	Infyna Chic™	40
One-Piece Pouching Systems	12	Touch Free Hydrophilic	
One-Piece Urostomy Pouches Systems	13	Intermittent Catheters	42
Two-Piece Pouches	15	VaPro™	42
Two-Piece Pouches.....	16	VaPro™ Pocket	43
High Output Pouches.....	17	Advance Plus™	44
Ostomy Kits	17	Intermittent Catheters	45
Urostomy Kits.....	17	Apogee™	45
Urostomy Pouches.....	18	Male External Catheters	46
Two-Piece Barriers.....	21	Hollister Latex.....	46
Flextend Barriers	22	InView™ Silicone	46
FlexWear Barriers	23	Urinary Collection Systems.....	47
One-Piece Pouching Systems	24	Leg Bags.....	47
Flextend Soft Convex	26	m9 Cleaner/Decrystallizer	47
Flextend Convex	26	External Collection Devices	
Flextend Flat.....	27	Urinary Pouches.....	48
FlexWear Flat.....	27	Fecal Collectors.....	48
SoftFlex Flat	28	Critical Care.....	50
Flextend High Output	29	Oral Endotracheal Tube Fasteners.....	51
Flextend Urostomy	30	Tube Attachment Devices.....	51
Accessories	33	Numerical Index	52
CeraPlus™ Accessories	34		
Adapt™	35		
m9 Odour Eliminators.....	36		
Pouch Accessories.....	36		
Skin Sheets.....	36		
High Output Collection Bag	36		
Irrigation Products	36		
Wound Drainage Collectors.....	36		

Hollister Customer Service Department

Business Hours: Monday to Friday from 8:30 a.m. to 5:00 p.m. (Eastern Standard Time)

Telephone: 1-800-263-7400

Fax: 1-800-432-8846

Email: Canadacustomerservice@hollister.com




Helping your patients live life on **their own terms**

You spoke, we listened. CeraPlus™ Pouches just got even better^{1,2,3}. Because we're passionate about helping you improve the quality of life for people living with a stoma.**

With the upgraded pouches, we have developed our best product combination. A wide choice of configurations to achieve a secure and comfortable **Fit** from day one, protecting from leakage. A moisturizing ceramide **Formulation** that helps maintain healthy skin around the stoma, minimizing transepidermal water loss, dryness and itching⁴. And a high-quality **Pouch** with improvements made to the closure, fabric and aesthetics.

We love that feeling when everything comes together. It's what makes CeraPlus™ Products' different. The outcome? A better quality of life for your patients, with more comfort and dignity.

Upgraded pouch™ features include:



Fabric

- Softer fabric that's comfortable to wear and smooth on the skin
- Enhanced water repellent fabric that dries fast

Aesthetics

- Improved beige color that is discreet under clothes and visually appealing
- The fabric isn't see through to ensure patient privacy

Closure

- Secure fastener provides reassuring feedback to help indicate a secure seal
- The closure has soft rounded corners

www.hollister.com

*CeraPlus™ Skin Barriers contain the Remois Technology of Alcare Co.,Ltd.

**compared to our previous pouches. 1 Data on File, 2023; n=43. 2 Data on File, 2024. 3 Data on File, 2023; n=110. 4 Data on File, 2023; n=68 OSO Registry.

Prior to the use, be sure to read the Instructions for Use for information regarding Intended Use, Contraindications, Warnings, Precautions, and Instructions.

Not all products are CE marked.

Hollister, the Hollister logo, CeraPlus are trademarks of Hollister Incorporated. © 2024 Hollister Incorporated.



CeraPlus™
Ostomy Products

Ostomy Care

Collin Jarvis

Living with an ostomy since 2013

"I have found that CeraPlus™ Pouches and Skin Barriers are easy to use and to remove, and they are also flexible while maintaining dependable adhesion."













Scan Me
To Learn More About
Collin's Journey

Healthy skin. Positive outcomes.

Healthy Skin. Positive Outcomes.

Choosing the Right Skin Barrier FIT

FIT is critical to get a good seal around the stoma to prevent leakage and promote comfort.

Possible Clinical Case Scenarios	Barrier Examples
 <ul style="list-style-type: none">• Flat peristomal plane• Adequate stoma protrusion• Centered stoma opening	 <p>Flat Skin Barrier</p>
 <ul style="list-style-type: none">• Gentle pressure required to peristomal plane• Flush stoma• Peristomal injury	 <p>Integrated Soft Convex Barrier</p>
 <ul style="list-style-type: none">• Stoma flush with skin• Uneven peristomal plane• Off-centered stoma opening	 <p>Flat Skin Barrier with Round or Oval Convex Barrier Ring</p>
 <ul style="list-style-type: none">• Creases or folds in peristomal plane• Off-centered stoma opening• Stoma retracted below skin surface	 <p>Integrated Convex Barrier</p>
 <ul style="list-style-type: none">• Recessed or flushed stoma• Uneven abdominal contours• Off-centered stoma opening	 <p>Integrated Convex Barrier with Round or Oval Convex Barrier Ring</p>

Prior to use, be sure to read the Instructions for Use for information regarding intended use, contraindications, warnings, precautions and instructions.

1. Colwell JC, Goldberg MT & Carmel JE (Eds). (2015). WOCN Society Core Curriculum Ostomy Management. Chapter 10. Philadelphia, PA: Wolters Kluwer
2. Hoefloek J et al. Use of Convexity in Ostomy Care: Results of an International Consensus Meeting. JWOCN, January/February 2017

Characteristics of Convex Skin Barriers and Clinical Application: Results of an International Consensus Panel

McNichol, L, Cobb, T, Depaifve, Y, Quigley, M, Smitka, K & Gray, M 2021, 'Characteristics of Convex Skin Barriers and Clinical Application - Results of an International Consensus Panel' - *Journal Wound Ostomy Continence Nursing*, 48(6): 524-532.

Method Consensus of 12 nurses from 11 countries with extensive experience in ostomy care

Results Five convex skin barrier characteristics and 16 clinical application statements reached

Conclusion Standard definitions to help clinicians select the right convex barrier based on assessment criteria for better patient outcomes




Scan the QR Code to read the full consensus results!

The Five Characteristics of Convexity

1. Depth

The **depth** of the convex skin barrier is defined as the measurement from the apex of the dome to the base.



2. Compressibility

Compressibility is the capacity of the convex dome to be displaced or flattened.




3. Flexibility

Flexibility is how easily the convex skin barrier can bend.



4. Tension Location

Tension location is the position in which the Convex dome exerts downward and outward forces on the peristomal topography.



5. Slope

The **Slope** is the angle from the base of the convex skin barrier to the apex of the dome.



Hollister CeraPlus™

Ostomy Products* have received official dermatological accreditation from the Skin Health Alliance




Who is the Skin Health Alliance?

The **Skin Health Alliance** is a leading international skin health organization that works closely with dermatologists, scientists and the skin care industry with a common interest in the promotion of healthy skin. The **Skin Health Alliance** awards professional accreditation and endorsement to companies, services and brands seeking independent dermatological recognition for their product research.



- **Hollister CeraPlus™** Ostomy Products* have received official dermatological accreditation from the Skin Health Alliance
- In addition to the current rigorous safety testing our products go through, the SHA accreditation signifies that the **Hollister CeraPlus™** Ostomy Products have been independently verified by a specialist panel.
- **Hollister CeraPlus™** Ostomy Products may help prevent skin problems and help keep skin healthy from Day 1.
- By choosing **Hollister CeraPlus™** Ostomy Products, people using stoma pouches can have the confidence and pride of knowing their barrier helps protect their skin, empowering them to live their life to the fullest.



Prior to use, be sure to read the Instructions for Use for information regarding Intended Use, Contraindications, Warnings, Precautions, and Instructions. Hollister, the Hollister logo and CeraPlus are trademarks of Hollister Incorporated. All other trademarks and copyright are the property of their respective owners. © 2023 Hollister Incorporated. 923539

*CeraPlus skin barriers contain the Remois technology of Alcare Co., Ltd.



Ostomy Care
Healthy skin. Positive outcomes.

CeraPlus™ Ostomy Products*

Patrick Belton
Hollister End-User,
Canmore, AB

"Using the CeraPlus™ Ostomy products* has been so gentle on the skin that I have not had any type of ongoing peristomal issues."

Protection Where it Matters Most.



Two-Piece Barriers

Soft Convex Skin Barriers



Cut-to-fit * **NEW!**

Product #	Flange Size	For Stomas Up To	Box of
11702	1-3/4" (44 mm)	1" (25 mm)	5 With Tape Border
11703	2-1/4" (57 mm)	1-1/2" (38 mm)	
11704	2-3/4" (70 mm)	2" (51 mm)	

Product #	Flange Size	For Stomas Up To	Box of
* 15702	1-3/4" (44 mm)	1" (25 mm)	5 Without Tape Border
* 15703	2-1/4" (57 mm)	1-1/2" (38 mm)	
* 15704	2-3/4" (70 mm)	2" (51 mm)	

Pre-sized * **NEW!**

Product #	Flange Size	Barrier Opening	Box of
* 11902	1-3/4" (44 mm)	3/4" (19 mm)	5 With Tape Border
11903	1-3/4" (44 mm)	7/8" (22 mm)	
11904	1-3/4" (44 mm)	1" (25 mm)	
11905	2-1/4" (57 mm)	1-1/8" (29 mm)	
11906	2-1/4" (57 mm)	1-1/4" (32 mm)	
* 11907	2-1/4" (57 mm)	1-3/8" (35 mm)	
* 11908	2-1/4" (57 mm)	1-1/2" (38 mm)	
* 11909	2-3/4" (70 mm)	1-5/8" (41 mm)	
* 11910	2-3/4" (70 mm)	1-3/4" (44 mm)	
* 11911	2-3/4" (70 mm)	2" (51 mm)	

Product #	Flange Size	Barrier Opening	Box of
* 157102	1-3/4" (44 mm)	3/4" (19 mm)	5 Without Tape Border
* 157103	1-3/4" (44 mm)	7/8" (22 mm)	
* 157104	1-3/4" (44 mm)	1" (25 mm)	
* 157105	2-1/4" (57 mm)	1-1/8" (29 mm)	
* 157106	2-1/4" (57 mm)	1-1/4" (32 mm)	
* 157107	2-1/4" (57 mm)	1-3/8" (35 mm)	
* 157108	2-1/4" (57 mm)	1-1/2" (38 mm)	
* 157109	2-3/4" (70 mm)	1-5/8" (41 mm)	
* 157110	2-3/4" (70 mm)	1-3/4" (44 mm)	
* 157111	2-3/4" (70 mm)	2" (51 mm)	



Please go to page 16 for corresponding 2 piece pouch options.

Convex Skin Barriers



Cut-to-fit

Product #	Flange Size	For Stomas Up To	Box of
11402	1-3/4" (44 mm)	1" (25 mm)	5 With Tape Border
11403	2-1/4" (57 mm)	1-1/2" (38 mm)	
11404	2-3/4" (70 mm)	2" (51 mm)	
15302	1-3/4" (44 mm)	1" (25 mm)	5 Without Tape Border
15303	2-1/4" (57 mm)	1-1/2" (38 mm)	
15304	2-3/4" (70 mm)	2" (51 mm)	

Pre-sized

Product #	Flange Size	Barrier Opening	Box of
11503	1-3/4" (44 mm)	7/8" (22 mm)	5 With Tape Border
11504	1-3/4" (44 mm)	1" (25 mm)	
11505	2-1/4" (57 mm)	1-1/8" (29 mm)	
11506	2-1/4" (57 mm)	1-1/4" (32 mm)	5 Without Tape Border
153103	1-3/4" (44 mm)	7/8" (22 mm)	
153104	1-3/4" (44 mm)	1" (25 mm)	
153105	2-1/4" (57 mm)	1-1/8" (29 mm)	
153106	2-1/4" (57 mm)	1-1/4" (32 mm)	

Colour Matching System

To help you select a compatible two-piece pouch, every box label is colour-coded.

1. Select the skin barrier flange size with the appropriate barrier stoma opening and note its colour code.
2. Select a suitable pouch with the desired features in the matching colour.

- 1-3/4" (44 mm)
- 2-1/4" (57 mm)
- 2-3/4" (70 mm)
- 4" (102 mm)



Two-Piece Barriers + Kits

Flat Skin Barriers



Cut-to-fit

Product #	Flange Size	For Stomas Up To	Box of	
11202	1-3/4" (44 mm)	1-1/4" (32 mm)	5	With Tape Border
11203	2-1/4" (57 mm)	1-3/4" (44 mm)		
11204	2-3/4" (70 mm)	2-1/4" (57 mm)		
15102	1-3/4" (44 mm)	1-1/4" (32 mm)	5	Without Tape Border
15103	2-1/4" (57 mm)	1-3/4" (44 mm)		
15104	2-3/4" (70 mm)	2-1/4" (57 mm)		

Pre-sized

Product #	Flange Size	Barrier Opening	Box of	
112013	1-3/4" (44 mm)	7/8" (22 mm)	5	With Tape Border
112014	1-3/4" (44 mm)	1" (25 mm)		
112015	1-3/4" (44 mm)	1-1/8" (29 mm)		
112016	1-3/4" (44 mm)	1-1/4" (32 mm)		
112017	2-1/4" (57 mm)	1-3/8" (35 mm)		
112018	2-1/4" (57 mm)	1-1/2" (38 mm)	5	Without Tape Border
151113	1-3/4" (44 mm)	7/8" (22 mm)		
151114	1-3/4" (44 mm)	1" (25 mm)		
151115	1-3/4" (44 mm)	1-1/8" (29 mm)		
151116	1-3/4" (44 mm)	1-1/4" (32 mm)		
151117	2-1/4" (57 mm)	1-3/8" (35 mm)		
151118	2-1/4" (57 mm)	1-1/2" (38 mm)		

Soft Convex Kits

COMING SOON Ostomy Kits

Product #	Flange Size	For Stomas Up To	Box of	
19002C	1-3/4" (44 mm)	1" (25 mm)	5	With Tape Border
19003C	2-1/4" (57 mm)	1-1/2" (38 mm)		
19004C	2-3/4" (70 mm)	2" (51 mm)		
19612	1-3/4" (44 mm)	1" (25 mm)	5	Without Tape Border
19613	2-1/4" (57 mm)	1-1/2" (38 mm)		
19614	2-3/4" (70 mm)	2" (51 mm)		

Includes: Cut-to-fit Soft Convex Skin Barrier & transparent drainable pouch

COMING SOON Urostomy Kits

Product #	Flange Size	For Stomas Up To	Box of	
19202C	1-3/4" (44 mm)	1" (25 mm)	5	With Tape Border
19203C	2-1/4" (57 mm)	1-1/2" (38 mm)		
19204C	2-3/4" (70 mm)	2" (51 mm)		
19912	1-3/4" (44 mm)	1" (25 mm)	5	Without Tape Border
19913	2-1/4" (57 mm)	1-1/2" (38 mm)		
19914	2-3/4" (70 mm)	2" (51 mm)		

Includes: Cut-to-fit Soft Convex Skin Barrier, transparent urostomy pouch & drain tube adapter

Flat Kits

Ostomy Kits

Product #	Flange Size	For Stomas Up To	Box of
19603	2-1/4" (57 mm)	1-3/4" (44 mm)	5
19604	2-3/4" (70 mm)	2-1/4" (57 mm)	

Includes: Cut-to-fit Flat Skin Barrier with tape border & transparent drainable pouch

Urostomy Kits

Product #	Flange Size	For Stomas Up To	Box of
19903	2-1/4" (57 mm)	1-3/4" (44 mm)	5

Includes: Cut-to-fit Flat Skin Barrier with tape border, transparent urostomy pouch & drain tube adapter

One-Piece Pouching Systems

Soft Convex Skin Barriers

Drainable Pouch



Cut-to-fit With Tape Border					
Product #	Colour	For Stomas Up To	Pouch Length	Filter	Box of
8964	Transparent	1" (25 mm)	12" (30 cm)	AF300	5
8968		1-1/2" (38 mm)	12" (30 cm)		
89611		2-1/8" (55 mm)	12" (30 cm)		
8954	Beige with Viewing Option	1" (25 mm)	12" (30 cm)	AF300	5
8958		1-1/2" (38 mm)	12" (30 cm)		
89511		2-1/8" (55 mm)	12" (30 cm)		

Cut-to-fit Without Tape Border					
Product #	Colour	For Stomas Up To	Pouch Length	Filter	Box of
89518	Beige with Viewing Option	1-1/2" (38 mm)	12" (30 cm)	AF300	5
8195111		2-1/8" (55 mm)	12" (30 cm)		

Pre-sized With Tape Border					
Product #	Colour	Barrier Opening	Pouch Length	Filter	Box of
8960	Beige with Viewing Option	3/4" (19 mm)	12" (30 cm)	AF300	5
8961		7/8" (22 mm)	12" (30 cm)		
8962		1" (25 mm)	12" (30 cm)		
8963		1-1/8" (29 mm)	12" (30 cm)		

Flat Skin Barriers

Drainable Pouch



Cut-to-fit With Tape Border					
Product #	Colour	For Stomas Up To	Pouch Length	Filter	Box of
8931	Ultra-clear	2-1/2" (64 mm)	12" (30 cm)	No	10
8901	Beige with Viewing Option	2-1/2" (64 mm)	12" (30 cm)	AF300	

Cut-to-fit Without Tape Border					
Product #	Colour	For Stomas Up To	Pouch Length	Filter	Box of
89011	Beige with Viewing Option	2-1/2" (64 mm)	12" (30 cm)	AF300	10

Pre-sized With Tape Border					
Product #	Colour	Barrier Opening	Pouch Length	Filter	Box of
8925	Beige with Viewing Option	1" (25 mm)	12" (30 cm)	AF300	10
8930		1-3/16" (30 mm)	12" (30 cm)		
8935		1-3/8" (35 mm)	12" (30 cm)		
8940		1-9/16" (40 mm)	12" (30 cm)		

Convex Skin Barriers

Drainable Pouch



Cut-to-fit With Tape Border					
Product #	Colour	For Stomas Up To	Pouch Length	Filter	Box of
8984	Ultra-clear	1" (25 mm)	12" (30 cm)	No	5
8988		1-1/2" (38 mm)	12" (30 cm)		
89811		2" (51 mm)	12" (30 cm)		
8914	Beige with Viewing Option	1" (25 mm)	12" (30 cm)	AF300	5
8918		1-1/2" (38 mm)	12" (30 cm)		
89111		2" (51 mm)	12" (30 cm)		

Pre-sized With Tape Border					
Product #	Colour	Barrier Opening	Pouch Length	Filter	Box of
8990	Beige with Viewing Option	3/4" (19 mm)	12" (30 cm)	AF300	5
8991		7/8" (22 mm)	12" (30 cm)		
8992		1" (25 mm)	12" (30 cm)		
8993		1-1/8" (29 mm)	12" (30 cm)		

The CeraPlus™ Skin Barrier delivers both Fit and Formulation

Ceramide Infused Barrier

- **Protects** the skin's natural moisture barrier
- **Decreases** TEWL (transepidermal water loss) from damaged or eroded skin
- **Maintains** peristomal skin health

Start using CeraPlus™ Skin Barrier from day one to help ensure the healthy skin your patients deserve.

Contact your Hollister sales representative today for details.

SCAN ME







One-Piece Urostomy Pouching Systems

Soft Convex Skin Barriers Urostomy Pouch

Enhanced Design




Cut-to-fit				With Tape Border 
Product #	Colour	For Stomas Up To	Pouch Length	Box of
84134	Ultra-clear	1" (25 mm)	9" (23 cm)	5
84138		1-1/2" (38 mm)	9" (23 cm)	
841311		2-1/8" (55 mm)	9" (23 cm)	


Pre-sized				With Tape Border 
Product #	Colour	Barrier Opening	Pouch Length	Box of
8412	Beige with Viewing Option	3/4" (19 mm)	9" (23 cm)	5
8413		7/8" (22 mm)	9" (23 cm)	
8414		1" (25 mm)	9" (23 cm)	
8415		1-1/8" (29 mm)	9" (23 cm)	

Convex Skin Barriers Urostomy Pouch

Enhanced Design




Cut-to-fit				With Tape Border 
Product #	Colour	For Stomas Up To	Pouch Length	Box of
84404	Ultra-clear	1" (25 mm)	9" (23 cm)	5
84408		1-1/2" (38 mm)	9" (23 cm)	
844011		2" (51 mm)	9" (23 cm)	

Pre-sized				With Tape Border 
Product #	Colour	Barrier Opening	Pouch Length	Box of
84292	Beige with Viewing Option	3/4" (19 mm)	9" (23 cm)	5
84293		7/8" (22 mm)	9" (23 cm)	
84294		1" (25 mm)	9" (23 cm)	
84295		1-1/8" (29 mm)	9" (23 cm)	

Flat Skin Barriers Urostomy Pouch

Enhanced Design



Cut-to-fit				With Tape Border 
Product #	Colour	For Stomas Up To	Pouch Length	Box of
84200	Ultra-clear	2-1/2" (64 mm)	9" (23 cm)	10

A Great Pouch Even Better

Enhanced Flow-Control Valve

The new improved valve is made with softer material and a streamlined shape.

Multi-Chamber Design

The updated pouches have multiple built-in chambers, which are designed to facilitate even distribution of urine.



Multi Chamber
Front View



Multi Chamber,
Top View, Sliced

Upgraded pouches your patients can rely on^{1,2}

Our customer research helped us identify opportunities to make improvements to our pouches*:



Fabric

- Softer fabric that's comfortable to wear and smooth on the skin
- More breathable fabric to reduce sweating sensation
- Enhanced water repellent fabric that dries fast



Aesthetics

- Improved beige color that is discreet under clothes and visually appealing
- The fabric isn't see-through to ensure patient privacy



Closure

- Secure fastener provides reassuring feedback to help indicate a secure seal
- The closure forms a bigger and wider opening, designed for easier access to the inside and to clean the edges
- The closure has soft rounded corners and is durable, for a focus on patient comfort

We will upgrade all existing CeraPlus™ Pouches, so every patient using Hollister can benefit.

To find out more about the pouch upgrade, request a call from your Hollister representative who can take you through the positive changes and product availability.



www.hollister.com

Prior to the use, be sure to read the Instructions for Use for information regarding Intended Use, Contraindications, Warnings, Precautions, and Instructions. Hollister, the Hollister logo, CeraPlus, CeraRing are trademarks of Hollister Incorporated. Not all products are CE marked. Person received compensation from Hollister Incorporated. *CeraPlus skin barriers contain the Remois Technology of Alcare Co., Ltd. ©2024 Hollister Incorporated.

*compared to our previous pouch. 1 Data on File, 2023; n=43. 2 Data on File, 2024.

Find out more at **Hollister.com/update** or call: **1-800-263-7400**



CeraPlus™
Ostomy Products

Two-Piece Pouches



Scan Me
To Learn More About
LeeAnne's Journey

LeeAnne Hayden

Living with an ostomy since 2012

"I am passionate about yoga, and to do what I love I need an ostomy product that will protect my skin and that I can rely on."

Two-Piece Pouches

New Image

Drainable Pouches



Colour: Ultra-clear

Product #	Flange Size	Pouch Length	Filter	Box of
18002	1-3/4" (44 mm)	12" (30 cm)	No	10
18003	2-1/4" (57 mm)	12" (30 cm)		
18004	2-3/4" (70 mm)	12" (30 cm)		
18006	4" (102 mm)	12" (30 cm)		

Colour: Transparent

Product #	Flange Size	Pouch Length	Filter	Box of
18132	1-3/4" (44 mm)	12" (30 cm)	No	10
18133	2-1/4" (57 mm)	12" (30 cm)		
18134	2-3/4" (70 mm)	12" (30 cm)		
18192	1-3/4" (44 mm)	12" (30 cm)	AF300	10
18193	2-1/4" (57 mm)	12" (30 cm)		
18194	2-3/4" (70 mm)	12" (30 cm)		
18292	1-3/4" (44 mm)	7" (18 cm)	AF300	20
18293	2-1/4" (57 mm)	7" (18 cm)		
18294	2-3/4" (70 mm)	7" (18 cm)		

Colour: Beige

Product #	Flange Size	Pouch Length	Filter	Box of
18112	1-3/4" (44 mm)	12" (30 cm)	No	10
18113	2-1/4" (57 mm)	12" (30 cm)		
18114	2-3/4" (70 mm)	12" (30 cm)		
18182	1-3/4" (44 mm)	12" (30 cm)	AF300	10
18183	2-1/4" (57 mm)	12" (30 cm)		
18184	2-3/4" (70 mm)	12" (30 cm)		
18282	1-3/4" (44 mm)	7" (18 cm)	AF300	20
18283	2-1/4" (57 mm)	7" (18 cm)		
18284	2-3/4" (70 mm)	7" (18 cm)		

Clamp Closure



Colour: Ultra-clear

Product #	Flange Size	Pouch Length	Filter	Box of
18172	1-3/4" (44 mm)	12" (30 cm)	No	10
18173	2-1/4" (57 mm)	12" (30 cm)		
18174	2-3/4" (70 mm)	12" (30 cm)		
18176	4" (102 mm)	12" (30 cm)		

Colour: Transparent

Product #	Flange Size	Pouch Length	Filter	Box of
18102	1-3/4" (44 mm)	12" (30 cm)	No	10
18103	2-1/4" (57 mm)	12" (30 cm)		
18104	2-3/4" (70 mm)	12" (30 cm)		

Colour: Beige

Product #	Flange Size	Pouch Length	Filter	Box of
18122	1-3/4" (44 mm)	12" (30 cm)	No	10
18123	2-1/4" (57 mm)	12" (30 cm)		
18124	2-3/4" (70 mm)	12" (30 cm)		
18142	1-3/4" (44 mm)	12" (30 cm)	AF300	10
18143	2-1/4" (57 mm)	12" (30 cm)		
18144	2-3/4" (70 mm)	12" (30 cm)		
18202	1-3/4" (44 mm)	9" (23 cm)	No	10
18203	2-1/4" (57 mm)	9" (23 cm)		
18204	2-3/4" (70 mm)	9" (23 cm)		

Closed Pouches



Colour: Transparent

Product #	Flange Size	Pouch Length	Filter	Box of
18362	1-3/4" (44 mm)	9" (23 cm)	AF300	60
18363	2-1/4" (57 mm)	9" (23 cm)		
18364	2-3/4" (70 mm)	9" (23 cm)		

Colour: Beige

Product #	Flange Size	Pouch Length	Filter	Box of
18352	1-3/4" (44 mm)	7" (18 cm)	No	30
18353	2-1/4" (57 mm)	7" (18 cm)		
18354	2-3/4" (70 mm)	7" (18 cm)		
18372-30	1-3/4" (44 mm)	9" (23 cm)	AF300	30
18373-30	2-1/4" (57 mm)	9" (23 cm)		
18374-30	2-3/4" (70 mm)	9" (23 cm)		
18392	1-3/4" (44 mm)	7" (18 cm)	AF300	60
18393	2-1/4" (57 mm)	7" (18 cm)		
18394	2-3/4" (70 mm)	7" (18 cm)		

Two-Piece Pouches + Kits

New Image

High Output Pouches

Enhanced Design



Colour: Transparent

Product #	Flange Size	Pouch Length	Filter	Box of
18022	1-3/4" (44 mm)	12" (30 cm)	AF300	10
18023	2-1/4" (57 mm)	12" (30 cm)		
18024	2-3/4" (70 mm)	12" (30 cm)		

* 5551 Collection bag can be found on page 36


Original Design



Colour: Ultra-clear


Product #	Flange Size	Pouch Length	Filter	Box of
18012	1-3/4" (44 mm)	12" (30 cm)	No	10
18013	2-1/4" (57 mm)	12" (30 cm)		
18014	2-3/4" (70 mm)	12" (30 cm)		
18016	4" (102 mm)	12" (30 cm)		

Enhanced Design Features Include:




A reimagined wide-bore soft tap minimizes the burden of blockages

- A large inner diameter
- Soft, flexible, and durable material
- Plug-style closure with a large knob



Large collection bag and long tubing offers patients more freedom — day and night

- 2 litre capacity collection bag – Less frequent draining may improve quality of sleep
- Corrugated tubing reduces kinking during movement
- 152.4cm tube offers freedom of movement



Intuitive locking mechanism connects patients to an extra sense of security

- New twist and lock feature for a secure pouch drain tap and collection bag connection

Ostomy Kits

FlexWear
Drainable Pouching System

Non-sterile				
Product #	Flange Size	For Stomas Up To	Pouch Length	Box of
19006	4" (102 mm)	3-1/2" (89 mm)	12" (30 mm)	5
Sterile				
Product #	Flange Size	For Stomas Up To	Pouch Length	Box of
19054	2-3/4" (70 mm)	2-1/4" (57 mm)	12" (30 mm)	5

Urostomy Kits

Flexlend

Sterile			
Product #	Flange Size	For Stomas Up To	Box of
19254	2-3/4" (70 mm)	2-1/4" (57 mm)	5

Includes:

- *Transparent pouch*
- *Cut-to-fit Flexlend Skin Barrier with tape border*
- *Drain tube adapter*

Two-Piece Pouches

New Image

Urostomy Pouches

Enhanced Design



Colour: Ultra-clear

Product #	Flange Size	Pouch Length	Box of
18922	1-3/4" (44 mm)	9" (23 cm)	10
18923	2-1/4" (57 mm)	9" (23 cm)	
18924	2-3/4" (70 mm)	9" (23 cm)	

Colour: Transparent

Product #	Flange Size	Pouch Length	Box of
18902	1-3/4" (44 mm)	9" (23 cm)	10
18903	2-1/4" (57 mm)	9" (23 cm)	
18904	2-3/4" (70 mm)	9" (23 cm)	

Colour: Beige

Product #	Flange Size	Pouch Length	Box of
18912	1-3/4" (44 mm)	9" (23 cm)	10
18913	2-1/4" (57 mm)	9" (23 cm)	
18914	2-3/4" (70 mm)	9" (23 cm)	

A Great Pouch, Even Better

Enhanced Flow-Control Valve

The new improved valve is made with softer material and a streamlined shape.

Multi-Chamber Design

The updated pouches have multiple built-in chambers, which are designed to facilitate even distribution of urine.



Multi Chamber
Front View



Multi Chamber,
Top View, Sliced



Original Design

Colour: Transparent

Product #	Flange Size	Pouch Length	Box of
18402	1-3/4" (44 mm)	9" (23 cm)	10
18403	2-1/4" (57 mm)	9" (23 cm)	
18404	2-3/4" (70 mm)	9" (23 cm)	

Selecting the Right Ostomy Product Can Make All the Difference




The skin around your stoma deserves advanced protection. Let us help you find the CeraPlus™ Product* that gives you the security you need to enable healthy skin in your daily life.

- By answering seven simple questions, let us help you find the CeraPlus™ Ostomy Product* you need.



Seven questions. One result.
See what works for you.



**Contains the Remoils Technology of Alcare Co., Ltd.*

Prior to use, be sure to read the Instructions for Use for information regarding Intended Use, Contraindications, Warnings, Precautions, and Instructions.

Hollister, the Hollister logo and CeraPlus are trademarks of Hollister Incorporated.
©2025 Hollister Incorporated.



Healthy Skin. Positive Outcomes.

Choosing the Right Skin Barrier FORMULATION

FORMULATION is the customized ingredients responsible for adhesion, erosion resistance and absorption that must be carefully formulated to create an environment where healthy skin can thrive.

Extended Wear Skin Barriers

CeraPlus™ Skin Barriers*



- Contains skin barrier characteristics (adhesion, absorption and erosion resistance) and is **infused with ceramide** — a naturally occurring lipid in the skin which protects against dryness
- Formulation designed to help protect the skin's natural moisture barrier and help maintain good peristomal skin health

Flextend

- Durable with resistance to erosion
- Swells when it comes in contact with liquid to help form a seal between the skin and the stoma
- Suitable for those who perspire heavily
- Strong adhesion and flexibility

Standard Wear Skin Barriers

FlexWear

- Good balance of durability and gentleness providing moderate adhesion, and can handle frequent or infrequent pouch changes
- Resistant to fluid erosion
- Suitable for all skin types, depending on change frequency
- Designed to be worn more than one day

SoftFlex

- Appropriate choice for:
 - Those who have thin, fragile, or sensitive skin (such as those who have undergone chemo/radiation therapy)
 - People who desire frequent removal
- Designed to be gentle to the skin, which allows for frequent pouch changes

Did you know?

Hollister CeraPlus™ Ostomy Products* have received official dermatological accreditation from the Skin Health Alliance.

Who is the Skin Health Alliance? It is a leading international skin health organization that works closely with dermatologists, scientists and skin care industry with a common interest in the promotion of healthy skin. The Skin Health Alliance awards professional accreditation and endorsement to companies, services and brands seeking independent dermatological recognition for their product research.



*CeraPlus products contain the Remois Technology of Alcare Co., Ltd.



Two-Piece Barriers

Healthy skin. Positive outcomes.

Two-Piece Barriers

New Image

Convex Skin Barriers

Flextend



Cut-to-fit

Product #	Flange Size	For Stomas Up To	Box of	
14802	1-3/4" (44 mm)	1" (25 mm)	5	With Tape Border
14803	2-1/4" (57 mm)	1-1/2" (38 mm)		
14804	2-3/4" (70 mm)	2" (51 mm)		
15802	1-3/4" (44 mm)	1" (25 mm)	5	Without Tape Border
15803	2-1/4" (57 mm)	1-1/2" (38 mm)		
15804	2-3/4" (70 mm)	2" (51 mm)		

Pre-sized

Product #	Flange Size	Barrier Opening	Box of	
13901	1-3/4" (44 mm)	5/8" (16 mm)	5	With Tape Border
13902	1-3/4" (44 mm)	3/4" (19 mm)		
13903	1-3/4" (44 mm)	7/8" (22 mm)		
13904	1-3/4" (44 mm)	1" (25 mm)		
13905	2-1/4" (57 mm)	1-1/8" (29 mm)		
13906	2-1/4" (57 mm)	1-1/4" (32 mm)		
13907	2-1/4" (57 mm)	1-3/8" (35 mm)		
13908	2-1/4" (57 mm)	1-1/2" (38 mm)		
13909	2-3/4" (70 mm)	1-5/8" (41 mm)		
13910	2-3/4" (70 mm)	1-3/4" (44 mm)		
13911	2-3/4" (70 mm)	2" (51 mm)		
15903	1-3/4" (44 mm)	7/8" (22 mm)	5	Without Tape Border
15904	1-3/4" (44 mm)	1" (25 mm)		
15905	2-1/4" (57 mm)	1-1/8" (29 mm)		
15906	2-1/4" (57 mm)	1-1/4" (32 mm)		

Flat Skin Barriers

Flextend



Cut-to-fit

Product #	Flange Size	For Stomas Up To	Box of	
14602	1-3/4" (44 mm)	1-1/4" (32 mm)	5	With Tape Border
14603	2-1/4" (57 mm)	1-3/4" (44 mm)		
14604	2-3/4" (70 mm)	2-1/4" (57 mm)		
14606	4" (102 mm)	3-1/2" (89 mm)		
15602	1-3/4" (44 mm)	1-1/4" (32 mm)	5	Without Tape Border
15603	2-1/4" (57 mm)	1-3/4" (44 mm)		
15604	2-3/4" (70 mm)	2-1/4" (57 mm)		





Pre-sized

Product #	Flange Size	Barrier Opening	Box of	
14701	1-3/4" (44 mm)	5/8" (16 mm)	5	With Tape Border
14702	1-3/4" (44 mm)	3/4" (19 mm)		
14703	1-3/4" (44 mm)	7/8" (22 mm)		
14704	1-3/4" (44 mm)	1" (25 mm)		
14705	1-3/4" (44 mm)	1-1/8" (29 mm)		
14706	1-3/4" (44 mm)	1-1/4" (32 mm)		
14707	2-1/4" (57 mm)	1-3/8" (35 mm)		
14708	2-1/4" (57 mm)	1-1/2" (38 mm)		
14709	2-1/4" (57 mm)	1-3/4" (44 mm)		
16101	1-3/4" (44 mm)	5/8" (16 mm)	5	Without Tape Border
16102	1-3/4" (44 mm)	3/4" (19 mm)		
16103	1-3/4" (44 mm)	7/8" (22 mm)		
16104	1-3/4" (44 mm)	1" (25 mm)		
16105	1-3/4" (44 mm)	1-1/8" (29 mm)		
16106	1-3/4" (44 mm)	1-1/4" (32 mm)		
16107	2-1/4" (57 mm)	1-3/8" (35 mm)		
16108	2-1/4" (57 mm)	1-1/2" (38 mm)		

Colour Matching System

To help you select a compatible two-piece pouch, every box label is colour-coded.

1. Select the skin barrier flange size with the appropriate barrier stoma opening and note its colour code.
2. Select a suitable pouch with the desired features in the matching colour.


	1-3/4" (44 mm)
	2-1/4" (57 mm)
	2-3/4" (70 mm)
	4" (102 mm)

Two-Piece Barriers

New Image

Convex Skin Barriers

FlexWear



Cut-to-fit


Product #	Flange Size	For Stomas Up To	Box of	
13402	1-3/4" (44 mm)	1" (25 mm)	5	With Tape Border
13403	2-1/4" (57 mm)	1-1/2" (38 mm)		
13404	2-3/4" (70 mm)	2" (51 mm)		
15402	1-3/4" (44 mm)	1" (25 mm)	5	Without Tape Border
15403	2-1/4" (57 mm)	1-1/2" (38 mm)		
15404	2-3/4" (70 mm)	2" (51 mm)		

Pre-sized

Product #	Flange Size	Barrier Opening	Box of	
13503	1-3/4" (44 mm)	7/8" (22 mm)	5	With Tape Border
13504	1-3/4" (44 mm)	1" (25 mm)		
13505	2-1/4" (57 mm)	1-1/8" (29 mm)		
13506	2-1/4" (57 mm)	1-1/4" (32 mm)	5	Without Tape Border
15503	1-3/4" (44 mm)	7/8" (22 mm)		
15504	1-3/4" (44 mm)	1" (25 mm)		
15505	2-1/4" (57 mm)	1-1/8" (29 mm)		
15506	2-1/4" (57 mm)	1-1/4" (32 mm)		

Adhesive Border

The adhesive border may help provide peace of mind and increase confidence in patients with skin creases and varying topography.



Flat Skin Barriers

FlexWear



Cut-to-fit

Product #	Flange Size	For Stomas Up To	Box of	
14202	1-3/4" (44 mm)	1-1/4" (32 mm)	5	With Tape Border
14203	2-1/4" (57 mm)	1-3/4" (44 mm)		
14204	2-3/4" (70 mm)	2-1/4" (57 mm)		
14206	4" (102 mm)	3-1/2" (89 mm)		
15202	1-3/4" (44 mm)	1-1/4" (32 mm)	5	Without Tape Border
15203	2-1/4" (57 mm)	1-3/4" (44 mm)		
15204	2-3/4" (70 mm)	2-1/4" (57 mm)		

Pre-sized

Product #	Flange Size	Barrier Opening	Box of	
14301	1-3/4" (44 mm)	5/8" (16 mm)	5	With Tape Border
14302	1-3/4" (44 mm)	3/4" (19 mm)		
14303	1-3/4" (44 mm)	7/8" (22 mm)		
14304	1-3/4" (44 mm)	1" (25 mm)		
14305	1-3/4" (44 mm)	1-1/8" (29 mm)		
14306	1-3/4" (44 mm)	1-1/4" (32 mm)		
14307	2-1/4" (57 mm)	1-3/8" (35 mm)		
14308	2-1/4" (57 mm)	1-1/2" (38 mm)		
14309	2-1/4" (57 mm)	1-3/4" (44 mm)		
16401	1-3/4" (44 mm)	5/8" (16 mm)	5	Without Tape Border
16402	1-3/4" (44 mm)	3/4" (19 mm)		
16403	1-3/4" (44 mm)	7/8" (22 mm)		
16404	1-3/4" (44 mm)	1" (25 mm)		
16405	1-3/4" (44 mm)	1-1/8" (29 mm)		
16406	1-3/4" (44 mm)	1-1/4" (32 mm)		
16407	2-1/4" (57 mm)	1-3/8" (35 mm)		
16408	2-1/4" (57 mm)	1-1/2" (38 mm)		

Floating Flange


The floating flange may help allow a secure pouch and barrier coupling without pressing on the abdomen.



Start With Soft Convexity

A flat skin barrier works well with a flat peristomal plane and adequate stoma protrusion. However, for many patients, the presence of skin folds, recessed stomas, or extra weight can make getting the right fit challenging.

CeraPlus™ Soft Convex allows for gentle pressure for a comfortable and secure fit as it conforms to naturally uneven body contours.



DID YOU KNOW?

People with creases in the peristomal area were **3X more likely to experience a skin complication** than those who did not have creases?²

1. Based on evaluations completed in France, Germany, Canada, US, Australia, UK, Belgium, New Zealand, and Switzerland; Hollister Data on File, ref-02530. May 2021. Types of ostomies evaluated: Ileostomy (42%), Colostomy (29%), Urostomy (25%), Other/Multiple (4%). Patients evaluated ranged from post-op to 20+ years living with an ostomy; 50% had been living with an ostomy for less than 4 weeks.
2. Salvadalena G, Colwell JC, Skountrianos G, Pittman J. Lessons Learned About Peristomal Skin Complications: Secondary Analysis of the ADVOCATE Trial. J Wound Ostomy Continence Nurs. 2020;47(4):357-363.

Clinicians Agree:

CeraPlus™ Soft Convex is an Ideal Choice

70%

would recommend it for post-operative use.¹

Based on 300 evaluations

91%

were satisfied or very satisfied with the ability to get the best fit.¹

Based on 322 evaluations



One-Piece Pouching Systems

Healthy skin. Positive outcomes.

One-Piece Pouching Systems

Premier

Soft Convex Skin Barriers

Flextend
Drainable Pouch





Cut-to-fit			With Tape Border		
Product #	Colour	For Stomas Up To	Pouch Length	Filter	Box of
8674	Beige with Viewing Option	5/8"– 1" (13 - 25 mm)	12" (30 cm)	AF300	5
8678		5/8"– 1-1/2" (13 - 38 mm)	12" (30 cm)		
86711		5/8"– 2-1/8" (13 - 55 mm)	12" (30 cm)		

Pre-sized			With Tape Border		
Product #	Colour	Barrier Opening	Pouch Length	Filter	Box of
8661	Beige with Viewing Option	7/8" (22 mm)	12" (30 cm)	AF300	5
8662		1" (25 mm)	12" (30 cm)		
8663		1-1/8" (29 mm)	12" (30 cm)		
8664		1-1/4" (32 mm)	12" (30 cm)		
8665		1-1/2" (38 mm)	12" (30 cm)		
8666		1-3/4" (44 mm)	12" (30 cm)		

Convex Skin Barriers

Flextend
Clamp Closure





Cut-to-fit			With Tape Border		
Product #	Colour	For Stomas Up To	Pouch Length	Box of	
82624	Transparent	1" (25 mm)	12" (30 cm)	5	
82628		1-1/2" (38 mm)	12" (30 cm)		
826211		2" (51 mm)	12" (30 cm)		
Pre-sized			With Tape Border		
Product #	Colour	Barrier Opening	Pouch Length	Box of	
82610	Transparent	3/4" (19 mm)	12" (30 cm)	5	
82611		7/8" (22 mm)	12" (30 cm)		
82612		1" (25 mm)	12" (30 cm)		
82613		1-1/8" (29 mm)	12" (30 cm)		
82614		1-1/4" (32 mm)	12" (30 cm)		
82615		1-3/8" (35 mm)	12" (30 cm)		
82616		1-1/2" (38 mm)	12" (30 cm)		
82617		1-5/8" (41 mm)	12" (30 cm)		
82618		1-3/4" (44 mm)	12" (30 cm)		
82619		2" (51 mm)	12" (30 cm)		

Convex Skin Barriers

Flextend
Drainable Pouch



Cut-to-fit			With Tape Border 		
Product #	Colour	For Stomas Up To	Pouch Length	Filter	Box of
8244	Transparent	1" (25 mm)	12" (30 cm)	No	5
8248		1-1/2" (38 mm)	12" (30 cm)		
82411		2" (51 mm)	12" (30 cm)		
8274	Transparent	1" (25 mm)	12" (30 cm)	AF300	
8278		1-1/2" (38 mm)	12" (30 cm)		
82711		2" (51 mm)	12" (30 cm)		
8214	Beige	1" (25 mm)	12" (30 cm)	No	5
82111		2" (51 mm)	12" (30 cm)		

Pre-sized			With Tape Border 		
Product #	Colour	Barrier Opening	Pouch Length	Box of	
8280	Transparent	3/4" (19 mm)	12" (30 cm)	5	
8281		7/8" (22 mm)	12" (30 cm)		
8282		1" (25 mm)	12" (30 cm)		
8283		1-1/8" (29 mm)	12" (30 cm)		
8284		1-1/4" (32 mm)	12" (30 cm)		
8285		1-3/8" (35 mm)	12" (30 cm)		
8286		1-1/2" (38 mm)	12" (30 cm)		
8288		1-3/4" (44 mm)	12" (30 cm)		
8289		2" (51 mm)	12" (30 cm)		
8290	Beige	3/4" (19 mm)	12" (30 cm)	5	
8291		7/8" (22 mm)	12" (30 cm)		
8292		1" (25 mm)	12" (30 cm)		
8293		1-1/8" (29 mm)	12" (30 cm)		
8294		1-1/4" (32 mm)	12" (30 cm)		
8295		1-3/8" (35 mm)	12" (30 cm)		
8296		1-1/2" (38 mm)	12" (30 cm)		

One-Piece Pouching Systems

Premier

Flat Skin Barriers

Flextend Drainable Pouch



Cut-to-fit With Tape Border					
Product #	Colour	For Stomas Up To	Pouch Length	Filter	Box of
8531	Ultra-clear	2-1/2" (64 mm)	12" (30 cm)	No	10
8331	Transparent	2-1/2" (64 mm)	12" (30 cm)	AF300	
8551	Beige	2-1/2" (64 mm)	12" (30 cm)	No	

Pre-sized With Tape Border				
Product #	Colour	Barrier Opening	Pouch Length	Box of
8538	Transparent	1-1/4" (32 mm)	12" (30 cm)	10
8552	Beige	1" (25 mm)	12" (30 cm)	10
8558		1-1/4" (32 mm)	12" (30 cm)	

FlexWear Drainable Pouch



Cut-to-fit With Tape Border					
Product #	Colour	For Stomas Up To	Pouch Length	Filter	Box of
8031	Ultra-clear	2-1/2" (64 mm)	12" (30 cm)	No	10
8081	Transparent	2-1/2" (64 mm)	12" (30 cm)	AF300	
8181	Beige	2-1/2" (64 mm)	12" (30 cm)	AF300	

Flextend Clamp Closure



Cut-to-fit With Tape Border				
Product #	Colour	For Stomas Up To	Pouch Length	Box of
8651	Ultra-clear	2-1/2" (64 mm)	12" (30 cm)	10
8631	Transparent	2-1/2" (64 mm)	12" (30 cm)	10
8641		2-1/2" (64 mm)	9" (23 cm)	

Pre-sized With Tape Border				
Product #	Colour	Barrier Opening	Pouch Length	Box of
8637	Transparent	3/4" (19 mm)	12" (30 cm)	10
8632		1" (25 mm)	12" (30 cm)	
8638		1-1/4" (32 mm)	12" (30 cm)	
8633		1-1/2" (38 mm)	12" (30 cm)	
8639		1-3/4" (44 mm)	12" (30 cm)	
8635	Beige	2-1/2" (64 mm)	12" (30 cm)	10
8647		3/4" (19 mm)	9" (23 cm)	
8642		1" (25 mm)	9" (23 cm)	
8648		1-1/4" (32 mm)	9" (23 cm)	
8643		1-1/2" (38 mm)	9" (23 cm)	
8649		1-3/4" (44 mm)	9" (23 cm)	

One-Piece Pouching Systems

Premier

Flat Skin Barriers

SoftFlex
Drainable Pouch



Cut-to-fit Without Tape Border					
Product #	Colour	For Stomas Up To	Pouch Length	Filter	Box of
88300	Beige	2-1/8" (55 mm)	12" (30 cm)	AF300	10
88700		2-1/8" (55 mm)	7" (18 cm)		20

Cut-to-fit — Oval Skin Barrier Without Tape Border					
Product #	Colour	For Stomas Up To	Pouch Length	Filter	Box of
88402	Transparent	2-1/2" x 3" (66 x 77 mm)	12" (30 cm)	AF300	10
88302	Beige	2-1/2" x 3" (66 x 77 mm)	12" (30 cm)		

Pre-sized Without Tape Border					
Product #	Colour	Barrier Opening	Pouch Length	Filter	Box of
88325	Beige	1" (25 mm)	12" (30 cm)	AF300	10
88330		1-3/16" (30 mm)	12" (30 cm)		
88335		1-3/8" (35 mm)	12" (30 cm)		
88340		1-9/16" (40 mm)	12" (30 cm)		
88725		1" (25 mm)	7" (18 cm)	AF300	20
88730		1-3/16" (30 mm)	7" (18 cm)		
88735		1-3/8" (35 mm)	7" (18 cm)		
88740		1-9/16" (40 mm)	7" (18 cm)		

Flat Skin Barriers

SoftFlex
Closed Pouch



Cut-to-fit Without Tape Border					
Product #	Colour	For Stomas Up To	Pouch Length	Filter	Box of
82300	Beige	2-1/8" (55 mm)	9" (23 cm)	AF300	30
82100		2-1/8" (55 mm)	7" (18 cm)		

Cut-to-fit — Oval Skin Barrier Without Tape Border					
Product #	Colour	For Stomas Up To	Pouch Length	Filter	Box of
82402	Transparent	2-1/2" x 3" (66 mm x 77 mm)	9" (23 cm)	AF300	30
82302	Beige	2-1/2" x 3" (66 mm x 77 mm)	9" (23 cm)		

Pre-sized Without Tape Border					
Product #	Colour	Barrier Opening	Pouch Length	Filter	Box of
82325	Beige	1" (25 mm)	9" (23 cm)	AF300	30
82330		1-3/16" (30 mm)	9" (23 cm)		
82335		1-3/8" (35 mm)	9" (23 cm)		
82125		1" (25 mm)	7" (18 cm)		
82130		1-3/16" (30 mm)	7" (18 cm)		
82135		1-3/8" (35 mm)	7" (18 cm)		

Stoma Caps
with Microporous Tape



Pre-sized With Tape Border					
Product #	Colour	For Stomas Up To	Pouch Length	Filter	Box of
3184	White	2" (51 mm)	4-1/4" (11 cm)	Yes	30
3186		3" (76 mm)	4-1/4" (11 cm)		

SoftFlex Skin Barrier



Pre-sized Without Tape Border					
Product #	Colour	For Stomas Up To	Pouch Length	Filter	Box of
1796	Beige	1-5/16" (50 mm)	4" (10 cm)	Yes	30

One-Piece Pouching Systems

Premier

High Output Pouches

Flextend
Enhanced Design

Soft Convex



Cut-to-fit

Product #	Colour	For Stomas Up to	Pouch Length	Box of	Box of
83025		1" (25 mm)	12" (30 cm)		
83038	Transparent	1-1/2" (38 mm)	12" (30 cm)	AF300	5
83055		2-1/8" (55 mm)	12" (30 cm)		

Pre-Sized

Product #	Colour	Barrier Opening	Pouch Length	Box of	Box of
83125		1" (25 mm)	12" (30 cm)		
83130	Transparent	1-3/16" (30 mm)	12" (30 cm)	AF300	5
83135		1-3/8" (35 mm)	12" (30 cm)		

Flat

Cut-to-fit

Product #	Colour	For Stomas Up to	Pouch Length	Filter	Box of
26064b	Transparent	2-1/2" (64 mm)	12" (30 cm)	AF300	10
26364	Beige	2-1/2" (64 mm)	12" (30 cm)		

* 5551 Collection bag can be found on page 36

Flextend
Original Design




Non-sterile

Without Tape Border


Product #	Colour	Barrier Opening	Pouch Length	Box of
80070	Ultra-clear	2-3/4" (70 mm)	12" (30 cm)	10
80110		4-1/3" (110 mm)	12" (30 cm)	

Enhanced Design features include:




Unique combination of fit and formulation

- Flat and Soft Convex skin barrier options
- Highly durable formulation maximises wear time and withstands erosion
- Adhesive Border Barrier helps provide peace of mind with an extra sense of security




A reimagined wide-bore soft tap

- A large inner diameter
- Soft, flexible, and durable material
- Plug-style closure with a large knob



Large collection bag and long tubing

- 2 litre capacity collection bag – Less frequent draining may improve quality of sleep
- Corrugated tubing reduces kinking during movement
- 152.4cm tube offers freedom of movement



Intuitive locking mechanism

- New twist and lock feature for a secure pouch drain tap and collection bag connection

One-Piece Pouching Systems

Premier


Convex Skin Barriers


Urostomy Pouch

Flextend

Enhanced Design




Cut-to-fit				With Tape Border 
Product #	Colour	For Stomas Up To	Pouch Length	Box of
84794	Ultra-clear	1" (25 mm)	9" (23 cm)	5
84798		1-1/2" (38 mm)	9" (23 cm)	
8479111		2" (51 mm)	9" (23 cm)	
84398	Beige	1-1/2" (38 mm)	9" (23 cm)	
8439111		2" (51 mm)	9" (23 cm)	


Pre-sized				With Tape Border 
Product #	Colour	Barrier Opening	Pouch Length	Box of
84890	Ultra-clear	1/2" (13 mm)	9" (23 cm)	5
84891		5/8" (16 mm)	9" (23 cm)	
84892		3/4" (19 mm)	9" (23 cm)	
84893		7/8" (22 mm)	9" (23 cm)	
84894		1" (25 mm)	9" (23 cm)	
84895		1-1/8" (29 mm)	9" (23 cm)	
84896		1-1/4" (32 mm)	9" (23 cm)	
84897		1-3/8" (35 mm)	9" (23 cm)	
84898		1-1/2" (38 mm)	9" (23 cm)	
84899		1-3/4" (44 mm)	9" (23 cm)	
84990	Beige	1/2" (13 mm)	9" (23 cm)	5
84991		5/8" (16 mm)	9" (23 cm)	
84992		3/4" (19 mm)	9" (23 cm)	
84993		7/8" (22 mm)	9" (23 cm)	
84994		1" (25 mm)	9" (23 cm)	
84995		1-1/8" (29 mm)	9" (23 cm)	
84996		1-1/4" (32 mm)	9" (23 cm)	
84997		1-3/8" (35 mm)	9" (23 cm)	
84998		1-1/2" (38 mm)	9" (23 cm)	
84999		1-3/4" (44 mm)	9" (23 cm)	

Flextend

Original Design



Cut-to-fit				With Tape Border 
Product #	Colour	For Stomas Up To	Pouch Length	Box of
8474	Transparent	1" (25 mm)	9" (23 cm)	5
8478		1-1/2" (38 mm)	9" (23 cm)	
84711		2" (51 mm)	9" (23 cm)	

Pre-sized				With Tape Border 
Product #	Colour	Barrier Opening	Pouch Length	Box of
8480	Transparent	1/2" (13 mm)	9" (23 cm)	5
8481		5/8" (16 mm)	9" (23 cm)	
8482		3/4" (19 mm)	9" (23 cm)	
8483		7/8" (22 mm)	9" (23 cm)	
8484		1" (25 mm)	9" (23 cm)	
8485		1-1/8" (29 mm)	9" (23 cm)	
8486		1-1/4" (32 mm)	9" (23 cm)	
8487		1-3/8" (35 mm)	9" (23 cm)	
8488		1-1/2" (38 mm)	9" (23 cm)	
8489		1-3/4" (44 mm)	9" (23 cm)	
8490	Beige	1/2" (13 mm)	9" (23 cm)	5
8491		5/8" (16 mm)	9" (23 cm)	
8492		3/4" (19 mm)	9" (23 cm)	
8493		7/8" (22 mm)	9" (23 cm)	
8494		1" (25 mm)	9" (23 cm)	
8495		1-1/8" (29 mm)	9" (23 cm)	
8496		1-1/4" (32 mm)	9" (23 cm)	
8497		1-3/8" (35 mm)	9" (23 cm)	
8498		1-1/2" (38 mm)	9" (23 cm)	

One-Piece Pouching Systems

Premier



Flat Skin Barriers

Urostomy Pouch

Flextend

Enhanced Design





Cut-to-fit			With Tape Border 	
Product #	Colour	For Stomas Up To	Pouch Length	Box of
84590	Ultra-clear	2-1/2" (64 mm)	9" (23 cm)	10
84490	Beige	2-1/2" (64 mm)	9" (23 cm)	
Pre-sized			With Tape Border 	
Product #	Colour	Barrier Opening	Pouch Length	Box of
84697	Transparent	3/4" (19 mm)	9" (23 cm)	10
84692		1" (25 mm)	9" (23 cm)	
84698		1-1/4" (32 mm)	9" (23 cm)	
84693		1-1/2" (38 mm)	9" (23 cm)	
84699		1-3/4" (44 mm)	9" (23 cm)	
84694		2" (51 mm)	9" (23 cm)	
84497	Beige	3/4" (19 mm)	9" (23 cm)	10
84492		1" (25 mm)	9" (23 cm)	
84498		1-1/4" (32 mm)	9" (23 cm)	
84493		1-1/2" (38 mm)	9" (23 cm)	
84499		1-3/4" (44 mm)	9" (23 cm)	
84494		2" (51 mm)	9" (23 cm)	

Flextend

Original Design



Cut-to-fit			With Tape Border 	
Product #	Colour	For Stomas Up To	Pouch Length	Box of
8450	Ultra-clear	2-1/2" (64 mm)	9" (23 cm)	10
8460	Transparent	2-1/2" (64 mm)	9" (23 cm)	
8440	Beige	2-1/2" (64 mm)	9" (23 cm)	
Pre-sized			With Tape Border 	
Product #	Colour	Barrier Opening	Pouch Length	Box of
8467	Transparent	3/4" (19 mm)	9" (23 cm)	10
8462		1" (25 mm)	9" (23 cm)	
8468		1-1/4" (32 mm)	9" (23 cm)	
8463		1-1/2" (38 mm)	9" (23 cm)	
8469		1-3/4" (44 mm)	9" (23 cm)	
8464		2" (51 mm)	9" (23 cm)	
8447	Beige	3/4" (19 mm)	9" (23 cm)	10
8442		1" (25 mm)	9" (23 cm)	
8448		1-1/4" (32 mm)	9" (23 cm)	
8443		1-1/2" (38 mm)	9" (23 cm)	
8449		1-3/4" (44 mm)	9" (23 cm)	

New!



RESEALABLE

for when Ring needs to be
used more than once



DURABLE

for product protection



60%

**LESS* PLASTIC
PACKAGING**
per barrier ring

>25%

SMALLER* BOX
for improved discretion
and efficient storage

Same trusted barrier rings, **less plastic packaging waste.**

Hollister barrier ring packaging is changing to be more environmentally conscious. The new design reduces plastic waste, furthering Hollister's commitment to sustainability practices.

The plastic savings is equivalent to removing ~132 million plastic water bottles from circulation.**



* Data on File 2025, packaging material analysis. Compared to previous Hollister packaging design. Average calculated across Hollister barrier ring portfolio.

** Data on File, 2025. Compared to previous Hollister packaging design. Average calculated across Hollister flat barrier ring portfolio annual sales 2024-2025.

Accessories

LaTesha Harrison-Thompson

Living with an ostomy since she was 14 years old

"Prior to using the CeraPlus™ Barrier Extenders, I was too scared to go in the water with fear that my pouch would come off. Now I can spend hours in the pool without any worry."



Scan Me
To Learn More About
LaTesha's Journey

Security and Skin Health

Sealed with the confidence of ceramide

The security + confidence of Adapt Accessories meet the science of ceramide



- Infused with ceramide — a natural component of the skin
- Helps protect the skin's natural moisture barrier
- Helps prevent leakage and protect skin—and may even extend barrier wear time
- Helps decrease transepidermal water loss (TEWL) from damaged or eroded skin stists erosion—even when wet or exposed to discharge from the stoma

CeraPlus™ Rings

Flat



Product #	Type	Outer Diameter	Inner Diameter	Thickness	Box of
8805	Standard	2" (48 mm)	3/4" (18 mm)	4.5 mm	10
8815	Slim	2" (48 mm)	3/4" (18 mm)	2.3 mm	
881527	Slim	2-1/4" (57 mm)	1-1/16" (27 mm)	2.3 mm	

Convex



Product #	Inner Diameter	Stretch Up To	Box of
89520	13/16" (20 mm)	1" (25 mm)	10
89530	1-3/16" (30 mm)	1-3/8" (35 mm)	
89540	1-9/16" (40 mm)	1-3/4" (45 mm)	

Oval Convex



Product #	Inner Diameter	Stretch Up To	Box of
89601	7/8" x 1-1/2" (22 x 38 mm)	1-1/8" x 1-3/4" (27 x 43 mm)	10
89602	1-3/16" x 1-7/8" (30 x 48 mm)	1-3/8" x 2-1/8" (35 x 53 mm)	
89603	1-1/2" x 2-3/16" (38 x 56 mm)	1-3/4" x 2-3/8" (43 x 61 mm)	

CeraPlus™ Barrier Extenders



Product #	Usage	Size	Box of
79402	Helps provide extra security around the edges of skin barrier, reduces skin barrier edge lifting and may help increase skin barrier wear time. Provides extra adhesion to accommodate a peristomal hernia.	Each barrier extender covers a surface area of 48 cm ²	30

CeraPlus™ Skin Barrier



Product #	Size	Box of
8900	4 x 4" (10 x 10 cm)	5

Accessories








Adapt™ Barrier Rings

Flat
FlexTend





Product #	Type	Outer Diameter	Inner Diameter	Thickness	Box of
7805	Standard	2" (48 mm)	3/4" (18 mm)	4.5 mm	10
7806	Slim	4" (102 mm)	3/4" (18 mm)	2.3 mm	
7815	Slim	2" (48 mm)	3/4" (18 mm)	2.3 mm	

Adapt™ Accessories




Product #		Description	Usage		Size	Box of
7737		Adapt Medical Adhesive Remover Spray <i>*featuring no-sting, silicone-based, alcohol free formulation</i>	360° spray technology to deliver controlled, quiet spraying from any angle, with minimal waste	Spray Can	1.7 fl oz (50 ml)	1
7760		Universal Remover Wipes <i>*featuring no-sting silicone-based, alcohol free formulation</i>	This skin friendly product helps remove adhesive residues.	50 Wipes per Carton		50
7917		Adapt Skin Protective Wipe <i>*featuring an alcohol-free formulation improved with the addition of aloe vera</i>	Protects skin during adhesive removal Not recommended for use with extended wear skin barriers	50 Wipes per Carton		50
7906		Adapt Stoma Powder	Use to absorb moisture from broken skin around the stoma. Dust on. Brush off excess powder. Discontinue use once the condition is resolved.	Puff Bottle	1 oz (28 g)	1
79301		Adapt Paste	Used to fill or caulk uneven skin contours to create a flatter surface. It can help to prevent drainage from getting under the ostomy skin barrier.	Tube	.5 oz (14 g)	20
79300					2 oz (60 g)	1
79400		Adapt Barrier Strips with FlexTend Barrier (Alcohol Free)	Used to fill in or caulk uneven skin contours near the stoma to create a flatter surface and help prevent ostomy drainage from getting under the ostomy barrier.	Triangular Shape	5-1/8" (13 cm)	10
78500		Adapt Lubricating Deodorant	Neutralizes odours, eases emptying, and helps prevent pouch static and sticking. It is not a perfume or masking agent. Should be used in drainable and closed pouching systems only.	Bottle	8 oz (236 ml)	1
78501				Packets	.27 oz (8 ml)	50

Accessories

m9 Odour Eliminators

Product #		Description	Size	Box of
7715		Drops	1 oz (30 ml)	12
7717			8 oz (240 ml)	6
7732		Spray (Unscented)	2 oz (60 ml)	12
7733			8 oz (240 ml)	6
7735		Spray (Apple Scent)	8 oz (240 ml)	

Pouch Accessories

Product #		Description	Size	Box of
8770		Drainable Pouch Clamp, Beige		20
7300		Beige Colour Ostomy Belt	Medium 23" - 43" (58 - 109 cm)	10
7299			Large 34" - 65" (86 - 165 cm)	
7331		Urostomy Drain Tube Adapter		10

Skin Sheets



Product #	Description	Size	Box of
8800	Flexextend	4" x 4" (10 cm x 10 cm)	5
8801		8" x 8" (20 cm x 20 cm)	3
7800	Premium	4" x 4" (10 cm x 10 cm)	5
7801		8" x 8" (20 cm x 20 cm)	3
7700	HolliHesive	4" x 4" (10 cm x 10 cm)	5

High Output Collection Bag



Product #	Description	Box of
5551	2L Collection Bag	10

* to be used with the following high output pouches:
18022, 18023, 18024, 83025, 83038, 83055, 83125, 83130, 83135, 26064b and 26364

Irrigation Products



Product #	Description	Information		Box of
7721	Cone Irrigator Kit	Includes: (1) Irrigator (1) Stoma Cone with Connector (1) Cleaning Brush		1
7723	Stoma Cone	With Connector Replacement Unit		10
7728	Irrigator Sleeve	With Attached Closure and Belt Tabs	2" (51 mm)	20
7724			3" (76 mm)	
18152	Two-Piece Irrigator Sleeve	Flange: 1-3/4" (44 mm)	Pouch: 30" (76 cm)	5
18153		Flange: 2-1/4" (57 mm)	Pouch: 30" (76 cm)	
18154		Flange: 2-3/4" (70 mm)	Pouch: 30" (76 cm)	

Wound Drainage Collectors



Product #	For Wounds Up To		Box of
9773	3" (7.6 cm)	FlexWear (Non-sterile)	10
9776	3.75" (9.5 cm)		10
9778	4" x 8" (10 cm x 20 cm)		3

A photograph of a man and a young child with curly hair, both smiling and engaged in cooking. The man is wearing a white t-shirt and a grey apron, and the child is wearing a yellow t-shirt and a white apron. They are standing in a kitchen with white tiled walls and wooden shelves. On the shelves are several glass jars filled with various ingredients. In the foreground, there is a wooden cutting board with vegetables, a bowl of tomatoes, and a pot of spaghetti. The man is holding the child's hand, and they are both looking at the food on the counter.

Pouchkins™


With the Gentleness and Flexibility of SoftFlex

Pouchkins

Pediatric Pouching System

One-Piece SoftFlex

Premie/Newborn Pouches



Premie Newborn

Cut-to-fit Without Tape Border

Product #	Description	For Stomas Up To	Pouch Length	Box of
Premie				
3777	Ultra-clear, no starter hole	5/8" (15 mm)	4-1/2" (11.5 cm)	15
Newborn				
3778	Ultra-clear, no starter hole	1-3/8" x 7/8" (35 x 22 mm)	6-1/2" (16 cm)	

Drainable Pouches



Cut-to-fit Without Tape Border

Product #	Description	For Stomas Up To	Pouch Length	Box of
3795	Transparent, no starter hole	2" (51 mm)	7" (18 cm)	10
3796		1-1/2" (38 mm)	6-1/2" (16 cm)	

Urostomy Pouch



Cut-to-fit Without Tape Border

Product #	Description	For Stomas Up To	Pouch Length	Box of
3797	Transparent, no starter hole	1-1/2" (38 mm)	6-1/2" (16 cm)	10

Pouchkins SoftFlex Skin Barrier:

- Designed to be gentle on sensitive infant skin
- No starter hole, which enables the opening to be placed where needed
- No tape by design
- Flexible and sized for babies and children
- None of the *Pouchkins* pediatric pouching systems are made with natural rubber latex
- Any 1-3/4" (44mm) two-piece pouch will fit this skin barrier

Two-Piece SoftFlex

Flat Skin Barrier



Cut-to-fit Without Tape Border

Product #	Description	Flange Size	For Stomas Up To	Box of
3761	No starter hole	1-3/4" (44 mm)	1-1/4" (32 mm)	5

Colour: Transparent

Product #	Flange Size	Pouch Length	Filter	Box of
3799	1-3/4" (44 mm)	6-1/2" (16 cm)	No	10

Colour: Beige

18282	1-3/4" (44 mm)	7" (18 cm)	AF300	20
--------------	----------------	------------	-------	----

Urostomy Pouches



Enhanced Design

Product #	Colour	Flange Size	Pouch Length	Box of
18922	Ultra-clear	1-3/4" (44 mm)	9" (23 cm)	10
18912	Beige	1-3/4" (44 mm)	9" (23 cm)	10

Ostomy Belt



Product #	Length	Box of
3774	10 - 17" (25 - 43 cm)	1

Shadow Buddies dolls

- Exclusive to **Hollister**, *Shadow Buddies* dolls are special friends for children who have had or will have ostomy surgery
- Great teaching tool to help with education and care giving
- Huggable dolls with heart-shaped eyes and a stoma provide comfort and support



Shadow Buddies is a trademark of the Shadow Buddies Foundation

Continence Care

A photograph of two women, one older with grey hair and glasses, and one younger with brown hair, both smiling and looking at a laptop screen. The older woman is wearing a brown collared shirt, and the younger woman is wearing a beige knit sweater. The background is a soft-focus indoor setting with a plant visible on the left.

People are at the center of everything we do.

The improved Infyna Chic™ Intermittent Catheter

New, clear tube

to aid visualisation.

Longer catheter tube

for confident, full bladder drainage.

More lubricious coating

for more comfortable insertion.



HydraBalance™
Lubricating Technology



Recyclable casing



Made with 11%
less plastic







PVC free

A brighter choice

We know that even the small details can make a difference to how someone feels, which is why Infyna Chic™ catheter is available in French sizes 8-14 in 4 pastel colour funnels.



Size	Catheter length	Tip	SKU	Quantity per box
 8 Fr	14cm	Straight	7008	30
 10 Fr	14cm	Straight	7010	30
 12 Fr	14cm	Straight	7012	30
 14 Fr	14cm	Straight	7014	30



To learn more, visit www.Hollister.com/en-ca/InfynaChic

A photograph of a man and a woman in a lush forest. The man, wearing a yellow jacket and blue jeans, is sitting on a mossy rock. The woman, wearing a red jacket and dark pants, is crouching in a stream, splashing water. A cluster of white, translucent bubbles is floating in the air between them. A red rectangular box with white text is overlaid on the lower left of the image.

**The catheter glides effortlessly for ease of use
giving users confidence and more time for living.¹**

Jonny uses VaPro™ catheters
Katie uses Infyna Chic™ catheters

Experience the ease of use² and comfort³ of our HydraBalance™ Lubricating Technology

Made with naturally sourced ingredients, its unique combination of hydrophilic coating and hydration fluid creates a smoothly gliding hydrogel cushion that not only helps protect the urethra, but also provides smooth and effortless insertion and withdrawal for a more confident and natural catheter experience.^{4, 5}

Key features of our HydraBalance™ Lubricating Technology:

- **Smooth and comfortable on catheter insertion and withdrawal.⁶**
- **Helps protect the urethra from potential micro trauma and Urinary Tract Infections.⁷**
- **Designed to capture the right balance of naturally sourced ingredients for as natural of an intermittent catheter experience as possible.⁸**





Touch Free Hydrophilic Intermittent Catheters

VaPro Catheters Provide 100% No Touch Protection with a Protective Tip and Sleeve:

VaPro catheters are hygienic by design to help protect patients from germs throughout the entire catheterization process.

Protective tip:

- helps shield the sterile catheter during insertion from bacteria located within the first 15 mm of the distal urethra
- helps reduce the risk of carrying bacteria into the urinary tract

Protective sleeve:

- allows for catheter to be gripped anywhere
- provides a barrier to the catheter to help keep germs away



VaPro™



Product #	Description	Size	Length	Box of
72082	Straight tip	8 Fr	8" (20 cm)	30
72102	Straight tip	10 Fr	8" (20 cm)	30
72122	Straight tip	12 Fr	8" (20 cm)	30
72084	Straight tip	8 Fr	16" (40 cm)	30
72104	Straight tip	10 Fr	16" (40 cm)	30
72124	Straight tip	12 Fr	16" (40 cm)	30
73124	Coudé tip		16" (40 cm)	
72142	Straight tip		8" (20 cm)	
72144	Straight tip	14 Fr	16" (40 cm)	30
73144	Coudé tip		16" (40 cm)	
72164	Straight tip	16 Fr	16" (40 cm)	30
73164	Coudé tip		16" (40 cm)	

VaPro Plus™



Product #	Description	Size	Length	Box of
74122	Straight tip	12 Fr	8" (20 cm)	30
74142	Straight tip	14 Fr	8" (20 cm)	30
74124	Straight tip	12 Fr	16" (40 cm)	30
74144	Straight tip	14 Fr	16" (40 cm)	30

VaPro Plus is a closed system with an integrated 1000 ml collection bag.

VaPro F-Style is a **FIRMER** Touch Free Hydrophilic Catheter.

VaPro™ F-Style



Product #	Description	Size	Length	Box of
7600084	Straight tip	8 Fr	16" (40 cm)	30
7600104	Straight tip	10 Fr	16" (40 cm)	30
7600124	Straight tip	12 Fr	16" (40 cm)	30
7600144	Straight tip	14 Fr	16" (40 cm)	30
7600164	Straight tip	16 Fr	16" (40 cm)	30


VaPro Plus™ F-Style



Product #	Description	Size	Length	Box of
7700124	Straight tip	12 Fr	16" (40 cm)	30
7700144	Straight tip	14 Fr	16" (40 cm)	30

VaPro Plus F-Style is a closed system with an integrated 1000 ml collection bag.

VaPro™ Pocket



VaPro Pocket™ catheters offer the features that patients want.¹

Protection. Ease of Use. Discretion.

VaPro Pocket catheters provide 100% No Touch Protection with a protective tip and sleeve, and they are easy to use, with an improved protective tip that offers a more comfortable insertion. The pocket-size package design is more discreet—a quality users prefer.

Consider VaPro Pocket catheters right from the start.



VaPro™ Pocket

Product #	Description	Size	Length	Box of
70104	Straight tip	10 Fr	16" (40 cm)	30
70124	Straight tip	12 Fr	16" (40 cm)	30
70144	Straight tip	14 Fr	16" (40 cm)	30
70164	Straight tip	16 Fr	16" (40 cm)	30
77124	Coudé tip	12 Fr	16" (40 cm)	30
77144	Coudé tip	14 Fr	16" (40 cm)	30
77164	Coudé tip	16 Fr	16" (40 cm)	30

VaPro™ Plus Pocket

Product #	Description	Size	Length	Box of
71084	Straight tip	8 Fr	16" (40 cm)	30
71104	Straight tip	10 Fr	16" (40 cm)	30
71122	Straight tip	12 Fr	8" (20 cm)	30
71124	Straight tip		16" (40 cm)	
71142	Straight tip	14 Fr	8" (20 cm)	30
71144	Straight tip		16" (40 cm)	
71164	Straight tip	16 Fr	16" (40 cm)	30

VaPro Plus Pocket is a closed system with an integrated 1000 ml collection bag.

Touch Free Intermittent Catheters

Closed System



Product #	Description	Size	Length	Box of
94064IC	Straight tip	6 Fr	16" (40 cm)	100
94084IC	Straight tip	8 Fr	16" (40 cm)	
94104IC	Straight tip	10 Fr	16" (40 cm)	
94124IC	Straight tip	12 Fr	16" (40 cm)	
95124IC	Coudé tip		16" (40 cm)	
94144IC	Straight tip	14 Fr	16" (40 cm)	
95144IC	Coudé tip		16" (40 cm)	
94164IC	Straight tip	16 Fr	16" (40 cm)	
95164IC	Coudé tip		16" (40 cm)	
94184IC	Straight tip	18 Fr	16" (40 cm)	

Closed System Kit



Product #	Description	Size	Length	Box of
96064IC	Straight tip	6 Fr	16" (40 cm)	100
96084IC	Straight tip	8 Fr	16" (40 cm)	
96104IC	Straight tip	10 Fr	16" (40 cm)	
96124IC	Straight tip	12 Fr	16" (40 cm)	
97124	Coudé tip		16" (40 cm)	
96144IC	Straight tip	14 Fr	16" (40 cm)	
97144	Coudé tip		16" (40 cm)	
96164IC	Straight tip	16 Fr	16" (40 cm)	
97164	Coudé tip		16" (40 cm)	
96184IC	Straight tip	18 Fr	16" (40 cm)	


Kit Includes:

- Advance Plus Closed System
- Large Gloves
- Benzalkonium Chloride antiseptic swab (BZK)
- Waterproof underpad

Advance Plus Touch Free Intermittent Catheters:

- Have a protective tip that helps shield the catheter from bacteria located in the first 15 mm of the distal urethra
- Have a user-regulated gel reservoir which allows for control of the catheter lubrication
- Smooth eyelets help provide trouble-free insertion and withdrawal





Intermittent Catheters



Product #	Description	Size	Length	Box of
10610	Straight tip	6 Fr	10" (25 cm)	30
10616	Straight tip	6 Fr	15" (38 cm)	
10806	Straight tip	8 Fr	6" (15 cm)	30
10810	Straight tip		10" (25 cm)	
10816	Straight tip		16" (40 cm)	
10826	Coudé tip	8 Fr	16" (40 cm)	
11006	Straight tip	10 Fr	6" (15 cm)	30
11010	Straight tip		10" (25 cm)	
11026	Coudé tip		16" (40 cm)	
11016	Straight tip	12 Fr	16" (40 cm)	30
11206	Straight tip		6" (15 cm)	
11216	Straight tip		16" (40 cm)	
11226	Coudé tip	14 Fr	16" (40 cm)	30
11406	Straight tip		6" (15 cm)	
11446	Straight tip, No funnel		6" (15 cm)	
11416	Straight tip		16" (40 cm)	
B1061	Soft tip		16" (40 cm)	
11426	Coudé tip	16 Fr	16" (40 cm)	30
B1065	Straight tip with curved packaging		16" (40 cm)	
11616	Straight tip	18 Fr	16" (40 cm)	30
11626	Coudé tip		16" (40 cm)	
11816	Straight tip	18 Fr	16" (40 cm)	30
11826	Coudé tip	18 Fr	16" (40 cm)	

Features of the Apogee Family of Intermittent Catheters:

- Two staggered smooth eyelets help provide trouble-free insertion
- Not made with natural rubber latex
- Designed for comfort
- Medical grade Polyvinyl Chloride (PVC)
- Non lubricated, clear catheter
- Straight and coudé tip options available

Colour	Size French	Size Millimeter
	6	2.0
	8	2.7
	10	3.3
	12	4.0
	14	4.7
	16	5.3
	18	6.0

Male External Catheters

Hollister

Latex




Extended Wear with Unique Inner Flap			Sheath length 2-3/4" (7 cm), adhesive length 1-3/5" (4 cm)
Product #	Description	Size	Box of
9206	Small/Pediatric	7/8"- 1" (22 - 25 mm)	30
9207	Medium	1"- 1-1/8" (26 - 30 mm)	
9208	Intermediate	1-1/8"- 1-1/4" (31 - 35 mm)	
9209	Large	1-3/8"- 1-1/2" (36 - 39 mm)	
Everyday Wear			Sheath length 2-3/4" (7 cm), adhesive length 1-3/5" (4 cm)
Product #	Description	Size	Box of
9106	Small/Pediatric	7/8"- 1" (22 - 25 mm)	30
9107	Medium	1"- 1-1/8" (26 - 30 mm)	
9108	Intermediate	1-1/8"- 1-1/4" (31 - 35 mm)	
9109	Large	1-3/8"- 1-1/2" (36 - 39 mm)	

Contains Natural Rubber Latex

Size Selection Guide

Use the Hollister Sizing Guide to determine the most appropriate size catheter

- Place sizing guide just behind the glans of the penis
- Select the size that most closely fits



InView™

Silicone



Standard			Sheath length 2-3/4" (7 cm), adhesive length 1-3/5" (4 cm)
Product #	Description	Size	Box of
97525	Small	1" (25 mm)	30
97529	Medium	1-1/8" (29 mm)	
97532	Intermediate	1-1/4" (32 mm)	
97536	Large	1-1/2" (36 mm)	
97541	Extra Large	1-5/8" (41 mm)	100
97525-100	Small	1" (25 mm)	
97529-100	Medium	1-1/8" (29 mm)	
97532-100	Intermediate	1-1/4" (32 mm)	
97536-100	Large	1-1/2" (36 mm)	
97541-100	Extra Large	1-5/8" (41 mm)	

Special			Sheath length 1-3/5" (4 cm), adhesive length 1-3/5" (4 cm)
Product #	Description	Size	Box of
97425	Small	1" (25 mm)	30
97429	Medium	1-1/8" (29 mm)	
97432	Intermediate	1-1/4" (32 mm)	
97436	Large	1-1/2" (36 mm)	
97441	Extra Large	1-5/8" (41 mm)	100
97425-100	Small	1" (25 mm)	
97429-100	Medium	1-1/8" (29 mm)	
97432-100	Intermediate	1-1/4" (32 mm)	
97436-100	Large	1-1/2" (36 mm)	
97441-100	Extra Large	1-5/8" (41 mm)	

Extra			Sheath length 2-3/4" (7 cm), adhesive length 2-3/4" (7 cm)
Product #	Description	Size	Box of
97625	Small	1" (25 mm)	30
97629	Medium	1-1/8" (29 mm)	
97632	Intermediate	1-1/4" (32 mm)	
97636	Large	1-1/2" (36 mm)	
97641	Extra Large	1-5/8" (41 mm)	100
97625-100	Small	1" (25 mm)	
97629-100	Medium	1-1/8" (29 mm)	
97632-100	Intermediate	1-1/4" (32 mm)	
97636-100	Large	1-1/2" (36 mm)	
97641-100	Extra Large	1-5/8" (41 mm)	

Not Made with Natural Rubber Latex

Urinary Collection Systems

Leg Bags

Leg Bag



Product #	Description	Capacity	Box of
9814	Medium	18 oz (540 ml)	10
9805	Large	30 oz (900 ml)	

Combination Pack



Product #	Description	Capacity	Box of
9348	Medium	18 oz (540 ml)	12 packs
9349	Large	30 oz (900 ml)	

Extension Tubing with Connector



Product #	Description	Length	Box of
9345	Non-sterile	18" (45 cm)	10
9346	Sterile	18" (45 cm)	

Vented Combination Pack



Product #	Description	Capacity	Box of
9645	Medium	18 oz (540 ml)	10 packs
9655	Large	30 oz (900 ml)	

Fabric Strap



Product #	Description	Capacity	Box of
9342	Small	15" (38 cm)	10 pairs
9343	Medium	23" (58 cm)	

Bedside Drainage



Product #	Description	Capacity	Box of
9839	2000 ml Collection System	68 oz (2000 ml)	20

Hollister Leg Bags Are:

- Sterile
- Not Made with Natural Rubber Latex
- Designed with an anti-reflux valve

m9 Cleaner/Decrystallizer



Concentrate & Wash Bottle

Product #	Description	Size	Box of
7736	1 Bottle of Concentrate 1 Washing Bottle	16 oz (480 ml)	12

External Collection Devices

Urinary Pouches

Female

SoftFlex Skin Barrier



Product #	Description	Box of
9840	Length 7-1/2" (19 cm) Width 4-1/2" (11cm)	10

Male

Flextend Skin Barrier

Thin barrier for flexibility, no additional tape needed



Product #	Description	Box of
9873	Length 7-1/2" (19 cm) Width 4-1/2" (11cm)	10

SoftFlex Skin Barrier



Product #	Description	Box of
9811	Length 7-1/2" (19 cm) Width 4-1/2" (11cm)	10

Fecal Collectors

Pouches

Flextend Skin Barrier



Product #	Description	Pouch Length	Box of
9823	Medium, Ultra-clear	10" (25 cm)	10
9880	Large, Transparent	12" (30 cm)	

Transparent Pouches


SoftFlex Skin Barrier



Product #	Description	Pouch Length	Box of
9822	Medium	10" (25 cm)	10
9821	Large	12" (30 cm)	

Hollister's Urinary Pouches:

- Not Made with Natural Rubber Latex
- Have a cut-to-fit skin barrier
- Have a odour-barrier pouch film



VaPro™ No Touch
Intermittent Catheters

Study¹ Confirms Superior Protection*

Hollister asked an independent laboratory to conduct a study to see if the **superior protection** of VaPro™ No Touch Intermittent Catheters could be measured. The study compared VaPro catheters to three other brands.

*In comparison to other tested catheters.

Reference: 1. Hollister data on file, TR-00453, 2021.

Find Out What
Happened



Continence Care

A photograph of two healthcare professionals, a woman with dark hair and a woman with short curly hair, both wearing blue scrubs. They are standing in a brightly lit room with large windows in the background. The woman on the left is holding a tablet, and the woman on the right is pointing at the screen. The image has a soft, professional feel with a slight blur in the background.

Critical Care

Dedicated to excellence.

Oral Endotracheal Tube Fasteners

AnchorFast™



Product #	Description	Box of
9799	<ul style="list-style-type: none"> ET tube wrap with one-click security clamp locks the ET tube in place Gliding tube shuttle allows the ET tube to be easily repositioned from side to side while the tube remains secure Flextend Skin Barrier For tubes sized 5.0-10.0mm inner diameter 	12

AnchorFast SlimFit™



Product #	Description	Box of
9787	<ul style="list-style-type: none"> ET tube wrap with one-click security Gliding tube shuttle allows for side-to-side repositioning of the ET tube Thumb pad with increased surface area for finger or thumb placement when opening Designed to accommodate smaller faces Secures ET tubes ranging in size 5.0-8.0 mm 	12

AnchorFast Guard Select™



Product #	Description	Box of
9700	<ul style="list-style-type: none"> ET tube wrap with one-click security clamp locks the ET tube in place Gliding tube shuttle allows the ET tube to be easily repositioned from side to side while the tube remains secure Integrated tube protection sleeve designed to help prevent occlusion Secures oral ET tubes ranging in size from 5.0 to 9.0 mm inner diameter and ET tubes with integrated subglottic suction lumen sizes 6.0 to 8.0 mm inner diameter 	12

AnchorFast Replacement Neck Band Straps



Product #	Description	Box of
9788	Replacement Neck Band Straps	4

Tube Attachment Devices

Horizontal Tube (HTAD)



Product #	Description	Box of
9781	Single use, sterile	5

Vertical Tube (VTAD)



Product #	Description	Box of
9782	Single use, sterile	5

Drain Tube (DTAD)



Product #	Description	Box of
9780	Single use, sterile	5

Feeding Tube (FTAD)



Product #	Description	Box of
9786	<ul style="list-style-type: none"> FlexWear skin barrier conforms to patient's nose Single use, sterile 	20

Suction Tube (STAD)



Product #	Description	Box of
9785	<ul style="list-style-type: none"> FlexWear skin barrier conforms to patient's nose Single use, sterile 	20

Numerical Index

Product #	Page	Product #	Page	Product #	Page	Product #	Page	Product #	Page	Product #	Page
1796	28	8248	26	8484	30	8801	36	9776	36	11504	10
3184	28	8274	26	8485	30	8805	34	9778	36	11505	10
3186	28	8278	26	8486	30	8815	34	9780	51	11506	10
3761	38	8280	26	8487	30	8900	34	9781	51	11616	45
3774	38	8281	26	8488	30	8901	12	9782	51	11626	45
3777	38	8282	26	8489	30	8914	12	9785	51	11702	10
3778	38	8283	26	8490	30	8918	12	9786	51	11703	10
3795	38	8284	26	8491	30	8925	12	9787	51	11704	10
3796	38	8285	26	8492	30	8930	12	9788	51	11816	45
3797	38	8286	26	8493	30	8931	12	9799	51	11826	45
3799	38	8288	26	8494	30	8935	12	9805	47	11902	10
5551	36	8289	26	8495	30	8940	12	9811	48	11903	10
7008	40	8290	26	8496	30	8954	12	9814	47	11904	10
7010	40	8291	26	8497	30	8958	12	9821	48	11905	10
7012	40	8292	26	8498	30	8960	12	9822	48	11906	10
7014	40	8293	26	8531	27	8961	12	9823	48	11907	10
7299	36	8294	26	8538	27	8962	12	9839	47	11908	10
7300	36	8295	26	8551	27	8963	12	9840	48	11909	10
7331	36	8296	26	8552	27	8964	12	9873	48	11910	10
7700	36	8331	27	8558	27	8968	12	9880	48	11911	10
7715	36	8412	13	8631	27	8984	12	10610	45	13402	23
7717	36	8413	13	8632	27	8988	12	10616	45	13403	23
7721	36	8414	13	8633	27	8990	12	10806	45	13404	23
7723	36	8415	13	8635	27	8991	12	10810	45	13503	23
7724	36	8440	31	8637	27	8992	12	10816	45	13504	23
7728	36	8442	31	8638	27	8993	12	10826	45	13505	23
7732	36	8443	31	8639	27	9106	46	11006	45	13506	23
7733	36	8447	31	8641	27	9107	46	11010	45	13901	22
7735	36	8448	31	8642	27	9108	46	11016	45	13902	22
7736	47	8449	31	8643	27	9109	46	11026	45	13903	22
7737	35	8450	31	8647	27	9206	46	11202	11	13904	22
7760	35	8460	31	8648	27	9207	46	11203	11	13905	22
7800	36	8462	31	8649	27	9208	46	11204	11	13906	22
7801	36	8463	31	8651	27	9209	46	11206	45	13907	22
7805	35	8464	31	8661	26	9342	47	11216	45	13908	22
7806	35	8467	31	8662	26	9343	47	11226	45	13909	22
7815	35	8468	31	8663	26	9345	47	11402	10	13910	22
7906	35	8469	31	8664	26	9346	47	11403	10	13911	22
7917	35	8474	30	8665	26	9348	47	11404	10	14202	23
8031	27	8478	30	8666	26	9349	47	11406	45	14203	23
8081	27	8480	30	8674	26	9645	47	11416	45	14204	23
8181	27	8481	30	8678	26	9655	47	11426	45	14206	23
8214	26	8482	30	8770	36	9700	51	11446	45	14301	23
8244	26	8483	30	8800	36	9773	36	11503	10	14302	23

Numerical Index

Product #	Page	Product #	Page	Product #	Page	Product #	Page	Product #	Page	Product #	Page
14303	23	15704	10	18132	16	18902	18	72102	42	82617	26
14304	23	15802	22	18133	16	18903	18	72104	42	82618	26
14305	23	15803	22	18134	16	18904	18	72122	42	82619	26
14306	23	15804	22	18142	16	18912	18	72124	42	82624	26
14307	23	15903	22	18143	16	18912	38	72142	42	82628	26
14308	23	15904	22	18144	16	18913	18	72144	42	82711	26
14309	23	15905	22	18152	36	18914	18	72164	42	83025	29
14602	22	15906	22	18153	36	18922	18	73124	42	83038	29
14603	22	16101	22	18154	36	18922	38	73144	42	83055	29
14604	22	16102	22	18172	16	18923	18	73164	42	83125	29
14606	22	16103	22	18173	16	18924	18	74122	42	83130	29
14701	22	16104	22	18174	16	19002C	11	74124	42	83135	29
14702	22	16105	22	18176	16	19003C	11	74142	42	84134	13
14703	22	16106	22	18182	16	19004C	11	74144	42	84138	13
14704	22	16107	22	18183	16	19006	17	77124	43	84200	13
14705	22	16108	22	18184	16	19054	17	77144	43	84292	13
14706	22	16401	23	18192	16	19202C	11	77164	43	84293	13
14707	22	16402	23	18193	16	19203C	11	78500	35	84294	13
14708	22	16403	23	18194	16	19204C	11	78501	35	84295	13
14709	22	16404	23	18202	16	19254	17	79300	35	84398	30
14802	22	16405	23	18203	16	19603	11	79301	35	84404	13
14803	22	16406	23	18204	16	19604	11	79400	35	84408	13
14804	22	16407	23	18282	16	19612	11	79402	34	84490	31
15102	11	16408	23	18282	38	19613	11	80070	29	84492	31
15103	11	18002	16	18283	16	19614	11	80110	29	84493	31
15104	11	18003	16	18284	16	19903	11	82100	28	84494	31
15202	23	18004	16	18292	16	19912	11	82111	26	84497	31
15203	23	18006	16	18293	16	19913	11	82125	28	84498	31
15204	23	18012	17	18294	16	19914	11	82130	28	84499	31
15302	10	18013	17	18352	16	26064b	29	82135	28	84590	31
15303	10	18014	17	18353	16	26364	29	82300	28	84692	31
15304	10	18016	17	18354	16	70104	43	82302	28	84693	31
15402	23	18022	17	18362	16	70124	43	82325	28	84694	31
15403	23	18023	17	18363	16	70144	43	82330	28	84697	31
15404	23	18024	17	18364	16	70164	43	82335	28	84698	31
15503	23	18102	16	18372-30	16	71084	43	82402	28	84699	31
15504	23	18103	16	18373-30	16	71104	43	82411	26	84711	30
15505	23	18104	16	18374-30	16	71122	43	82610	26	84794	30
15506	23	18112	16	18392	16	71124	43	82611	26	84798	30
15602	22	18113	16	18393	16	71142	43	82612	26	84890	30
15603	22	18114	16	18394	16	71144	43	82613	26	84891	30
15604	22	18122	16	18402	18	71164	43	82614	26	84892	30
15702	10	18123	16	18403	18	72082	42	82615	26	84893	30
15703	10	18124	16	18404	18	72084	42	82616	26	84894	30

Numerical Index

Product #	Page	Product #	Page	Product #	Page
84895	30	94164IC	44	112013	11
84896	30	94184IC	44	112014	11
84897	30	95124IC	44	112015	11
84898	30	95144IC	44	112016	11
84899	30	95164IC	44	112017	11
84990	30	96064IC	44	112018	11
84991	30	96084IC	44	151113	11
84992	30	96104IC	44	151114	11
84993	30	96124IC	44	151115	11
84994	30	96144IC	44	151116	11
84995	30	96164IC	44	151117	11
84996	30	96184IC	44	151118	11
84997	30	97124	44	153103	10
84998	30	97144	44	153104	10
84999	30	97164	44	153105	10
86711	26	97425	46	153106	10
88300	28	97425-100	46	157102	10
88302	28	97429	46	157103	10
88325	28	97429-100	46	157104	10
88330	28	97432	46	157105	10
88335	28	97432-100	46	157106	10
88340	28	97436	46	157107	10
88402	28	97436-100	46	157108	10
88700	28	97441	46	157109	10
88725	28	97441-100	46	157110	10
88730	28	97525	46	157111	10
88735	28	97525-100	46	826211	26
88740	28	97529	46	841311	13
89011	12	97529-100	46	844011	13
89111	12	97532	46	881527	34
89511	12	97532-100	46	7600084	42
89518	12	97536	46	7600104	42
89520	34	97536-100	46	7600124	42
89530	34	97541	46	7600144	42
89540	34	97541-100	46	7600164	42
89601	34	97625	46	7700124	42
89602	34	97625-100	46	7700144	42
89603	34	97629	46	8195111	12
89611	12	97629-100	46	8439111	30
89811	12	97632	46	8479111	30
94064IC	44	97632-100	46	B1061	45
94084IC	44	97636	46	B1065	45
94104IC	44	97636-100	46		
94124IC	44	97641	46		
94144IC	44	97641-100	46		

General Information

Terms of Payment

Net 30 days from date of invoice. Electronic Fund Transfers (EFT) are possible. Please contact: CanadaContact@Hollister.com

Minimum Order

The minimum order requirement is \$250.00.

A charge of \$10.00 will be applied to any order below the minimum requirement.

Shipping

Prices for shipments to destinations in Canada are FOB buyer's address.

If the method and route of shipment are specified by the buyer, the buyer is required to pay the delivery charges, upon delivery of products to the carrier, all risk of loss and other incidents of ownership pass to the buyer. These terms may not be altered.

Pricing & Terms of Sale

GST/HST will be charged where applicable. Buyer agrees that they will not at anytime, directly or indirectly, sell or distribute, to customers which are located outside of Canada.

Hollister Return Goods Policy

All returns for credit require written authorization from Hollister Limited. A Return Goods Authorization can be obtained by contacting a Customer Service representative at Hollister Limited via email at canadacustomerservice@hollister.com. Only unopened, full boxes of items sold and fulfilled by Hollister Limited may be returned. Items must be in re-sellable condition, with no external labelling or markings. Customer has **14 days** from the date of receipt to return for a full refund.

- A 15% restock fee will be applied to all returned goods over **14 days** from the date of receipt.
- No returns allowed on any goods over **30 days** from date of receipt.
- Copy of invoice or purchase order number must be provided for the return request.
- Product SKU and Lot# must match PO.
- Returned goods must match what is on the approved Return Goods Authorization.
- Any goods included in the return not part of the approved Return Goods Authorization will be refused and returned to sender at the customer's expense.
- All returns must be shipped to address listed on the Return Goods Authorization.
- Returns shipped to any other address will not be credited and will be returned at the customer's expense.
- Unauthorized returns will not be accepted. Unauthorized returns will be refused and returned at the customer's expense.

www.hollister.com

The Hollister logo, Adapt, Advance, Advance Plus, Anchorfast, Anchorfast Slim Fit, Anchorfast Guard Select, Apogee, CeraPlus, Extended Wear, Everyday, Flexend, FlexWear, Infyna Chic, InView, m9, New Image, Pouchkins, Premier, Secure Start, SoftFlex, VaPro, VaPro F-Style, VaPro Plus, VaPro Pocket and VaPro Plus Pocket are trademarks of Hollister Incorporated.

All other trademarks and copyrights are the property of their respective owners.

©2025 Hollister Incorporated

MRP15000251

Damaged Product

- Customer must contact Hollister Limited to receive a Return Goods Authorization within 10 days of receipt.
- **Customer must receive written approval and has 30 days to return Damaged Product.**
- Receiving dock must show in writing any package received damaged for Hollister to claim from Freight Company.
- Customer must provide copy of invoice or purchase order number for the return request.
- Credit will not be issued if customer disposes the damaged product without approval by Customer Service.
- If product cannot be returned due to extreme damage (leakage/broken) customer must provide pictures to Customer Service within **10 days** of receipt to get approval to dispose of and thereby obtain credit for the same.

Shortages

- Customer must notify Hollister Limited Customer Service within **10 days** of receipt on shortages. Receiving dock must show in writing any package received short.

Ostomy & Continence Retail Sample Policy

All sample requests at the Retail level are to be processed through Hollister Customer Service by calling 1.800.263.7400 or by email canadasamples@hollister.com.

Retailer/client is to contact Customer Service and provide consumer information for sample processing. Samples will be shipped directly to consumers from Hollister Limited. Samples will **no longer** be available to ship directly to a retail setting.

Prior to Use

Prior to use, be sure to read the Instructions for Use for information regarding intended use, contraindications, warnings, precautions and instructions.

Latex Statement

The specifications for the Hollister Continence Care products listed below include natural rubber latex as a component; our suppliers add natural rubber latex as part of the production process. If a Hollister product(s) found in this catalog is not listed below, the specification for that product(s) and packaging does not include natural rubber latex as a component, and our component suppliers do not add natural rubber as part of their production process. Product # 9206, 9207, 9208, 9209, 9106, 9107, 9108, 9109



*CeraPlus products contain the Remois technology of Alcare Co., Ltd.



Hollister Secure Start Services

Enroll on the Go! It's easy and free!

Hollister Secure Start™ Services provides you with a lifetime of personalized support, so you can live your life, your way!

Enroll today to receive:

- One-to-one personalized support
- Product-related care tips
- Educational Resources
- Product samples that best fit your need
- Finding reimbursement for your products

Start enjoying the benefits of Secure Start services today!



1-866-789-7574



1-800-263-7421

