

# Making a Difference in the Journey of Life

2020 Product Catalogue

## Contact Information

### Customer Service

Toll Free: 1.800.263.7400

Local: 905.727.4344

### Fax

Toll Free: 1.800.432.8846

Local: 905.727.1614

### Email

canadacustomerservice@hollister.com

### Business Hours

Monday – Friday; 8:30 am - 5:00 pm (EST)

Excluding Holidays

### Hollister Limited

Aurora, ON



# Hollister Secure Start Services

Hollister Secure Start Services provides a lifetime of customized support for people in the ostomy and continence care communities.

We are proud to offer dedicated support for each and every patient along the continuum of care. Once enrolled, a Hollister Secure Start Services Coordinator will call to explain our services, answer any questions, and provide ongoing support.

A Hollister Secure Start Services Coordinator can assist with the following:

- Finding the right products
- Providing product information and condition-specific education
- Finding local resources like retailers and clinicians
- Product usage care tips
- Find reimbursement for your catheter products

#### Ostomy Care:

**1.866.789.7574**, or email us at  
**securestartcanada@hollister.com**

#### Continence Care:

**1.800.263.7421**, or email us at  
**canadassccccontact@hollister.com**



# Table of Contents

General Information	2	<b>Accessories</b>	
<b>Ostomy Care</b>		Adapt Barrier Rings	25
Choosing the Right Skin Barrier FIT	4	Adapt Accessories	26
Choosing the Right Skin Barrier FORMULATION	5	m9 Odour Eliminators	27
<b>New Image</b>		Pouch Accessories	27
<b>Two-Piece Pouching System</b>		Karaya Accessories	27
Flat Skin Barriers	6	Skin Sheets	27
Soft Convex Skin Barriers	8	Loop Ostomy Bridge	27
Convex Skin Barriers	8	Irrigation Products	27
Drainable Pouches	10	Wound Drainage Collectors	28
Closed Pouches	11	Access Port	28
High Output Pouches	11	Stoma Caps	28
Urostomy Pouches	12	<b>Continence Care</b>	
Ostomy Kits	13	Infyna Chic	30
Urostomy Kits	13	<b>VaPro</b>	
<b>Premier</b>		Touch Free Hydrophilic Intermittent Catheters	32
<b>One-Piece Pouching System</b>		<b>Advance Plus</b>	
Flat Skin Barriers - Drainable Pouch	16	Touch Free Intermittent Catheters	34
Flat Skin Barriers - Closed Pouch	17	<b>Apogee</b>	
High Output Pouches	17	Apogee Family of Intermittent Catheters	35
Soft Convex Skin Barriers - Drainable Pouch	18	<b>Male External Catheters</b>	
Convex Skin Barriers - Drainable Pouch	18	Hollister Extended and Everyday Wear	36
Flat Skin Barriers - Urostomy Pouch	20	InView Silicone	36
Soft Convex Skin Barriers - Urostomy Pouch	21	<b>Urinary Collection Systems</b>	
Convex Skin Barriers - Urostomy Pouch	21	Leg Bags	37
<b>Karaya</b>		Bedside Drainage	37
<b>One-Piece Pouching System—Natural Skin Barriers</b>		m9 Cleaner/Decrystallizer	37
Drainable Pouches	23	<b>External Collection Devices</b>	
Closed Pouches	23	Urinary Pouches	38
Urostomy Pouches	23	Fecal Collectors	38
<b>Pouchkins</b>		<b>Critical Care</b>	
<b>Pediatric Pouching System</b>		Oral Endotracheal Tube Fasteners	41
One-Piece	24	Tube Attachment Devices	41
Two-Piece	24	Numerical Index	42



# General Information

## Terms of Payment

Net 30 days from date of invoice. Electronic Fund Transfers (ETF) are possible. Please contact the customer order department.

## Minimum Order

The minimum order requirement for a dealer is \$250.00. A charge of \$10.00 will be applied to any order below the minimum requirement.

## Shipping

Prices for shipments to destinations in Canada are FOB buyer's address.

If the method and route of shipment are specified by the buyer, the buyer is required to pay the delivery charges, upon delivery of products to the carrier, all risk of loss and other incidents of ownership pass to the buyer. These terms may not be altered.

## Pricing & Terms of Sale

The prices and specifications in this catalogue are subject to change without notice. GST/HST will be charged where applicable. Buyer agrees that they will not at anytime, directly or indirectly, sell or distribute, to customers which are located outside of Canada.

## Return Goods Policy

All returns for credit require an authorization number from Hollister Limited. A Returns Goods Authorization (RGA) can be obtained by contacting a representative at Customer Service 1.800.263.7400 or by email [canadacustomerservice@hollister.com](mailto:canadacustomerservice@hollister.com). Written request for authorization is required.

**Hollister will authorize returns under the following conditions:**

1. A 15% restock fee will be applied to all returned goods over 30 days from the date of invoice.
2. Copy of invoice or purchase order number must be provided for the return request.
3. **Product returned by customers without an authorized RGA number will be disposed of by Hollister Limited and no credit will be issued.**

**Hollister will NOT authorize the return of the following:**

1. Goods over (1) year from the date of invoice.
2. Partial boxes of merchandise.
3. Boxes in unsaleable condition (i.e markings, sticker for sale, sun fading, dented boxes due to improper handling/return packaging, sticker glue, etc.)

**Damaged Product:**

1. Customer must call Hollister Limited to receive an RGA number within 5 business days of receipt.
2. Receiving dock is responsible for documenting, on the bill of lading, when a shipment is received in damaged condition.
3. Customer must provide copy of invoice, invoice number or purchase order number for the return request.

4. Credit will **NOT** be issued if customer disposes of damaged/ defective product without approval by Hollister Limited.
5. If product cannot be returned due to extreme damage (leakage/broken) customer must provide pictures to Hollister Limited within 5 business days of receipt to gain approval to dispose and thereby obtain credit for the same.

## Ostomy & Continence Retail Sample Policy

All sample requests at the Retail level are to be processed through Hollister Customer Service by calling 1.800.263.7400 or by email [canadacustomerservice@hollister.com](mailto:canadacustomerservice@hollister.com).

Retailer/client is to contact Customer Service and provide consumer information for sample processing. Samples will be shipped directly to consumers from Hollister Limited. Samples will **no longer** be available to ship directly to a retail setting.

## Shortages

Customer must notify Hollister Limited within 5 business days of receipt on shortages. A shortage is a difference between what was physically received and the master pack list.

## Customer Satisfaction Return Policy

All Customer Satisfaction Returns require authorization by your Hollister Territory Manager or a Customer Service Representative.

Unauthorized Customer Satisfaction Returns will not be processed.

Hollister will authorize Customer Satisfaction Returns under the following conditions:

1. The customer has been given a no-charge, full box replacement of Hollister product for their partial box return, and
2. A Customer Satisfaction Return form has been completed, and
3. The returned box has been given to, and the completed Customer Satisfaction Return form has been authorized by your Hollister Territory Manager or a Customer Service Representative.

## Prior to Use

Prior to use, be sure to read the Instructions for Use for information regarding intended use, contraindications, warnings, precautions and instructions.

## Latex Statement

The specifications for the Hollister Continence Care products listed below include natural rubber latex as a component; our suppliers add natural rubber latex as part of the production process. If a Hollister product(s) found in this catalog is not listed below, the specification for that product(s) and packaging does not include natural rubber latex as a component, and our component suppliers do not add natural rubber as part of their production process. Product # 9206, 9207, 9208, 9209, 9106, 9107, 9108, 9109

# Ostomy Care



*Healthy skin. Positive outcomes.*

# Healthy Skin. Positive Outcomes.

## Choosing the Right Skin Barrier FIT

FIT is critical to get a good seal around the stoma to prevent leakage and promote comfort.

### Barrier Examples



Flat Skin Barrier



Flat Skin Barrier  
with Round or Oval  
Convex Barrier Ring



Integrated Soft  
Convex Barrier



Integrated Convex  
Barrier



Integrated Convex  
Barrier with Round  
or Oval Convex  
Barrier Ring

### Possible Clinical Case Scenarios



Flat Skin Barrier

- Flat peristomal plane
- Adequate stoma protrusion
- Centered stoma opening



Flexible Convexity

- Stoma flush with skin
- Uneven peristomal plane
- Off-centered stoma opening



Soft Convexity

- Gentle pressure required to  
peristomal plane
- Flush stoma
- Peristomal injury



Firm Convexity

- Creases or folds in peristomal plane
- Off-centered stoma opening
- Stoma retracted below  
skin surface



Firm Convexity with  
More Depth

- Recessed or flushed stoma
- Uneven abdominal contours
- Off-centered stoma opening

Prior to use, be sure to read the Instructions for Use for information regarding intended use, contraindications, warnings, precautions and instructions.

1. Colwell JC, Goldberg MT & Carmel JE (Eds). (2015). WOCN Society Core Curriculum Ostomy Management. Chapter 10. Philadelphia, PA: Wolters Kluwer
2. Hoeflok J et al. Use of Convexity in Ostomy Care: Results of an International Consensus Meeting. JWOCN, January/February 2017

# Healthy Skin. Positive Outcomes.

## Choosing the Right Skin Barrier FORMULATION

FORMULATION is the customized ingredients responsible for adhesion, erosion resistance and absorption that must be carefully formulated to create an environment where healthy skin can thrive.

### Extended Wear Skin Barriers

---

#### CeraPlus

- Contains skin barrier characteristics (adhesion, absorption and erosion resistance) and is infused with ceramide—a naturally occurring lipid in the skin which protects against dryness
- Formulation designed to help protect the skin's natural moisture barrier and help maintain good peristomal skin health

#### Flexextend

- Durable with resistance to erosion
- Swells when it comes in contact with liquid to help form a seal between the skin and the stoma
- Suitable for those who perspire heavily
- Strong adhesion and flexibility

#### FormaFlex

- Maintains its shape after being stretched or shaped
- Provides an easy-to-customize opening due to barrier material
- Easy to shape to individual stomas without the use of scissors

### Standard Wear Skin Barriers

---

#### FlexWear

- Good balance of durability and gentleness providing moderate adhesion, and can handle frequent or infrequent pouch changes
- Resistant to fluid erosion
- Suitable for all skin types, depending on change frequency
- Designed to be worn more than one day

#### SoftFlex

- Appropriate choice for:
  - Those who have thin, fragile, or sensitive skin (such as those who have undergone chemo/radiation therapy)
  - People who desire frequent removal
- Designed to be gentle to the skin, which allows for frequent pouch changes

#### Karaya 5

- Well-suited for sensitive skin or in some cases where allergic reactions to synthetic barriers has occurred
- Mildly adhesive—needs a belt to secure pouching system



# New Image

## Two-Piece Pouching System

### Flat Skin Barriers

CeraPlus with Remois Technology\*




Cut-to-fit				
Product #	Flange Size	For Stomas Up To	Box of	
<b>11202</b>	1-3/4" (44 mm)	1-1/4" (32 mm)	5	With Tape Border
<b>11203</b>	2-1/4" (57 mm)	1-3/4" (44 mm)		
<b>11204</b>	2-3/4" (70 mm)	2-1/4" (57 mm)		
<b>15102</b>	1-3/4" (44 mm)	1-1/4" (32 mm)	5	Without Tape Border
<b>15103</b>	2-1/4" (57 mm)	1-3/4" (44 mm)		
<b>15104</b>	2-3/4" (70 mm)	2-1/4" (57 mm)		

\*Remois is a technology of Alcare Co., Ltd.

#### Designed for Diverse Bodies

Adhesive (Tape) Border Barriers are an appropriate option for most patients, especially those with:

- Skin irregularities such as hernias, wrinkles, or scars
- Anatomical challenges due to obesity or cachexia
- Added need for a sense of security, or extra assurances to accommodate an active lifestyle
- Increased desire for comfort and flexibility



Flextend



Cut-to-fit				
Product #	Flange Size	For Stomas Up To	Box of	
<b>14602</b>	1-3/4" (44 mm)	1-1/4" (32 mm)	5	With Tape Border
<b>14603</b>	2-1/4" (57 mm)	1-3/4" (44 mm)		
<b>14604</b>	2-3/4" (70 mm)	2-1/4" (57 mm)		
<b>14606</b>	4" (102 mm)	3-1/2" (89 mm)	5	Without Tape Border
<b>15602</b>	1-3/4" (44 mm)	1-1/4" (32 mm)		
<b>15603</b>	2-1/4" (57 mm)	1-3/4" (44 mm)		
<b>15604</b>	2-3/4" (70 mm)	2-1/4" (57 mm)		

#### Pre-sized

Product #	Flange Size	Barrier Opening	Box of	
<b>14701</b>	1-3/4" (44 mm)	5/8" (16 mm)	5	With Tape Border
<b>14702</b>	1-3/4" (44 mm)	3/4" (19 mm)		
<b>14703</b>	1-3/4" (44 mm)	7/8" (22 mm)		
<b>14704</b>	1-3/4" (44 mm)	1" (25 mm)		
<b>14705</b>	1-3/4" (44 mm)	1-1/8" (29 mm)		
<b>14706</b>	1-3/4" (44 mm)	1-1/4" (32 mm)		
<b>14707</b>	2-1/4" (57 mm)	1-3/8" (35 mm)		
<b>14708</b>	2-1/4" (57 mm)	1-1/2" (38 mm)		
<b>14709</b>	2-1/4" (57 mm)	1-3/4" (44 mm)		
<b>16101</b>	1-3/4" (44 mm)	5/8" (16 mm)	5	Without Tape Border
<b>16102</b>	1-3/4" (44 mm)	3/4" (19 mm)		
<b>16103</b>	1-3/4" (44 mm)	7/8" (22 mm)		
<b>16104</b>	1-3/4" (44 mm)	1" (25 mm)		
<b>16105</b>	1-3/4" (44 mm)	1-1/8" (29 mm)		
<b>16106</b>	1-3/4" (44 mm)	1-1/4" (32 mm)		
<b>16107</b>	2-1/4" (57 mm)	1-3/8" (35 mm)		
<b>16108</b>	2-1/4" (57 mm)	1-1/2" (38 mm)		




# New Image

## Two-Piece Pouching System

### Flat Skin Barriers

#### FormaFlex






Product #	Flange Size	For Stomas Up To	Box of
<b>14102</b>	1-3/4" (44 mm)	1-1/4" (32 mm)	5 
<b>14103</b>	2-1/4" (57 mm)	1-11/16" (43 mm)	
<b>14104</b>	2-3/4" (70 mm)	2-1/4" (57 mm)	

#### Colour Matching System

##### Two-Piece Pouching Systems

To help you select compatible *New Image* pouching system components, every box label is colour-coded.



1. Select the skin barrier flange size with the appropriate barrier stoma opening and note its colour code.
2. Select a suitable pouch with the desired features in the matching colour.

	1-3/4" (44 mm)
	2-1/4" (57 mm)
	2-3/4" (70 mm)
	4" (102 mm)



#### FlexWear



##### Cut-to-fit

Product #	Flange Size	For Stomas Up To	Box of
<b>14202</b>	1-3/4" (44 mm)	1-1/4" (32 mm)	5 
<b>14203</b>	2-1/4" (57 mm)	1-3/4" (44 mm)	
<b>14204</b>	2-3/4" (70 mm)	2-1/4" (57 mm)	
<b>14206</b>	4" (102 mm)	3-1/2" (89 mm)	
<b>15202</b>	1-3/4" (44 mm)	1-1/4" (32 mm)	5 
<b>15203</b>	2-1/4" (57 mm)	1-3/4" (44 mm)	
<b>15204</b>	2-3/4" (70 mm)	2-1/4" (57 mm)	

##### Pre-sized

Product #	Flange Size	Barrier Opening	Box of
<b>14301</b>	1-3/4" (44 mm)	5/8" (16 mm)	5 
<b>14302</b>	1-3/4" (44 mm)	3/4" (19 mm)	
<b>14303</b>	1-3/4" (44 mm)	7/8" (22 mm)	
<b>14304</b>	1-3/4" (44 mm)	1" (25 mm)	
<b>14305</b>	1-3/4" (44 mm)	1-1/8" (29 mm)	
<b>14306</b>	1-3/4" (44 mm)	1-1/4" (32 mm)	
<b>14307</b>	2-1/4" (57 mm)	1-3/8" (35 mm)	
<b>14308</b>	2-1/4" (57 mm)	1-1/2" (38 mm)	
<b>14309</b>	2-1/4" (57 mm)	1-3/4" (44 mm)	
<b>16401</b>	1-3/4" (44 mm)	5/8" (16 mm)	5 
<b>16402</b>	1-3/4" (44 mm)	3/4" (19 mm)	
<b>16403</b>	1-3/4" (44 mm)	7/8" (22 mm)	
<b>16404</b>	1-3/4" (44 mm)	1" (25 mm)	
<b>16405</b>	1-3/4" (44 mm)	1-1/8" (29 mm)	
<b>16406</b>	1-3/4" (44 mm)	1-1/4" (32 mm)	
<b>16407</b>	2-1/4" (57 mm)	1-3/8" (35 mm)	
<b>16408</b>	2-1/4" (57 mm)	1-1/2" (38 mm)	

# New Image

## Two-Piece Pouching System

### Soft Convex Skin Barriers **New!**

CeraPlus with Remois Technology\*



#### Cut-to-fit

Product #	Flange Size	For Stomas Up To	Box of
<b>11702</b>	1-3/4" (44 mm)	1" (25 mm)	5 With Tape Border
<b>11703</b>	2-1/4" (57 mm)	1-1/2" (38 mm)	
<b>11704</b>	2-3/4" (70 mm)	2" (51 mm)	

#### Pre-sized

Product #	Flange Size	Barrier Opening	Box of
<b>11903</b>	1-3/4" (44 mm)	7/8" (22 mm)	5 With Tape Border
<b>11904</b>	1-3/4" (44 mm)	1" (25 mm)	
<b>11905</b>	2-1/4" (57 mm)	1-1/8" (29 mm)	
<b>11906</b>	2-1/4" (57 mm)	1-1/4" (32 mm)	

#### Soft Convexity



Flexible design may help provide gentle pressure around the stoma to help obtain the right fit.

### Convex Skin Barriers

CeraPlus with Remois Technology\*



#### Cut-to-fit

Product #	Flange Size	For Stomas Up To	Box of
<b>11402</b>	1-3/4" (44 mm)	1" (25 mm)	5 With Tape Border
<b>11403</b>	2-1/4" (57 mm)	1-1/2" (38 mm)	
<b>11404</b>	2-3/4" (70 mm)	2" (51 mm)	
<b>15302</b>	1-3/4" (44 mm)	1" (25 mm)	5 Without Tape Border
<b>15303</b>	2-1/4" (57 mm)	1-1/2" (38 mm)	
<b>15304</b>	2-3/4" (70 mm)	2" (51 mm)	

#### Pre-sized

Product #	Flange Size	Barrier Opening	Box of
<b>11503</b>	1-3/4" (44 mm)	7/8" (22 mm)	5 With Tape Border
<b>11504</b>	1-3/4" (44 mm)	1" (25 mm)	
<b>11505</b>	2-1/4" (57 mm)	1-1/8" (29 mm)	
<b>11506</b>	2-1/4" (57 mm)	1-1/4" (32 mm)	

\*Remois is a technology of Alcare Co., Ltd.

### The CeraPlus Skin Barrier delivers both Fit and Formulation

#### Ceramide Infused Barrier

- **Protects** the skin's natural moisture barrier
- **Decreases** TEWL (transepidermal water loss) from damaged or eroded skin
- **Maintains** peristomal skin health

# New Image

## Two-Piece Pouching System

### Convex Skin Barriers

#### Flextend



##### Cut-to-fit

Product #	Flange Size	For Stomas Up To	Box of	
<b>14802</b>	1-3/4" (44 mm)	1" (25 mm)	5	With Tape Border
<b>14803</b>	2-1/4" (57 mm)	1-1/2" (38 mm)		
<b>14804</b>	2-3/4" (70 mm)	2" (51 mm)		
<b>15802</b>	1-3/4" (44 mm)	1" (25 mm)	5	Without Tape Border
<b>15803</b>	2-1/4" (57 mm)	1-1/2" (38 mm)		
<b>15804</b>	2-3/4" (70 mm)	2" (51 mm)		

##### Pre-sized

Product #	Flange Size	Barrier Opening	Box of	
<b>13901</b>	1-3/4" (44 mm)	5/8" (16 mm)	5	With Tape Border (Tapered)
<b>13902</b>	1-3/4" (44 mm)	3/4" (19 mm)		
<b>13903</b>	1-3/4" (44 mm)	7/8" (22 mm)		
<b>13904</b>	1-3/4" (44 mm)	1" (25 mm)		
<b>13905</b>	2-1/4" (57 mm)	1-1/8" (29 mm)		
<b>13906</b>	2-1/4" (57 mm)	1-1/4" (32 mm)		
<b>13907</b>	2-1/4" (57 mm)	1-3/8" (35 mm)		
<b>13908</b>	2-1/4" (57 mm)	1-1/2" (38 mm)		
<b>13909</b>	2-3/4" (70 mm)	1-5/8" (41 mm)		
<b>13910</b>	2-3/4" (70 mm)	1-3/4" (44 mm)		
<b>13911</b>	2-3/4" (70 mm)	2" (51 mm)		
<b>14901</b>	1-3/4" (44 mm)	5/8" (16 mm)	5	With Tape Border
<b>14902</b>	1-3/4" (44 mm)	3/4" (19 mm)		
<b>14903</b>	1-3/4" (44 mm)	7/8" (22 mm)		
<b>14904</b>	1-3/4" (44 mm)	1" (25 mm)		
<b>14905</b>	2-1/4" (57 mm)	1-1/8" (29 mm)		
<b>14906</b>	2-1/4" (57 mm)	1-1/4" (32 mm)		
<b>14907</b>	2-1/4" (57 mm)	1-3/8" (35 mm)		
<b>14908</b>	2-1/4" (57 mm)	1-1/2" (38 mm)		
<b>14909</b>	2-3/4" (70 mm)	1-5/8" (41 mm)		
<b>14910</b>	2-3/4" (70 mm)	1-3/4" (44 mm)		
<b>14911</b>	2-3/4" (70 mm)	2" (51 mm)		
<b>15903</b>	1-3/4" (44 mm)	7/8" (22 mm)	5	Without Tape Border
<b>15904</b>	1-3/4" (44 mm)	1" (25 mm)		
<b>15905</b>	2-1/4" (57 mm)	1-1/8" (29 mm)		
<b>15906</b>	2-1/4" (57 mm)	1-1/4" (32 mm)		

#### FlexWear



##### Cut-to-fit

Product #	Flange Size	For Stomas Up To	Box of	
<b>13402</b>	1-3/4" (44 mm)	1" (25 mm)	5	With Tape Border (Tapered)
<b>13403</b>	2-1/4" (57 mm)	1-1/2" (38 mm)		
<b>13404</b>	2-3/4" (70 mm)	2" (51 mm)		
<b>14402</b>	1-3/4" (44 mm)	1" (25 mm)	5	With Tape Border
<b>14403</b>	2-1/4" (57 mm)	1-1/2" (38 mm)		
<b>14404</b>	2-3/4" (70 mm)	2" (51 mm)		
<b>15402</b>	1-3/4" (44 mm)	1" (25 mm)	5	Without Tape Border
<b>15403</b>	2-1/4" (57 mm)	1-1/2" (38 mm)		
<b>15404</b>	2-3/4" (70 mm)	2" (51 mm)		

##### Pre-sized

Product #	Flange Size	Barrier Opening	Box of	
<b>13503</b>	1-3/4" (44 mm)	7/8" (22 mm)	5	With Tape Border (Tapered)
<b>13504</b>	1-3/4" (44 mm)	1" (25 mm)		
<b>13505</b>	2-1/4" (57 mm)	1-1/8" (29 mm)		
<b>13506</b>	2-1/4" (57 mm)	1-1/4" (32 mm)		
<b>14503</b>	1-3/4" (44 mm)	7/8" (22 mm)	5	With Tape Border
<b>14504</b>	1-3/4" (44 mm)	1" (25 mm)		
<b>14505</b>	2-1/4" (57 mm)	1-1/8" (29 mm)		
<b>14506</b>	2-1/4" (57 mm)	1-1/4" (32 mm)		
<b>15503</b>	1-3/4" (44 mm)	7/8" (22 mm)	5	Without Tape Border
<b>15504</b>	1-3/4" (44 mm)	1" (25 mm)		
<b>15505</b>	2-1/4" (57 mm)	1-1/8" (29 mm)		
<b>15506</b>	2-1/4" (57 mm)	1-1/4" (32 mm)		

#### Choosing the Right Skin Barrier FIT

- FIT is critical to get a good seal around the stoma to prevent leakage and promote comfort.
- Convex barrier products can be utilized to help provide a secure and predictable seal.
- Product selection is made by matching stomal protrusion and peristomal contours to the appropriate barriers.
- The outward curving of the convex skin barrier is designed to interface with the immediate peristomal skin.

Refer to page 4 for additional information on The Continuum of Convex Barrier Options.



# New Image

## Two-Piece Pouching System

### Drainable Pouches

#### Lock 'n Roll Closure System



##### Colour: Ultra-clear

Product #	Flange Size	Pouch Length	Filter	Box of
<b>18002</b>	1-3/4" (44 mm)	12" (30 cm)	No	10
<b>18003</b>	2-1/4" (57 mm)	12" (30 cm)		
<b>18004</b>	2-3/4" (70 mm)	12" (30 cm)		
<b>18006</b>	4" (102 mm)	12" (30 cm)		

##### Colour: Transparent

Product #	Flange Size	Pouch Length	Filter	Box of
<b>18132</b>	1-3/4" (44 mm)	12" (30 cm)	No	10
<b>18133</b>	2-1/4" (57 mm)	12" (30 cm)		
<b>18134</b>	2-3/4" (70 mm)	12" (30 cm)		
<b>18192</b>	1-3/4" (44 mm)	12" (30 cm)	AF300	10
<b>18193</b>	2-1/4" (57 mm)	12" (30 cm)		
<b>18194</b>	2-3/4" (70 mm)	12" (30 cm)		
<b>18292</b>	1-3/4" (44 mm)	7" (18 cm)	AF300	20
<b>18293</b>	2-1/4" (57 mm)	7" (18 cm)		
<b>18294</b>	2-3/4" (70 mm)	7" (18 cm)		

##### Colour: Beige

Product #	Flange Size	Pouch Length	Filter	Box of
<b>18112</b>	1-3/4" (44 mm)	12" (30 cm)	No	10
<b>18113</b>	2-1/4" (57 mm)	12" (30 cm)		
<b>18114</b>	2-3/4" (70 mm)	12" (30 cm)		
<b>18182</b>	1-3/4" (44 mm)	12" (30 cm)	AF300	10
<b>18183</b>	2-1/4" (57 mm)	12" (30 cm)		
<b>18184</b>	2-3/4" (70 mm)	12" (30 cm)		
<b>18282</b>	1-3/4" (44 mm)	7" (18 cm)	AF300	20
<b>18283</b>	2-1/4" (57 mm)	7" (18 cm)		
<b>18284</b>	2-3/4" (70 mm)	7" (18 cm)		

#### Lock 'n Roll Microseal Closure Features:

- A secure closure with tiny microseals
- An easy-to-open-and-close seal
- The material used on the closure is water-resistant and can be securely closed even when it is wet, but will work best when dry.

#### Clamp Closure



##### Colour: Ultra-clear

Product #	Flange Size	Pouch Length	Filter	Box of
<b>18172</b>	1-3/4" (44 mm)	12" (30 cm)	No	10
<b>18173</b>	2-1/4" (57 mm)	12" (30 cm)		
<b>18174</b>	2-3/4" (70 mm)	12" (30 cm)		
<b>18176</b>	4" (102 mm)	12" (30 cm)		

##### Colour: Transparent

Product #	Flange Size	Pouch Length	Filter	Box of
<b>18102</b>	1-3/4" (44 mm)	12" (30 cm)	No	10
<b>18103</b>	2-1/4" (57 mm)	12" (30 cm)		
<b>18104</b>	2-3/4" (70 mm)	12" (30 cm)		

##### Colour: Beige

Product #	Flange Size	Pouch Length	Filter	Box of
<b>18122</b>	1-3/4" (44 mm)	12" (30 cm)	No	10
<b>18123</b>	2-1/4" (57 mm)	12" (30 cm)		
<b>18124</b>	2-3/4" (70 mm)	12" (30 cm)		
<b>18142</b>	1-3/4" (44 mm)	12" (30 cm)	AF300	10
<b>18143</b>	2-1/4" (57 mm)	12" (30 cm)		
<b>18144</b>	2-3/4" (70 mm)	12" (30 cm)		
<b>18202</b>	1-3/4" (44 mm)	9" (23 cm)	No	10
<b>18203</b>	2-1/4" (57 mm)	9" (23 cm)		
<b>18204</b>	2-3/4" (70 mm)	9" (23 cm)		

# New Image

## Two-Piece Pouching System

### Closed Pouches



#### Colour: Transparent

Product #	Flange Size	Pouch Length	Filter	Box of
<b>18362</b>	1-3/4" (44 mm)	9" (23 cm)	AF300	60
<b>18363</b>	2-1/4" (57 mm)	9" (23 cm)		
<b>18364</b>	2-3/4" (70 mm)	9" (23 cm)		

#### Colour: Beige QuietWear Fabric

Product #	Flange Size	Pouch Length	Filter	Box of
<b>18332</b>	1-3/4" (44 mm)	9" (23 cm)	No	30
<b>18333</b>	2-1/4" (57 mm)	9" (23 cm)		
<b>18334</b>	2-3/4" (70 mm)	9" (23 cm)		
<b>18322</b>	1-3/4" (44 mm)	9" (23 cm)	AF300	30
<b>18323</b>	2-1/4" (57 mm)	9" (23 cm)		
<b>18324</b>	2-3/4" (70 mm)	9" (23 cm)		

#### Colour: Beige

Product #	Flange Size	Pouch Length	Filter	Box of
<b>18352</b>	1-3/4" (44 mm)	7" (18 cm)	No	30
<b>18353</b>	2-1/4" (57 mm)	7" (18 cm)		
<b>18354</b>	2-3/4" (70 mm)	7" (18 cm)		
<b>18372-30</b>	1-3/4" (44 mm)	9" (23 cm)	AF300	30
<b>18373-30</b>	2-1/4" (57 mm)	9" (23 cm)		
<b>18374-30</b>	2-3/4" (70 mm)	9" (23 cm)		
<b>18392</b>	1-3/4" (44 mm)	7" (18 cm)	AF300	60
<b>18393</b>	2-1/4" (57 mm)	7" (18 cm)		
<b>18394</b>	2-3/4" (70 mm)	7" (18 cm)		

### High Output Pouches



#### Soft Tap Closure System

#### Colour: Ultra-Clear

Product #	Flange Size	Pouch Length	Box of
<b>18012</b>	1-3/4" (44 mm)	12" (30 cm)	10
<b>18013</b>	2-1/4" (57 mm)	12" (30 cm)	
<b>18014</b>	2-3/4" (70 mm)	12" (30 cm)	
<b>18016</b>	4" (102 mm)	12" (30 cm)	

### Colour Matching System

#### Two-Piece Pouching Systems

To help you select compatible *New Image* pouching system components, every box label is colour-coded.

1. Select the skin barrier flange size with the appropriate barrier stoma opening and note its colour code.
2. Select a suitable pouch with the desired features in the matching colour.

	1-3/4" (44 mm)
	2-1/4" (57 mm)
	2-3/4" (70 mm)
	4" (102 mm)

# New Image

## Two-Piece Pouching System

### Urostomy Pouches

#### Enhanced Design



##### Colour: Ultra-Clear

Product #	Flange Size	Pouch Length	Box of
<b>18922</b>	1-3/4" (44 mm)	9" (23 cm)	10
<b>18923</b>	2-1/4" (57 mm)	9" (23 cm)	
<b>18924</b>	2-3/4" (70 mm)	9" (23 cm)	

##### Colour: Transparent

Product #	Flange Size	Pouch Length	Box of
<b>18902</b>	1-3/4" (44 mm)	9" (23 cm)	10
<b>18903</b>	2-1/4" (57 mm)	9" (23 cm)	
<b>18904</b>	2-3/4" (70 mm)	9" (23 cm)	

##### Colour: Beige

Product #	Flange Size	Pouch Length	Box of
<b>18912</b>	1-3/4" (44 mm)	9" (23 cm)	10
<b>18913</b>	2-1/4" (57 mm)	9" (23 cm)	
<b>18914</b>	2-3/4" (70 mm)	9" (23 cm)	

#### Original Design



##### Colour: Ultra-Clear

Product #	Flange Size	Pouch Length	Box of
<b>18422</b>	1-3/4" (44 mm)	9" (23 cm)	10
<b>18423</b>	2-1/4" (57 mm)	9" (23 cm)	

##### Colour: Transparent

Product #	Flange Size	Pouch Length	Box of
<b>18402</b>	1-3/4" (44 mm)	9" (23 cm)	10
<b>18403</b>	2-1/4" (57 mm)	9" (23 cm)	
<b>18404</b>	2-3/4" (70 mm)	9" (23 cm)	

##### Colour: Beige

Product #	Flange Size	Pouch Length	Box of
<b>18412</b>	1-3/4" (44 mm)	9" (23 cm)	10
<b>18413</b>	2-1/4" (57 mm)	9" (23 cm)	
<b>18414</b>	2-3/4" (70 mm)	9" (23 cm)	


### A Great Pouch Becomes Even Better

#### Enhanced Flow-Control Valve

The new valve is made with softer material and a streamlined shape.

#### New Multi-Chamber Design

The updated pouches now have multiple built-in chambers, which are designed to facilitate even distribution of urine.



Multi Chamber  
Front View



Multi Chamber,  
Top View, Sliced



# New Image

## Two-Piece Pouching System

### Ostomy Kits

#### Drainable









#### CeraPlus Barrier

##### Lock 'n Roll Closure System

Non-sterile			
Product #	Flange Size	For Stomas Up To	Box of
<b>19602</b>	1-3/4" (44 mm)	1-1/4" (32 mm)	5
<b>19603</b>	2-1/4" (57 mm)	1-3/4" (44 mm)	
<b>19604</b>	2-3/4" (70 mm)	2-1/4" (57 mm)	

#### FlexWear Barrier

##### Lock 'n Roll Closure System

Non-sterile				
Product #		Flange Size	For Stomas Up To	Box of
<b>19002</b>		1-3/4" (44 mm)	1-1/4" (32 mm)	5
<b>19003</b>		2-1/4" (57 mm)	1-3/4" (44 mm)	
<b>19004</b>		2-3/4" (70 mm)	2-1/4" (57 mm)	
<b>19006</b>		4" (102 mm)	3-1/2" (89 mm)	
Sterile				
Product #		Flange Size	For Stomas Up To	Box of
<b>19052</b>		1-3/4" (44 mm)	1-1/4" (32 mm)	5
<b>19053</b>		2-1/4" (57 mm)	1-3/4" (44 mm)	
<b>19054</b>		2-3/4" (70 mm)	2-1/4" (57 mm)	
<b>19056</b>		4" (102 mm)	3-1/2" (89 mm)	









#### CeraPlus Barrier

##### Clamp Closure System

Non-sterile			
Product #	Flange Size	For Stomas Up To	Box of
<b>19802</b>	1-3/4" (44 mm)	1-1/4" (32 mm)	5
<b>19803</b>	2-1/4" (57 mm)	1-3/4" (44 mm)	
<b>19804</b>	2-3/4" (70 mm)	2-1/4" (57 mm)	

#### FlexWear Barrier

##### Clamp Closure System

Non-sterile				
Product #		Flange Size	For Stomas Up To	Box of
<b>19102</b>		1-3/4" (44 mm)	1-1/4" (32 mm)	5
<b>19103</b>		2-1/4" (57 mm)	1-3/4" (44 mm)	
<b>19104</b>		2-3/4" (70 mm)	2-1/4" (57 mm)	
<b>19106</b>		4" (102 mm)	3-1/2" (89 mm)	
Sterile				
Product #		Flange Size	For Stomas Up To	Box of
<b>19152</b>		1-3/4" (44 mm)	1-1/4" (32 mm)	5
<b>19153</b>		2-1/4" (57 mm)	1-3/4" (44 mm)	
<b>19154</b>		2-3/4" (70 mm)	2-1/4" (57 mm)	
<b>19156</b>		4" (102 mm)	3-1/2" (89 mm)	

Ostomy Kits Include:

- Skin barrier with tape border
- Ultra-clear drainable pouch
- Loop ostomy bridge with 4" flange

### Urostomy Kits







#### CeraPlus Barrier

Non-sterile			
Product #	Flange Size	For Stomas Up To	Box of
<b>19902</b>	1-3/4" (44 mm)	1-1/4" (32 mm)	
<b>19903</b>	2-1/4" (57 mm)	1-3/4" (44 mm)	
<b>19904</b>	2-3/4" (70 mm)	2-1/4" (57 mm)	

Urostomy Kits Include:

- Ultra-clear pouch
- Cut-to-fit Flexextend Skin Barrier
- Drain tube adapter

#### Flexextend Barrier

Non-sterile				
Product #		Flange Size	For Stomas Up To	Box of
<b>19202</b>		1-3/4" (44 mm)	1-1/4" (32 mm)	5
<b>19203</b>		2-1/4" (57 mm)	1-3/4" (44 mm)	
<b>19204</b>		2-3/4" (70 mm)	2-1/4" (57 mm)	
Sterile				
Product #		Flange Size	For Stomas Up To	Box of
<b>19252</b>		1-3/4" (44 mm)	1-1/4" (32 mm)	5
<b>19253</b>		2-1/4" (57 mm)	1-3/4" (44 mm)	
<b>19254</b>		2-3/4" (70 mm)	2-1/4" (57 mm)	



# Follow the evidence to help ensure healthy skin

After surgery: Peristomal skin starts healthy



Top layers of the epidermis are removed with *every barrier change*.


It's unavoidable, and the damage isn't always visible.



Epidermal stripping is linked to increased levels  
of transepidermal water loss (TEWL)<sup>1,2</sup>

## Healthy Skin:

**Keeps moisture in and irritants out**




Low TEWL  
value = higher  
skin integrity<sup>3</sup>

Ceramides are abundant in healthy skin.  
They form a waterproof barrier to help keep  
moisture in and bacteria out.

— = Ceramides

## Damaged Skin:

**Allows moisture out and irritants in**



Higher TEWL  
value = lower  
skin integrity<sup>3</sup>

Damaged skin contains fewer ceramides.

— = Water droplet illustrates moisture loss from skin.

### References:

1. McNichol L, Lund C, Rosen T, Gray M. Medical adhesives and patient safety: state of the science: consensus statements for the assessment, prevention, and treatment of adhesive-related skin injuries. J Wound Ostomy Continence Nurs. 2013;40(4):365-380.
2. Grove G, Houser T, Sibbald G, Salvadalena G. Measuring epidermal effects of ostomy skin barriers. Skin Res Technol. 2019;25:179-186. <https://doi.org/10.1111/srt.12630>
3. Nybaek H, Lophagen S, Karlsmark T, Bang Knudsen D, Jemec GBE. Stratum corneum integrity as a predictor for peristomal skin problems in ostomates. Br J Dermatol. 2010;162:357-361.

# Healthy Skin. Positive Outcomes.

## What Happens When Skin is Compromised

Some common signs of a Peristomal Skin Complication (PSC) include:  
**ITCHING • REDNESS • PAIN • WEEPING SKIN  
BLEEDING • BURNING**



PSCs are often classified as **Peristomal Medical Adhesive-Related Skin Injury (PMARSI)<sup>9</sup>** and **Peristomal Moisture-Associated Skin Damage (PMASD)<sup>1</sup>**



PSCs associated with  
**INCREASED COST  
OF CARE<sup>6</sup>**



PSCs associated with  
**LOWER QUALITY  
OF LIFE<sup>7</sup>**

## What Happens When You Choose CeraPlus Skin Barrier



**RELATIVE REDUCTION**  
in the number of PSCs<sup>5</sup>



**INCREASED SATISFACTION**  
with prevention of itching<sup>5</sup>



Significantly  
**LOWER STOMA-RELATED  
COST OF CARE<sup>5</sup>**



Significantly  
**MORE LIKELY TO  
RESOLVE PSCs<sup>5</sup>**



**TEWL IS  
MINIMIZED<sup>2,8</sup>**




Healthier peristomal skin =  
**HIGHER QUALITY OF LIFE<sup>7</sup>**

**You have a choice. Choose CeraPlus skin barrier and help ensure healthy peristomal skin for your patients from **DAY ONE**.**

### References:

- Colwell J, Pittman J, Raizman R, Salvadalena G. A Randomized Controlled Trial Determining Variances in Ostomy Skin Conditions (ADVOCATE). J Wound Ostomy Continence Nurs. 2018;45(1):37-42.
- Taneja, C., et al., Risk and Economic Burden of Peristomal Skin Complications Following Ostomy Surgery. J Wound Ostomy Continence Nurs. 2019. 46(2): 143-149
- Nichols T., Inglese G. The Burden of Peristomal Skin Complications on an Ostomy Population as Assessed by Health Utility and the Physical Component Summary of the SF-36v2 International Society for Pharmacoeconomics and Outcomes Research (ISPOR), Science Direct Value in Health; 2018; p. 92-93, Tables 3-4
- Nichols T, Houser T, Grove G. Comparing the skin stripping effects of three ostomy skin barriers infused with ceramide, honey or aloe; Journal of Stomal Therapy Australia 2019; 39(2):14-18
- Gray M, Colwel J, Doughty D, Goldberg M, Hoek J, Manson A, McNichol L, Rao S. (2013) Peristomal Moisture-Associated Skin Damage in Adults with Fecal Ostomies: A Comprehensive Review and Consensus. J Wound Ostomy Continence Nurs: 38(5): 541-53





# Premier


## One-Piece Pouching System


### Flat Skin Barriers

#### Drainable Pouch

CeraPlus with Remois Technology\*  
Lock 'n Roll Closure System



Cut-to-fit With Tape Border 					
Product #	Colour	For Stomas Up To	Pouch Length	Filter	Box of
<b>8931</b>	Ultra-clear	2-1/2" (64 mm)	12" (30 cm)	No	10
<b>8901</b>	Beige with Viewing Option	2-1/2" (64 mm)	12" (30 cm)	AF300	


Pre-sized With Tape Border 					
Product #	Colour	Barrier Opening	Pouch Length	Filter	Box of
<b>8925</b>	Beige with Viewing Option	1" (25 mm)	12" (30 cm)	AF300	10
<b>8930</b>		1-3/16" (30 mm)	12" (30 cm)		
<b>8935</b>		1-3/8" (35 mm)	12" (30 cm)		
<b>8940</b>		1-9/16" (40 mm)	12" (30 cm)		


\*Remois is a technology of Alcare Co., Ltd.

### FlexTend

#### Lock 'n Roll Closure System




Cut-to-fit With Tape Border 					
Product #	Colour	For Stomas Up To	Pouch Length	Filter	Box of
<b>8531</b>	Ultra-clear	2-1/2" (64 mm)	12" (30 cm)	No	10
<b>8331</b>	Transparent	2-1/2" (64 mm)	12" (30 cm)	AF300	
<b>8551</b>	Beige	2-1/2" (64 mm)	12" (30 cm)	No	


Pre-sized With Tape Border 				
Product #	Colour	Barrier Opening	Pouch Length	Box of
<b>8538</b>	Transparent	1-1/4" (32 mm)	12" (30 cm)	10
<b>8552</b>	Beige	1" (25 mm)	12" (30 cm)	10
<b>8558</b>		1-1/4" (32 mm)	12" (30 cm)	
<b>8559</b>		1-3/4" (44 mm)	12" (30 cm)	

### FlexTend

#### Clamp Closure System




Cut-to-fit With Tape Border 				
Product #	Colour	For Stomas Up To	Pouch Length	Box of
<b>8651</b>	Ultra-clear	2-1/2" (64 mm)	12" (30 cm)	10
<b>8631</b>	Transparent	2-1/2" (64 mm)	12" (30 cm)	10
<b>8641</b>		2-1/2" (64 mm)	9" (23 cm)	

Pre-sized With Tape Border 				
Product #	Colour	Barrier Opening	Pouch Length	Box of
<b>8637</b>	Transparent	3/4" (19 mm)	12" (30 cm)	10
<b>8632</b>		1" (25 mm)	12" (30 cm)	
<b>8638</b>		1-1/4" (32 mm)	12" (30 cm)	
<b>8633</b>		1-1/2" (38 mm)	12" (30 cm)	
<b>8639</b>		1-3/4" (44 mm)	12" (30 cm)	
<b>8634</b>		2" (51 mm)	12" (30 cm)	
<b>8635</b>	Beige	2-1/2" (64 mm)	12" (30 cm)	10
<b>8647</b>		3/4" (19 mm)	9" (23 cm)	
<b>8642</b>		1" (25 mm)	9" (23 cm)	
<b>8648</b>		1-1/4" (32 mm)	9" (23 cm)	
<b>8643</b>		1-1/2" (38 mm)	9" (23 cm)	
<b>8649</b>		1-3/4" (44 mm)	9" (23 cm)	
<b>8644</b>		2" (51 mm)	9" (23 cm)	

### FlexWear

#### Lock 'n Roll Closure System




Cut-to-fit With Tape Border 					
Product #	Colour	For Stomas Up To	Pouch Length	Filter	Box of
<b>8031</b>	Ultra-clear	2-1/2" (64 mm)	12" (30 cm)	No	10
<b>8081</b>	Transparent	2-1/2" (64 mm)	12" (30 cm)	AF300	
<b>8181</b>	Beige	2-1/2" (64 mm)	12" (30 cm)	AF300	

### FlexWear

#### Clamp Closure System



Cut-to-fit With Tape Border 					
Product #	Colour	For Stomas Up To	Pouch Length	Filter	Box of
<b>8131</b>	Ultra-clear	2-1/2" (64 mm)	12" (30 cm)	No	10

# Premier


## One-Piece Pouching System


### Flat Skin Barriers


#### Drainable Pouch

SoftFlex  
Lock 'n Roll Closure System



Cut-to-fit Without Tape Border 					
Product #	Colour	For Stomas Up To	Pouch Length	Filter	Box of
<b>88300</b>	Beige	2-1/8" (55 mm)	12" (30 cm)	AF300	10
<b>88700</b>		2-1/8" (55 mm)	7" (18 cm)		20

Cut-to-fit—Oval Skin Barrier Without Tape Border 					
Product #	Colour	For Stomas Up To	Pouch Length	Filter	Box of
<b>88402</b>	Transparent	2-1/2" x 3" (66 x 77 mm)	12" (30 cm)	AF300	10
<b>88302</b>	Beige	2-1/2" x 3" (66 x 77 mm)	12" (30 cm)		

Pre-sized Without Tape Border 					
Product #	Colour	Barrier Opening	Pouch Length	Filter	Box of
<b>88325</b>	Beige	1" (25 mm)	12" (30 cm)	AF300	10
<b>88330</b>		1-3/16" (30 mm)	12" (30 cm)		
<b>88335</b>		1-3/8" (35 mm)	12" (30 cm)		
<b>88340</b>		1-9/16" (40 mm)	12" (30 cm)		
<b>88725</b>		1" (25 mm)	7" (18 cm)	AF300	20
<b>88730</b>		1-3/16" (30 mm)	7" (18 cm)		
<b>88735</b>		1-3/8" (35 mm)	7" (18 cm)		
<b>88740</b>		1-9/16" (40 mm)	7" (18 cm)		

#### SoftFlex


- Appropriate choice for:
  - Those who have thin, fragile, or sensitive skin (such as those who have undergone chemo/radiation therapy)
  - People who desire frequent removal
- Designed to be gentle to the skin, which allows for frequent pouch changes


### Flat Skin Barriers


#### Closed Pouch

SoftFlex



Cut-to-fit Without Tape Border 					
Product #	Colour	For Stomas Up To	Pouch Length	Filter	Box of
<b>82300</b>	Beige	2-1/8" (55 mm)	9" (23 cm)	AF300	30
<b>82100</b>		2-1/8" (55 mm)	7" (18 cm)		

Cut-to-fit—Oval Skin Barrier Without Tape Border 					
Product #	Colour	For Stomas Up To	Pouch Length	Filter	Box of
<b>82402</b>	Transparent	2-1/2" x 3" (66 mm x 77 mm)	9" (23 cm)	AF300	30
<b>82302</b>	Beige	2-1/2" x 3" (66 mm x 77 mm)	9" (23 cm)		

Pre-sized Without Tape Border 					
Product #	Colour	Barrier Opening	Pouch Length	Filter	Box of
<b>82325</b>	Beige	1" (25 mm)	9" (23 cm)	AF300	30
<b>82330</b>		1-3/16" (30 mm)	9" (23 cm)		
<b>82335</b>		1-3/8" (35 mm)	9" (23 cm)		
<b>82125</b>		1" (25 mm)	7" (18 cm)		
<b>82130</b>		1-3/16" (30 mm)	7" (18 cm)	AF300	30
<b>82135</b>		1-3/8" (35 mm)	7" (18 cm)		

### High Output Pouches

Flextend



Non-sterile Without Tape Border 			
Product #	Barrier Opening	Pouch Length	Box of
<b>80070</b>	2-3/4" (70 mm)	12" (30 cm)	10
<b>80110</b>	4-1/3" (110 mm)	12" (30 cm)	


# Premier


## One-Piece Pouching System

### Soft Convex Skin Barriers Drainable Pouch

CeraPlus with Remois Technology\*  
Lock 'n Roll Closure System




Cut-to-fit			With Tape Border 		
Product #		For Stomas Up To	Pouch Length	Filter	Box of
8954	Beige with Viewing Option	1" (25 mm)	12" (30 cm)	AF300	5
8958		1-1/2" (38 mm)	12" (30 cm)		
89511		2-1/8" (55 mm)	12" (30 cm)		


Pre-sized			With Tape Border 		
Product #		Barrier Opening	Pouch Length	Filter	Box of
8960	Beige with Viewing Option	3/4" (19 mm)	12" (30 cm)	AF300	5
8961		7/8" (22 mm)	12" (30 cm)		
8962		1" (25 mm)	12" (30 cm)		
8963		1-1/8" (29 mm)	12" (30 cm)		

\*Remois is a technology of Alcare Co., Ltd.

### Flextend Lock 'n Roll Closure System




Cut-to-fit				With Tape Border	
Product #	Colour	Stoma Size	Filter	Box of	
8674	Beige with Viewing Option	5/8" – 1" (13 - 25 mm)	AF300	5	
8678		5/8" – 1-1/2" (13 - 38 mm)			
86711		5/8" – 2-1/8" (13 - 55 mm)			


Pre-sized				With Tape Border	
Product #	Colour	Barrier Opening	Filter	Box of	
8661	Beige with Viewing Option	7/8" (22 mm)	AF300	5	
8662		1" (25 mm)			
8663		1-1/8" (29 mm)			

### Convex Skin Barriers Drainable Pouch

CeraPlus with Remois Technology\*  
Lock 'n Roll Closure System



Cut-to-fit				With Tape Border	
Product #		For Stomas Up To	Pouch Length	Filter	Box of
8984	Ultra-clear	1" (25 mm)	12" (30 cm)	No	5
8988		1-1/2" (38 mm)	12" (30 cm)		
89811		2" (51 mm)	12" (30 cm)		
8914	Beige with Viewing Option	1" (25 mm)	12" (30 cm)	AF300	
8918		1-1/2" (38 mm)	12" (30 cm)		
89111		2" (51 mm)	12" (30 cm)		

Pre-sized				With Tape Border	
Product #		Barrier Opening	Pouch Length	Filter	Box of
8990	Beige with Viewing Option	3/4" (19 mm)	12" (30 cm)	AF300	5
8991		7/8" (22 mm)	12" (30 cm)		
8992		1" (25 mm)	12" (30 cm)		
8993		1-1/8" (29 mm)	12" (30 cm)		

\*Remois is a technology of Alcare Co., Ltd.

#### Soft Convex Skin Barriers

Soft convexity is designed to help provide gentle peristomal pressure for an individualized fit. It helps conform to accommodate uneven skin surfaces, and is a more flexible alternative to firm convexity.

- Lock 'n Roll closure
- Tape border
- Beige with viewing option
- AF300 filter
- Belt tabs



# Premier

## One-Piece Pouching System

### Convex Skin Barriers

#### Drainable Pouch

Flexextend  
Lock 'n Roll Closure System



**Cut-to-fit** With Tape Border

Colour: Transparent

Product #		For Stomas Up To	Pouch Length	Filter	Box of
Tapered Border	Non-Tapered Border				
8244	8524	1" (25 mm)	12" (30 cm)	No	5
8248	8528	1-1/2" (38 mm)	12" (30 cm)		
82411	85211	2" (51 mm)	12" (30 cm)		
8274	8584	1" (25 mm)	12" (30 cm)	AF300	
8278	8588	1-1/2" (38 mm)	12" (30 cm)		
82711	85811	2" (51 mm)	12" (30 cm)		

Colour: Beige

Product #		For Stomas Up To	Pouch Length	Filter	Box of
Tapered Border	Non-Tapered Border				
<b>8214</b>	<b>8504</b>	1" (25 mm)	12" (30 cm)	No	5
<b>8218</b>	<b>8508</b>	1-1/2" (38 mm)	12" (30 cm)		
<b>82111</b>	<b>85011</b>	2" (51 mm)	12" (30 cm)		

**Pre-sized** With Tape Border

Colour: Transparent

Product #		Barrier Opening	Pouch Length	Box of
Tapered Border	Non-Tapered Border			
<b>8280</b>	<b>8510</b>	3/4" (19 mm)	12" (30 cm)	5
<b>8281</b>	<b>8511</b>	7/8" (22 mm)	12" (30 cm)	
<b>8282</b>	<b>8512</b>	1" (25 mm)	12" (30 cm)	
<b>8283</b>	<b>8513</b>	1-1/8" (29 mm)	12" (30 cm)	
<b>8284</b>	<b>8514</b>	1-1/4" (32 mm)	12" (30 cm)	
<b>8285</b>	<b>8515</b>	1-3/8" (35 mm)	12" (30 cm)	
<b>8286</b>	<b>85116</b>	1-1/2" (38 mm)	12" (30 cm)	
<b>8288</b>	<b>85118</b>	1-3/4" (44 mm)	12" (30 cm)	
<b>8289</b>	<b>85119</b>	2" (51 mm)	12" (30 cm)	

Colour: Beige

Product #		Barrier Opening	Pouch Length	Box of
Tapered Border	Non-Tapered Border			
<b>8290</b>	<b>8590</b>	3/4" (19 mm)	12" (30 cm)	5
<b>8291</b>	<b>8591</b>	7/8" (22 mm)	12" (30 cm)	
<b>8292</b>	<b>8592</b>	1" (25 mm)	12" (30 cm)	
<b>8293</b>	<b>8593</b>	1-1/8" (29 mm)	12" (30 cm)	
<b>8294</b>	<b>8594</b>	1-1/4" (32 mm)	12" (30 cm)	
<b>8295</b>	<b>8595</b>	1-3/8" (35 mm)	12" (30 cm)	
<b>8296</b>	<b>8596</b>	1-1/2" (38 mm)	12" (30 cm)	
<b>8297</b>	<b>8597</b>	1-5/8" (41 mm)	12" (30 cm)	

Flexextend  
Clamp Closure System



**Cut-to-fit** With Tape Border

Colour: Transparent

Product #		For Stomas Up To	Pouch Length	Box of
Tapered Border	Non-Tapered Border			
<b>82624</b>	<b>8624</b>	1" (25 mm)	12" (30 cm)	5
<b>82628</b>	<b>8628</b>	1-1/2" (38 mm)	12" (30 cm)	
<b>826211</b>	<b>86211</b>	2" (51 mm)	12" (30 cm)	

**Pre-sized** With Tape Border

Colour: Transparent

Product #		Barrier Opening	Pouch Length	Box of
Tapered Border	Non-Tapered Border			
<b>82610</b>	<b>8610</b>	3/4" (19 mm)	12" (30 cm)	5
<b>82611</b>	<b>8611</b>	7/8" (22 mm)	12" (30 cm)	
<b>82612</b>	<b>8612</b>	1" (25 mm)	12" (30 cm)	
<b>82613</b>	<b>8613</b>	1-1/8" (29 mm)	12" (30 cm)	
<b>82614</b>	<b>8614</b>	1-1/4" (32 mm)	12" (30 cm)	
<b>82615</b>	<b>8615</b>	1-3/8" (35 mm)	12" (30 cm)	
<b>82616</b>	<b>8616</b>	1-1/2" (38 mm)	12" (30 cm)	
<b>82617</b>	<b>8617</b>	1-5/8" (41 mm)	12" (30 cm)	
<b>82618</b>	<b>8618</b>	1-3/4" (44 mm)	12" (30 cm)	
<b>82619</b>	<b>8619</b>	2" (51 mm)	12" (30 cm)	

#### Tapered Border is Designed for Enhanced Sense of Security

- Tapered barrier edges help prevent lift from moisture or motion, which can help achieve desired wear time
- Gentle on the skin while still giving the feeling of extra security
- Adhesive flexes and bends to move comfortably with the body
- A new flexible plastic release liner that covers both the skin barrier and the adhesive border to help improve ease-of-use



# Premier


## One-Piece Pouching System

### Flat Skin Barriers

#### Urostomy Pouch

CeraPlus with Remois Technology\*  
Enhanced Design



Cut-to-fit					With Tape Border 
Product #	Colour	For Stomas	Pouch Length	Filter	Box of
<b>84200</b>	Ultra-clear	2-1/2" (64 mm)	9" (23 cm)	No	10

\*Remois is a technology of Alcare Co., Ltd.


The CeraPlus Skin Barrier  
delivers both Fit and Formulation


#### Ceramide Infused Barrier

- **Protects** the skin's natural moisture barrier
- **Decreases** TEWL (transepidermal water loss) from damaged or eroded skin
- **Maintains** peristomal skin health

Flextend  
Original Design





Cut-to-fit					With Tape Border 
Product #	Colour	For Stomas Up To	Pouch Length	Box of	
<b>8450</b>	Ultra-clear	2-1/2" (64 mm)	9" (23 cm)	10	
<b>8460</b>	Transparent	2-1/2" (64 mm)	9" (23 cm)		
<b>8440</b>	Beige	2-1/2" (64 mm)	9" (23 cm)		

Pre-sized					With Tape Border 
Product #	Colour	Barrier Opening	Pouch Length	Box of	
<b>8467</b>	Transparent	3/4" (19 mm)	9" (23 cm)	10	
<b>8462</b>		1" (25 mm)	9" (23 cm)		
<b>8468</b>		1-1/4" (32 mm)	9" (23 cm)		
<b>8463</b>		1-1/2" (38 mm)	9" (23 cm)		
<b>8469</b>		1-3/4" (44 mm)	9" (23 cm)		
<b>8464</b>		2" (51 mm)	9" (23 cm)	10	
<b>8447</b>	Beige	3/4" (19 mm)	9" (23 cm)		
<b>8442</b>		1" (25 mm)	9" (23 cm)		
<b>8448</b>		1-1/4" (32 mm)	9" (23 cm)		
<b>8443</b>		1-1/2" (38 mm)	9" (23 cm)		
<b>8449</b>		1-3/4" (44 mm)	9" (23 cm)		

Flextend  
Enhanced Design



Cut-to-fit					With Tape Border 
Product #	Colour	For Stomas Up To	Pouch Length	Box of	
<b>84590</b>	Ultra-clear	2-1/2" (64 mm)	9" (23 cm)	10	
<b>84490</b>	Beige	2-1/2" (64 mm)	9" (23 cm)		

Pre-sized					With Tape Border 
Product #	Colour	Barrier Opening	Pouch Length	Box of	
<b>84697</b>	Transparent	3/4" (19 mm)	9" (23 cm)	10	
<b>84692</b>		1" (25 mm)	9" (23 cm)		
<b>84698</b>		1-1/4" (32 mm)	9" (23 cm)		
<b>84693</b>		1-1/2" (38 mm)	9" (23 cm)		
<b>84699</b>		1-3/4" (44 mm)	9" (23 cm)		
<b>84694</b>		2" (51 mm)	9" (23 cm)	10	
<b>84497</b>	Beige	3/4" (19 mm)	9" (23 cm)		
<b>84492</b>		1" (25 mm)	9" (23 cm)		
<b>84498</b>		1-1/4" (32 mm)	9" (23 cm)		
<b>84493</b>		1-1/2" (38 mm)	9" (23 cm)		
<b>84499</b>		1-3/4" (44 mm)	9" (23 cm)		
<b>84494</b>		2" (51 mm)	9" (23 cm)		




# Premier


## One-Piece Pouching System

### Soft Convex Skin Barriers Urostomy Pouch

CeraPlus with Remois Technology\*  
Enhanced Design



Cut-to-fit				With Tape Border 
Product #	Colour	For Stomas Up To	Pouch Length	Box of
<b>84134</b>	Ultra-clear	1" (25 mm)	9" (23 cm)	5
<b>84138</b>		1-1/2" (38 mm)	9" (23 cm)	
<b>841311</b>		2-1/8" (55 mm)	9" (23 cm)	


Pre-sized				With Tape Border 
Product #	Colour	Barrier Opening	Pouch Length	Box of
<b>8412</b>	Beige with Viewing Option	3/4" (19 mm)	9" (23 cm)	5
<b>8413</b>		7/8" (22 mm)	9" (23 cm)	
<b>8414</b>		1" (25 mm)	9" (23 cm)	
<b>8415</b>		1-1/8" (29 mm)	9" (23 cm)	

\*Remois is a technology of Alcare Co., Ltd.

### Convex Skin Barriers Urostomy Pouch

CeraPlus with Remois Technology\*  
Enhanced Design



Cut-to-fit				With Tape Border	
Product #	Colour	For Stomas Up To	Pouch Length	Filter	Box of
<b>84404</b>	Ultra-clear	1" (25 mm)	9" (23 cm)	No	5
<b>84408</b>		1-1/2" (38 mm)	9" (23 cm)		
<b>844011</b>		2" (51 mm)	9" (23 cm)		

Pre-sized				
Product #	Colour	Barrier Opening	Pouch Length	Box of
<b>84292</b>	Beige with Viewing Option	3/4" (19 mm)	9" (23 cm)	5
<b>84293</b>		7/8" (22 mm)	9" (23 cm)	
<b>84294</b>		1" (25 mm)	9" (23 cm)	
<b>84295</b>		1-1/8" (29 mm)	9" (23 cm)	

\*Remois is a technology of Alcare Co., Ltd.

#### Premier Soft Convex CeraPlus Pouching Systems:

- Are designed to conform to the body without compromising fit
- Deliver both convex support and flexibility
- Protect the skin's natural moisture barrier and help maintain good peristomal skin health



#### A Great Pouch Becomes Even Better

##### Enhanced Flow-Control Valve

The new valve is made with softer material and a streamlined shape.

##### New Multi-Chamber Design

The updated pouches now have multiple built-in chambers, which are designed to facilitate even distribution of urine.



Multi Chamber  
Front View



Multi Chamber,  
Top View, Sliced

# Premier


## One-Piece Pouching System


### Convex Skin Barriers

#### Urostomy Pouch

Flextend  
Enhanced Design





Cut-to-fit			With Tape Border 	
Product #	Colour	For Stomas Up To	Pouch Length	Box of
<b>84794</b>	Ultra-clear	1" (25 mm)	9" (23 cm)	5
<b>84798</b>		1-1/2" (38 mm)	9" (23 cm)	
<b>8479111</b>		2" (51 mm)	9" (23 cm)	
<b>84394</b>	Beige	1" (25 mm)	9" (23 cm)	
<b>84398</b>		1-1/2" (38 mm)	9" (23 cm)	
<b>8439111</b>		2" (51 mm)	9" (23 cm)	

Pre-sized			With Tape Border 	
Product #	Colour	Barrier Opening	Pouch Length	Box of
<b>84890</b>	Ultra-clear	1/2" (13 mm)	9" (23 cm)	5
<b>84891</b>		5/8" (16 mm)	9" (23 cm)	
<b>84892</b>		3/4" (19 mm)	9" (23 cm)	
<b>84893</b>		7/8" (22 mm)	9" (23 cm)	
<b>84894</b>		1" (25 mm)	9" (23 cm)	
<b>84895</b>		1-1/8" (29 mm)	9" (23 cm)	
<b>84896</b>		1-1/4" (32 mm)	9" (23 cm)	
<b>84897</b>		1-3/8" (35 mm)	9" (23 cm)	
<b>84898</b>		1-1/2" (38 mm)	9" (23 cm)	
<b>84899</b>		1-3/4" (44 mm)	9" (23 cm)	
<b>84990</b>	Beige	1/2" (13 mm)	9" (23 cm)	5
<b>84991</b>		5/8" (16 mm)	9" (23 cm)	
<b>84992</b>		3/4" (19 mm)	9" (23 cm)	
<b>84993</b>		7/8" (22 mm)	9" (23 cm)	
<b>84994</b>		1" (25 mm)	9" (23 cm)	
<b>84995</b>		1-1/8" (29 mm)	9" (23 cm)	
<b>84996</b>		1-1/4" (32 mm)	9" (23 cm)	
<b>84997</b>		1-3/8" (35 mm)	9" (23 cm)	
<b>84998</b>		1-1/2" (38 mm)	9" (23 cm)	
<b>84999</b>		1-3/4" (44 mm)	9" (23 cm)	

Flextend  
Original Design



Cut-to-fit			With Tape Border 	
Product #	Colour	For Stomas Up To	Pouch Length	Box of
<b>8474</b>	Transparent	1" (25 mm)	9" (23 cm)	5
<b>8478</b>		1-1/2" (38 mm)	9" (23 cm)	
<b>84711</b>		2" (51 mm)	9" (23 cm)	
Pre-sized				
With Tape Border 				
Product #	Colour	Barrier Opening	Pouch Length	Box of
<b>8480</b>	Transparent	1/2" (13 mm)	9" (23 cm)	5
<b>8481</b>		5/8" (16 mm)	9" (23 cm)	
<b>8482</b>		3/4" (19 mm)	9" (23 cm)	
<b>8483</b>		7/8" (22 mm)	9" (23 cm)	
<b>8484</b>		1" (25 mm)	9" (23 cm)	
<b>8485</b>		1-1/8" (29 mm)	9" (23 cm)	
<b>8486</b>		1-1/4" (32 mm)	9" (23 cm)	
<b>8487</b>		1-3/8" (35 mm)	9" (23 cm)	
<b>8488</b>		1-1/2" (38 mm)	9" (23 cm)	
<b>8489</b>		1-3/4" (44 mm)	9" (23 cm)	
<b>8490</b>	Beige	1/2" (13 mm)	9" (23 cm)	5
<b>8491</b>		5/8" (16 mm)	9" (23 cm)	
<b>8492</b>		3/4" (19 mm)	9" (23 cm)	
<b>8493</b>		7/8" (22 mm)	9" (23 cm)	
<b>8494</b>		1" (25 mm)	9" (23 cm)	
<b>8495</b>		1-1/8" (29 mm)	9" (23 cm)	
<b>8496</b>		1-1/4" (32 mm)	9" (23 cm)	
<b>8497</b>		1-3/8" (35 mm)	9" (23 cm)	
<b>8498</b>		1-1/2" (38 mm)	9" (23 cm)	


### Adhesive (Tape) Border Barriers


- Feature a thin layer of adhesive that surrounds the hydrocolloid barrier.
- Designed to be gentle on the skin, the adhesive moves with the body to easily conform to skin irregularities, creases and folds.
- Application is simple, with no additional accessories required to gain an added sense of protection. Go with the skin barrier that comes with built-in confidence.


### Drainable Pouches

#### Clamp Closure System




Colour: Transparent With Tape Border 			
Product #	Stoma Size	Pouch Length	Box of
<b>3232</b>	5/8" (15 mm)	9" (23 cm)	30
<b>3238</b>	7/8" (22 mm)	9" (23 cm)	
<b>3233</b>	1-1/8" (29mm)	9" (23 cm)	
<b>3239</b>	1-3/8" (35 mm)	9" (23 cm)	
<b>3234</b>	1-1/2" (38 mm)	9" (23 cm)	
<b>3222</b>	5/8" (15 mm)	12" (30 cm)	30
<b>3228</b>	7/8" (22 mm)	12" (30 cm)	
<b>3223</b>	1-1/8" (29mm)	12" (30 cm)	
<b>3229</b>	1-3/8" (35 mm)	12" (30 cm)	
<b>3224</b>	1-1/2" (38 mm)	12" (30 cm)	
<b>3225</b>	2" (51 mm)	12" (30 cm)	
<b>3226</b>	2-1/2" (64 mm)	12" (30 cm)	
<b>3278</b>	7/8" (22 mm)	16" (40 cm)	30
<b>3273</b>	1-1/8" (29mm)	16" (40 cm)	
<b>3279</b>	1-3/8" (35 mm)	16" (40 cm)	
<b>3274</b>	1-1/2" (38 mm)	16" (40 cm)	
<b>3275</b>	2" (51 mm)	16" (40 cm)	
<b>3276</b>	2-1/2" (64 mm)	16" (40 cm)	
<b>3608</b>	7/8" (22 mm)	12" (30 cm)	10
<b>3603</b>	1-1/8" (29mm)	12" (30 cm)	
<b>3609</b>	1-3/8" (35 mm)	12" (30 cm)	
<b>3604</b>	1-1/2" (38 mm)	12" (30 cm)	
<b>3605</b>	2" (51 mm)	12" (30 cm)	

Colour: Beige With Tape Border 			
Product #	Stoma Size	Pouch Length	Box of
<b>3112</b>	5/8" (15 mm)	12" (30 cm)	30
<b>3119</b>	1-3/8" (35 mm)	12" (30 cm)	


Colour: Transparent With Replacement Filter With Tape Border 			
Product #	Stoma Size	Pouch Length	Box of
<b>3662</b>	5/8" (15 mm)	12" (30 cm)	10
<b>3668</b>	7/8" (22 mm)	12" (30 cm)	

### Closed Pouches




Colour: Transparent Without Tape Border 			
Product #	Stoma Size	Pouch Length	Box of
<b>7164</b>	1-1/2" (38 mm)	9" (23 cm)	30

Colour: Transparent With Integrated Filter With Tape Border 			
Product #	Stoma Size	Pouch Length	Box of
<b>3328</b>	7/8" (22 mm)	9" (23 cm)	30
<b>3323</b>	1-1/8" (29mm)	9" (23 cm)	
<b>3329</b>	1-3/8" (35 mm)	9" (23 cm)	
<b>3324</b>	1-1/2" (38 mm)	9" (23 cm)	
<b>3325</b>	2" (51 mm)	9" (23 cm)	
<b>3326</b>	2-1/2" (64 mm)	9" (23 cm)	

### Urostomy Pouches



Colour: Transparent With Tape Border 			
Product #	Stoma Size	Pouch Length	Box of
<b>1432</b>	5/8" (15 mm)	9" (23 cm)	10
<b>1438</b>	7/8" (22 mm)	9" (23 cm)	
<b>1433</b>	1-1/8" (29mm)	9" (23 cm)	
<b>1439</b>	1-3/8" (35 mm)	9" (23 cm)	
<b>1434</b>	1-1/2" (38 mm)	9" (23 cm)	


- To select the appropriate *Karaya* pouch, determine the stoma size and select the matching Product Code
- The labeling on the *Karaya* box and pouch you have selected indicates the gasket size, NOT the stoma size
- The *Karaya* seal extends into the gasket opening, ensuring that the Product Code you have selected fits the stoma size you have measured


# Pouchkins

## Pediatric Pouching System

### One-Piece SoftFlex

#### Premie/Newborn Pouches



Cut-to-fit			Without Tape Border	
Product #	Description	For Stomas Up To	Box of	
Premie				
3777	Ultra-clear, no starter hole	5/8" (15 mm)	15	
Newborn				
3778	Ultra-clear, no starter hole	1-3/8" x 7/8" (35 x 22 mm)		

#### Drainable Pouches Lock 'n Roll Closure System



Cut-to-fit Without Tape Border				
Product #	Description	For Stomas Up To	Pouch Length	Box of
<b>3795</b>	Transparent, no starter hole	2" (51 mm)	7" (18 cm)	10
<b>3796</b>		1-1/2" (38 mm)	6-1/2" (16 cm)	

#### Urostomy Pouch



Cut-to-fit Without Tape Border				
Product #	Description	For Stomas Up To	Pouch Length	Box of
<b>3797</b>	Transparent, no starter hole	1-1/2" (38 mm)	6-1/2" (16 cm)	10

#### **Pouchkins SoftFlex Skin Barrier:**

- Designed to be gentle on sensitive infant skin
- No starter hole, which enables the opening to be placed where needed
- No tape by design
- Flexible and sized for babies and children
- None of the *Pouchkins* pediatric pouching systems are made with natural rubber latex
- Any 1-3/4" *New Image* pouch will fit this skin barrier

### Two-Piece SoftFlex

#### Flat Skin Barrier



Cut-to-fit Without Tape Border				
Product #	Description	Flange Size	For Stomas Up To	Box of
<b>3761</b>	No starter hole	1-3/4" (44 mm)	1-1/4" (32 mm)	5

#### Drainable Pouches Lock 'n Roll Closure System



Colour: Transparent Without Tape Border				
Product #	Flange Size	Pouch Length	Filter	Box of
<b>3799</b>	1-3/4" (44 mm)	6-1/2" (16 cm)	No	10
<b>18292</b>	1-3/4" (44 mm)	7" (18 cm)	AF300	20

#### Urostomy Pouches



Without Tape Border				
Product #	Flange Size	Pouch Length	Filter	Box of
<b>18402</b>	1-3/4" (44 mm)	9" (23 cm)	No	10

#### Enhanced Design

Product #	Colour	Flange Size	Pouch Length	Box of
<b>18922</b>	Ultra-clear	1-3/4" (44 mm)	9" (23 cm)	10
<b>18912</b>	Beige	1-3/4" (44 mm)	9" (23 cm)	10

#### Ostomy Belt



Without Tape Border		
Product #	Length	Box of
<b>3774</b>	10 - 17" (25 - 43 cm)	1

# Accessories

## Adapt Barrier Rings

### Flat CeraRing



Product #	Type	Outer Diameter	Thickness	Box of
<b>8805</b>	Standard	2" (48 mm)	4.5 mm	10
<b>8815</b>	Slim	2" (48 mm)	2.30 mm	

### Convex CeraRing



Product #	Inner Diameter	Stretch Up To	Box of
<b>89520</b>	13/16" (20 mm)	1" (25 mm)	10
<b>89530</b>	1-3/16" (30 mm)	1-3/8" (35 mm)	
<b>89540</b>	1-9/16" (40 mm)	1-3/4" (45 mm)	

### Oval Convex CeraRing



Product #	Inner Diameter	Stretch Up To	Box of
<b>89601</b>	7/8" x 1-1/2" (22 x 38 mm)	1-1/8" x 1-3/4" (27 x 43 mm)	10
<b>89602</b>	1-3/16" x 1-7/8" (30 x 48 mm)	1-3/8" x 2-1/8" (35 x 53 mm)	
<b>89603</b>	1-1/2" x 2-3/16" (38 x 56 mm)	1-3/4" x 2-3/8" (43 x 61 mm)	

### Flat Flextend



Product #	Type	Outer Diameter	Thickness	Box of
<b>7805</b>	Standard	2" (48 mm)	4.5 mm	10
<b>7806</b>	Standard	4" (102 mm)	4.5 mm	
<b>7815</b>	Slim	2" (48 mm)	2.30 mm	

### Convex Flextend



Product #	Inner Diameter	Stretch Up To	Box of
<b>79520</b>	13/16" (20 mm)	1" (25 mm)	10
<b>79530</b>	1-3/16" (30 mm)	1-3/8" (35 mm)	
<b>79540</b>	1-9/16" (40 mm)	1-3/4" (45 mm)	

### Oval Convex Flextend



Product #	Inner Diameter	Stretch Up To	Box of
<b>79601</b>	7/8" x 1-1/2" (22 x 38 mm)	1-1/8" x 1-3/4" (27 x 43 mm)	10
<b>79602</b>	1-3/16" x 1-7/8" (30 x 48 mm)	1-3/8" x 2-1/8" (35 x 53 mm)	
<b>79603</b>	1-1/2" x 2-3/16" (38 x 56 mm)	1-3/4" x 2-3/8" (43 x 61 mm)	

## Adapt CeraRing

**Sealed with the confidence of ceramide**

**The security + confidence of Adapt meets the science of ceramide**











- Available in Flat and Convex for each individual need
- Helps prevent leakage and protect skin—and may even extend barrier wear time
- Resists erosion—even when wet or exposed to discharge from the stoma
- Offers ultimate comfort and ease of use—alcohol-free rings are designed to be gentle for all skin types
- Extended wear
- Infused with ceramide — a natural component of the skin
- Help protect the skin's natural moisture barrier
- Help decrease transepidermal water loss (TEWL) from damaged or eroded skin







# Accessories

## Adapt Accessories





Product #		Description	Usage	Size		Box of
<b>7760</b>		Universal Remover Wipes <i>*featuring no-sting silicone-based, alcohol free formulation</i>	Silicone-based and alcohol-free, this skin friendly product helps remove adhesive residues.	50 Wipes per Carton		50
<b>7737</b>		Adapt Medical Adhesive Remover Spray <i>*featuring no-sting, silicone-based, alcohol free formulation</i>	360° spray technology to deliver controlled, quiet spraying from any angle, with minimal waste	Spray Can	1.7 fl oz (50 ml)	1
<b>7906</b>		Adapt Stoma Powder	Use to absorb moisture from broken skin around the stoma. Dust on. Brush off excess powder. Stop using when skin heals.	Puff bottle	1 oz (28 g)	1
<b>79402</b>		Adapt Barrier Extenders	Helps provide extra security around the edges of skin barrier, reduces skin barrier edge lifting and may help increase skin barrier wear time. Provides extra adhesion to accommodate a peristomal hernia.	Each barrier extender covers a surface area of 48 cm <sup>2</sup>		30
<b>78500</b>		Adapt Lubricating Deodorant	Neutralizes odours, eases emptying, and helps prevent pouch static and sticking. It is not a perfume or masking agent. Should be used in drainable and closed pouching systems only.	Bottle	8 oz (236 ml)	1
<b>78501</b>				Packets	.27 oz (8 ml)	50
<b>79301</b>		Adapt Paste	Used to fill or caulk uneven skin contours to create a flatter surface. It can help to prevent drainage from getting under the ostomy skin barrier.	Tube	.5 oz (14 g)	20
<b>79300</b>					2 oz (60 g)	1
<b>79400</b>		Adapt Barrier Strips with FlexTend Barrier (Alcohol Free)	Used to fill in or caulk uneven skin contours near the stoma to create a flatter surface and help prevent ostomy drainage from getting under the ostomy barrier.	Triangular Shape	5-1/8" (13 cm)	10
<b>7730</b>		Medical Adhesive Spray	Improves adhesive contact between skin and barrier Not recommended for use with extended wear barriers	Spray Can <i>Contains no CFCs</i>	3.8 oz (112 ml)	1
<b>7917</b>		Adapt Skin Protective Wipe <i>*featuring no-sting silicone-based, alcohol free formulation</i>	Protects skin during adhesive removal Not recommended for use with extended wear skin barriers	50 Wipes per Carton		50
<b>HBRF2650</b>		Hydrofera Blue Antibacterial Ostomy Ring	Antibacterial foam dressing that helps manage both wound exudate and bioburden. To be applied to skin under ostomy skin barrier.	Outer Diameter 2-1/2" (6.4 cm)		10

# Accessories

## m9 Odour Eliminators

Product #	Description	Size	Box of
<b>7715</b>	 Drops	1 oz (30 ml)	12
<b>7717</b>		8 oz (240 ml)	6
<b>7732</b>	 Spray Unscented	2 oz (60 ml)	12
<b>7733</b>		8 oz (240 ml)	6
<b>7735</b>	Spray Apple Scent	8 oz (240 ml)	

## Pouch Accessories

Product #	Description	Size	Box of
<b>8770</b>	 Drainable Pouch Clamp, Beige		20
<b>7300</b>	 Beige Colour Ostomy Belt	Medium 23" - 43" (58 - 109 cm)	10
<b>7299</b>		Large 34" - 65" (86 - 165 cm)	
<b>7766</b>	 Replacement Filter	20 filters and 2 plugs/bag (product series 366x, 383x, and 389x)	5 bags
<b>7331</b>	 Enhanced Design Urostomy Drain Tube Adapter		10

## Karaya Accessories

Product #	Description	Size	Box of
<b>7910</b>	 Paste	Tube 4.5 oz (128 g)	12
<b>7905</b>	 Powder	Puff Bottle 2.5 oz (71 g)	

## Skin Sheets



Product #	Description	Size	Box of
<b>8800</b>	Flexextend (Non-sterile)	4" x 4" (10 cm x 10 cm)	5
<b>8801</b>		8" x 8" (20 cm x 20 cm)	3
<b>7800</b>	Premium (Non-sterile)	4" x 4" (10 cm x 10 cm)	5
<b>7801</b>		8" x 8" (20 cm x 20 cm)	3
<b>7700</b>	HolliHesive (Non-sterile)	4" x 4" (10 cm x 10 cm)	5

Skin barrier sheets are designed to be the most resistant to discharge and is the barrier of choice for many urostomies and ileostomies. It may be cut for customized use in stoma care. Available in extended wear (Flexextend) and standard wear (Premium & HolliHesive).

## Loop Ostomy Bridge



Product #	Description	Length	Box of
<b>7767</b>	Sterile	3-1/2" (89 mm)	20

## Irrigation Products



Product #	Description	Information	Box of
<b>7721</b>	Cone Irrigator Kit	Includes: (1) Irrigator (1) Stoma Cone with Connector (1) Cleaning Brush	1
<b>7723</b>	Stoma Cone	With Connector Replacement Unit	10
<b>7728</b>	Irrigator Sleeve	With Attached Closure and Belt Tabs	2" (51 mm)
<b>7724</b>			3" (76 mm)
<b>18152</b>	New Image Two-Piece Irrigator Sleeve	Flange: 1-3/4" (44 mm)	Pouch: 30" (76 cm)
<b>18153</b>		Flange: 2-1/4" (57 mm)	Pouch: 30" (76 cm)
<b>18154</b>		Flange: 2-3/4" (70 mm)	Pouch: 30" (76 cm)
<b>7740</b>	Stoma Lubricant	Bottle	4 oz (118 ml)

# Accessories

## Wound Drainage Collectors

### With Skin Barrier



Product #	For Wounds Up To		Box of
<b>9701</b>	3" (7.6 cm)	Sterile	5
<b>9702</b>	3.75" (9.5 cm)		5
<b>9703</b>	4" x 8" (10 cm x 20 cm)		3
<b>9773</b>	3" (7.6 cm)	Non-sterile	10
<b>9776</b>	3.75" (9.5 cm)		10
<b>9778</b>	4" x 8" (10 cm x 20 cm)		3

### With Microporous Tape

Product #	For Wounds Up To		Box of
<b>9772</b>	3" (7.6 cm)	Sterile	10
<b>9775</b>	3.75" (9.5 cm)		10
<b>9777</b>	4" x 8" (10 cm x 20 cm)		3

### Wound Drainage Collectors

For management of draining wounds, fistulas, or draining tube sites

- Available with and without FlexWear skin barriers
- Sterile and non-sterile drainable pouches
- Odour-barrier film
- Resealable drainage cap and large access window
- Without starter hole

## Access Port



Product #	Description		Box of
<b>9779</b>	Graduated nipple and punch allows introduction of catheter or drain through Wound Drainage Collector or pouch system	Sterile	8

## Stoma Caps



Pre-sized					With Tape Border
Product #	Colour	For Stomas Up To	Pouch Length	Filter	Box of
<b>3184</b>	White	2" (51 mm)	4-1/4" (11 cm)	Yes	30
<b>3186</b>		3" (76 mm)	4-1/4" (11 cm)		

## SoftFlex Skin Barrier



Pre-sized					Without Tape Border
Product #	Colour	For Stomas Up To	Pouch Length	Filter	Box of
<b>1796</b>	Beige	1-5/16" (50 mm)	4" (10 cm)	Yes	30

# Continence Care

A man with a beard and glasses, wearing a blue denim shirt over a white t-shirt, is sitting in a wheelchair. He is smiling and looking off to the side. The background is a blurred outdoor setting with trees and a path. The bottom of the image has a red overlay.

People are at the center of everything we do.



## Discretion by design

Introducing a discreet catheter with design that is truly inspired.

The Infyna Chic catheter was designed with input from clinicians to help provide a high level of discretion for women who use catheters. It is beautiful and easy to use. The Infyna Chic catheter can help a woman feel better about having to use one.

### Infyna Chic **New!**

Product #	Description	Size	Length	Box of
<b>7008</b>	Straight tip	8 Fr	5" (13 cm)	30
<b>7010</b>	Straight tip	10 Fr	5" (13 cm)	30
<b>7012</b>	Straight tip	12 Fr	5" (13 cm)	30
<b>7014</b>	Straight tip	14 Fr	5" (13 cm)	30



## Features that help change her intermittent catheter experience


- Flip-top cap can be opened and closed with one hand
- Ready to use, hydrophilic catheter
- Just right stiffness helps support easy, touch-free insertion
- Catheter length helps give confidence that the bladder is fully drained
- Can attach a collection bag without removing the catheter from the case



Easy to open with just one hand



Easy to grip and remove from case



Does not leak when reclosed





# VaPro

Touch Free Intermittent Catheter

F-Style

## Another reason to choose VaPro Catheters

The VaPro F-Style catheter, a firmer-style catheter, offers patients the freedom to choose more firmness, increased control, and the extra protection that the portfolio of VaPro intermittent catheters delivers.

### VaPro F-Style Catheters

More firm.  
More control.  
More choice.



Continence Care



## Touch Free Hydrophilic Intermittent Catheters

### VaPro



Product #	Description	Size	Length	Box of
<b>72124</b>	Straight tip	12 Fr	16" (40 cm)	30
<b>73124</b>	Coudé tip		16" (40 cm)	
<b>72142</b>	Straight tip	14 Fr	8" (20 cm)	30
<b>72144</b>	Straight tip		16" (40 cm)	
<b>73144</b>	Coudé tip	16 Fr	16" (40 cm)	30
<b>72164</b>	Straight tip		16" (40 cm)	
<b>73164</b>	Coudé tip		16" (40 cm)	

### VaPro Plus



Product #	Description	Size	Length	Box of
<b>74124</b>	Straight tip	12 Fr	16" (40 cm)	30
<b>74144</b>	Straight tip	14 Fr	16" (40 cm)	30

*VaPro Plus is a closed system with an integrated 1000 ml collection bag.*

### VaPro F-Style

**New!**



Product #	Description	Size	Length	Box of
<b>7600124</b>	Straight tip	12 Fr	16" (40 cm)	30
<b>7600144</b>	Straight tip	14 Fr	16" (40 cm)	30
<b>7600164</b>	Straight tip	16 Fr	16" (40 cm)	30

### VaPro Plus F-Style



Product #	Description	Size	Length	Box of
<b>7700124</b>	Straight tip	12 Fr	16" (40 cm)	30
<b>7700144</b>	Straight tip	14 Fr	16" (40 cm)	30

*VaPro Plus F-Style is a closed system with an integrated 1000 ml collection bag.*

### VaPro Pocket



Product #	Description	Size	Length	Box of
<b>70124</b>	Straight tip	12 Fr	16" (40 cm)	30
<b>70144</b>	Straight tip	14 Fr	16" (40 cm)	30
<b>70164</b>	Straight tip	16 Fr	16" (40 cm)	30

### VaPro Plus Pocket



Product #	Description	Size	Length	Box of
<b>71084</b>	Straight tip	8 Fr	16" (40 cm)	30
<b>71104</b>	Straight tip	10 Fr	16" (40 cm)	30
<b>71122</b>	Straight tip	12 Fr	8" (20 cm)	30
<b>71124</b>	Straight tip		16" (40 cm)	
<b>71142</b>	Straight tip	14 Fr	8" (20 cm)	30
<b>71144</b>	Straight tip		16" (40 cm)	
<b>71164</b>	Straight tip	16 Fr	16" (40 cm)	30

*VaPro Plus Pocket is a closed system with an integrated 1000 ml collection bag.*

### VaPro Catheters Provide 100% No Touch Protection with a Protective Tip and Sleeve:

VaPro catheters are hygienic by design to help protect patients from germs throughout the entire catheterization process.

#### Protective tip:

- helps shield the sterile catheter during insertion from bacteria located within the first 15 mm of the distal urethra
- helps reduce the risk of carrying bacteria into the urinary tract

#### Protective sleeve:

- allows for catheter to be gripped anywhere
- provides a barrier to the catheter to help keep germs away

#### The protective tip is designed for male and female use

The protective tip helps reduce the risk of bacteria being carried into the urinary tract during the insertion process.





## The Hygienic Benefits of a No-Touch Intermittent Catheter with a Protective Tip and/or Sleeve

The 'No-Touch' Method of Intermittent Urinary Catheter Insertion:  
Can it Reduce the Risk of Bacteria Entering the Bladder?

The in vitro model showed that intermittent catheters with the hygienic feature of a no-touch sleeve helped reduce the potential for external contamination during preparation and insertion of an intermittent catheter. This reduction may help reduce the risk of bacteria entering the bladder.<sup>1</sup>

For more information on this study or to request your free sample, please visit us at [www.hollister.com](http://www.hollister.com) or call us at 1.800.263.7400.

1. Hudson E, and Murahata R. Spinal Cord. 2005; 43; 611-614.

Prior to use, be sure to read the Instructions for Use for information regarding Intended Use, Contraindications, Warnings, Precautions, and Instructions. The Hollister logo, VaPro, and "Continence Care. People First." are trademarks of Hollister Incorporated. ©2018 Hollister Incorporated.



## Touch Free Intermittent Catheters

### Closed System



Product #	Description	Size	Length	Box of
<b>94064IC</b>	Straight tip	6 Fr	16" (40 cm)	100
<b>94084IC</b>	Straight tip	8 Fr	16" (40 cm)	
<b>94104IC</b>	Straight tip	10 Fr	16" (40 cm)	
<b>94124IC</b>	Straight tip	12 Fr	16" (40 cm)	
<b>95124IC</b>	Coudé tip		16" (40 cm)	
<b>94144IC</b>	Straight tip	14 Fr	16" (40 cm)	
<b>95144IC</b>	Coudé tip		16" (40 cm)	
<b>94164IC</b>	Straight tip	16 Fr	16" (40 cm)	
<b>95164IC</b>	Coudé tip		16" (40 cm)	
<b>94184IC</b>	Straight tip	18 Fr	16" (40 cm)	

### Closed System Kit



Product #	Description	Size	Length	Box of
<b>96064IC</b>	Straight tip	6 Fr	16" (40 cm)	100
<b>96084IC</b>	Straight tip	8 Fr	16" (40 cm)	
<b>96104IC</b>	Straight tip	10 Fr	16" (40 cm)	
<b>96124IC</b>	Straight tip	12 Fr	16" (40 cm)	
<b>97124</b>	Coudé tip		16" (40 cm)	
<b>96144IC</b>	Straight tip	14 Fr	16" (40 cm)	
<b>97144</b>	Coudé tip		16" (40 cm)	
<b>96164IC</b>	Straight tip	16 Fr	16" (40 cm)	
<b>97164</b>	Coudé tip		16" (40 cm)	
<b>96184IC</b>	Straight tip	18 Fr	16" (40 cm)	

#### Kit Includes:

- Advance Plus Closed System
- Gloves
- Benzalkonium Chloride antiseptic swab (BZK)
- Waterproof underpad

#### Advance Plus Touch Free Intermittent Catheters:

- Have a protective tip that helps shield the catheter from bacteria located in the first 15 mm of the distal urethra
- Have a user-regulated gel reservoir which allows for control of the catheter lubrication
- Smooth eyelets help provide trouble-free insertion and withdrawal



## Apogee Family of Intermittent Catheters



Product #	Description	Size	Length	Box of
<b>B1007</b>	Straight tip	6 Fr	15" (38 cm)	30
<b>B1016</b>	Straight tip	8 Fr	6" (15 cm)	30
<b>B1018</b>	Straight tip		16" (40 cm)	
<b>B1033</b>	Straight tip	10 Fr	6" (15 cm)	30
<b>1035</b>	Coudé, Olive tip		16" (40 cm)	
<b>B1034</b>	Straight tip		16" (40 cm)	
<b>B1047</b>	Straight tip	12 Fr	6" (15 cm)	30
<b>B1048</b>	Straight tip		16" (40 cm)	
<b>B1050</b>	Coudé tip		16" (40 cm)	
<b>B1062</b>	Straight tip	14 Fr	6" (15 cm)	30
<b>B1063</b>	Straight tip, No funnel		6" (15 cm)	
<b>B1064</b>	Straight tip		16" (40 cm)	
<b>B1061</b>	Soft tip		16" (40 cm)	
<b>B1066</b>	Coudé tip		16" (40 cm)	
<b>B1065</b>	Straight tip with curved packaging	16 Fr	16" (40 cm)	50
<b>B1082</b>	Straight tip		16" (40 cm)	
<b>B1086</b>	Coudé tip		16" (40 cm)	
<b>B1096</b>	Straight tip	18 Fr	16" (40 cm)	30

### Features of the Apogee Family of Intermittent Catheters:

- Two staggered smooth eyelets help provide trouble-free insertion
- Not made with natural rubber latex
- Designed for comfort
- Medical grade Polyvinyl Chloride (PVC)
- Non lubricated, clear catheter
- Straight and coudé tip options available



# Male External Catheters

## Hollister

### Latex



<b>Extended Wear</b> with Unique Inner Flap			Sheath length 2-3/4" (7 cm), adhesive length 1-3/5" (4 cm)
Product #	Description	Size	Box of
<b>9206</b>	Small/Pediatric	7/8"- 1" (22 - 25 mm)	30
<b>9207</b>	Medium	1"- 1-1/8" (26 - 30 mm)	
<b>9208</b>	Intermediate	1-1/8"- 1-1/4" (31 - 35 mm)	
<b>9209</b>	Large	1-3/8"- 1-1/2" (36 - 39 mm)	
<b>Everyday Wear</b>			Sheath length 2-3/4" (7 cm), adhesive length 1-3/5" (4 cm)
Product #	Description	Size	Box of
<b>9106</b>	Small/Pediatric	7/8"- 1" (22 - 25 mm)	30
<b>9107</b>	Medium	1"- 1-1/8" (26 - 30 mm)	
<b>9108</b>	Intermediate	1-1/8"- 1-1/4" (31 - 35 mm)	
<b>9109</b>	Large	1-3/8"- 1-1/2" (36 - 39 mm)	

Contains Natural Rubber Latex

### Size Selection Guide

Use the Hollister Sizing Guide to determine the most appropriate size catheter

- Place sizing guide just behind the glans of the penis
- Select the size that most closely fits



## InView

### Silicone



<b>Standard</b>			Sheath length 2-3/4" (7 cm), adhesive length 1-3/5" (4 cm)
Product #	Description	Size	Box of
<b>97525</b>	Small	1" (25 mm)	30
<b>97529</b>	Medium	1-1/8" (29 mm)	
<b>97532</b>	Intermediate	1-1/4" (32 mm)	
<b>97536</b>	Large	1-1/2" (36 mm)	
<b>97541</b>	Extra Large	1-5/8" (41 mm)	
<b>97525-100</b>	Small	1" (25 mm)	100
<b>97529-100</b>	Medium	1-1/8" (29 mm)	
<b>97532-100</b>	Intermediate	1-1/4" (32 mm)	
<b>97536-100</b>	Large	1-1/2" (36 mm)	
<b>97541-100</b>	Extra Large	1-5/8" (41 mm)	

<b>Special</b>			Sheath length 1-3/5" (4 cm), adhesive length 1-3/5" (4 cm)
Product #	Description	Size	Box of
<b>97425</b>	Small	1" (25 mm)	30
<b>97429</b>	Medium	1-1/8" (29 mm)	
<b>97432</b>	Intermediate	1-1/4" (32 mm)	
<b>97436</b>	Large	1-1/2" (36 mm)	
<b>97441</b>	Extra Large	1-5/8" (41 mm)	
<b>97425-100</b>	Small	1" (25 mm)	100
<b>97429-100</b>	Medium	1-1/8" (29 mm)	
<b>97432-100</b>	Intermediate	1-1/4" (32 mm)	
<b>97436-100</b>	Large	1-1/2" (36 mm)	
<b>97441-100</b>	Extra Large	1-5/8" (41 mm)	

<b>Extra</b>			Sheath length 2-3/4" (7 cm), adhesive length 2-3/4" (7 cm)
Product #	Description	Size	Box of
<b>97625</b>	Small	1" (25 mm)	30
<b>97629</b>	Medium	1-1/8" (29 mm)	
<b>97632</b>	Intermediate	1-1/4" (32 mm)	
<b>97636</b>	Large	1-1/2" (36 mm)	
<b>97641</b>	Extra Large	1-5/8" (41 mm)	
<b>97625-100</b>	Small	1" (25 mm)	100
<b>97629-100</b>	Medium	1-1/8" (29 mm)	
<b>97632-100</b>	Intermediate	1-1/4" (32 mm)	
<b>97636-100</b>	Large	1-1/2" (36 mm)	
<b>97641-100</b>	Extra Large	1-5/8" (41 mm)	

Not Made with Natural Rubber Latex

# Urinary Collection Systems

## Leg Bags

### Leg Bag



Product #	Description	Capacity	Box of
<b>9814</b>	Medium	18 oz (540 ml)	10
<b>9805</b>	Large	30 oz (900 ml)	

### Extension Tubing with Connector



Product #	Description	Length	Box of
<b>9345</b>	Non-sterile	18" (45 cm)	10
<b>9346</b>	Sterile	18" (45 cm)	

### Fabric Strap



Product #	Description	Capacity	Box of
<b>9342</b>	Small	15" (38 cm)	10 pairs
<b>9343</b>	Medium	23" (58 cm)	

#### Hollister Leg Bags Are:

- Sterile
- Not Made with Natural Rubber Latex
- Designed with an anti-reflux valve

### Combination Pack



Product #	Description	Capacity	Box of
<b>9348</b>	Medium	18 oz (540 ml)	12 packs
<b>9349</b>	Large	30 oz (900 ml)	

### Vented Combination Pack



Product #	Description	Capacity	Box of
<b>9645</b>	Medium	18 oz (540 ml)	10 packs
<b>9655</b>	Large	30 oz (900 ml)	

### Bedside Drainage



Product #	Description	Capacity	Box of
<b>9839</b>	2000 ml Collection System	68 oz	20

## m9 Cleaner/Decrystallizer



### Concentrate & Wash Bottle

Product #	Description	Size	Box of
<b>7736</b>	1 Bottle of Concentrate 1 Washing Bottle	16 oz (480 ml)	12

# External Collection Devices

## Urinary Pouches

### Female

#### SoftFlex Skin Barrier



Product #	Information	Box of
<b>9840</b>	Adapt Paste, tape strips and Skin Gel Wipe included	10

### Male

#### Flextend Skin Barrier



Product #	Information	Box of
<b>9873</b>	Thin barrier for flexibility, no additional tape needed	10

#### SoftFlex Skin Barrier



Product #	Information	Box of
<b>9811</b>	Tape strips and Skin Gel Wipe included	10

## Fecal Collectors

### Pouches

#### Flextend Skin Barrier



Product #	Description	Pouch Length	Box of
<b>9823</b>	Medium, Ultra-clear	10" (25 cm)	10
<b>9880</b>	Large, Transparent	12" (30 cm)	

### Transparent Pouches

#### SoftFlex Skin Barrier



Product #	Description	Pouch Length	Box of
<b>9822</b>	Medium	10" (25 cm)	10
<b>9821</b>	Large	12" (30 cm)	

#### Hollister's Urinary Pouches:

- Not Made with Natural Rubber Latex
- Have a cut-to-fit skin barrier
- Have a odour-barrier pouch film



With 100% No Touch Protection,  
whatever you touch  
won't touch VaPro catheters.



Continence Care



# Critical Care

A woman with dark hair, wearing a pink sweater, is smiling and looking down at a young girl with brown hair in pigtails. The girl is wearing a white t-shirt and blue overalls, and is focused on drawing with a black marker on a piece of paper. They are sitting at a wooden table. Several other colored markers are scattered on the table. The background is a bright, out-of-focus room with a white brick wall and a bookshelf.

Dedicated to excellence.



## Oral Endotracheal Tube Fasteners

### AnchorFast SlimFit

**New!**



Product #	Description	Box of
<b>9787</b>	<ul style="list-style-type: none"> <li>• Secure without use of adhesive tape</li> <li>• Designed to accommodate smaller faces</li> <li>• Gliding tube shuttle</li> <li>• Upper lip foam designed to increase surface area contact</li> <li>• Adjustable neck strap</li> </ul>	12

### AnchorFast Guard



Product #	Description	Box of
<b>9700</b>	<ul style="list-style-type: none"> <li>• Integrated Tube Protection Sleeve</li> <li>• Allows side-to-side repositioning of the tube to help optimize oral care while the tube remains secure</li> <li>• For tubes sized 5.0-9.0mm inner diameter</li> </ul>	12

### AnchorFast



Product #	Description	Box of
<b>9799</b>	<ul style="list-style-type: none"> <li>• Flextend Skin Barrier</li> <li>• Allows side-to-side repositioning of the tube to help optimize oral care while the tube remains secure</li> </ul>	12

### AnchorFast Replacement Neck Band Straps



Product #	Description	Box of
<b>9788</b>	Replacement Neck Band Straps	4

## Tube Attachment Devices

### Horizontal Tube (HTAD)



Product #	Description	Box of
<b>9781</b>	Single use, sterile	5

### Vertical Tube (VTAD)



Product #	Description	Box of
<b>9782</b>	Single use, sterile	5

### Drain Tube (DTAD)



Product #	Description	Box of
<b>9780</b>	Sterile	5

### Feeding Tube (FTAD)



Product #	Description	Box of
<b>9786</b>	FlexWear skin barrier conforms to patient's nose	20

### Suction Tube (STAD)



Product #	Description	Box of
<b>9785</b>	FlexWear skin barrier conforms to patient's nose	20

# Numerical Index

Product #	Page	Product #	Page	Product #	Page	Product #	Page	Product #	Page	Product #	Page
1035	35	3778	24	8244	19	8483	22	8613	19	8960	18
1432	23	3795	24	8248	19	8484	22	8614	19	8961	18
1433	23	3796	24	8274	19	8485	22	8615	19	8962	18
1434	23	3797	24	8278	19	8486	22	8616	19	8963	18
1438	23	3799	24	8280	19	8487	22	8617	19	8984	18
1439	23	7008	30	8281	19	8488	22	8618	19	8988	18
1796	28	7010	30	8282	19	8489	22	8619	19	8990	18
3112	23	7012	30	8283	19	8490	22	8624	19	8991	18
3119	23	7014	30	8284	19	8491	22	8628	19	8992	18
3184	28	7164	23	8285	19	8492	22	8631	16	8993	18
3186	28	7299	27	8286	19	8493	22	8632	16	9106	36
3222	23	7300	27	8288	19	8494	22	8633	16	9107	36
3223	23	7331	27	8289	19	8495	22	8634	16	9108	36
3224	23	7700	27	8290	19	8496	22	8635	16	9109	36
3225	23	7715	27	8291	19	8497	22	8637	16	9206	36
3226	23	7717	27	8292	19	8498	22	8638	16	9207	36
3228	23	7721	27	8293	19	8504	19	8639	16	9208	36
3229	23	7723	27	8294	19	8508	19	8641	16	9209	36
3232	23	7724	27	8295	19	8510	19	8642	16	9342	37
3233	23	7728	27	8296	19	8511	19	8643	16	9343	37
3234	23	7730	26	8297	19	8512	19	8644	16	9345	37
3238	23	7732	27	8331	16	8513	19	8647	16	9346	37
3239	23	7733	27	8412	21	8514	19	8648	16	9348	37
3273	23	7735	27	8413	21	8515	19	8649	16	9349	37
3274	23	7736	37	8414	21	8524	19	8651	16	9645	37
3275	23	7737	26	8415	21	8528	19	8661	18	9655	37
3276	23	7740	27	8440	20	8531	16	8662	18	9700	41
3278	23	7760	26	8442	20	8538	16	8663	18	9701	28
3279	23	7766	27	8443	20	8551	16	8674	18	9702	28
3323	23	7767	27	8447	20	8552	16	8678	18	9703	28
3324	23	7800	27	8448	20	8558	16	8770	27	9772	28
3325	23	7801	27	8449	20	8559	16	8800	27	9773	28
3326	23	7805	25	8450	20	8584	19	8801	27	9775	28
3328	23	7806	25	8460	20	8588	19	8805	25	9776	28
3329	23	7815	25	8462	20	8590	19	8815	25	9777	28
3603	23	7905	27	8463	20	8591	19	8901	16	9778	28
3604	23	7906	26	8464	20	8592	19	8914	18	9779	28
3605	23	7910	27	8467	20	8593	19	8918	18	9780	41
3608	23	7917	26	8468	20	8594	19	8925	16	9781	41
3609	23	8031	16	8469	20	8595	19	8930	16	9782	41
3662	23	8081	16	8474	22	8596	19	8931	16	9785	41
3668	23	8131	16	8478	22	8597	19	8935	16	9786	41
3761	24	8181	16	8480	22	8610	19	8940	16	9787	41
3774	24	8214	19	8481	22	8611	19	8954	18	9788	41
3777	24	8218	19	8482	22	8612	19	8958	18	9799	41

Product #	Page	Product #	Page	Product #	Page	Product #	Page	Product #	Page	Product #	Page
9805	37	14102	7	14907	9	16407	7	18292	24	19052	13
9811	38	14103	7	14908	9	16408	7	18293	10	19053	13
9814	37	14104	7	14909	9	18002	10	18294	10	19054	13
9821	38	14202	7	14910	9	18003	10	18322	11	19056	13
9822	38	14203	7	14911	9	18004	10	18323	11	19102	13
9823	38	14204	7	15102	6	18006	10	18324	11	19103	13
9839	37	14206	7	15103	6	18012	11	18332	11	19104	13
9840	38	14301	7	15104	6	18013	11	18333	11	19106	13
9873	38	14302	7	15202	7	18014	11	18334	11	19152	13
9880	38	14303	7	15203	7	18016	11	18352	11	19153	13
11202	6	14304	7	15204	7	18102	10	18353	11	19154	13
11203	6	14305	7	15302	8	18103	10	18354	11	19156	13
11204	6	14306	7	15303	8	18104	10	18362	11	19202	13
11402	8	14307	7	15304	8	18112	10	18363	11	19203	13
11403	8	14308	7	15402	9	18113	10	18364	11	19204	13
11404	8	14309	7	15403	9	18114	10	18372-30	11	19252	13
11503	8	14402	9	15404	9	18122	10	18373-30	11	19253	13
11504	8	14403	9	15503	9	18123	10	18374-30	11	19254	13
11505	8	14404	9	15504	9	18124	10	18392	11	19602	13
11506	8	14503	9	15505	9	18132	10	18393	11	19603	13
11702	8	14504	9	15506	9	18133	10	18394	11	19604	13
11703	8	14505	9	15602	6	18134	10	18402	12	19802	13
11704	8	14506	9	15603	6	18142	10	18402	24	19803	13
11903	8	14602	6	15604	6	18143	10	18403	12	19804	13
11904	8	14603	6	15802	9	18144	10	18404	12	19902	13
11905	8	14604	6	15803	9	18152	27	18412	12	19903	13
11906	8	14606	6	15804	9	18153	27	18413	12	19904	13
13402	9	14701	6	15903	9	18154	27	18414	12	70124	32
13403	9	14702	6	15904	9	18172	10	18422	12	70144	32
13404	9	14703	6	15905	9	18173	10	18423	12	70164	32
13503	9	14704	6	15906	9	18174	10	18902	12	71084	32
13504	9	14705	6	16101	6	18176	10	18903	12	71104	32
13505	9	14706	6	16102	6	18182	10	18904	12	71122	32
13506	9	14707	6	16103	6	18183	10	18912	12	71124	32
13901	9	14708	6	16104	6	18184	10	18912	24	71142	32
13902	9	14709	6	16105	6	18192	10	18913	12	71144	32
13903	9	14802	9	16106	6	18193	10	18914	12	71164	32
13904	9	14803	9	16107	6	18194	10	18922	12	72124	32
13905	9	14804	9	16108	6	18202	10	18922	24	72142	32
13906	9	14901	9	16401	7	18203	10	18923	12	72144	32
13907	9	14902	9	16402	7	18204	10	18924	12	72164	32
13908	9	14903	9	16403	7	18282	10	19002	13	73124	32
13909	9	14904	9	16404	7	18283	10	19003	13	73144	32
13910	9	14905	9	16405	7	18284	10	19004	13	73164	32
13911	9	14906	9	16406	7	18292	10	19006	13	74124	32

# Numerical Index

Product #	Page	Product #	Page	Product #	Page	Product #	Page	Product #	Page
74144	32	84294	21	85118	19	97144	34	B1033	35
78500	26	84295	21	85119	19	97164	34	B1034	35
78501	26	84394	22	85211	19	97425	36	B1047	35
79300	26	84398	22	85811	19	97425-100	36	B1048	35
79301	26	84404	21	86211	19	97429	36	B1050	35
79400	26	84408	21	86711	18	97429-100	36	B1061	35
79402	26	84490	20	88300	17	97432	36	B1062	35
79520	25	84492	20	88302	17	97432-100	36	B1063	35
79530	25	84493	20	88325	17	97436	36	B1064	35
79540	25	84494	20	88330	17	97436-100	36	B1065	35
79601	25	84497	20	88335	17	97441	36	B1066	35
79602	25	84498	20	88340	17	97441-100	36	B1082	35
79603	25	84499	20	88402	17	97525	36	B1086	35
80070	17	84590	20	88700	17	97525-100	36	B1096	35
80110	17	84692	20	88725	17	97529	36	HBRF2650	26
82100	17	84693	20	88730	17	97529-100	36		
82111	19	84694	20	88735	17	97532	36		
82125	17	84697	20	88740	17	97532-100	36		
82130	17	84698	20	89111	18	97536	36		
82135	17	84699	20	89511	18	97536-100	36		
82300	17	84711	22	89520	25	97541	36		
82302	17	84794	22	89530	25	97541-100	36		
82325	17	84798	22	89540	25	97625	36		
82330	17	84890	22	89601	25	97625-100	36		
82335	17	84891	22	89602	25	97629	36		
82402	17	84892	22	89603	25	97629-100	36		
82411	19	84893	22	89811	18	97632	36		
82610	19	84894	22	94064IC	34	97632-100	36		
82611	19	84895	22	94084IC	34	97636	36		
82612	19	84896	22	94104IC	34	97636-100	36		
82613	19	84897	22	94124IC	34	97641	36		
82614	19	84898	22	94144IC	34	97641-100	36		
82615	19	84899	22	94164IC	34	826211	19		
82616	19	84990	22	94184IC	34	841311	21		
82617	19	84991	22	95124IC	34	844011	21		
82618	19	84992	22	95144IC	34	7600124	32		
82619	19	84993	22	95164IC	34	7600144	32		
82624	19	84994	22	96064IC	34	7600164	32		
82628	19	84995	22	96084IC	34	7700124	32		
82711	19	84996	22	96104IC	34	7700144	32		
84134	21	84997	22	96124IC	34	8439111	22		
84138	21	84998	22	96144IC	34	8479111	22		
84200	20	84999	22	96164IC	34	B1007	35		
84292	21	85011	19	96184IC	34	B1016	35		
84293	21	85116	19	97124	34	B1018	35		

A man in a dark suit, white shirt, patterned tie, and a fedora-style hat is sitting on a wooden chair outdoors. He is looking off to the side with a thoughtful expression, his hand resting on his chin. The background shows a rustic wooden building and other wooden chairs, suggesting a historical or heritage setting.

The spirit of Hollister  
remains the same today  
as it did in 1921.

#### Our Mission

To make life more rewarding and dignified for the people who use our products and services.

#### Our Vision

To grow and prosper as an independent, employee-owned company and in the process become better human beings.

## Our Heritage

The roots of Hollister Incorporated date back to 1921, when 23-year-old entrepreneur John Dickinson Schneider founded a small printing company in Chicago, which he named JDS Printer Craftsman. From the very beginning, John Schneider's business philosophy was "only first class is good enough," and the company became known for its quality products, quality service, and quality people.



**CeraPlus™ skin barrier**  
with Remois Technology\*

**vs. SenSura Mio**

## The Results are Clear: CeraPlus Skin Barrier Helps Ensure Healthy Skin

In the Measuring Epidermal Effects of Ostomy Skin Barriers study (Grove et al, 2019), the CeraPlus skin barrier was significantly more effective at maintaining skin health than the SenSura Mio barrier.



**CeraPlus  
Skin Barrier**

**SenSura  
Mio Barrier**

Infused with ceramide, the skin's natural protection against water loss	✓	✗
Observed to be better at minimizing skin irritation	✓	✗
Observed to be better at managing Transepidermal water loss (TEWL)	✓	✗
Observed to be better at protecting skin against damage from repeated barrier removal	✓	✗
Observed to be better at maintaining skin health	✓	✗



**CeraPlus Skin Barrier Minimizes Potential Impact  
of Peristomal Skin Complications Right from the Start!**

*For more information, please contact your Hollister Representative*

Grove G, Houser T, Sibbald G, Salvadlena G. Measuring epidermal effects of ostomy skin barriers. *Skin Res Technol*. 2019;25:179–186. <https://doi.org/10.1111/srt.12630>  
The Hollister logo, CeraPlus and "Healthy skin. Positive Outcomes." are trademarks of Hollister Incorporated.  
All other trademarks and copyrights are the property of their respective owners.  
©2019 Hollister Incorporated

CE Prior to use, be sure to read the Instructions for Use for information regarding Intended Use, Contraindications, Warnings, Precautions, and Instructions.  
Not all products are CE marked.



\*Remois is a technology of Alcare Co., Ltd.



**Ostomy Care**  
*Healthy skin. Positive outcomes.*