



# VaPro

Touch Free Intermittent Catheter

F-Style

## Another reason to choose VaPro Catheters

The VaPro F-Style catheter, a firmer-style catheter, offers patients the freedom to choose more firmness, increased control, and the extra protection that the portfolio of VaPro intermittent catheters delivers.

### VaPro F-Style Catheters

More firm.  
More control.  
More choice.



Continence Care

# VaPro F-Style Catheters and VaPro Plus F-Style Catheters

## Another choice in our portfolio



### Firmness patients want

- Designed to increase control during insertion\*
- Designed to help patients with dexterity issues\*
- Designed for patients who prefer a firmer catheter\*

\*Responses from 27 nurses from UK, US, and NL in Hollister Clinical Advisory Board meetings in October 2016.

### Extra protection patients need

- Protective introducer tip helps shield the sterile catheter from bacteria in the initial portion of the urethra
- Protective introducer tip helps reduce the risk of carrying bacteria into the urinary tract
- Protective sleeve allows for catheter to be gripped anywhere
- Protective sleeve provides a barrier to the catheter to help keep germs away

### Easy to use

The unique combination of tip and sleeve creates a truly touch-free system that is easy to use.

- Ready to use with no extra steps required
- Evenly lubricated for smooth insertion and easier withdrawal
- Easy-open packaging is designed even for those with limited dexterity

#### Ring cap

helps keep protective tip clean and protected after package is opened

#### Protective tip

helps bypass bacteria in the first 15 mm of the distal urethra

#### Protective sleeve

promotes easy handling while also helping to keep bacteria away

#### F-Style

firmer-style catheter tubing



### Available in a range of sizes

#### VaPro F-Style

Ready to use hydrophilic intermittent catheter

Product No	Length (cm/inches)	Size	Box Qty
7600124	40 cm/16 in	12	30
7600144	40 cm/16 in	14	30
7600164	40 cm/16 in	16	30

#### VaPro Plus F-Style

Ready to use hydrophilic intermittent catheter with collection bag

Product No	Length (cm/inches)	Size	Box Qty
7700124	40 cm/16 in	12	30
7700144	40 cm/16 in	14	30

This intermittent catheter is a flexible tubular device that is inserted through the urethra by male, female, and pediatric (paediatric) patients who need to drain urine from the bladder.

**WARNING: To help reduce the potential for infection and/or other complications, do not reuse.**

**WARNING:** If discomfort or any sign of trauma occurs, discontinue use immediately and consult your healthcare practitioner.

**PRECAUTION: Please consult your healthcare practitioner before using this product if any of the following conditions are present:**

- Severed urethra
- Pronounced stricture
- Urethritis – inflammation of the urethra
- Epididymitis – inflammation of the epididymis (testicle tube)
- Unexplained urethral bleeding
- False passage
- Prostatitis – inflammation of the prostate gland

**PRECAUTION:** Self-catheterization (catheterisation) should follow the plan of care and advice given by your healthcare practitioner and be carried out only in accordance with the instructions for use provided. Because catheterization (catheterisation) frequency varies by person, the recommended frequency of your catheterization (catheterisation) should be provided by your healthcare practitioner. For any other questions about your catheterization (catheterisation), please contact your healthcare practitioner.

**NOTE:** Store boxes in a flat position and at normal room temperature.

This product consists of: polyurethane sleeve, thermoplastic elastomer tip, and tubing not made with DEHP.