



## Maintain a low profile.

Hollister leg bags maintain a low profile as they rest discreetly against the leg with soft, pleated sides and are held comfortably in place with wide fabric straps.



## Urine Collection Bags and Accessories

## Leg Bag

### Urinary Leg Bag

Box of 10

9814/9805

- Includes:
  - (1) Sterile, Urinary Leg Bag — Medium or Large
- Not Made with Natural Rubber Latex

HCPCS: A9270

HCPCS: A4358 (if straps are added)

Stock No	Description	Length	Width
9814	Medium — 18 oz (540 mL)	10"	4-1/4"
9805	Large — 30 oz (900 mL)	11-1/8"	5"

## Leg Bag Packs

### Urinary Leg Bag Combination Pack

Box of 10

9348/9349

- Includes:
  - (1) Sterile, Urinary Leg Bag — Medium or Large
  - (1) Non-Sterile, Extension Tubing (18" Long, Oval, 11/32" Effective Inside Diameter) with Connector
  - (1) 1 Pair, Fabric Leg Bag Straps—Medium (23" Long, 1-1/8" Wide)
- Not Made with Natural Rubber Latex

HCPCS: A4358

Stock No	Description	Length	Width
9348	Medium — 18 oz (540 mL)	10"	4-1/4"
9349	Large — 30 oz (900 mL)	11-1/8"	5"

### Vented Urinary Leg Bag Combination Pack

Box of 10

9645/9655

- Includes:
  - (1) Sterile, Urinary Leg Bag — Medium or Large
  - (1) Sterile, Extension Tubing with Vacuum Relief Valve, (18" Long, Oval, 11/32" Effective Inside Diameter) with Connector
  - (1) 1 Pair, Fabric Leg Bag Straps — Medium (23" Long, 1-1/8" Wide)
- Not Made with Natural Rubber Latex

HCPCS: A4358

Stock No	Description	Length	Width
9645	Medium — 18 oz (540 mL)	10"	4-1/4"
9655	Large — (30 oz / 900 mL)	11-1/8"	5"

## Accessories

### Fabric Leg Bag Straps

10 Pair

9342/9343

- Not Made with Natural Rubber Latex

HCPCS: A5114

Stock No	Description	Length	Width
9342	Small	15"	1-1/8"
9343	Medium	23"	1-1/8"

### Extension Tubing with Connector

Box of 10

9342/9343

- Not Made with Natural Rubber Latex

HCPCS: A5114

Stock No	Description	Length
9345	Non-sterile	18"
9346	Sterile	18"

## Bedside Drainage Collection System with Anti-Reflux Valve

### Bedside Drainage Collection System with Anti-Reflux Valve

Box of 10

9839


- Sample Port
- Sterile Fluid Pathway
- Not Made with Natural Rubber Latex

HCPCS: A4357

Stock No	Description	Length	Width
9839	68 oz (2000 mL)	12"	9-3/4"

# Urine Collection Bags and Accessories

The Hollister vented and non-vented leg bags use an anti-reflux valve to help prevent urine from back flowing into the tubing. Our urinary leg bags are not made with natural rubber latex and are designed to be discreet and secure, promoting independence and quality of life.



#### Anti-Reflux Valve

- Helps prevent urine from back flowing into the tubing

#### Buttons

- Low profile buttons help secure the bag

#### Extension Tubing

- Our oval extension tubing helps keep the fluid pathway open
- Sterile and non-sterile options are available

#### Sterile

#### Soft Vinyl Film

- Helps ensure comfort against the skin and limits unwanted odor

#### Convenience Drain

- Easy to operate with one hand, even for those with limited dexterity

#### Fabric Straps

- Wide straps help hold the bag comfortably and securely in place

#### Optional Vented System

#### Vent

- Helps prevent air build-up in the leg bag, providing a lower profile.

#### Bedside drainage collection system with anti-reflux valve also available

Prior to use, be sure to read the Instructions for Use for information regarding Intended Use, Contraindications, Warnings, Precautions, and Instructions.

Hollister and the Hollister logo are trademarks of Hollister Incorporated.

© 2021 Hollister Incorporated.

922093-0121 US-01338

